

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UptiBond™ is Interchim's proprietary GC capillary column technology that has been developed around rigorous, innovative manufacturing processes from the highest quality raw material. Every step of column manufacturing is strictly controlled, establishing a new standard in GC chromatography. UptiBond™ exhibits enhanced column to column reproducibility and the utilization of Interchim's in-house chemistry resources establish major product benefits such as ultra low levels of column bleeding.

Stationary phase selection guide

Applications	Composition	Polarity Isothermal/ Programmed	Temp. Range	Phases with Similar Polarity
Bonded Phases Amines, hydrocarbons, pesticides, PCBs, phenols, sulfur compounds,	100% dimethylpolysiloxane	Nonpolar	a : -60 to 325/350 b : -60 to 300/320 c : -60 to 260/280	UBP1, DB-1, BP-1, HP-1 CP-Sil5, Rtx-1, OV-1,007-1
Alkaloids, drugs, FAMES Halogenated compounds, Semi-volatiles, pesticides	5% phenyl - 95 % dimethyl polysiloxane	Nonpolar	a : -60 to 325/350 b : -60 to 300/320 c : -60 to 260/280	UB5P, DB-5, DB-5ms, DB-5ht, XTI-5, CP-Sil8CB, 007-2, OV-5, Rtx-5, Rtx-5, BPX-5, BP-5, SPB-5
Aroclors, alcohols, pesticides, VOAs	6% cyanopropyl-phenyl 94% dimethyl polysiloxane	Mid-polar	a : -20 to 280/300 b : -20 to 260/280	UB1301, DB-1301, DB-624, Rtx-624, Mtx-1301, CP-624, RTX-1301
Aroclors, Amines, pesticides pharmaceuticals	35% dimethyl polysiloxane	Mid-polar	a : 40 to 300/320	DB-35, RTX-35, SPB-35, AT-35
Volatile compounds		Mid-polar	a : 260/260	UB624, DB-624, OV624, ZB-624 HP-VOC, CP-Select624CB
Pesticides, herbicides, TMS Sugars, aroclors	14 % cyanopropyl-phenyl 86 % dimethyl polysiloxane	Mid-polar	a : -20 to 280/300 b : -20 to 260/280	UB1701, DB-1701, Rtx-1701, CP-Sil 19 CB, OV1701
Drugs, glycols, pesticides, steroids	50% phenyls, 50% methyl- polysiloxane	Mid-polar	a : 40 to 280/300 b : 40 to 260/280	UB17, UB17P, DB-17, DB-17ht, CP-Sil24CB, BPX200, 007-17, OV-17
Aldehydes, pesticides, Herbicides, organo-chlorine Organo-phosphorous		Polar	a : -45 to 240/260 b : -45 to 220/240	UB210, DB-210, Rtx-200
FAMES, alditol acetates Neutral sterols	50 % cyanopropyl-phenyl 50 % methyl-polysiloxane	Polar Polar	a : 40 to 220/240 b : 40 to 200/220	UB225, DB-225, SP-2330, OV-225, Rtx-225, BP-225, 007-225
Solvents, glycols, alcohols	Polyethylene	Polar	a : 20 to 250/260 b : 20 to 230/240	UBWAX, DB-WAX, Carbowax 20M Rtx-WAX, ZB-WAX, HP-20M
Acids, alcohols, aldehydes Acrylates, ketones, nitriles	Polyethylene glycol-acid	Polar	a : 60 to 240/250 b : 60 to 230/240	UBFFAP, DB-FFAP, Stabilwax DA 007-FFAP, BP21,

a : d.i. < 530 µm - b : d.i. : 530 µm - c : d.i. > 2 µm

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UB1P - [100% Dimethyl Polysiloxane]

- Low bleed MS grade column
- Suitable upgrade from : DB-1, HP-1, Ultra-1, SPB-1, CP-Sil5CB, ZB-1, OV1
- Excellent general purpose column
- Hydrocarbon, pesticides, phenols, amines, flavors and fragrances

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.25	325/350	UB1P152525
0.25	30	0.25	325/350	UB1P302525
0.25	60	0.25	325/350	UB1P602525
0.32	15	0.25	325/350	UB1P153225
0.32	30	0.25	325/350	UB1P303225
0.32	60	0.25	325/350	UB1P603225

UB5P - [5% phenyl - 95% Dimethyl Polysiloxane]

- Low bleed MS grade column
- Suitable Upgrade from : DB-5, DB-5ms, HP-5, HP-5ms, SPB-5, VF-5ms, ZB-5, OV5, Rtx-5, Rtx-5ms, Rtx-5SiIMS
- Excellent general purpose column
- Alkaloids, drugs, FAMES, halogenated compounds, aromatic compounds, pesticides

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.10	325/350	UB5P152501
0.25	15	0.25	325/350	UB5P152525
0.25	15	0.50	325/350	UB5P152550
0.25	30	0.10	325/350	UB5P302501
0.25	30	0.25	325/350	UB5P302525
0.25	30	0.50	325/350	UB5P302550
0.25	30	1.0	325/350	UB5P302510
0.25	60	0.10	325/350	UB5P602501
0.25	60	0.25	325/350	UB5P602525
0.32	15	0.10	325/350	UB5P153201
0.32	15	0.25	325/350	UB5P153225
0.32	15	0.50	325/350	UB5P153205
0.32	30	0.10	325/350	UB5P303201
0.32	30	0.25	325/350	UB5P303225
0.32	30	0.50	325/350	UB5P303205
0.32	60	0.10	325/350	UB5P603201
0.32	60	0.25	325/350	UB5P603225



Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UB17P - [50% phenyl - 50% Dimethyl Polysiloxane]

- Ideal for pesticide analysis
- Polar compounds, FAMES, free phenols, chlorinated pesticides, basic drugs, phthalate esters
- Suitable upgrade from : DB-17, DB-17ms, Rtx-50, ZB-50, CP-Sil19 CB, BPX-50, 007-17

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.25	320/350	UB17P152525
0.25	30	0.25	320/350	UB17P302525
0.25	60	0.25	320/350	UB17P602525
0.32	15	0.25	320/350	UB17P153225
0.32	30	0.25	320/350	UB17P303225
0.32	60	0.25	320/350	UB17P603225



Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UB1301 - [6% Cyanopropylphenyl - 94% Dimethylpolysiloxane]

- Ideal for pesticide analysis
- Polar compounds, FAMES, free phenols, chlorinated pesticides, basic drugs, phatolate esters
- Suitable Upgrade from : DB-1301, Rtx-1301, CPSil1301 etc..

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.25	280/300	UB1301152525
0.25	15	0.5	280/300	UB1301152550
0.25	15	1	260/280	UB1301152510
0.25	30	0.25	280/300	UB1301302525
0.25	30	0.5	280/300	UB1301302505
0.25	30	1	260/280	UB1301302510
0.25	60	0.25	280/300	UB1301602525
0.25	60	0.5	280/300	UB1301602550
0.25	60	1	260/280	UB1301602510
0.32	15	0.25	280/300	UB1301153225
0.32	15	0.5	280/300	UB1301153205
0.32	15	1	260/280	UB1301153210
0.32	30	0.25	280/300	UB1301303225
0.32	30	0.5	280/300	UB1301303205
0.32	30	1	260/280	UB1301303210
0.32	60	0.25	280/300	UB1301603225
0.32	60	0.5	280/300	UB1301603205
0.32	60	1	260/280	UB1301603210
0.53	15	1	260/280	UB1301155310
0.53	30	1	260/280	UB1301305310

UB624

- UB624 is tested with volatile halogenated compounds, making it an excellent choice for Method 624, 501, 503.1, 524.2, 601-603 and 8260.
- Suitable Upgrade from : DB-624, Rtx-624, CP-Select 624CB, SPB-624, BP-624, AT-624, ZB-624

Volatile organic compounds (VOCs), residual solvents, volatile aromatic and unsaturated hydrocarbons.

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	30	1.4	260/260	UB624302514
0.25	60	1.4	260/260	UB624602514
0.32	30	1.8	260/260	UB624303218
0.32	30	3	260/260	UB624303230
0.32	60	1.8	260/260	UB624603218
0.53	30	3	260/260	UB624305330
0.53	75	3	260/260	UB624755330



Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UB1701 - [14% Cyanopropylphenyl - 86% Methylpolysiloxane]

- Suitable upgrade from : DB-1701, Rtx-1701, SPB-1701, CP-Sil19CB, OV-1701, 007-1701, AT-1701
- Pesticide organochlorides, Sugars, TMS, Triazines, Herbicides, Drugs

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.25	280/300	UB1701152525
0.25	15	0.50	280/300	UB1701152505
0.25	15	1.0	260/280	UB1701152510
0.25	30	0.25	280/300	UB1701302525
0.25	30	0.50	280/300	UB1701302505
0.25	30	1.0	260/280	UB1701302510
0.25	60	0.25	280/300	UB1701602525
0.25	60	0.50	280/300	UB1701602505
0.25	60	1.0	260/280	UB1701602510
0.32	15	0.25	280/300	UB1701153225
0.32	15	0.50	280/300	UB1701153205
0.32	15	1.0	260/280	UB1701153210
0.32	30	0.25	280/300	UB1701303225
0.32	30	0.25	280/300	UB1701303225
0.32	30	0.50	280/300	UB1701303205
0.32	30	1.0	260/280	UB1701303210
0.32	60	0.25	280/300	UB1701603225
0.32	60	0.50	280/300	UB1701603205
0.32	60	1.0	260/280	UB1701603210
0.53	15	1.0	260/280	UB1701155310
0.53	30	1.0	260/280	UB1701305310



Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UB17P - [50% Diphenyl – 50% Dimethylpolysiloxane]

- Ideal for pesticide analysis
- Polar compounds, FAMES, free phenols, chlorinated pesticides, basic drugs, phatolate esters
- Suitable upgrade from : DB-17, DB-17ms, Rtx-50, ZB-50, CP-Sil19 CB, BPX-50, 007-17

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	15	0.25	320/340	UB17152525
0.25	30	0.25	320/340	UB17302525
0.25	60	0.25	320/340	UB17602525
0.32	30	0.25	320/340	UB17303225
0.32	60	0.25	320/340	UB17603225
0.53	15	1.0	300/320	UB17155310
0.53	30	1.0	300/320	UB17305310

UB210

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	30	0.25	240/260	UB210302525

UB225

I.D mm	Length M	Thickness μm	Temp. Range ($^{\circ}\text{C}$)	P/N
0.25	30	0.25	240/260	UB225302525



Analysis - GC Capillary Columns

UptiBond™ GC capillary columns

UBWAX - [Polyethylene Glycol]


I.D mm	Length M	Thickness µm	Temp. Range (°C)	P/N
0.25	15	0.25	250/260	UBWAX152525
0.25	30	0.25	250/260	UBWAX302525
0.25	30	0.50	250/260	UBWAX302505
0.25	60	0.25	250/260	UBWAX602525
0.25	60	0.50	250/260	UBWAX602505
0.32	15	0.25	250/260	UBWAX153225
0.32	30	0.25	250/260	UBWAX303225
0.32	30	0.50	250/260	UBWAX303205
0.32	60	0.25	250/260	UBWAX603225
0.32	60	0.50	250/260	UBWAX603205
0.53	15	1.0	230/240	UBWAX155310
0.53	15	2.0	230/240	UBWAX155320
0.53	30	1.0	230/240	UBWAX305310
0.53	30	2.0	230/240	UBWAX305320

UBAmines

I.D mm	Length M	Thickness µm	Temp. Range (°C)	P/N
0.32	15		265/300	UBAMIN1532
0.32	30		265/300	UBAMIN3032
0.32	60		265/275	UBAMIN6032

UBFFAP - [Polyethylene Glycol]

I.D mm	Length M	Thickness µm	Temp. Range (°C)	P/N
0.25	15	0.25	230/240	UBFFAP152525
0.25	30	0.25	230/240	UBFFAP302525
0.25	30	0.50	230/240	UBFFAP302505
0.25	60	0.25	230/240	UBFFAP602525
0.25	60	0.50	230/240	UBFFAP602505
0.32	30	0.25	230/240	UBFFAP303225
0.32	60	0.25	230/240	UBFFAP603225
0.53	15	1.0	220/230	UBFFAP155310
0.53	30	1.0	220/230	UBFFAP305310



Analysis - GC Capillary Columns


Uptpacked™ GC columns

Interchim can supply packed columns made with any combination of phases

Premium grade stainless tubing (only manufactured for GC only)
Including brass Swagelock fittings
Delivery 2 or 3 weeks

For order enquiries please indicate :

1. Type of tubing required i.e SS or Glass
2. Length and OD of the column
3. Loading % of liquid phase
4. Nature of liquid phase
5. Solid support and particle size (mesh)
6. Name and model of GC chromatograph



For GC glass columns please specify your column configuration according to the drawing highlighted above.

Order enquiry example

6'	x	1/8"	Stainless steel	10%	OV1	Chromosorb WMP	80/100	for HP 6890
Length		OD	Tubing type	Percentage	Liquid phase	Solid support	Mesh size	GC model

Standard SS Uptpacked

- Available from stock
- 1.8 m x 1/8"



Description	P/N	Description	P/N
10% Apiezon L on Chromosorb W-HP 80/100	793980	10% OV-1 on Chromosorb W-HP 80/100	867980
10% Carbowax 1500 on Chrom. W-HP 80/100	816520	3% OV-17 on Chromosorb W-HP 80/100	393090
10% Carbowax 20 M on Chromosorb W-HP 80/100	245691	10% OV-17 on Chromosorb W-HP 80/100	222540
10% Carbowax 20 M on Chromosorb W-AW 80/100	184630	3% OV-101 on Chromosorb W-HP 80/100	743800
10% Carbowax 20 M-TPA on Chrom. W-AW 80/100	742060	10% OV-101 on Chromosorb W-HP 80/100	802070
Chromosorb 101, 80/100	834260	3% OV-225 on Chromosorb W-HP 80/100	848440
Chromosorb 102, 80/100	249630	10% OV-225 on Chromosorb W-HP 80/100	882870
10% DEGS on Chromosorb W-AW 80/100	787230	3% SE-30 GC Grade on Chromosorb W-HP 80/100	388810
5% FFAP on Chromosorb W-AW 80/100	835850	Porapak Q, 80/100	520800
3% OV-1 on Chromosorb W-HP 80/100	765440	Porapak N, 80/100	978770


Analysis - GC Capillary Columns

UptiBond™ GC capillary columns - applications

Interchim columns are tested according Grob test process. Each compound of this test mixture measures, efficiency, film thickness, acid & basic characteristics and column activity.

UptiBond™ UB1P


Compound	Rt (min)	Compound	Rt (min)		
1	2,3-Butanediol	4.159 et 4.293	7	2-Ethylhexanoic acid	11.599
2	n-Decane	9.019	8	2,6-Dimethylaniline	12.078
3	1-Octanol	10.254	9	Methyl decanoate	15.984
4	2,6-Dimethylphenol	10.829	10	Dicyclohexylamine	18.074
5	1-Nonanal	10.920	11	Methyl undecanoate	18.074
6	n-Undecane	11.357	12	Methyl dodecanoate	20.058



UB1P303225 / N° Series : 1P-53542
Temp. : 40°C to 190°C at 6°C /min
Carrier Gas : Helium, 20 cm/sec (FID)

UptiBond™ UB5P


Compound	Rt (min)	Compound	Rt (min)		
1	2,3-Butanediol	4.742 et 4.901	5	1-Nonanal	11.812
2	n-Decane	9.332	6	2,6-Dimethylphenol	11.914
3	1-Octanol	10.996	7	2-Ethylhexanoic acid	12.203
4	n-Undecane	11.724	8	2,6-Dimethylaniline	13.314
			9	Methyl decanoate	16.805
			10	Dicyclohexylamine	18.926
			11	Methyl undecanoate	18.926
			12	Methyl dodecanoate	20.938



UB5P303225 / N° Series : 5P-61700
Temp. : 40°C to 190°C at 6°C /min
Carrier Gas : Helium, 20 cm/sec (FID)

UptiBond™ UB1301

Compound	Rt (min)	Compound	Rt (min)		
1	2,3-Butanediol	6.096 et 6.249	5	1-Nonanal	12.148
2	n-Decane	8.922	6	2,6-Dimethylphenol	13.185
3	n-Undecane	11.254	7	2,6-Dimethylaniline	14.028
4	1-Octanol	11.627	8	2-Ethylhexanoic acid	14.609
			9	Methyl decanoate	16.802
			10	Dicyclohexylamine	18.521
			11	Methyl undecanoate	18.861
			12	Methyl dodecanoate	20.816




UB1301303225 / N° Series : 1301-62915
Temp. : 40°C to 190°C at 6°C /min
Carrier Gas : Helium, 20 cm /sec (FID)

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns - applications

UptiBond™ UB624


Compound	Rt (min)	Compound	Rt (min)		
1	2,3-Butanediol	10.245 et 10.469	5	1-Nonanal	17.847
2	n-Decane	14.090	6	2,6-Dimethylphenol	18.973
3	n-Undecane	16.762	7	2-Ethylhexanoic acid	19.343
4	1-Octanol	17.103	8	2,6-Dimethylaniline	20.139
			9	Methyl decanoate	22.771
			10	Methyl undecanoate	24.956
			11	Dicyclohexylamine	25.090
			12	Methyl dodecanoate	27.023



UB624303218 / N° Series : 624-62911
 Temp. : 40°C to 230°C at 6°C /min
 Carrier Gas : Helium, 20 cm /sec (FID)

UptiBond™ UBWAX


Compound	Rt (min)	Compound	Rt (min)		
1	n-Decane	4.553	5	1-Octanol	14.999
2	n-Undecane	5.973	6	Methyl decanoate	15.702
3	1-Nonanal	11.771	7	Dicyclohexylamine	16.710
4	2,3-Butanediol	14.765 et 15.477	8	Methyl undecanoate	17.612
			9	Methyl dodecanoate	19.445
			10	2,6-Dimethylaniline	20.403
			11	2,6-Dimethylphenol	21.447
			12	2-Ethylhexanoic acid	22.090



UBWAX303225 / N° Series : WAX-62918
 Temp. : 40°C to 190°C at 6°C /min
 Carrier Gas : Helium, 20 cm /sec (FID)

UptiBond™ UB1701


Compound	Rt (min)	Compound	Rt (min)		
1	2,3-Butanediol	7.842 et 8.079	5	1-Nonanal	13.238
2	n-Decane	8.740	6	2,6-Dimethylphenol	15.200
3	n-Undecane	11.081	7	2-Ethylhexanoic acid	15.650
4	1-Octanol	12.866	8	2,6-Dimethylaniline	15.914
			9	Methyl decanoate	17.606
			10	Dicyclohexylamine	19.189
			11	Methyl undecanoate	19.655
			12	Methyl dodecanoate	21.601



UB1701303225 / N° Series : 170162914
 Temp. : 40°C to 190°C at 6°C /min
 Carrier Gas : Helium, 20 cm /sec (FID)

UptiBond™ UBFFAP

Compound	Rt (min)	Compound	Rt (min)		
1	n-Decane	4.562	5	1-Octanol	15.062
2	n-Undecane	5.974	6	Methyl decanoate	15.745
3	1-Nonanal	11.846	7	Dicyclohexylamine	17.650
4	2,3-Butanediol	14.905 et 15.617	8	Methyl undecanoate	17.650
			9	Methyl dodecanoate	19.477
			10	2,6-Dimethylaniline	20.592
			11	2,6-Dimethylphenol	21.608
			12	2-Ethylhexanoic acid	22.179




UBFFAP303225 / N° Series : FFAP-62916
 Temp. : 40°C to 190°C at 6°C /min
 Carrier Gas : Helium, 20 cm /sec (FID)

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns - applications

Fatty acid methyl ester

- 1 - 3.129 C8:0
- 2 - 3.623 C9:0
- 3 - 4.182 C10:0
- 4 - 4.789 C11:0
- 5 - 5.487 C12:0




UB5-Premium 0.25 µm, 30 m x 0.25 mm

Carrier Gas : He, 25 cm /sec, constant flow
 Oven : 130°C (hold 1 min)
 200°C at 20°C /min (hold 2 min)
 Injection : 1 µl, split 1:100, 270°C
 Detection : FID, 280°C

Alcohols

- 1 - 2.671 Methanol
- 2 - 2.884 Ethanol
- 3 - 3.025 Acetone
- 4 - 3.106 Isopropanol
- 5 - 4.040 Butanone
- 6 - 7.071 3-Methylbutan-1-ol
- 7 - 7.163 2-Methylbutan-1-ol
- 8 - 7.816 Pentan-1-ol



UB1-Premium 0.25µm, 30 m x 0.25 mm

Carrier Gas : He, 25 cm /sec, constant flow
 Oven : 30°C (hold 4 min)
 100°C at 10°C /min
 Injection : 1 µl, split 1:100, 200°C
 Detection : FID, 200°C

Solvents

- 1 - 2.555 Acetone
- 2 - 3.018 Butanone
- 3 - 3.628 Benzene
- 4 - 5.074 Pentanol
- 5 - 6.847 Ethylbenzene
- 6 - 7.039 p-Xylene
- 7 - 7.430 Styrene
- 8 - 7.480 Cyclohexanone



UB5-Premium 0.25µm, 30 m x 0.25 mm


Carrier Gas : He, 25 cm /sec, constant flow
 Oven : 30°C (hold 3 min)
 100°C at 10°C /min
 (hold 2 min)
 Injection : 1 µl,
 split 1:100, 250°C
 Detection : FID, 300°C

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns - applications

Hydrocarbons

- 1 - 2.511 unknown
- 2 - 3.680 Hexane
- 3 - 4.812 Benzene
- 4 - 5.944 Heptane
- 5 - 9.721 Octane
- 6 - 11.183 Ethylbenzene
- 7 - 11.988 Nonane
- 8 - 13.562 Decane
- 9 - 14.802 Undecane
- 10 - 15.870 Dodecane
- 11 - 16.839 Tridecane




UB5-Premium 0.25 µm, 30 m x 0.25 mm

Carrier Gas : He, 25 cm /sec, constant flow
Oven : 35°C (hold 7.5 min)
200°C at 15°C /min (hold 5 min)
Injection : 1 µl, split 1:100, 250°C
Detection : FID, 250°C

HAPs

- | | | | |
|------------------------|------------------|--------------------|------------------|
| 1- Naphthalene | 7- Fluorene | anthracene | fluoranthene |
| 2- 2-Methylnaphthalene | 8- Phenanthrene | 14- Benzo[b] | 19- Dibenz[a,h] |
| 3- 1-Methylnaphthalene | 9- Anthracene | anthracene | anthracene |
| 4- Acenaphthylene | 10- Fluoranthene | 15- Benzo[k] | 20- Benzo[g,h,i] |
| 5- Acenaphthene | 11- Pyrene | anthracene | perylene |
| 6- Biphenyl | 12- Chrysene | 17- Benzo[e]pyrene | 21- Indeno |
| | 13- Benzo[a] | 18- Benzo[e] | [1,2,3-cd]pyrene |




UB17P 0.25 mm I.D. x 30 m df = 0.25 µm

Oven : 80°C (3 min hold) ; 5°C/min 340°C (5 min hold)
Detection : MS
Sample Conc. : each 10 µg/L

Mixture

- | | | | |
|--------------------------------|---------------------------|---------------------------|-------------------------------|
| 1- Methanol | 16- Chloroform | ketone | 46- Cyclohexanone |
| 2- Ethanol | 17- i-Butanol | 32- n-Amyl alcohol | 47- o-Xylene |
| 3- Acetone | 18- Tetrahydrofuran | 33- Toluene | 48- Cellosolve acetate |
| 4- i-Propanol | 19- Methyl cellosolve | 34- i-Butyl acetate | 49- Butyl cellosolve |
| 5- Ethyl ether | 20- 1,1,1-Trichloroethane | 35- N,N-Dimethylformamide | 50- n-Amyl acetate |
| 6- Methyl acetate | 21- 1,2-Dichloroethane | 36- Methyl n-butyl ketone | 51- 1,1,2,2-tetrachloroethane |
| 7- Dichloromethane | 22- i-Propyl acetate | 37- n-Butyl acetate | 52- Methylcyclohexanol |
| 8- Carbon disulfide | 23- Benzene | 38- Tetrachloroethylene | 53- Methylcyclohexanone |
| 9- n-Propanol | 24- n-Butanol | 39- Chlorobenzene | 54- PPhenol |
| 10- trans-1,2-Dichloroethylene | 25- Carbon tetrachloride | 40- Ethylbenzene | 55- o-Dichlorobenzene |
| 11- Methyl ethyl ketone | 26- Trichloroethylene | 41- m-Xylene | 56- o-Cresol |
| 12- Hexane | 27- 1,4-Dioxane | 42- p-Xylene | 57- m-Cresol |
| 13- 2-Butanol | 28- Ethyl cellosolve | 43- i-Amyl Acetate | 58- p-Cresol |
| 14- cis-1,1-Dichloroethylene | 29- n-Propyl acetate | 44- Cyclohexanol | |
| 15- Ethyl acetate | 30- i-Amyl alcohol | 45- Styrene | |
| | 31- Methyl i-butyl | | |



UBWAX 0.25 mm I.D. x 60 m df = 0.25 µm


Oven : 40°C (5 min hold) 4°C/min 230°C
Injection : Split
Detection : FID

Analysis - GC Capillary Columns

UptiBond™ GC capillary columns - applications

Mixture

- 1- Methanol
- 2- Dimethylamine
- 3- Ethanol
- 4- Acetonitrile
- 5- Acetone
- 6- Iso-Propanol
- 7- Acetic acid
- 8- Diethylamine
- 9- Ethylacetate
- 10- Triethylamine
- 11- Pyridine
- 12- Dimethylformamide
- 13- Toluene




UB Amines 0.32 mm I.D. x 60 m

Oven : 50°C (hold 3 min) 10°C /min 220°C (5 min. Hold)
Detection : FID 250°C
Injecteur : 260°C
Split ratio : 1:50
Samples : solvents & amines

Xylene compound


- 1- p-Xylene
- 2- m-Xylene
- 3- o-Xylene
- 4- Styrene



UBBFAP 0.25 mm I.D. x 60 m df = 0.5 µm

Oven : 150°C
Injection : 250°C
Detection : FID 250°C

Fragrances



UB5-Premium 0.25 µm, 30 m x 0.25 mm


Carrier Gas : He, 0.1 cm /sec, constant flow
Oven : 40°C ; 85°C at 100°C /min ; 180°C at 3°C /min
Injection : 1 µl, split 1:100, 250°C
Detection : FID, 250°C

Analysis - GC Capillary Columns


UptiBond™ GC capillary columns - applications

Amino acids analysis


Amino acids standard



Mushroom tissue with nitrogen



Mushroom tissue with nitrogen



UB5-Premium 0.25 μ m, 30 m x 0.25 mm

Carrier Gas : Helium 1.2 ml /min

Oven : 110°C for 3 min

110°C - 260°C at 6°C/min

260°C for 4 min

Injection : 250°C

1/10 μ l injection

Detection : MSD

Sample preparation :

- Lyophilization of mushrooms
- Ethanol extraction
- Dry evaporation
- Dissolution in HCl
- Dry evaporation
- Derivatization with MTBSTFA/
DMF (5/1)
3 min 83°C