

Biochromatography

Introduction

Interchim Biochromatography encompasses a wide range of technologies, techniques and combination processes to provide the research and developer with a diverse array of tools for the purification of complex biological samples such as proteins, peptides and polypeptides.

Whilst the technologies highlighted have been developed to predominately serve these applications, the industry has also demonstrated that these techniques can often be useful for purification applications of other complex molecules.

The table below provides a general guide to product selection. It cross references the product versus the technology it can be utilized within.

Product	Reverse phase	HIC	Ion Exchange	Gel Filtration	Affinity	IMAC
Affarose™					X	X
AmiR gel™					X	
Atoll™	RP					
Cellufine®		HIC	IE	GF		
Hydrocell®	RP	HIC	IE			
UptiBeads™					X	
Uptisphere® 300Å	RP					
Uptisphere® Protein A					X	
Uptispin™					X	X

Protein background

Proteins are long chain covalently linked amino acids. Protein structures are folded rather than linear due to varying types of binding forces such as hydrogen, van der Waals, and disulfide bonds. Proteins vary in nature relative to chain length and amino acid sequence.

The primary structure of proteins and peptides refers to the linear number and sequence of amino acids present.

Protein Secondary Structure

The amino acid sequence in a protein determines the form within a particular protein. This constitutes the secondary structures of a protein. In general, proteins fall into two main classes of structure, namely globular proteins or fibrous proteins. Globular proteins are compactly folded and coiled, whereas, fibrous proteins are more filamentous and elongated.

Several kinds of secondary structure exist :

- Alpha-Helix
- 3.10-Helix
- Sheets
- Turns



Super-secondary Structure

Some proteins contain an order within secondary structures that forms a distinct sub-structural classification. These are referred as super-secondary structures.

Tertiary Structure

Tertiary structures are complete three-dimensional polypeptide units of a given protein. This incorporates the spatial relationship of different secondary structures from one another within a polypeptide chain and how secondary structures fold into the three-dimensional form of the protein.


Secondary structures of proteins often constitute distinct classifications, therefore a tertiary structure also refers to the relationship of different classifications to one another within a protein. Interactions of different classifications is governed by a range of binding forces, these include hydrogen bonding, hydrophobic interactions, electrostatic interactions and van der Waals forces.

Quaternary Structure

Many proteins contain two or more different polypeptide chains that are held in association by the same non-covalent forces that stabilize the tertiary structures of proteins. Proteins with multiple polypeptide chains are termed oligomeric proteins. The structure formed by monomer-monomer interaction in an oligomeric protein is known as a quaternary structure.

Complex Protein Structures

Proteins are also covalently conjugated with carbohydrates. Proteins covalently associated with carbohydrates are termed glycoproteins. Structural complexes involving protein associated with lipid via non-covalent interactions are termed lipoproteins. Such modifications occur following the synthesis (translation) of proteins and are termed post-translational modifications.



Protein & Peptide Standards

Interchim protein standards and peptides are detailed within the reverse phase section of this chapter (p.C20). Whilst this is considered by Interchim as the primary application for these standards they are also suited to other applications such as Gel Filtration.

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Method development kits

Oligopeptides and peptides

Three same dimension columns packed with the following chemistry :

- Uptisphere® WOD, 300 Å C18 mono functional
- Uptisphere® PXP, 300 Å C18 poly functional type I
- Uptisphere® WTF, 300 Å C18 tri functional

Item	Dimension	P/N
Oligopeptide & peptide kit	250 x 4.6 mm	CI4220
Oligopeptide & peptide kit	250 x 2.1 mm	CI4230

Polypeptides and proteins

Three same dimension columns packed with the following chemistry :

- Uptisphere® WC4, 300Å C4 mono functional
- Uptisphere® WD4, 300Å C4 poly functional type I
- Uptisphere® WT4, 300Å C4 tri functional

Item	Dimension	P/N
Polypeptide & protein kit	250 x 4.6 mm	CI4240
Polypeptide & protein kit	250 x 2.1 mm	CI4250

Client specification

Three same dimension columns packed with three chemistries from our range of seven Uptisphere® 300Å.

Item	Dimension	P/N
Method development Kit 5 µm	250 x 4.6 mm	DEV052546
Method development Kit 5 µm	250 x 3.0 mm	DEV052530
Method development Kit 5 µm	250 x 2.0 mm	DEV052520

Validation method kit

Three columns packed with the same chemistry with three different batches

Item	Dimension	P/N
Validation method Kit 5 µm	250 x 4.6 mm	PP052546
Validation method Kit 5 µm	250 x 2.1 mm	PP052521

Kit specifications

Column capacities :

Analytical, micro, nano, and preparative.

Particle Size :

3, 5, 10, & 15 µm.

Method Development Kits :

Set of three columns with three bonding technologies in C4 or C18.

Prep Method Development Kits :


contain one analytical column and one preparative column packed with the same batch of silica.

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Selection Guide & Column Characteristics



This chart provides a basic guide to column selection based on the molecular weight and hydrophobicity of the peptides and proteins.



Separation of lysozyme

Uptisphere 5 µm WTF, 250 x 4.6 mm vs competitor Column C18

Uptisphere 5 µm WTF, 250 x 4.6 mm
Flow rate : 1 µl /min
Temp : 50°C
Volume : 30.00 µl
Solvent : H₂O/TFA 0.05%/CH₃CN 80 %



Competitor 5 µm C18, 250 x 4.6 mm
Flow rate : 1 µl /min
Volume : 30.00 µl
Solvent : H₂O/TFA 0.05%/CH₃CN 80 %

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns - Hardware

Modulo-Cart QS and QK

Modulo-Cart QS is Interchim's adopted hardware for the majority of analytical columns. Internal studies have demonstrated a 15% gain in efficiency over competitor hardware. The QK system includes a guard holder and cartridge.

Availability :

Lengths : 33 mm -to- 300 mm

i.d. : 1.0 mm -to- 4.6 mm

- Hand tight connection
- Direct connection guard cartridge -to- analytical column
- Frits : easily replaceable
- No drying during storage
- The PEEK - Carbon constitution guarantees a long life time

Column i.d. (mm)	Connection tube i.d. (inches)	Flow (ml /min)
2.1 - 3.0	0.005 - 0.007	0.20 - 0.5
3.0, 3.9, 4.0, 4.6	0.007 - 0.010	0.50 - 1.0

Modulo-Cart HS

Reduce your analysis time without sacrificing resolution

Modulo-Cart HS hardware allow you to reduce analysis time without sacrificing resolution. It is specifically dedicated to ultra efficient and fast analysis in the smallest geometry.

Availability :

Lengths : 20, 33, 50 & 75 mm

i.d. : 2.0 & 4.6 mm i.d. and length.

- Fast analysis
- Perfect for fast gradient and LC-MS
- Zero dead volume



Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns - Hardware

Modulo-Cart the MC & MCK

1.0 mm internal diameter for high sensitivity

MC columns display an approximate 20 fold increase in sensitivity when compared to a 4.6 mm column. They are particularly suited to trace analysis & research of impurities or metabolites in very low concentrations and are a valuable tool in biochromatography. In addition its high sensitivity leads to low solvent consumption, perfect for LC/MS and LC/RMN applications.

Availability : 20, 50, 100 & 150 mm lengths.

Modulo-Cart MCK columns have a 1 mm pre-installed guard cartridge packed with the stationary phase of your choice.



Column i.d. (mm)	Connection tube i.d. (inches)	Flow (ml /min)
1.0 - 2.1	0.005	0.05 - 0.2

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns - Hardware

Nano LC - Stainless Steel hardware w/500 µm i.d.

High value or very low quantity biomolecule analysis

For analysis requiring extreme sensitivity. Displays a 50 fold increase in sensitivity compared to a 4.6 mm column.

Frit design dramatically reduces hold up volume and facilitates sample to bed interaction.
Nano LC is ideal for biochemistry where typically samples are valuable and in micro quantities.

Availability :

Lengths : 50, 100 & 150 mm

i.d. 500 µm

i.d.	Flow	Optimal flow	Capacity
500 µm	4 - 12 µl	9 µl	0.1 µg

Nano LC - PEEK hardware w/300 µm i.d.

Extreme sensitivity applications

For analysis requiring extreme sensitivity. Displays a 100 fold increase in sensitivity compared to a 4.6 mm column.

Frit design dramatically reduces hold up volume and facilitates sample to bed interaction.
Nano LC is ideal for biochemistry where typically samples are valuable and in micro quantities.

Availability :

Lengths : 50, 100 & 150 mm

i.d. 300 µm

i.d.	Flow	Optimal flow	Capacity
300 µm	2 - 4 µl	3 µl	0.03 µg



Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart QS	Modulo-Cart QK
Octadecyl 300 Å - 100 m ² /g %C : 10 mono functional + "one step" end-capping	WOD	3 µm	100 x 2.0 mm	UP3WOD#10QS	UP3WOD#10QK
	WOD	3 µm	125 x 2.0 mm	UP3WOD#12QS	UP3WOD#12QK
	WOD	3 µm	150 x 2.0 mm	UP3WOD#15QS	UP3WOD#15QK
	WOD	5 µm	100 x 2.0 mm	UP5WOD#10QS	UP5WOD#10QK
	WOD	5 µm	125 x 2.0 mm	UP5WOD#12QS	UP5WOD#12QK
	WOD	5 µm	150 x 2.0 mm	UP5WOD#15QS	UP5WOD#15QK
	WOD	5 µm	250 x 2.0 mm	UP5WOD#25QS	UP5WOD#25QK
	WOD	3 µm	100 x 3.0 mm	UP3WOD\$10QS	UP3WOD\$10QK
	WOD	3 µm	125 x 3.0 mm	UP3WOD\$12QS	UP3WOD\$12QK
	WOD	3 µm	150 x 3.0 mm	UP3WOD\$15QS	UP3WOD\$15QK
	WOD	5 µm	100 x 3.0 mm	UP5WOD\$10QS	UP5WOD\$10QK
	WOD	5 µm	125 x 3.0 mm	UP5WOD\$12QS	UP5WOD\$12QK
	WOD	5 µm	150 x 3.0 mm	UP5WOD\$15QS	UP5WOD\$15QK
	WOD	5 µm	250 x 3.0 mm	UP5WOD\$25QS	UP5WOD\$25QK
	WOD	3 µm	100 x 4.6 mm	UP3WOD-10QS	UP3WOD-10QK
	WOD	3 µm	125 x 4.6 mm	UP3WOD-12QS	UP3WOD-12QK
	WOD	3 µm	150 x 4.6 mm	UP3WOD-15QS	UP3WOD-15QK
	WOD	5 µm	100 x 4.6 mm	UP5WOD-10QS	UP5WOD-10QK
	WOD	5 µm	125 x 4.6 mm	UP5WOD-12QS	UP5WOD-12QK
	WOD	5 µm	150 x 4.6 mm	UP5WOD-15QS	UP5WOD-15QK
WOD	5 µm	250 x 4.6 mm	UP5WOD-25QS	UP5WOD-25QK	
WOD	10 µm	250 x 4.6 mm	UP10WOD-25QS	UP10WOD-25QK	
Octadecyl 300 Å - 100 m ² /g %C : 8 poly functional type I	PXP	5 µm	100 x 2.0 mm	UP5PXP#10QS	UP5PXP#10QK
	PXP	5 µm	125 x 2.0 mm	UP5PXP#12QS	UP5PXP#12QK
	PXP	5 µm	150 x 2.0 mm	UP5PXP#15QS	UP5PXP#15QK
	PXP	5 µm	250 x 2.0 mm	UP5PXP#25QS	UP5PXP#25QK
	PXP	5 µm	100 x 3.0 mm	UP5PXP\$10QS	UP5PXP\$10QK
	PXP	5 µm	125 x 3.0 mm	UP5PXP\$12QS	UP5PXP\$12QK
	PXP	5 µm	150 x 3.0 mm	UP5PXP\$15QS	UP5PXP\$15QK
	PXP	5 µm	250 x 3.0 mm	UP5PXP\$25QS	UP5PXP\$25QK
	PXP	5 µm	100 x 4.6 mm	UP5PXP-10QS	UP5PXP-10QK
	PXP	5 µm	125 x 4.6 mm	UP5PXP-12QS	UP5PXP-12QK
	PXP	5 µm	150 x 4.6 mm	UP5PXP-15QS	UP5PXP-15QK
	PXP	5 µm	250 x 4.6 mm	UP5PXP-25QS	UP5PXP-25QK
	PXP	10 µm	250 x 4.6 mm	UP10PXP-25QS	UP10PXP-25QK

Modulo-cart information
p.C7 to C.9

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300 Å	Phase	Particle size	Dimension	Modulo-Cart QS	Modulo-Cart QK	
Octadecyl 300 Å - 100 m ² /g %C : 14 tri functional+ "one step" end-capping	WTF	5 µm	100 x 2.0 mm	UP5WTF#10QS	UP5WTF#10QK	
	WTF	5 µm	125 x 2.0 mm	UP5WTF#12QS	UP5WTF#12QK	
	WTF	5 µm	150 x 2.0 mm	UP5WTF#15QS	UP5WTF#15QK	
	WTF	5 µm	250 x 2.0 mm	UP5WTF#25QS	UP5WTF#25QK	
	WTF	5 µm	100 x 3.0 mm	UP5WTF\$10QS	UP5WTF\$10QK	
	WTF	5 µm	125 x 3.0 mm	UP5WTF\$12QS	UP5WTF\$12QK	
	WTF	5 µm	150 x 3.0 mm	UP5WTF\$15QS	UP5WTF\$15QK	
	WTF	5 µm	250 x 3.0 mm	UP5WTF\$25QS	UP5WTF\$25QK	
	WTF	5 µm	100 x 4.6 mm	UP5WTF-10QS	UP5WTF-10QK	
	WTF	5 µm	125 x 4.6 mm	UP5WTF-12QS	UP5WTF-12QK	
	WTF	5 µm	150 x 4.6 mm	UP5WTF-15QS	UP5WTF-15QK	
	WTF	5 µm	250 x 4.6 mm	UP5WTF-25QS	UP5WTF-25QK	
	Octyl 300 Å - 100 m ² /g %C : 8 mono functional + "one step" end-capping	WC8	5 µm	100 x 2.0 mm	UP5WC8#10QS	UP5WC8#10QK
		WC8	5 µm	125 x 2.0 mm	UP5WC8#12QS	UP5WC8#12QK
WC8		5 µm	150 x 2.0 mm	UP5WC8#15QS	UP5WC8#15QK	
WC8		5 µm	250 x 2.0 mm	UP5WC8#25QS	UP5WC8#25QK	
WC8		5 µm	100 x 3.0 mm	UP5WC8\$10QS	UP5WC8\$10QK	
WC8		5 µm	125 x 3.0 mm	UP5WC8\$12QS	UP5WC8\$12QK	
WC8		5 µm	150 x 3.0 mm	UP5WC8\$15QS	UP5WC8\$15QK	
WC8		5 µm	250 x 3.0 mm	UP5WC8\$25QS	UP5WC8\$25QK	
WC8		5 µm	100 x 4.6 mm	UP5WC8-10QS	UP5WC8-10QK	
WC8		5 µm	125 x 4.6 mm	UP5WC8-12QS	UP5WC8-12QK	
WC8		5 µm	150 x 4.6 mm	UP5WC8-15QS	UP5WC8-15QK	
WC8		5 µm	250 x 4.6 mm	UP5WC8-25QS	UP5WC8-25QK	
Butyl 300 Å - 100 m ² /g %C : 4 mono functional + "one step" end-capping		WC4	3 µm	100 x 2.0 mm	UP3WC4#10QS	UP3WC4#10QK
		WC4	3 µm	125 x 2.0 mm	UP3WC4#12QS	UP3WC4#12QK
	WC4	3 µm	150 x 2.0 mm	UP3WC4#15QS	UP3WC4#15QK	
	WC4	5 µm	100 x 2.0 mm	UP5WC4#10QS	UP5WC4#10QK	
	WC4	5 µm	125 x 2.0 mm	UP5WC4#12QS	UP5WC4#12QK	
	WC4	5 µm	150 x 2.0 mm	UP5WC4#15QS	UP5WC4#15QK	
	WC4	5 µm	250 x 2.0 mm	UP5WC4#25QS	UP5WC4#25QK	
	WC4	3 µm	100 x 3.0 mm	UP3WC4\$10QS	UP3WC4\$10QK	
	WC4	3 µm	125 x 3.0 mm	UP3WC4\$12QS	UP3WC4\$12QK	
	WC4	3 µm	150 x 3.0 mm	UP3WC4\$15QS	UP3WC4\$15QK	
	WC4	5 µm	100 x 3.0 mm	UP5WC4\$10QS	UP5WC4\$10QK	
	WC4	5 µm	125 x 3.0 mm	UP5WC4\$12QS	UP5WC4\$12QK	
	WC4	5 µm	150 x 3.0 mm	UP5WC4\$15QS	UP5WC4\$15QK	
	WC4	5 µm	250 x 3.0 mm	UP5WC4\$25QS	UP5WC4\$25QK	
	WC4	3 µm	100 x 4.6 mm	UP3WC4-10QS	UP3WC4-10QK	
	WC4	3 µm	125 x 4.6 mm	UP3WC4-12QS	UP3WC4-12QK	
	WC4	3 µm	150 x 4.6 mm	UP3WC4-15QS	UP3WC4-15QK	
	WC4	5 µm	100 x 4.6 mm	UP5WC4-10QS	UP5WC4-10QK	
	WC4	5 µm	125 x 4.6 mm	UP5WC4-12QS	UP5WC4-12QK	
	WC4	5 µm	150 x 4.6 mm	UP5WC4-15QS	UP5WC4-15QK	
	WC4	5 µm	250 x 4.6 mm	UP5WC4-25QS	UP5WC4-25QK	
	WC4	10 µm	250 x 4.6 mm	UP10WC4-25QS	UP10WC4-25QK	

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart QS	Modulo-Cart QK
Butyl 300 Å - 100 m ² /g %C : 4 poly functional type I	WD4	3 µm	100 x 2.0 mm	UP3WD4#10QS	UP3WD4#10QK
	WD4	3 µm	125 x 2.0 mm	UP3WD4#12QS	UP3WD4#12QK
	WD4	3 µm	150 x 2.0 mm	UP3WD4#15QS	UP3WD4#15QK
	WD4	5 µm	100 x 2.0 mm	UP5WD4#10QS	UP5WD4#10QK
	WD4	5 µm	125 x 2.0 mm	UP5WD4#12QS	UP5WD4#12QK
	WD4	5 µm	150 x 2.0 mm	UP5WD4#15QS	UP5WD4#15QK
	WD4	5 µm	250 x 2.0 mm	UP5WD4#25QS	UP5WD4#25QK
	WD4	3 µm	100 x 3.0 mm	UP3WD4\$10QS	UP3WD4\$10QK
	WD4	3 µm	125 x 3.0 mm	UP3WD4\$12QS	UP3WD4\$12QK
	WD4	3 µm	150 x 3.0 mm	UP3WD4\$15QS	UP3WD4\$15QK
	WD4	5 µm	100 x 3.0 mm	UP5WD4\$10QS	UP5WD4\$10QK
	WD4	5 µm	125 x 3.0 mm	UP5WD4\$12QS	UP5WD4\$12QK
	WD4	5 µm	150 x 3.0 mm	UP5WD4\$15QS	UP5WD4\$15QK
	WD4	5 µm	250 x 3.0 mm	UP5WD4\$25QS	UP5WD4\$25QK
	WD4	3 µm	100 x 4.6 mm	UP3WD4-10QS	UP3WD4-10QK
	WD4	3 µm	125 x 4.6 mm	UP3WD4-12QS	UP3WD4-12QK
	WD4	3 µm	150 x 4.6 mm	UP3WD4-15QS	UP3WD4-15QK
	WD4	5 µm	100 x 4.6 mm	UP5WD4-10QS	UP5WD4-10QK
	WD4	5 µm	125 x 4.6 mm	UP5WD4-12QS	UP5WD4-12QK
	WD4	5 µm	150 x 4.6 mm	UP5WD4-15QS	UP5WD4-15QK
	WD4	5 µm	250 x 4.6 mm	UP5WD4-25QS	UP5WD4-25QK
	WD4	10 µm	250 x 4.6 mm	UP10WD4-25QS	UP10WD4-25QK
	Butyl 300 Å - 100 m ² /g %C : 3 tri functional + end-capping "one step"	WT4	5 µm	100 x 2.0 mm	UP5WT4#10QS
WT4		5 µm	125 x 2.0 mm	UP5WT4#12QS	UP5WT4#12QK
WT4		5 µm	150 x 2.0 mm	UP5WT4#15QS	UP5WT4#15QK
WT4		5 µm	250 x 2.0 mm	UP5WT4#25QS	UP5WT4#25QK
WT4		5 µm	100 x 3.0 mm	UP5WT4\$10QS	UP5WT4\$10QK
WT4		5 µm	125 x 3.0 mm	UP5WT4\$12QS	UP5WT4\$12QK
WT4		5 µm	150 x 3.0 mm	UP5WT4\$15QS	UP5WT4\$15QK
WT4		5 µm	250 x 3.0 mm	UP5WT4\$25QS	UP5WT4\$25QK
WT4		5 µm	100 x 4.6 mm	UP5WT4-10QS	UP5WT4-10QK
WT4		5 µm	125 x 4.6 mm	UP5WT4-12QS	UP5WT4-12QK
WT4		5 µm	150 x 4.6 mm	UP5WT4-15QS	UP5WT4-15QK
WT4		5 µm	250 x 4.6 mm	UP5WT4-25QS	UP5WT4-25QK
WT4		10 µm	250 x 4.6 mm	UP10WT4-25QS	UP10WT4-25QK
Octadecyl 500 Å - 60 m ² /g %C : n.a. poly functional type I	XTF	5 µm	100 x 2.0 mm	UP5XTF#10QS	UP5XTF#10QK
	XTF	5 µm	150 x 2.0 mm	UP5XTF#15QS	UP5XTF#15QK
	XTF	5 µm	250 x 2.0 mm	UP5XTF#25QS	UP5XTF#25QK
	XTF	5 µm	100 x 3.0 mm	UP5XTF\$10QS	UP5XTF\$10QK
	XTF	5 µm	150 x 3.0 mm	UP5XTF\$15QS	UP5XTF\$15QK
	XTF	5 µm	250 x 3.0 mm	UP5XTF\$25QS	UP5XTF\$25QK
	XTF	5 µm	100 x 4.6 mm	UP5XTF-10QS	UP5XTF-10QK
	XTF	5 µm	150 x 4.6 mm	UP5XTF-15QS	UP5XTF-15QK
	XTF	5 µm	250 x 4.6 mm	UP5XTF-25QS	UP5XTF-25QK

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart HS	
Octadecyl 300 Å - 100 m ² /g %C : 10 mono functional + "one step" end-capping	WOD	3 µm	20 x 2.0 mm	UP3WOD#2HS	
	WOD	3 µm	33 x 2.0 mm	UP3WOD#3HS	
	WOD	3 µm	50 x 2.0 mm	UP3WOD#5HS	
	WOD	5 µm	33 x 2.0 mm	UP5WOD#3HS	
	WOD	5 µm	50 x 2.0 mm	UP5WOD#5HS	
	WOD	3 µm	50 x 3.0 mm	UP3WOD\$5HS	
	WOD	5 µm	50 x 3.0 mm	UP5WOD\$5HS	
	WOD	3 µm	20 x 4.6 mm	UP3WOD-2HS	
	WOD	3 µm	33 x 4.6 mm	UP3WOD-3HS	
	WOD	3 µm	50 x 4.6 mm	UP3WOD-5HS	
	WOD	3 µm	75 x 4.6 mm	UP3WOD-7HS	
	WOD	5 µm	33 x 4.6 mm	UP5WOD-3HS	
	WOD	5 µm	50 x 4.6 mm	UP5WOD-5HS	
	WOD	5 µm	75 x 4.6 mm	UP5WOD-7HS	
	Octadecyl 300 Å - 100 m ² /g %C : 14 tri functional + "one step" end-capping	WTF	5 µm	33 x 2.0 mm	UP5WTF#3HS
		WTF	5 µm	50 x 2.0 mm	UP5WTF#5HS
WTF		5 µm	50 x 3.0 mm	UP5WTF\$5HS	
WTF		5 µm	33 x 4.6 mm	UP5WTF-3HS	
WTF		5 µm	50 x 4.6 mm	UP5WTF-5HS	
WTF		5 µm	75 x 4.6 mm	UP5WTF-7HS	

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart HS	
Butyl 300 Å - 100 m ² /g %C : 4 mono functional + "one step" end-capping	WC4	3 µm	20 x 2.0 mm	UP3WC4#2HS	
	WC4	3 µm	33 x 2.0 mm	UP3WC4#3HS	
	WC4	3 µm	50 x 2.0 mm	UP3WC4#5HS	
	WC4	5 µm	33 x 2.0 mm	UP5WC4#3HS	
	WC4	5 µm	50 x 2.0 mm	UP5WC4#5HS	
	WC4	3 µm	50 x 3.0 mm	UP3WC4\$5HS	
	WC4	5 µm	50 x 3.0 mm	UP5WC4\$5HS	
	WC4	3 µm	20 x 4.6 mm	UP3WC4-2HS	
	WC4	3 µm	33 x 4.6 mm	UP3WC4-3HS	
	WC4	3 µm	50 x 4.6 mm	UP3WC4-5HS	
	WC4	3 µm	75 x 4.6 mm	UP3WC4-7HS	
	WC4	5 µm	33 x 4.6 mm	UP5WC4-3HS	
	WC4	5 µm	50 x 4.6 mm	UP5WC4-5HS	
	WC4	5 µm	75 x 4.6 mm	UP5WC4-7HS	
	Butyl 300 Å - 100 m ² /g %C : 4 poly functional type I	WD4	3 µm	20 x 2.0 mm	UP3WD4#2HS
		WD4	3 µm	33 x 2.0 mm	UP3WD4#3HS
WD4		3 µm	50 x 2.0 mm	UP3WD4#5HS	
WD4		5 µm	33 x 2.0 mm	UP5WD4#3HS	
WD4		5 µm	50 x 2.0 mm	UP5WD4#5HS	
WD4		3 µm	50 x 3.0 mm	UP3WD4\$5HS	
WD4		5 µm	50 x 3.0 mm	UP5WD4\$5HS	
WD4		3 µm	20 x 4.6 mm	UP3WD4-2HS	
WD4		3 µm	33 x 4.6 mm	UP3WD4-3HS	
WD4		3 µm	50 x 4.6 mm	UP3WD4-5HS	
WD4		3 µm	75 x 4.6 mm	UP3WD4-7HS	
WD4		5 µm	33 x 4.6 mm	UP5WD4-3HS	
WD4		5 µm	50 x 4.6 mm	UP5WD4-5HS	
WD4		5 µm	75 x 4.6 mm	UP5WD4-7HS	
Butyl 300 Å - 100 m ² /g %C : 3 tri functional + "one step" end-capping		WT4	5 µm	33 x 2.0 mm	UP5WT4#3HS
		WT4	5 µm	50 x 2.0 mm	UP5WT4#5HS
	WT4	5 µm	50 x 3.0 mm	UP5WT4\$5HS	
	WT4	5 µm	33 x 4.6 mm	UP5WT4-3HS	
	WT4	5 µm	50 x 4.6 mm	UP5WT4-5HS	
	WT4	5 µm	75 x 4.6 mm	UP5WT4-7HS	

Modulo-cart information
p.C7 to C.9

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart MC	Modulo-Cart MCK
Octadecyl 300 Å - 100 m ² /g %C : 10 mono functional + "one step" end-capping	WOD	3 µm	20 x 1.0 mm	UP3WODD2MC	
	WOD	3 µm	50 x 1.0 mm	UP3WODD5MC	UP3WODD5MCK
	WOD	3 µm	150 x 1.0 mm	UP3WODD15MC	UP3WODD15MCK
	WOD	5 µm	50 x 1.0 mm	UP5WODD5MC	UP5WODD5MCK
	WOD	5 µm	150 x 1.0 mm	UP5WODD15MC	UP5WODD15MCK
Octadecyl 300 Å - 100 m ² /g %C : 14 tri functional + "one step" end-capping	WTF	5 µm	50 x 1.0 mm	UP5WTFD5MC	UP5WTFD5MCK
	WTF	5 µm	150 x 1.0 mm	UP5WTFD15MC	UP5WTFD15MCK
Butyl 300 Å - 100 m ² /g %C : 4 mono functional + "one step" end-capping	WC4	3 µm	20 x 1.0 mm	UP3WC4D2MC	
	WC4	3 µm	50 x 1.0 mm	UP3WC4D5MC	UP3WC4D5MCK
	WC4	3 µm	150 x 1.0 mm	UP3WC4D15MC	UP3WC4D15MCK
	WC4	5 µm	50 x 1.0 mm	UP5WC4D5MC	UP5WC4D5MCK
	WC4	5 µm	150 x 1.0 mm	UP5WC4D15MC	UP5WC4D15MCK
	WD4	3 µm	20 x 1.0 mm	UP3WD4D2MC	
	WD4	3 µm	50 x 1.0 mm	UP3WD4D5MC	UP3WD4D5MCK
	WD4	3 µm	150 x 1.0 mm	UP3WD4D15MC	UP3WD4D15MCK
	WD4	5 µm	50 x 1.0 mm	UP5WD4D5MC	UP5WD4D5MCK
	WD4	5 µm	150 x 1.0 mm	UP5WD4D15MC	UP5WD4D15MCK
	Butyl 300 Å - 100 m ² /g %C: 3 tri functional + "one step" end-capping	WT4	5 µm	50 x 1.0 mm	UP5WT4D5MC
WT4		5 µm	150 x 1.0 mm	UP5WT4D15MC	UP5WT4D15MCK

Biochromatography - Reverse Phase

Uptisphere® 300 Å Columns

Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart SS	Uptisphere® 300Å	Phase	Particle size	Dimension	Modulo-Cart PKG
Octadecyl 300 Å - 100 m ² /g %C : 10 mono functional + "one step" end-capping "	WOD	3 µm	50 x 0.5 mm	UP3WODB5	Octadecyl 300 Å - 100 m ² /g %C : 10 mono functional + "one step" end-capping	WOD	3 µm	50 x 0.3 mm	UP3WODA5
	WOD	3 µm	100 x 0.5 mm	UP3WODB10		WOD	3 µm	100 x 0.3 mm	UP3WODA10
	WOD	3 µm	150 x 0.5 mm	UP3WODB15		WOD	3 µm	150 x 0.3 mm	UP3WODA15
	WOD	5 µm	50 x 0.5 mm	UP5WODB5		WOD	5 µm	50 x 0.3 mm	UP5WODA5
	WOD	5 µm	100 x 0.5 mm	UP5WODB10		WOD	5 µm	100 x 0.3 mm	UP5WODA10
	WOD	5 µm	150 x 0.5 mm	UP5WODB15		WOD	5 µm	150 x 0.3 mm	UP5WODA15
300 Å - 100 m ² /g %C : 14 tri functional + "one step" end-capping	WTF	5 µm	50 x 0.5 mm	UP5WTFB5	300 Å - 100 m ² /g %C : 14 tri functional + "one step" end-capping	WTF	5 µm	50 x 0.3 mm	UP5WTFA5
	WTF	5 µm	100 x 0.5 mm	UP5WTFB10		WTF	5 µm	100 x 0.3 mm	UP5WTFA10
	WTF	5 µm	150 x 0.5 mm	UP5WTFB15		WTF	5 µm	150 x 0.3 mm	UP5WTFA15
Butyl 300 Å - 100 m ² /g %C : 4 mono functional + "one step" end-capping	WC4	3 µm	50 x 0.5 mm	UP3WC4B5	Butyl 300 Å - 100 m ² /g %C : 4 mono functional + "one step" end-capping	WC4	3 µm	50 x 0.3 mm	UP3WC4A5
	WC4	3 µm	100 x 0.5 mm	UP3WC4B10		WC4	3 µm	100 x 0.3 mm	UP3WC4A10
	WC4	3 µm	150 x 0.5 mm	UP3WC4B15		WC4	3 µm	150 x 0.3 mm	UP3WC4A15
	WC4	5 µm	50 x 0.5 mm	UP5WC4B5		WC4	5 µm	50 x 0.3 mm	UP5WC4A5
	WC4	5 µm	100 x 0.5 mm	UP5WC4B10		WC4	5 µm	100 x 0.3 mm	UP5WC4A10
	WC4	5 µm	150 x 0.5 mm	UP5WC4B15		WC4	5 µm	150 x 0.3 mm	UP5WC4A15
300 Å - 100 m ² /g %C : 4 poly functional type I	WD4	3 µm	50 x 0.5 mm	UP3WD4B5	300 Å - 100 m ² /g %C : 4 poly functional type I	WD4	3 µm	50 x 0.3 mm	UP3WD4A5
	WD4	3 µm	100 x 0.5 mm	UP3WD4B10		WD4	3 µm	100 x 0.3 mm	UP3WD4A10
	WD4	3 µm	150 x 0.5 mm	UP3WD4B15		WD4	3 µm	150 x 0.3 mm	UP3WD4A15
	WD4	5 µm	50 x 0.5 mm	UP5WD4B5		WD4	5 µm	50 x 0.3 mm	UP5WD4A5
	WD4	5 µm	100 x 0.5 mm	UP5WD4B10		WD4	5 µm	100 x 0.3 mm	UP5WD4A10
	WD4	5 µm	150 x 0.5 mm	UP5WD4B15		WD4	5 µm	150 x 0.3 mm	UP5WD4A15
300 Å - 100 m ² /g %C : 3 tri functional + "one step" end-capping	WT4	5 µm	50 x 0.5 mm	UP5WT4B5	300 Å - 100 m ² /g %C : 3 tri functional + "one step" end-capping	WT4	5 µm	50 x 0.3 mm	UP5WT4A5
	WT4	5 µm	100 x 0.5 mm	UP5WT4B10		WT4	5 µm	100 x 0.3 mm	UP5WT4A10
	WT4	5 µm	150 x 0.5 mm	UP5WT4B15		WT4	5 µm	150 x 0.3 mm	UP5WT4A15

Modulo-cart information
p.C7 to C.9

Biochromatography - Reverse Phase

Atoll™ Polymeric columns

Atoll™ is a family of highly cross-linked, spherical, macro-porous co-polymers which exhibit extreme purity. This support has been designed for applications such as solid phase sample preparation and HPLC analysis of biomolecules. It does not swell and can be used with everyday solvents.

- Ultra pure
- High mechanical resistance
- pH range : 0 - 14

Unlike traditional silica based supports, Atoll™ does not suffer from irreversible adsorption. Sample bio-activity is therefore maintained and sample recovery optimized.

The Atoll Biochromatography family consists of :


Atoll™ NPR : non-porous resin with a 2 µm particle size dedicated to ultra-fast analysis of biomolecules in LC/MS, MS/MS...

Atoll™ MP : porous resin recommended for the purification of pharmaceutical molecules, small basic molecules and peptides

Atoll™ MP3 : macro-porous resin adapted for the separation of hydrophilic peptides, proteins, recombinant proteins and bio-purification.

Atoll™ is available in a variety of housing formats relative to application. i.e. columns, cartridges and 96 well plates for sample preparation

	Atoll NPR	Atoll MP	Atoll MP3
Pore size	0	n.a.	300 Å
Surface area	2.5 m ² /g	500 m ² /g	500 m ² /g
2 µm	X		
3 µm		on request	
5 µm		X	X
10 µm		on request	
15 µm		X	X
30 µm		on request	



Biochromatography - Reverse Phase

Atoll™ Polymeric columns

i.d.	Particle size	Dimension	Atoll NPR	Atoll MP	Atoll MP3
2.0 mm	2 µm	33 x 2.0 mm	ATL2NPR#3QS ATL2NPR#5QS	ATL5MP#3QS ATL5MP#5QS ATL5MP#7QS ATL5MP#10QS ATL5MP#15QS ATL5MP#25QS	ATL5MP3#3QS ATL5MP3#5QS ATL5MP3#7QS ATL5MP3#10QS ATL5MP3#15QS ATL5MP3#25QS
		50 x 2.0 mm			
	5 µm	33 x 2.0 mm			
		50 x 2.0 mm			
		75 x 2.0 mm			
		100 x 2.0 mm			
150 x 2.0 mm					
250 x 2.0 mm					
4.6 mm	2 µm	33 x 4.6 mm	ATL2NPR-3QS ATL2NPR-5QS	ATL5MP-3QS ATL5MP-5QS ATL5MP-7QS ATL5MP-10QS ATL5MP-15QS ATL5MP-25QS	ATL5MP3-3QS ATL5MP3-5QS ATL5MP3-7QS ATL5MP3-10QS ATL5MP3-15QS ATL5MP3-25QS
		50 x 4.6 mm			
	5 µm	33 x 4.6 mm			
		50 x 4.6 mm			
		75 x 4.6 mm			
		100 x 4.6 mm			
150 x 4.6 mm					
250 x 4.6 mm					
2.0 mm	2 µm	33 x 2.0 mm	ATL2NPR#3QK ATL2NPR#5QK	ATL5MP#3QK ATL5MP#5QK ATL5MP#7QK ATL5MP#10QK ATL5MP#15QK ATL5MP#25QK	ATL5MP3#3QK ATL5MP3#5QK ATL5MP3#7QK ATL5MP3#10QK ATL5MP3#15QK ATL5MP3#25QK
		50 x 2.0 mm			
	5 µm	33 x 2.0 mm			
		50 x 2.0 mm			
		75 x 2.0 mm			
		100 x 2.0 mm			
150 x 2.0 mm					
250 x 2.0 mm					
4.6 mm	2 µm	33 x 4.6 mm	ATL2NPR-3QK ATL2NPR-5QK	ATL5MP-3QK ATL5MP-5QK ATL5MP-7QK ATL5MP-10QK ATL5MP-15QK ATL5MP-25QK	ATL5MP3-3QK ATL5MP3-5QK ATL5MP3-7QK ATL5MP3-10QK ATL5MP3-15QK ATL5MP3-25QK
		50 x 4.6 mm			
	5 µm	33 x 4.6 mm			
		50 x 4.6 mm			
		75 x 4.6 mm			
		100 x 4.6 mm			
150 x 4.6 mm					
250 x 4.6 mm					
10.0 mm	15 µm	250 x 10.0 mm		ATL15MP.25M	ATL15MP3.25M
21.2 mm	15 µm	250 x 21.2 mm		ATL15MP/25M	AT15MP3/25M
28.0 mm	15 µm	250 x 28.0 mm		ATL15MP&25M	ATL15MP3&25M
10.0 mm	15 µm	250 x 10.0 mm		ATL15MP.25M	ATL15MP3.25M
21.2 mm	15 µm	250 x 21.2 mm		ATL15MP/25M	ATL15MP3/25M
28.0 mm	15 µm	250 x 28.0 mm		ATL15MP&25M	ATL15MP3&25M



Modulo-cart information
p.C7 to C.9