

PEO spacer containing linkers

MONOFUNCTIONNAL (mPEGs, blocked PEGs)

Products Description


The linkers displayed in this sheet contain:

- A **PEO chain** (PolyEthyleneOxy, a PEG chain that is monodisperse)
- a **Methyl/MethylOxy group** at one end of PEO chain (not-functional - blocked terminus)
- and one **functional group**: from a large variety: COOH, NH₂, SH, OH, blocking termini (t-but, t-boc, acetyl), as well as with other more reactive groups (NHS, TFP, Maleimide, Azide, Hydrazide).

mPEO _x	spacer	Group 2	-NH ₂	-NHS	-	-OH	-SH, Thiol	-N ₃	-Hyd
Group 1	PEO	COOH	Amine	and other NH ₂ -reacti	Maleimide and other SH-reactive	Hydroxyl	Sulfhydryl	Azide	Hydrazide
Methyl		See	See	See	See	See	See	See	


Blocking groups.

• **Methyl group** (CH₃-) and **Methoxy group** (CH₃-O-; Methyloxy, MeO) are equivalent groups when bound to a PEO unit :



Ex:

A methyl group coupled to the Oxygène of a PEO unit (PEO: [-O-]) = a MethOxy- group coupled to EthylOxy units



They are blocking group that end of the PEO chain of 'mPEO' and 'mPEG' linkers'.

• ask for reversibly-blocking groups are available such as **t-boc**, **t-**, **Fmoc**, **CBZ** that can be released in proper conditions. They are useful in organic synthesis to introduce protected functional groups. After deprotection, the functional group can be used for further synthesis or further applications.

ex.:

mPEO _x -amidoPEO _x -acid	➔ see FT-BH9511 or ask
mPEO _x -TFP ester	“
mPEO _x -amidoPEO _x -TFP ester	“
mPEO _x -NH-CO(CH ₂) ₃ CO-TFP ester	“
Carboxyl-dPEG@ ₄ -(m-PEO _x) ₃	“
NHS/TFP-PEO ₄ -(mPEO _x) ₃ -ester	“

Functional groups.

Carboxyl group (COOH) of the linkers can be coupled to Amine and other groups by standard chemistry.

Amine group (NH₂) of the linkers can be coupled to carboxyl using standard chemistry, i.e. using carbodiimides mediated amidation (EDC #UP52005). Guidelines for reaction of amino group with carboxyls:

Sulfhydryl group (SH) of the linkers can be coupled to free Sulfhydryls by addition, or to reduced sulfhydryl (dissulfides) by replacement, or to Maleimide or PyridylVinyl sulfone (rapid reaction).

Hydroxyl group (OH) of the linkers can be coupled by standard chemistry.

Succinimidyl group (NHSuc: N-Hydroxy-Succinimidyl) of the linker can be reacted with a amines containing compounds, e.g., a protein, a peptide, to yield a stable amide bound.

Maleimide group (MAL) of the linker can be reacted with a thiol/sulfhydryl modified compound, e.g., a peptide

For any question,
contact you local distributor
uptima@interchim.com

Uptima,
powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

FT-0B6601

drug, oligo, small molecule or thiol containing protein, to yield a stable bound.

See more information about the functional group, notably reactivity, in the technical sheet NT-PEO000

Azide group (**HYD**) of the linker is aldehyde/ketone reactive.

Hydrazide group (**HYD**) of the linker is aldehyde/ketone reactive.

t-boc-HYD is a protected hydrazide group, which can be deprotected using TFA or HCl.

Aldehyde group (**ALD**) present several interesting reactions (with amine, hydrazide,... groups). It exist à propionaldehyde (ALP) an butyraldehyde (ALB), but also¹ amide-propionaldehyde, urethane-propionaldehyde, amide-butyraldehyde, urethane-butyraldehyde.

mPEO_x (Methyl-PEO linkers, MeO-PEO' linkers)

mPEO_x-COOH

These pegylation reagents react with amines

mPEO₂-COOH **0B6581, 100mg 0B6582, 1g**

MW: 148.16 (single compound) – CAS : 149577-05-9

Spacer length: 10 atoms and 10.9 Å

mPEO₄-COOH **0B6591, 100mg 0B6592, 1g**

MW: 236.26 (single compound) - CAS: 67319-28-2

Spacer is 14 atoms and 15.6 Å

mPEO₈-COOH **0B6601, 100mg 0B6602, 1g**

MW: 412.47 (single compound)

Spacer is 26 atoms and 29.8 Å


mPEO₁₂-COOH **0B6611, 100mg 0B6612, 1g**

soluble in methylene chloride or water; other organic solvents like ethyl acetate and THF can be considered

MW: 568.68 (single compound)

Spacer is 38 atoms and 44.0 Å

mPEO₂₄-COOH **0B6621, 100mg 0B6622, 1g**




MW: 1117.31 (single compound) - CAS: 117786-94-4

Spacer is 74 atoms and 86.2 Å


mPEO₂₅-COOH **0B6631, 100mg 0B6632, 1g**

MW: 1161.37 (single compound)

Spacer is 77 atoms and 89.9 Å




mPEO₃₇-COOH **0B6641, 100mg 0B6642, 1g**



MW: 1690.0 (single compound) – CAS : 117786-94-4

Spacer is 112 atoms and 133.9 Å

mPEO₄₈-CO(CH₂)₃COOH **0B6661, 100mg 0B6662, 1g**



MW: 2259.70 (single compound)

Spacer is 151 atoms and 178.5 Å

Ask for polydisperse version mPEG-COOH [\[FT-W1\]](#)


FT-0B6601

mPEO_x-NH₂ (Amine)

These pegylation reagents react with acids, active esters and aldehydes (Carbonyl/carboxyl reactive PEOs)


mPEO₂-NH₂ 0B6671, 100mg 0B6673, 1g

MW: 119.16 (single compound)
Spacer length: 10 atoms and 10.9 Å



mPEO₄-NH₂ RJ2151, 100mg RJ2153, 1g

MW: 207.27 (single compound) - CAS #: 85030-56-4;
MethylOxy-PEO₃-Amine; Methyl-PEG₄-Amine; MA(PEG₄)
Spacer length: 14 atoms and 15.5 Å




mPEO₈-NH₂ RJ2161, 100mg RJ2163, 1g

MW: 383.48 (single compound) - CAS #: 869718-81-0
Spacer length: 26 atoms and 29.7 Å

mPEO₁₂-NH₂ RJ2171, 100mg RJ2173, 1g


MW: 559.69 (single compound) - CAS #: 32130-27-1
Spacer length: 38 atoms and 43.9 Å



soluble in methylene chloride or water; other organic solvents like ethyl acetate and THF can be considered

mPEO₂₄-NH₂ RJ2181, 100mg RJ2183, 1g

MW: 1088.32 (single compound) - CAS: 32130-27-1
Spacer length: 74 atoms and 86.1 Å




mPEO₃₆-NH₂ 0B6681, 100mg 0B6683, 1g

MW: 1616.95 (single compound) - CAS: 32130-27-1
Spacer length: 109 atoms and 1300 Å

mPEO₄₈-CO(CH₂)₃-NH₂ 0B6691, 100mg 0B6693, 1g

MW: 2145.56 (single compound) - CAS: 32130-27-1
Spacer length: 146 atoms and 173.1 Å



Ask for polydisperse version mPEG-Amine [\[FT-W1\]](#)

mPEO_x-NHSuc (N-Hydroxy-Succinimidyl ester)

These pegylation reagents react with Amines (NH₂)

mPEO₂-NHS A1LYS1, 100mg A1LYS3, 1g

MW: 245.23 (single compound) - Spacer is 8 atoms and 8.5 Å - CAS:1127247-34-0

mPEO₄-NHS BH9061, 100mg BH9063, 1g

MW: 333.33 (single compound) - Spacer is 14 atoms and 15.6 Å - CAS:622405-78-1

mPEO₈-NHS BH9131, 100mg BH9133, 1g

MW: 509.54 (single compound) - Spacer is 26 atoms and 29.8 Å - CAS:756525-90-3

mPEO₁₂-NHS BH9501, 100mg BH9503, 1g

MW: 685.75 (single compound) - Spacer is 38 atoms and 44 Å - CAS:756525-94-7

mPEO₂₄-NHS RJ2001, 100mg RJ2002, 1g

MW: 1214.39 (single compound) - Spacer is 74 atoms and 86.2 Å - CAS:756525-94-7

mPEO₂₅-NHS A1LYT1, 100mg A1LYT3, 1g

MW: 1258.44 (single compound) - Spacer is 77 atoms and 89.9 Å - CAS:N/A

mPEO₃₇-NHS A1LYU1, 100mg A1LYU3, 1g

MW: 1787.07 (single compound) - Spacer is 112 atoms and 133.9 Å - CAS:756525-94-7

Ask for polydisperse version mPEG-NHSester [\[FT-W1\]](#)

For any question,
contact you local distributor
uptima@interchim.com

Uptima,
powered by



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

FT-0B6601

mPEO_x-TFP ester

These pegylation reagents react with Amines (NH₂)


mPEO₁₂-TFP ester	A1LYP0, 100mg	A1LYP2, 1g
MW: 736.75(single compound) - Spacer is 38 atoms and 44.0 Å - CAS: N/A		
mPEO₂₄- TFP ester	A1LYQ0, 100mg	A1LYQ2, 1g
MW: 1265.38(single compound) - Spacer is 74 atoms and 86.2 Å - CAS: N/A		
mPEO₂₅- TFP ester	A1LYR0, 100mg	A1LYR2, 1g
MW: 1309.42(single compound) - Spacer is 77 atoms and 89.9 Å - CAS: N/A		

Ask for polydisperse version mPEG-TFPester^[FT-W]

mPEO_x-Maleimide (MAL)


These pegylation reagents react with sulfhydryls (thiols, -SH)

mPEO₄-Maleimide	A1LYG0, 100mg	A1LYG2, 1g
MW: 358.39 (single compound) – Spacer is 20 atoms and 22.1 Å - CAS #: 1263044-81-0		
mPEO₈-Maleimide	A1LYH0, 100mg	A1LYH2, 1g
MW: 534.60 (single compound) – Spacer is 32 atoms and 36.4 Å - CAS: 1334169-90-2		
mPEO₁₂-Maleimide	HH5771, 100mg	HH5772, 1g
MW: 710.81(single compound) - Spacer is 44 atoms and 50.7 Å - CAS: 88504-24-9		
mPEO₂₄-Maleimide	RJ2021, 100mg	RJ2023, 1g
MW: 1239.44(single compound) - Spacer is 80 atoms and 50.7 Å - CAS: 88504-24-9		
mPEO₃₆-Maleimide	A1LYI0, 100mg	A1LYI2, 1g
MW: 1768.07(single compound) - Spacer is 116 atoms and 137.2 Å - CAS: 88504-24-9		
mPEO₄₈-Maleimide	A1LYJ0, 100mg	A1LYJ2, 1g
MW: 2296.70(single compound) - Spacer is 152 atoms and 178.6 Å - CAS: 88504-24-9		



mPEO4-Maleimide

mPEO48-Maleimide



Other maleimide-PEOs and PEGs:

See also (purified) **mPEG_x-Maleimide** #DZ5761 (FT-[AYPMB0](#))

mPEG-Maleimide

(Methyl-PEO₁₂)₃-PEO₄-Maleimide RJ1911

see [Branched Maleimide PEOs \(LV5640\)](#)

Maleimide-PEO₁₁-Biotin BR4031

see Biotinylated PEOs

mPEO_x-Azide (-N=N₂)

mPEO₂-Azide **06581, 100mg**


MW: 1 (single compound) –

Spacer length: 10 atoms and 10.9 Å

See also (purified) **mPEG_x-Azide** #WU0001 (FT-[AYPMB0](#))

FT-0B6601

mPEO_x-Hydroxyl (OH)



BH911: mPEO₁₁-OH; Methyl-PEO₁₁-Hydroxyl; CH₃-PEO₁₁-OH; CH₃-O-PEO'₁₀-OH; MethylOxy-mPEO'₁₀-Hydroxyl; OH-PEO'₁₀-O-CH₃; m-dPEG₁₁-OH

mPEO₇-Hydroxyl

MW: 340.41 - Spacer is 23atoms and 26.1Å -CAS:4437-01-8; Hepta(ethylene glycol) methylether; Hepta(ethylene glycol) methylether

mPEO₁₁-Hydroxyl

MW: 516.62(single compound) - CAS: 114740-40-8(9004-74-4) Undecae(ethylene glycol) methylether- Spacer is 35 atoms and 40.3 Å

mPEO₁₅-Hydroxyl

MW: 692.83(single compound) - CAS: 114740-40-8(9004-74-4) - Spacer is 46 atoms and 54.7 Å

mPEO₂₃-Hydroxyl

MW: 1045.25(single compound) - CAS: 114740-40-8(9004-74-4) - Spacer is 70 atoms and 83.1 Å

mPEO₃₆-Hydroxyl

MW: 1617.93(single compound) - CAS: 114740-40-8(9004-74-4) - Spacer is 109 atoms and 130.1 Å

mPEO₄₈-Hydroxyl

MW: 2160.59(single compound) - CAS: N/A (9004-74-4) - Spacer is 146 atoms and 172.4 Å

AXCJZ0, 100mg

AXCJZ2, 1g

BH9110, 100mg

BH9112, 1g

A1LYL0, 100mg

A1LYL2, 1g

A1LYM0, 100mg

A1LYM2, 1g

A1LYN0, 100mg

A1LYN2, 1g

A1LYO0, 100mg

A1LYO2, 1g

Ask for purified versions of mPEG-OH (oligodisperse - 160 to 40000Da) #IO5021 (FT-[AYPMB0](#))

See also hydrolyzed PEGs (PEG-OH, 10-26KDa ; CAS: 9004-74-4)^[B36US0]

mPEO_x-Thiol (SH)

mPEO₂-Thiol

Inquire, 100mg , 1g

See also mPEG_x-Thiol MW:350-30000Da #B36EB1 (FT-[AYPMB0](#))

mPEO_x-Aldehyde (CHO)

mPEO₂-Aldehyde

Inquire

Methyl-PEO₈-propionaldehyde

BS4H31

AWIMM


CAS: 1234369-95-9 ; (mPEG8-pALD, dPEG8-Aldehyde, MeO-PEG7-CHO ; MW :396.5

Methoxy-PEO₂₄-propionaldehyde

RPX89

(mPEG24-pALD, dPEG24-Aldehyde, MeO-PEG24-CHO)

AWIMM



See also mPEG_x-Aldehydes MW:350-30000Da (e.g. #WT8911 in FT-[AYPMB0](#))

350D # B36BL, 550D #B36BM, 750D #B36BN; 1K # WT893, 2K #WT892, 5K #WT891, 10K #WT894, 20K #WT895, 30K #WT896, 40K #KV851

But also propionaldehyde (ALP) an butyraldehyde (ALB), but also[] amide-propionaldehyde, urethane-propionaldehyde, amide-butyraldehyde, urethane-butyraldehyde. E.g. mPEG_x-Propionaldehydes MW:350-30000Da (FT-[BS4GX1](#))

References:

Hermanson, Greg T, "Bioconjugate Techniques", Academic Press, Inc., San Diego, CA, 1996.

Associated document:

NT-PEO000

[Homobifunctional PEO linkers \(BH9511\)](#)

[mPEG reagents \(DZ3531\)](#): mono-functional PEGs (mPEG-NHS, -MAL,-SH, -OH,...)

For any question,
contact you local distributor
uptima@interchim.com

Uptima,
powered by



R10E
213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60