


FT-18405



DNase I (RNase-free)

Storage: at -20°C for one year

Description

Deoxyribonuclease I (DNase I) is an endonuclease that degrades double and single-strand DNA and chromatin. It functions by hydrolyzing phosphodiester linkages, producing mono and oligonucleotides with a 5'-phosphate and a 3'-hydroxyl group. Ribonuclease has been reduced to non-detectable levels. Its activity depends on Mg^{2+} or Mn^{2+} ion. DNase I with Mg^{2+} cuts randomly double strand DNA at any site, DNase I with Mn^{2+} cuts double strand DNA at the same site to form sticky-end with 1-2 nucleotide or blunt-end.

Source

purified from bovine pancreas

Molecular Weight

32 kDa (monomer)

Unit Definition

One unit is the amount of enzyme required to completely degrade 1 μg pBR322 plasmid DNA in 10 minutes at 37°C .

Activity Test Condition

40 mM Tris-HCl (pH 8.0), 10 mM MgSO_4 , 1 mM CaCl_2 , 1 μg of pBR322 DNA

Purity

Free of other DNA endonucleases and exonucleases, free of RNase

Storage Buffer

50 mM Tris-acetate (pH 7.5), 10 mM CaCl_2 , 50% (v/v) glycerol

10 \times Reaction Buffer

100 mM Tris-HCl (pH 7.5 at 25°C), 100 mM MgCl_2 , 1 mM CaCl_2

Components

DNase I (3 units/ μl)	1500 units
10 \times DNase I Reaction Buffer	2 \times 1 ml
200 mM EDTA	1 ml