

## Important recommendations

**Viromer® YELLOW** is fully compatible with cell culture media, sera or antibiotics. Seed cells in complete medium the day before transfection and replace with fresh medium before starting the experiment.

**Cell culture and plating:** Grow cells to reach 60...80% confluency at the day of transfection. Use the volume of complete medium as mentioned in the table below.

Multiwell plate type	96	24	6
<b>Adherent cells</b>			
Cells seeded per well	12,000	60,000	250,000
*Range	± 3,000	± 20,000	± 80,000
<b>Suspension cells</b>			
Cells seeded per well	48,000	240,000	1,000,000
*Range	± 12,000	± 80,000	± 320,000
Medium per well	0.1 ml	0.5 ml	2 ml

\* in reverse transfection protocols, cell numbers should be on the higher end

**Suspension cells:** These cells need more DNA, please start using the 1.5x transfection scale and go to 2.0x or 2.5x. See our complete online manual for further recommendations.

**Forward/reverse transfection:** In forward transfection protocols, cells are seeded the day before transfection and the transfection complexes are freshly prepared at the transfection day. Instructions for reverse transfection and use in high-throughput screening (HTS) are provided in the detailed manual on the website.



For further optimization, information and trouble-shooting, go to [www.viromer-transfection.com](http://www.viromer-transfection.com) and visit our support pages. ...do not hesitate to contact us!

Version Dec / 15

# VIROMER® YELLOW Quick Guide

## pDNA/mRNA transfection

Detailed manual on Lipocalyx website  
[www.viromer-transfection.com](http://www.viromer-transfection.com)

### Important!

Avoid contact of Viromer® YELLOW with dry ice. Always close the vial and tighten the cap immediately after use. Do not vortex. Storage: +2-8°C within the provided aluminum bag.

### Protocol steps:

- Day 0: Plating of cells
- Day 1: Transfection (preparation 10 min, incubation 15 min)
- Day 1-3: Final incubation before analysis (6-24h after transfection)

### Conditions of use and required materials:

Warm all reagents to room temperature. Complexes should be prepared freshly. Use sterile, DNase/RNase free and apyrogenic tips and tubes.

## Transfection Protocol: 3-condition optimization


Volumes given here support 24 or 96 well formats. For 6 well, scale up 4 fold.

- 1 Dilute your pDNA/mRNA to 11 ng/ $\mu$ l using Buffer YELLOW. Provide a volume of 135  $\mu$ l. >> Tube 1
- 2 Place a 0.6  $\mu$ l droplet of Viomer® YELLOW onto the wall of a fresh tube. Immediately add 14.4  $\mu$ l of Buffer YELLOW and vortex for 3-5 s. >> Tube 2  
**Always add Buffer YELLOW to Viomer® YELLOW, not vice versa!**
- 3 Pipette 135  $\mu$ l of the pDNA/mRNA solution from Tube 1 onto the 15  $\mu$ l of the Viomer® YELLOW solution in Tube 2. Mix swiftly and incubate for about 15 min at room temperature.
- 4 Add transfection complexes from step 3 to your cells. Titrate as per the table below to identify optimal conditions.

Transfection Scale	96 well		24 well		6 well	
	Transfer Volume per well	pDNA/ mRNA per well	Transfer Volume per well	pDNA/ mRNA per well	Transfer Volume per well	pDNA/ mRNA per well
low 0.5x	5 $\mu$ l	50 ng	25 $\mu$ l	250 ng	100 $\mu$ l	1000 ng
standard 1.0x	10 $\mu$ l	100 ng	50 $\mu$ l	500 ng	200 $\mu$ l	2000 ng
high 1.5x	15 $\mu$ l	150 ng	75 $\mu$ l	750 ng	300 $\mu$ l	3000 ng
	5x replicates		1x replicate		1x replicate	

- 5 Incubate cells as usual. Monitor pDNA/mRNA effects 6-24 h after transfection. Expression from mRNA can begin as early as 2 h.

## Workflow



## Final Transfection Protocol

During optimization you identified a specific transfer volume and scale. Please proceed with your optimal settings.

The table below is a protocol using the 1.0x transfection scale (standard). Please adjust all volumes according to your optimal transfection scale.

- 1 Start with diluting pDNA/mRNA to 11 ng/ $\mu$ l using Buffer YELLOW.
- |                  | 96 well      | 24 well      | 6 well       | comments   |
|------------------|--------------|--------------|--------------|--|
| 2 Viomer® YELLOW | 0.6 $\mu$ l  | 0.6 $\mu$ l  | 2.4 $\mu$ l  | Buffer YELLOW onto Viomer® YELLOW<br>Vortex immediately (3-5s) |
| Buffer YELLOW    | 14.4 $\mu$ l | 14.4 $\mu$ l | 57.6 $\mu$ l |  |
- 3 pDNA/mRNA in solution from Step 1  
Viomer® YELLOW solution from Step 2
- |                                     | 96 well     | 24 well     | 6 well      | comments                            |
|-------------------------------------|-------------|-------------|-------------|-------------------------------------|
| pDNA/mRNA in solution from Step 1   | 135 $\mu$ l | 135 $\mu$ l | 540 $\mu$ l | Mix swiftly and incubate for 15 min |
| Viomer® YELLOW solution from Step 2 | 15 $\mu$ l  | 15 $\mu$ l  | 60 $\mu$ l  |                                     |
- 4 Transfer volume replicates
- |                 | 96 well    | 24 well    | 6 well      | comments                |
|-----------------|------------|------------|-------------|-------------------------|
| Transfer volume | 10 $\mu$ l | 50 $\mu$ l | 200 $\mu$ l | Incubate cells as usual |
| replicates      | 15x        | 3x         | 3x          |                         |
- 5 Expression from pDNA/mRNA can be monitored from 6 to 24 h after transfection. Expression from mRNA can begin as early as 2 h.