

# Product Information

## Steady-Luc™ Firefly HTS Assay Kit

### Steady-Luc Firefly HTS Assay Kit (Lyophilized)

#### Kit Contents:

#### Steady-Luc Firefly HTS Assay

Component	30028-T (40 assays)	30028-1 (120 assays)	30028-2 (1000 assays)	30028-3 (10,000 assays)
D-Luciferin	99907 1 mg	99907 3 x 1 mg	30028A2 25 mg	30028A2 10 x 25 mg
Steady-Luc Assay Buffer	30028B-T 4 mL	30028B 12 mL	30028B2 100 mL	30028B2 10 x 100 mL

#### Steady Luc Firefly HTS Assay Kit (Lyophilized)

Component	30028L-1 (120 assays)	30028L-2 (1000 assays)	30028L-3 (10,000 assays)
D-Luciferin	99907 3 x 1 mg	30028A2 25 mg	30028A2 10 x 25 mg
Steady-Luc Assay Buffer (Lyophilized)	30028L-B1 1 bottle	30028L-B2 1 bottle	30028L-B2 10 bottles
Steady-Luc Reconstitution Buffer	30028L-C1 12 mL	30028L-C2 100 mL	30028L-C2 10 x 100 mL

Number of assays is based on 96-well plate format.

#### Storage and Handling

Store Steady-Luc Firefly HTS Assay Kit (30028 series) at  $-70^{\circ}\text{C}$ . Kit components are stable for at least six months from date of receipt when stored as recommended. Avoid repeated freeze-thaw cycles.

Store Steady-Luc Firefly HTS Assay Kit (Lyophilized) (30028L series) at  $-20^{\circ}\text{C}$ . Kit components are stable for at least six months from date of receipt when stored as recommended.

#### Product Description

Firefly luciferase is widely used as a reporter for studying gene regulation and function, and for pharmaceutical screening (1, 2). It is a very sensitive genetic reporter due to the lack of any endogenous activity in mammalian cells or tissues (3, 4). Firefly luciferase is a 62,000 Dalton protein, which is active as a monomer and does not require subsequent processing for its activity. The enzyme catalyzes ATP-dependent D-luciferin oxidation by oxygen into oxyluciferin with emission of light centered at 560 nm (Figure 1).

However, the light production resulting from the reaction leads to formation of suicidal adenylyl-oxyluciferin at the enzyme surface. It results in very short half-life of the light emission with a flash-type kinetics. Several substances have been described to prolong light production by regenerating enzyme through removing inhibitory oxyluciferin from the enzyme surface (5, 6), but the signal duration (10-15 minutes) is still too short for batch process screening.

Biotium's Steady-Luc™ HTS assay system is a proprietary mixture of substances that modify the enzymatic reaction to produce a long lasting signal (steady glow) by preventing the formation of adenylyl-oxyluciferin at the enzyme surface. It is a homogeneous high sensitivity firefly luciferase reporter gene assay kit for the quantification of firefly luciferase expression in mammalian cells with signal half life of about 3 hours (Figures 2). Glow-type luciferase assays like Steady-Luc have lower luminescence signal compared to flash-type assays. The sensitivity and limit of detection of the assay will depend on luciferase expression levels in your experimental system as well as luminometer sensitivity.

Biotium's original Steady-Luc™ Firefly HTS Assay Kit (30028 series) contains assay buffer in liquid format. Steady-Luc™ Firefly HTS Assay Kit (Lyophilized) (30028L series) is a newer packaging format that includes lyophilized assay buffer for convenient room temperature shipping and storage at  $-20^{\circ}\text{C}$ .

#### References

1. Alam, J. and J.L. Cook. 1990. *Anal. Biochem.* 188:245-254.
2. Bronstein, I., et al. 1994. *Anal. Biochem.* 219:169-181.
3. Gould, S.J. and S. Subramani. 1988. *Anal. Biochem.* 175:5-13.
4. Brasier, A.R., et al. 1989. *BioTechniques.* 7:1116-1122.
5. Wood, K.V. 1990. *Proceedings of the Vth International Symposium on Bioluminescence and Chemiluminescence*, Cambridge, Ed. By P. Stanley and L. J. Kricka. p 543.
6. Airth, R.L., et al. 1958. *Biochemica and Biophysica Acta* 27:519-532.




Figure 1. Bioluminescent reaction catalyzed by firefly luciferase.




Figure 2. Steady-Luc (30028) and Steady-Luc (Lyophilized) (30028L) assays in transfected CHO-K1 cells. CHO-K1 cells were grown in F12-K medium containing 10% FBS in a white 96-well plate. On the day after plating, cells were transiently transfected with varying amounts of pGL3 firefly luciferase expression vector (Promega) using Viafect transfection reagent (Promega). On the day after transfection, the medium was replaced with fresh growth medium (100  $\mu\text{L}$  per well). The plate was allowed to equilibrate to room temperature, and 100  $\mu\text{L}$  Steady-Luc assay buffer or reconstituted Steady-Luc Assay Buffer (Lyophilized) containing D-Luciferin was added to each well. The plate was placed in a Bio-Tek Synergy H1 microplate reader and mixed with fast orbital shaking for five minutes. Luminescence was read every five minutes for five hours, with three seconds of orbital shaking before each read. Background luminescence values from untransfected cells were subtracted from luminescence values for each time point. Averaged luminescence values for duplicate wells are shown for various amounts of pGL transfected (ng/per well).

Continued next page

## Assay Procedure

**Note:** Steady-Luc luminescence signal has a half-life of about 3 hours, but may fluctuate over time or with temperature variation, and may vary depending on culture medium used. Therefore, raw luminescence values should be directly compared only for samples in the same medium. For comparison of luminescence signal between plates that are read at different times, each plate should include the same common internal control. The luminescence signals from each plate can be normalized to the internal control from the same plate.

**Note:** Steady-Luc assay should be carried out on cells or samples in cell culture medium containing magnesium. Luminescence signal will be low in the absence of magnesium.

### Protocol

1. Equilibrate the kit components to room temperature (22°C).
2. For lyophilized format (30028L series) only: prepare Steady-Luc Assay Buffer by adding reconstitution buffer to the bottle containing the lyophilized buffer. For 30028L-B1, add 12 mL reconstitution buffer; for 30028L-B2, add 100 mL reconstitution buffer. Mix by rocking until the buffer is a homogenous solution. Reconstitution buffer contains detergent; mix gently to avoid excessive foaming. Reconstituted assay buffer is stable at -20 °C for at least 3 months or -70 °C for at least 6 months. Avoid repeated freeze-thaw cycles.
3. To prepare Steady-Luc working solution, mix D-luciferin substrate and Steady-Luc Assay Buffer in 1 mg to 4 mL ratio. For each 1 mg vial of D-luciferin, mix with 4 mL Assay Buffer. For each 25 mg vial of D-luciferin, mix with 100 mL Assay Buffer. Add a small volume of Assay Buffer to the D-luciferin vial and mix by inversion until the substrate is completely dissolved, then transfer the D-luciferin solution to the full volume of Assay Buffer required. Only prepare working solution as needed for one day.  
  
**Note:** D-luciferin in Assay Buffer has limited stability. Instead of dissolving the entire contents of the D-luciferin vial in Assay Buffer, you may prepare a D-luciferin stock solution at 10 mg/mL in dH<sub>2</sub>O, and store it at -20°C or below for repeated use. The D-luciferin stock solution should be stable for at least one month, depending on the frequency of freeze-thaw cycles. The required volume of working solution can be prepared by diluting D-luciferin in Assay Buffer to a final concentration of 0.25 mg/mL (2.5 uL of 10 mg/mL D-luciferin stock solution per 100 uL assay buffer).
4. Remove plates containing luciferase-expressing cells from the incubator. If plates will be read in luminescence microplate reader, make sure plates are compatible with the instrument.
5. Add a volume of assay solution equal to that of the culture medium in each well and mix well. For example, for 96-well plates, add 100 uL assay solution to each well containing 100 uL of cells in medium, for a final volume of 200 uL per well.
6. Wait at least 5 minutes for complete lysis of the cells. Mixing on an orbital shaker during cell lysis is recommended.
7. Immediately before reading luminescence, mix samples thoroughly. Measure luminescence with a microplate luminometer. Alternatively, cell lysates can be transferred to tubes to be measured in a single sample luminometer.

## Related Products

Catalog number	Product
30003	Firefly Luciferase Assay Kit
30075	Firefly Luciferase Assay Kit (Lyophilized)
30004	Renilla Luciferase Assay Kit
30005	Firefly & Renilla Dual Luciferase Assay Kit
30020	ATP-Glo Bioluminometric Cell Viability Assay
22003	Mini Cell Scrapers, pack of 200

Please visit our website at [www.biotium.com](http://www.biotium.com) for information on our life science research products, including environmentally friendly EvaGreen® qPCR master mixes, fluorescent CF™ dye antibody conjugates and reactive dyes, apoptosis reagents, fluorescent probes, and kits for cell biology research.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.