



~The Kit for a new type smear preparation~

# *Smear Gell*<sup>TM</sup>

《 Manual ~ brief version ~ 》

~ Please read this manual before use ~

《 Any inquiries about this product 》

E-mail : [ipgell@genostaff.com](mailto:ipgell@genostaff.com)

Genostaff Co., Ltd.  
Kawauchi Bldg. 6F,  
1-4-4, Nezu, Bunkyo-ku, Tokyo  
113-0031 JAPAN  
81-3-3824-5856: TEL  
81-3-3824-5851: FAX  
HP: <http://www.genostaff.com>

## § 1. Preparation

### 1. Sample

1) To prepare cell suspension, re-suspend cells in small amount of cell culture medium or neutral buffer after centrifugation.

\* Fig. 1 is rough standard for concentrations of cell suspension. If the suspension looks to have enough conc. (approx.  $4 \times 10^3$  cells/ $\mu\text{L}$ ), it is unnecessary to carry out centrifugation and re-suspension.

\* When the cell sample is fixated, centrifuge it and discard the fixative. After washing the cells with PBS, re-suspend with mediums or buffers.

Re-fix the smear preparation with an appropriate fixative, after solidifying the sample on slideglass.

\* When the cell suspension contains cryoprotectant like  $\alpha$ Cell Bankerö or Glycerol, wash the cells and re-suspend with mediums or buffers.

2) Prepare 3 $\mu\text{L}$  of cell suspension per one slideglass.

< Table-1 Volume of mixture (cell suspension and solution I ) required for each number of slide >

		Number of slideglass				
		one	two	three	five	ten
Use as mixture	Cell suspension	3 $\mu\text{L}$	6 $\mu\text{L}$	9 $\mu\text{L}$	15 $\mu\text{L}$	30 $\mu\text{L}$
	Solution I	2 $\mu\text{L}$	4 $\mu\text{L}$	6 $\mu\text{L}$	10 $\mu\text{L}$	20 $\mu\text{L}$

\*Use 5 $\mu\text{L}$  of solution II per one slideglass

### 2. Reagents and Others

1) Thaw reagents (solution I , II ) of the Smear Gell kit in hands, just before to use and keep them at room temperature.

\* Solution I and II are able to repeat only 2 ~ 3 times FREEZE-THAW cycles.

2) Prepare the APS coated slideglasses (MATSUNAMI).

\* If you use slideglasses of kit attachment, return them to room temperature before use.

3) Prepare an appropriate fixative.


\* When you fix smear preparations using this kit on slideglasses, soak them whole in the fixative.

\* Choose appropriate fixative suitable for your staining method. DO NOT USE ACETONE for this kit.

## § 2. Method of preparing smear preparation using Smear Gell

### 【Preparing only one slideglass at once】

#### < Schematic >



#### 【Protocol】

1. Put 3  $\mu\text{L}$  of cell suspension on a slideglass.
2. Add 2  $\mu\text{L}$  of solution I to the cell suspension and mix by pipetting on same slideglass.
3. Add 5  $\mu\text{L}$  of solution II to the mixture on same slideglass, and pipetting 5 times immediately after, and furthermore smear the total mixture thinly with same tip, in a series of actions.

It is easier to smear that the tip is slanted like Fig. 3.


\* [Do not pipetting more than 5 times when you mix mixture A and solution II, and start smearing no sooner after finish 5 times pipetting.](#)

\* If you feel difficulty in smearing with normal tips, cut the head of tips like Fig. 2 and try to use it in Step 3.

4. Leave the slideglass at room temperature for 2 minutes to solidify the gel. **AVOID DRYING.**  
\* Do not leave more than 2 minutes.
5. Soak the slideglass in an appropriate fixative which suitable for your objective staining method.  
\* About 15 ~ 30 minutes fixation is enough.

## 【Preparing some slideglasses at once】




### < Schematic >



### 【Protocol】


1. Add solution I to the cell suspension and pipetting well. ( mixture A)
  - \* See Table-1
2. Put 5  $\mu\text{L}$  of solution II on a slideglass.
3. Add 5  $\mu\text{L}$  of mixture A to solution II on the slideglass, and pipetting 5 times immediately after, and furthermore smear the total mixture thinly with same tip, in a series of actions.  
It is easier to smear that the tip is slanted like Fig. 3.
  - \* [Do not pipetting more than 5 times when you mix mixture A and solution II, and start smearing no sooner after finish 5 times pipetting.](#)
  - \* If you feel difficulty in smearing with normal tips, cut the head of tips like Fig. 2 and try to use it in Step 3.
4. Leave the slideglass at room temperature for 2 minutes to solidify the gel. **AVOID DRYING.**
  - \* Do not leave more than 2 minutes.
5. Soak the slideglass in an appropriate fixative which suitable for your objective staining method.
  - \* About 15 ~ 30 minutes fixation is enough.

< Fig. 1 Rough standard for concentration of cell suspension >

7.5 x 10 <sup>4</sup> cells / slide (= 2.5 x 10 <sup>4</sup> cells /μL)	1x10 <sup>4</sup> cells / slide (= 3.3 x 10 <sup>3</sup> cells /μL)	1x10 <sup>3</sup> cells /slide (= 3.3 x 10 <sup>2</sup> cells /μL)
		


Hematopoietic cells (almost leukocyte) mag. x 100

< Fig. 2 Tip >




< Fig. 3 Point of smearing >

《SIDE VIEW》





《OVERHEAD VIEW》



It is suitable for smearing at appropriate thickness to make the tip shuttle about 5 times in approx. 2cm × 1.5cm area.  
**DO NOT CARRY OUT RE-SMEARING.**  
 Finish smearing in one action.

< Fig. 4 Example of smear preparation using Smear Gell >

- 
- \* Uniformly spreaded and appropriate thick Smear Gell containing cells.
  - \* Cells were almost monolayer.
  - \* It is easy to observe by a microscope.



HE stain

- △
- \* Smear Gell is likely to unstick due to its inequality.
  - \* Cells were multi-layer.
  - \* It seems difficult to observe by a microscope.

