


n-Octyl-β-D-Glucopyranoside

Product Information

Reference:	UP263700, 1 g UP263701, 5 g
Chemical name:	n-Octyl-β-D-Glucopyranoside (OG)
Structure :	C ₁₄ H ₂₈ O ₆
CAS :	29836-26-8
Molecular Weight :	292.38



Test	Specification – typical values
Appearance:	White to Off-White Powder
Purity	(HPLC)>99%
Absorbance (Solvent : D.I. water, Concen. 1%)	340 nm <0.01 280 nm <0.01 260 nm <0.015 225 nm <0.04
Percent Alpha	<0.6%
Percent Alcohol	<0.006%
% Fluorescence due to detergent (exc. 280 nm/em. 345 nm of a solution at 0.1% in H ₂ O+0.01%BSA)	0.2%
pH (1% water, 20°C)	6.8
Solubility in water / Solubilité dans l'eau	>20%
Conductance (10% solution)	<18μS

Technical and Scientific information

n-Octyl-β-D-Glucopyranoside (OG) is biological non-ionic detergent (neutral) widely utilized for solubilizing membrane bound proteins in their native state (used typically at 1%).

OG has been used at 100mM to incorporate Tissue factor (TF, Thromboplastin) into phospholipid vesicles ("relipidation")¹.

OG is used also in crystallography², in Mass Spectrometry (0.1%, improves Trypsin digestion)³.

OG removal can be achieved easily by dialysis because of the low CMC (20-25mM). Use a MWCO preferably above the OG micelles size (~25kD), even OG single molecule diffuse out from large micelle at smaller cut-off.

FT-263703

Related products and documents

Other detergents: see BioSciences catalogue

NT-Protein Extraction

Interchim, Division Biologie Interbiotech; Hotline : 04 70 03 73 06
interbiotech@interchim.com ; <http://www.interchim.com>

J08E