


Oxonol V

Sensitive slow-response membrane potential probe

Product Description

Name :	Oxonol V
Catalog Number :	Bis-(3-phenyl-5-oxoisoxazol-4-yl)pentamethine oxonol
	FP-352021 25 mg
	FP-352023 50 mg
	FP-352022 100 mg
Molecular Weight :	MW= 384.39
	C ₂₃ H ₁₆ N ₂ O ₄
Solubility:	DMSO, EtOH
Fluorescence:	$\lambda_{exc}/\lambda_{em}$ (MeOH) = 610/639 nm



Storage: Room temperature. Protect from light and moisture

Introduction

Oxonol V is a sensitive slow-response membrane potential probe that is widely used for measuring membrane potentials of many biological systems. The fluorescence of Oxonol V decreases upon membrane hyperpolarization. In general, slow-response probes exhibit potential-dependent changes in their transmembrane distribution that are accompanied by a fluorescence change. The magnitude of their optical responses is much larger than that of fast-response probes (typically a 1% fluorescence change per mV). Slow-response probes, which include cationic carbocyanines, rhodamines and anionic oxonols, are suitable for detecting changes in average membrane potentials of nonexcitable cells caused by respiratory activity, ion-channel permeability, drug binding and other factors.

Directions for use

Directions for Use

- Prepare a stock solution at 0.7 mM in ethanol and store at 4°C
- Dilute freshly to a final concentration of 3 μ M before each experiment
- Incubate coverslips containing wounded monolayers for 30 min at room temperature in the appropriate solution containing 3 μ M Oxonol V.
- Mount the coverslips and placed under a fluorescence microscope.

Other protocol may be found in the litterature.

References

- **Chifflet S. et al.**, A possible role for membrane depolarization in epithelial wound healing, *Am J Physiol Cell Physiol* 288: C1420-C1430 (2005) [Article](#)
- **Das T. et al.**, Mechanism of response of potential-sensitive dyes studied by time-resolved fluorescence, *Biophys. J.*, 64: 1122 - 1132 (1993) [Article](#)
- **Freedman J. et al.**, Quantitative analysis of oxonol V fluorescence in submitochondrial particles, *Ann. N.Y. Acad. Sci.*, 671: 493 (1992)

FT-352022

- **Morimoto T. et al.**, Voltage-Sensitive Oxonol Dyes Are Novel Large-Conductance Ca²⁺-Activated K⁺ Channel Activators Selective for beta1 and beta4 but Not for beta2 Subunits, *Mol. Pharmacol.*, 71: 1075 - 1088 (2007) [Article](#)
- **Moriyama Y. et al.**, One-step Purification of *Escherichia coli* H⁺-ATPase (F0F1) and Its Reconstitution into Liposomes with Neurotransmitter Transporters, *J. Biol. Chem.*, 266, 33: 22141-22146 (1991) [Article](#)

Related / associated products and documents

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

- DiBAC₄(3), FP-46600A
- DiBAC₄(5), FP-465992
- DiSBAC₂(3), FP-46603A
- DiSBAC₂(5), FP-HH6390
- DiOC₂(5), FP-BZ9350
- Oxonol VI, FP393141
- Acridine orange, FP-21092A
- Propidium iodide, FP-36774A

Ordering information

[Catalog size quantities and prices may be found at www.interchim.com/](#)

Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.