

# Nuclease S1

## Product Description

<b>Name:</b>	<b>Nuclease S1</b>
<b>Reference:</b>	35601A, 10000 Units                      35601B, 50000 Units
<b>Source:</b>	<i>Aspergillus oryzae</i> Chromatographically purified. Specific for single-stranded DNA (ssDNA) degradation. Activity on native (ds) DNA undetectable under the assay conditions. A frozen solution in 30mM sodium acetate, pH 4.6, 50mM NaCl, 1mM ZnCl <sub>2</sub> , and 50% glycerol.
<b>Unit Definition:</b>	One Unit hydrolyzes one microgram of denatured calf thymus DNA per minute at 37°C, pH 4.6.
<b>Molecular Weight:</b>	Approximately 32,000 - 36,000 daltons, exists as a monomer
<b>Optimum pH :</b>	4.0 - 4.6
<b>Activators:</b>	Zn <sup>2+</sup> and/or Ca <sup>2+</sup>
<b>Inhibitors:</b>	EDTA, citrate and a high concentration of SDS
<b>Storage:</b>	-20°C

## Introduction

Nuclease S1, isolated from certain *Neurospora* and *Aspergillus* species, specifically hydrolyzes both terminal and internal phosphodiester bonds of single-stranded DNA and RNA. It is used to eliminate non-annealed polynucleotide tails and hair-pin loops in DNA-RNA or DNA-DNA duplexes in hybridization studies and in genetic recombination experiments.

## Directions for use

### Handling and Storage

For long term storage in solution, for up to six months, dilute NUCSI to ≥6000 u/ml in water and freeze in aliquots. Dilute solutions can be stabilized by adding 0.1% albumin (Worthington Code: BSANF) and 10% glycerol.

### Protocol

One unit is the amount of enzyme liberating 1μg (0.033 A260) of acid-soluble nucleotides from heat-denatured DNA per minute at 37°C and at pH 4.6

### Reagents

Buffer: 0.2M NaCl, 0.002M ZnCl<sub>2</sub>, 0.06M CH<sub>3</sub>COONa, pH 4.6: dissolve 5.844 gms NaCl (MW 58.44), 136 mg ZnCl<sub>2</sub> (MW 136.29) and 1.85 ml concentrated glacial acetic acid in 450 ml reagent grade water. Adjust pH to 4.6 with 10M NaOH. Bring to a final volume of 500 ml with reagent grade water.

Enzyme Diluent: Dissolve 40 mg BSA in 200 ml Buffer

Substrate: Shred into small fibers, 60 mg calf thymus DNA and dissolve in 50ml reagent grade water by standing at room temperature for at least 18 hours. Additional stirring may be necessary to effect solution. Remove 10ml DNA solution to 10 ml of buffer. This is native calf thymus DNA solution (Substrate B)

Heat the remaining DNA solution in a large Pyrex test tube with a stir bar in boiling water on a heater/stirrer while stirring for 20 minutes. Immediately pour into **PRE-FROZEN** 1 liter beaker on ice. Mix equal volumes of the DNA

Contact your local distributor

Uptima, powered by

FT-35601A

solution and **cold** buffer. This is heat-denatured calf thymus DNA solution (Substrate A). Use as soon as possible to prevent the blank from elevating.

15% Perchloric Acid: Add 21.5 ml concentrated perchloric acid (70%) to 78.5 ml deionized water.

## Procedure

1. To clean glass tubes (two for each point) add 2ml Substrate B for tests. Include 2 tubes with 2ml substrate B for blanks (no enzyme added) and 2 tubes with 2ml Substrate A (native DNA test).
2. Incubate for 5 minutes before adding the enzyme.
3. Add 0.1ml enzyme dilution
4. Incubate at 37°C for 10 minutes.
5. Stop reaction by adding 2ml 15% perchloric acid.
6. Leave on ice for 10 minutes.
7. Centrifuge on a bench-top centrifuge for 15 minutes at 2000 rpm.
8. Withdraw 3ml supernatant and read A260.

## Calculation

$$\text{units / ml} = \frac{\text{A260 of sample} - \text{A260 of blank} \times \text{dilution} \times 1242}{10}$$

Where 1242 is a factor derived by dividing the reaction volume, 4.1ml, by the A260 of 1µg, which is 0.033, and dividing by the enzyme sample volume used, which is 0.1ml.

## References

- **Ando, T.:** A Nuclease Specific for Heat-Denatured DNA is Isolated from a Product of *Aspergillus oryzae*, *Biochim Biophys Acta* 114, 158, 1966
- **Britten, R., et al.:** Methods in Enzymology Vol. 29, L. Grossman and K. Moldave, *Academic Press, NY*, 363, 1974
- **Brookes, A., et al.:** Evaluation of the Use of S1 Nuclease to Detect Small Length Variations in Genomic DNA, *Eur J Biochem* 183, 291, 1989
- **Esteban, J., et al.:** Activation of S1 Nuclease at Neutral pH, *Nucleic Acids Res* 20, 4932, 1992
- **Geigl, E., et al.:** Chromosome-specific Identification and Quantification of S1 Nuclease-Sensitive Sites in Yeast Chromatin by Pulsed-Field Gel Electrophoresis, *Mol Microbiol* 4, 801, 1990
- **Gite, S., et al.:** Single-Strand-Specific Nucleases, *CRC Crit Rev Microbiol* 21, 101, 1995
- **Grafi, G., et al.:** Characterization of S1/Mung Bean Type Nuclease Activity in Plant Cell Suspensions, *Plant Sci* 74, 107, 1991
- **Manale, A., et al.:** S1 Nuclease as a Probe for the Conformation of a Dimeric tRNA Precursor, *Biochemistry* 18, 77, 1979
- **Maniatis, T., et al.:** Amplification and Characterization of a beta-Globin Gene Synthesized in vitro, *Cell* 8, 163, 1976
- **Marmur, J., et al.:** Progress in Nucleic Acids Res Vol. 1, J. Davidson and W. Cohn, *Academic Press, NY*, 280, 1963
- **Nedospasov, S., et al.:** Control Region of SV40 Minichromosomes Is Preferentially Cleaved by Single-Strand Specific S1 Nuclease, *J Biomolec Struct Dynamics* 6, 907, 1989
- **Ostrem, J., et al.:** Purification of S1 Nuclease from *Aspergillus-Oryzae* by Recycling Isoelectric Focusing, *Electrophoresis* 11, 953, 1990
- **Rudert, W., et al.:** Double-Labeled Fluorescent Probes for 5' Nuclease Assays: Purification and Performance Evaluation, *Biotechniques* 22, 1140, 1997
- **Vogt, V.:** Purification and Properties of S1 Nuclease from *Aspergillus*, *Methods in Enzymology* Vol. 65, L. Grossman and K. Moldave, *Academic Press, NY*, 248, 1980

Contact your local distributor

[uptima@interchim.com](mailto:uptima@interchim.com)

Uptima, powered by  
 213 Avenue J.F. Kennedy - BP 1140  
 03103 Montluçon Cedex - France  
 Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

FT-35601A

## Related / associated products and documents

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

[#39375](#) IminoBiotin ([UP39375](#)) and NHS-IminoBiotin ([UP35329](#))  
[page A352+](#) [other saturating agents for immunoassays \(i.e. SeaBlock \)](#)

- Deoxyribonuclease I, [243540](#)
- Proteinase K, [858700](#)

## Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Uptima is not liable for any damage resulting from handling or contact with this product.

Contact your local distributor

[uptima@interchim.com](mailto:uptima@interchim.com)

Uptima, powered by  
 **interchim**  
213 Avenue J.F. Kennedy - BP 1140  
03103 Montluçon Cedex - France  
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60