



CFDA-SE

Product Information

| Product name cat.number | MW (g·mol ⁻¹) | $\lambda_{exc} \backslash \lambda_{em}$ max. Coupled (nm) | |
|-------------------------------------|------------------------------|---|--|
| CFDA-SE FP-52493A, 25mg | 557.47 | 495 / 519 | 5-(and-6)-carboxyfluorescein Diacetate Succinimidyl ester « mixed isomers » Also called CFSE in the litterature, but this abbreviation should be used for the non diacetate form. |
| 5-CFDA-SE FP-AM496A, 10mg | 557.47 | 495 / 519 | 5-carboxyfluorescein Diacetate Succinimidyl ester « single isomers » |
| 6-CFDA-SE FP-AM497A, 10mg | 557.47 | 495 / 519 | 6-carboxyfluorescein Diacetate Succinimidyl ester « single isomers » |

| | |
|--------------------------------|--|
| Soluble: | In DMF or DMSO |
| Absorption / Emission : | $\lambda_{exc} \backslash \lambda_{em}$ (coupled hydrophyle) = 495 / 519 nm; pKa=6.4 |
| Storage: | Store at -20°C >1 year. (M) Protect from light and moisture |

Introduction

CFDA-SE is the most popularized long term cell tracer in tissues and organisms for survival or redistribution purposes. It allows further tracing of cell becoming during cell division /i.e. in embryogenesis studies, cell migration or moving / cancerogenesis, cell transplantation.... It is also used for labeling amine containing probes (antibodies, aminoallyl nucleic acids...).

CFDA-SE is colorless and passively diffuses into cells. After its acetate groups are cleaved by intracellular esterases, it becomes a highly fluorescent amine-reactive fluorophore, labeling covalently intracellular proteins, that keeps the fluorophore inside cells.

The fluorescence is compatible with subsequent fixation with formaldehyde or glutaraldehyde.

CFDA-SE overcomes alternative methods to reduce extracellular fluorescence quenching, including acidification of the extracellular medium, addition of tryptan blue or an antifluorescein antibody. It also achieves a more stable cell staining than obtained with isothiocyanate derivatives, because the amide linkage formed by the coupling reaction is much more stable than a thiourea linkage.

The litterature reports the advantage of use of CFDA-SE in previous bacterial transport studies. It is efficiently stained bacteria without affect the transport or viability of bacteria.[\(1\)](#)

Directions for use

Handling and Storage

Prepare a stock solution of CFDA-SE at a concentration 1000-fold higher than the final usage concentration (for example, 2 mM if the final concentration is 2 μ M) in DMSO.

CFDA-SE is EXTREMELY moisture sensitive!

Use dry DMSO. CFDA-SE will hydrolyze quickly at room temperature in the presence of water, and much more slowly at -20°C under desiccating conditions. You can aliquot and store over desiccant at -20°C or -70°C but should not be for more than 2 months.

Guidelines for use – on primary cells and cell lines ⁰

If your cells show decreased labeling with the same stock of CFDA-SE, hydrolysis is the likely cause.

Generally cells are labeled at a final CFDA-SE concentration of 0.5 to 5 μ M.

For best results, do a titration and find the lowest concentration of CFDA-SE that will give effective cell labeling – this will vary from cell type to cell type, and also with the application.

CFDA-SE labeling is somewhat toxic and can induce growth arrest and apoptosis in some cell types – therefore, it is important to find the lowest acceptable labeling concentration and check the viability after labeling.

As a rough guide, 0.5 to 2 μ M is usually enough for in vitro experiments – cell tracking and generational analysis in transplanted cells may require 2 to 5 μ M.

1. Prepare cells in PBS or HBSS containing 0.1 BSA%. Cell concentrations can range widely from 1×10^6 cells/ml (for in vitro experiments) up to 5×10^7 cells/ml.
2. Prepare a solution of CFDA-SE from stock solution in PBS/0.1% BSA at 2X the final labeling concentration. For example, if you are labeling at 5 μ M, prepare a 10 μ M solution. Prepare a volume of CFDA-SE equal to the cell volume above (no more than 2 ml per labeling reaction).
3. Add an equal volume of CFDA-SE solution to the cell suspension. Mix gently and incubate for 5 to 10 minutes at $+37^{\circ}\text{C}$.
Note : titrate to find the minimal effective conditions
4. Immediately fill the labeling tube to the top with the tissue culture media intended for culture (such as RPMI/10% FBS) and centrifuge. Wash the cells three times with tissue culture media at room temperature. To reduce the amount of unbound CFDA-SE in cells, incubate the cells at $+37^{\circ}\text{C}$ for 5 minutes after the second wash and prior to the third. This allows free unreacted CFDA-SE to diffuse out of the cells and be removed in the final wash.

Guidelines for use – for protein labeling

This method uses carboxyfluorescein succinimidyl ester (CFDA-SE) rather than fluorescein isothiocyanate, resulting in more reliable labeling. Succinimidyl esters are excellent reagents for amine modification since the amide products formed are very stable. CFDA-SE has high reactivity with aliphatic amines, low reactivity with aromatic amines, including tyrosine.

1. Prepare or otherwise obtain pure protein; make sure it is free of other contaminating proteins (e.g. albumin).
2. Ensure protein to be labeled is in a suitable buffer. Buffers containing TRIS are NOT acceptable since the TRIS interferes with labeling. A reasonable buffer to use is PBS, pH 7.4. If necessary, exchange current buffer for PBS using one of three methods:
 - a. Gel filtration. Do not use if you have less than 1 mg protein.
 - b. Dialysis. Microdialysis is probably the best method if you do not have very much protein.
 - c. Centricon microconcentrators.
3. Adjust protein concentration to ~ 1 mg/ml.
4. Prepare CFDA-SE, 1.5 mg/ml in anhydrous DMSO. Refer to "handling" section. Dilute CFDA-SE at desired concentration for labeling in anhydrous DMSO immediately before use. A starting labeling molar ratio is 3-20 Fluorescein to protein, but should be optimized for each molecule.
5. Add 100 μ l CFDA-SE per 1 ml of protein solution.

FT-52493A

6. Incubate for 90 minutes in the dark at room temperature with continuous gentle agitation. Alternatively, incubate overnight in the dark at 4 °C with continuous gentle agitation.
7. Exact conjugation efficiency depends on temperature, length of incubation, concentration of protein, concentration of CFDA-SE, and nature of protein. For critical applications, conjugate and test a small amount first to verify the conditions.
8. Separate labeled protein from free Fluorescein compounds by extensive dialysis versus PBS or by gel filtration. Concentrate with Centriprep and/or Centricon concentrators as necessary.
9. Microfuge on high for 10 minutes and filter through 0.22 µm filter. Optionally add 0.5% azide as preservative.
10. Determine concentration and F/P ratio by measuring absorbance at 280 nm and at 495 nm. For **antibodies**, use the following formulae to get approximate values:

$$\text{Protein concentration (mg/ml)} = (\text{OD}_{280} - 0.35 \times \text{OD}_{495}) / 1.4$$
$$\text{F/P ratio} = (3.3 \times \text{OD}_{495}) / (\text{OD}_{280} - 0.35 \times \text{OD}_{495})$$

Guidelines for use – on bacteria using microplate spectrofluorometry ⁰

Protocol may found in the literature, and in our technical notice [NT-Amine reactive](#)

References

- Bryan B.**, *et al.*, Coordinated Vascular Endothelial Growth Factor Expression and Signaling During Skeletal Myogenic Differentiation, *Mol. Biol. Cell*, 19: 994 - 1006 (2008) [Article](#)
- Fuller M.E.**, *et al.*, Application of a Vital Fluorescent Staining Method for Simultaneous, Near-Real-Time Concentration Monitoring of Two Bacterial Strains in an Atlantic Coastal Plain Aquifer in Oyster, Virginia, *Applied and Environmental Microbiology*, **70**, 1680 (2004) [Article](#)
- Fuller M.E.**, *et al.*, « Development of a Vital Fluorescent Staining Method for Monitoring Bacterial Transport in Subsurface Environments », *Applied and Environmental Microbiology*, **66**, 4486 (2000) [Article](#)
- Lee Y. K.**, *et al.*, Permanent Colonization by *Lactobacillus casei* Is Hindered by the Low Rate of Cell Division in Mouse Gut, *Applied and Environmental Microbiology*, **70**, 670 (2004) [Article](#)
- Linke B.**, *et al.*, CXCL16/CXCR6-mediated adhesion of human peripheral blood mononuclear cells to inflamed endothelium, *Cytokine* (2017) <https://doi.org/10.1016/j.cyto.2017.06.008>
- Mintern J.**, *et al.*, The use of carboxyfluorescein diacetate succinimidyl ester to determine the site, duration and cell type responsible for antigen presentation *in vivo*, *Immunology and Cell Biology*, **77**, 539(1999) [Abstract](#)
- Monneaux F.**, *et al.*, Selective Modulation of CD4+ T Cells from Lupus Patients by a Promiscuous, Protective Peptide Analog., *The Journal of Immunology*, 175: 5839-5847 (2005) [Article](#)
- Ueno A.**, *et al.*, Transient Upregulation of Indoleamine 2,3-Dioxygenase in Dendritic Cells by Human Chorionic Gonadotropin Downregulates Autoimmune Diabetes, *Diabetes*, 56: 1686 - 1693 (2007) [Article](#)

Related products

- GebaFlex Dialysis tubes, 12000-14000 MWCO, Mini, [AZ390A](#)
- Dimethylsulfoxide, anhydrous, [FP-JW7390](#)
- PBS Buffer Tabs, 1tab/100ml, [UP307157](#)
- FluoProbes® 490 NHS-ester, [FP-JO2820](#)

Ordering information

Catalog size quantities and prices may be found at <http://www.fluoprobes.com>

Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

Rev.H11E-F01VB