

Reacti-Bind™ NeutrAvidin™ Coated 96-Well Plates

0611.2

Number	Description
15127	Reacti-Bind™ NeutrAvidin™ Coated Plates (clear, 8-well strips), 5 each
15129	Reacti-Bind™ NeutrAvidin™ Coated Plates (clear, 96-well), 5 each
15116	Reacti-Bind™ NeutrAvidin™ Coated Plates (white, 96-well), 5 each
15117	Reacti-Bind™ NeutrAvidin™ Coated Plates (black, 96-well), 5 each Activation Level: 100 µl Binding Capacity: ~15 pmol D-biotin/well Blocking Buffer: These plates are supplied blocked with SuperBlock® Blocking Buffer
15128	Reacti-Bind™ NeutrAvidin™ Coated Plates (clear, 8-well strips), 5 each
15123	Reacti-Bind™ NeutrAvidin™ Coated Plates (clear, 96-well), 5 each
15216	Reacti-Bind™ NeutrAvidin™ Coated Plates (white, 96-well), 5 each
15217	Reacti-Bind™ NeutrAvidin™ Coated Plates (black, 96-well), 5 each Activation Level: 200 µl Binding Capacity: > 15 pmol D-biotin/well Blocking Buffer: These plates are supplied blocked with Blocker™ BSA

Storage: Upon receipt store plates at 4°C in unopened pouches. Once opened, place unused plates in a resealable bag with desiccant and store at 4°C.

Introduction

The Reacti-Bind™ NeutrAvidin™ Coated Plates are ideal for capturing biotin-labeled molecules without interference from nonspecific binding. NeutrAvidin™ Protein is deglycosylated avidin, which reduces lectin binding to undetectable levels while retaining stability and biotin-binding affinity. NeutrAvidin™ Protein offers the advantages of a near-neutral pI (6.3), to minimize nonspecific adsorption, and the lack of the RYD sequence, which eliminates nonspecific binding to the RGD binding domain of adhesion receptors present in a variety of cells. NeutrAvidin™ Protein yields the lowest nonspecific binding among the known biotin-binding proteins. The clear, white and black plates can be used with colorimetric, chemiluminescent and fluorescent detection methods, respectively.

Example ELISA Protocol using NeutrAvidin™ Coated Plates

A. Materials Required

- Wash Buffer: Tris-buffered saline (25 mM Tris, 150 mM NaCl; pH 7.2; Product No. 28376), 0.1% BSA, 0.05% Tween® -20; alternatively, use Blocker™ BSA (Product No. 37520) supplemented with 0.05% Tween® -20
- Biotinylated capture antibody adjusted to 10 µg/ml, or other appropriate concentration, with Wash Buffer
- Antigen adjusted to appropriate concentration with Wash Buffer
- Primary antibody adjusted to appropriate concentration with Wash Buffer
- Enzyme-labeled secondary antibody adjusted to appropriate concentration with Wash Buffer
- Appropriate enzyme substrate: example substrates are the TMB Substrate Kit (Product No. 34021) for horseradish peroxidase and the Phosphatase Substrate Kit (Product No. 37620) for alkaline phosphatase

Warranty: Pierce Biotechnology products are warranted to meet stated product specifications and to conform to label descriptions when stored and used properly. Unless otherwise stated, this warranty is limited to one year from date of sale when used according to product instructions. Pierce Biotechnology's sole liability for the product is limited to replacement of the product or refund of the purchase price. Unless otherwise expressly authorized in writing by Pierce Biotechnology, Pierce products are supplied for Research Use Only and are intended to be used by a technically qualified individual. Pierce Biotechnology's quality system is certified to ISO 9001. Pierce Biotechnology products are not produced in accordance with FDA's current Good Manufacturing Practices. Pierce Biotechnology strives for 100% customer satisfaction. If you are not satisfied with the performance of a Pierce Biotechnology product, please contact Pierce Biotechnology or your local distributor.

B. Method

1. Wash each well three times with 200 µl of Wash Buffer. Add 100 µl of the biotinylated capture antibody to each well and incubate for 2 hours at room temperature.
2. Wash each well three times with 200 µl of Wash Buffer. Make a serial dilution of the antigen and add 100 µl to each well. Incubate plate for 30 minutes at room temperature.
3. Wash each well three times with 200 µl of Wash Buffer. Add 100 µl of the primary antibody to each well and incubate plate for 30 minutes at room temperature.
4. Wash each well three times with 200 µl of Wash Buffer. Add 100 µl of the enzyme-labeled secondary antibody to each well. Incubate plate for 30 minutes at room temperature.
5. Wash each well three times with 200 µl of Wash Buffer.
6. Follow the manufacturer’s instructions for the specific detection system.

Procedure for Determining Binding Activity of the NeutrAvidin™ Coated Plates

The binding activity of the plates may be tested using Biotinylated Alkaline Phosphatase (Product No. 29339) and PNPP (Product No. 37620) or Biotinylated Horseradish Peroxidase (Product No. 29139) and TMB (Product No. 34021).

1. Rinse each well with three times with 200 µl of wash buffer (e.g., TBS).
2. Prepare a 1 mg/ml solution of the biotinylated enzyme. Make 1:2 serial dilutions using a 1:1,000 dilution for the first well. Incubate the wells for 1 hour at room temperature.
3. Wash each well three times with 200 µl of TBS containing 0.05% Tween®-20.
4. Incubate with 100 µl of substrate solution for 15 minutes at room temperature.
5. Measure the absorbance of each well. Active plates result in an absorbance of 0.5 to 1.0 OD at 405 nm.

Related Pierce Products

37070	SuperSignal® ELISA Pico Chemiluminescent Substrate , 100 ml, peroxidase substrate
15169	QuantaBlu™ Fluorogenic Peroxidase Substrate Kit
34028	1-Step™ Ultra TMB-ELISA , 250 ml, colorimetric peroxidase substrate
37621	1-Step™ PNPP , 100 ml, colorimetric phosphatase substrate
29339	ImmunoPure® Biotinylated Alkaline Phosphatase , 1 mg
29139	ImmunoPure® Biotinylated Horseradish Peroxidase , 5 mg
15075	ImmunoWare™ Reagent Reservoirs , 200/pkg.
15082	ImmunoWare™ Microtube Racked System , 960 tubes
15036	Sealing Tape for 96-Well Plates , 100/pkg.

References

Denlinger, L.C., *et al.* (2001). Cutting Edge: The nucleotide receptor P2X7 contains multiple protein- and lipid-interaction motifs including a potential binding site for bacterial lipopolysaccharide. *J. Immunol.* **167**:1871-6.

Ferre-Aubineau, V., *et al.* (1995). Colorimetric microtiter plate hybridization assay using monoclonal antibody for detection of an amplified human immunodeficiency virus target. *J. Virol. Meth.* **55**:145-51.

Hiller, Y. *et al.* (1987). Biotin binding to avidin. Oligosaccharide side chain not required for ligand association. *Biochem. J.* **248**:167-71.

Holmstrom, K., *et al.* (1993). A highly sensitive and fast non-radioactive method for detection of polymerase chain reaction products. *Anal. Biochem.* **209**:278-283.

Simon, M.D., *et al.* (2004). A phage display selection of engrailed homeodomain mutants and the importance of residue Q50. *Nucl. Acid. Res.* **32(12)**:3623-31.

Tween® is a registered trademark of ICI Americas.

SuperSignal® Technology is protected by U.S. Patent # 6,432,662.

QuantaBlu™ Technology is protected by U.S. Patent # 6,040,150 and # 6,437,179.

Current versions of product instructions are available at www.piercenet.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

©Pierce Biotechnology, Inc., 3/2006. Printed in the USA.