

Rhod-2

Product Information

Product name cat.number	MW (g·mol ⁻¹)	$\lambda_{exc}/\lambda_{em}$ max. Low Ca ²⁺ (nm)	$\lambda_{exc}/\lambda_{em}$ max. High Ca ²⁺ (nm)	mol. abs. (M ⁻¹ cm ⁻¹)	Kd (nM)	Soluble In
Rhod-2, AM ester bromide FP-661582 , 1mg FP-661584, 1mg Pure grade FP-661585, 50mg Pure grade FP-661584, 20x50µg FP-HH6180, 1 ml, 1mM DMSO, Pure grade	1123.97	548 nm / negligible	552 / 578	low Ca ²⁺ 91 000 high[Ca ²⁺] 96 000 [		DMSO
					570	DMSO, water pH>6
Rhod-2, K salt FP-M2038A , 1mg	869.08					
Rhod-2, NH ₄ salt FP-66159A , 1 mg	805.08					Buffer pH>6
Rhod-2, Na salt FP-AM709A , 1 mg	820.8					DMSO Buffer pH>6

Storage: **Indicator salts** can be stored desiccated and protected from light at room temperature, +4°C or -20°C >1 year.



AM esters can be stored desiccated and protected from light at -20°C > 6 months.

Introduction

Rhod-2 is non fluorescent before Ca²⁺ binding but becomes fluorescent with increasing Ca²⁺ concentration. The fluorescent enhancement for rhod-2 from low [Ca²⁺] to high [Ca²⁺] is smaller than that for fluo-3, and also in general rhod-2 is somewhat less fluorescent than fluo-3. However, the wavelengths of absorption (556 nm) and emission (576 nm) maxima of Rhod-2, are longer than those of Fluo-3. This may make it useful for some applications where sample autofluorescence is a problem, where sample contains photoreceptors, when photoactivatable chelators or another fluorescent dye of shorter wavelength is used at the same time. Its Kd of 570 nM places Rhod-2 above Fluo-3 and Fura-2, thus for higher Ca²⁺ concentrations studies.

Rhod-2 salt forms are membrane-impermeant but can be loaded into cells via microinjection or scrape loading.

Membrane-permeant form and thus can be loaded into cells via incubation. Because of the relatively low water solubility of the AM ester, Pluronic F-127, a mild detergent, is often used as a dispersing agent to facilitate the loading. Rhod-2/AM itself does not bind Ca²⁺, but it is readily hydrolyzed to rhod-2 by endogenous esterases once the dye is inside the cells.



Absorption (dashed curve) and Ca²⁺-dependent fluorescence (solid curves) spectra of Rhod-2 (10 µM) in solutions with various Ca²⁺ concentrations (0.04-40 µM). (Congwu Du, *Biophys. J.*, Jan 2001; 80: 549)

Directions for use

Handling and Storage

For loading, a 1-10 mM stock solution is prepared using dry (anhydrous) dimethylsulfoxide and kept at -20°C. In order to prevent the deterioration of the AM esters that occurs with repeated thawing and freezing, aliquots should be prepared. Loading is usually performed in a serum-free culture medium at a concentration ranging from 1-10 µM. Pluronic may be added to the loading medium if there is difficulty in loading. (Not all cells require Pluronic.)

The cells can be incubated at 20°C (or 37°C), and the time of incubation may range from 30 to 60 minutes. They should be washed two to three times with fresh serum-free culture medium for clean results. Note that incubation conditions may vary from cell to cell and should be optimized for each cell type. Protocols abound in the literature for a variety of cells.

Guidelines for use – calcium activity in skeletal muscle [Q](#)

1. Isolate single muscle fibres
2. Load the mitochondria by incubating fibres with 5µM Rhod-3 AM for 2h. Then wash for a further 2h.
3. After loading, fibres were mounted in the perfusion channel of a muscle bath where they could be stimulated electrically to produce a tetanus.
4. Dyes were excited with 568 nm light and the emitted light collected through a 585 nm long-pass filter.

References

- Akao M., *et al.*, Mechanistically Distinct Steps in the Mitochondrial Death Pathway Triggered by Oxidative Stress in Cardiac Myocytes, *Circulation Research*, 92:186 (2003) [Article](#)
- Congwu Du, *et al.*, Calcium Measurements in Perfused Mouse Heart: Quantitating Fluorescence and Absorbance of Rhod-2 by Application of Photon Migration Theory, *Biophys. J.*, **80**, 549(2001) [Article](#)
- Guidarelli A. *et al.*, Inhibition of complex III promotes loss of Ca²⁺ dependence for mitochondrial superoxide formation and permeability transition evoked by peroxynitrite, *Journal of Cell Science* 120, 1908-1914 (2007) [Article](#)
- Meyer, T., *et al.*, *Biochemistry* **29**, 32 (1990)
- Minta, A., *et al.*, Fluorescent indicators for cytosolic calcium based on rhodamine and fluorescein chromophores, *J. Biol. Chem.* **264**, 8171 (1989)
- Parker, I., *et al.*, *Science* **250**, 977 (1990)
- Pedro J. Del Nido, *et al.*, Fluorescence measurement of calcium transients in perfused rabbit heart using rhod 2, *Am J Physiol Heart Circ Physiol*, **274**, 728 (1998) [Article](#)
- Suzuki K. *et al.*, Dynamic recruitment of phospholipase Cγ at transiently immobilized GPI-anchored receptor clusters induces IP3–Ca²⁺ signaling: single-molecule tracking study 2, *The Journal of Cell Biology*, Vol. 177, No. 4, 731-742 (2007) [Article](#)
- Vergara, J., *et al.*, Imaging of calcium transients in skeletal muscle fibers. *Biophys. J.* **59**, 12 (1991)
- Yan Y. *et al.*, Bidirectional regulation of Ca²⁺ sparks by mitochondria-derived reactive oxygen species in cardiac myocytes, *Cardiovasc Res*, 10.1093/cvr/cvm047 (2007) [Abstract](#)
- Zima A. *et al.*, Termination of Cardiac Ca²⁺ Sparks: Role of Intra-SR [Ca²⁺], Release Flux, and Intra-SR Ca²⁺ Diffusion, *Circ. Res.*, Oct 2008; 103: e105 – e115 [Article](#)

Related Products

- | | |
|-------------------------------------|--|
| - Rhod-4 AM, CQ6061 | - Fluo-8L AM, CP7551 |
| - Fluo-8 NW, CJ2560 | - Fluo-8H AM, CP7531 |
| - Fluo-8 AM, CP7501 | - Fura-2 AM, FP-42776A |

Ordering information

Catalog size quantities and prices may be found at <http://www.fluoprobes.com>
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

