

## Phosphate Buffered Saline (PBS)

### Products Description

*Phosphate Buffered Saline Buffer (PBS) is a classic solution with physiological buffering range*

Composition of PBS (1X) is 137mM Sodium Chloride, 10mM phosphate, 2.7mM Potassium Chloride; pH is 7.4

<b>Name:</b>	<b>Phosphate Buffered Saline Buffer (PBS 1X)</b>	
<b>Cat.number:</b>	<b>UP68723A, 1 pack (10L of PBS 1x)</b> <b>UP68723B, 10 packs</b>	<b>UP30715, 100 tabs</b> (each tab makes 100ml) other format in TechSheet <a href="#">AWJGB0</a> .
<b>Form:</b>	powder, pack 1 pack (9.88g), when dissolved with distilled water, makes 1L of 10X PBS, or 10L of 1X PBS	tabs 1 tablet, when dissolved with distilled water, will give 100 ml of 1X PBS
<b>Use:</b>	<i>Reconstitute PBS 1x with distilled water:</i> When reconstituted with water, PBS1x final concentrations are 137mM Sodium Chloride, 10mM phosphate, 2.7mM Potassium Chloride; pH is 7.4	

<b>Name:</b>	<b>Concentrated Phosphate Buffered PBS Buffer 20X</b>	<b>Saline Buffer solutions (conc.PBS) PBS Buffer 10X</b>
<b>Cat.number:</b>	<b>N1376A, 500 ml of 20X (10L of PBS 1x)</b> <b>N1376B, 1 L of 20X (10L of PBS 1x)</b> 0.2 um filtered, Sterilised by autoclaving at 121°C.	<b>N14012, 1 L of 10X (10L of PBS 1x)</b> <b>N14013, 4 L of 10X (40L of PBS 1x)</b>
<b>Form:</b>	Clear liquid	Clear liquid
<b>Use:</b>	<i>Dilute 1:20 with distilled water.</i> 500ml 20X makes 10L of PBS 1x	<i>Dilute 1:10 with distilled. 1L</i> 1L 10X makes 10L of PBS 1x

**Storage:** Room temperature

### Applications:

PBS is used in various laboratory techniques, including immunodetection, microarraying, biochemistry, purification, cell culture. It can be supplemented by a variety of additives (chelators such as EDTA, detergents such as Tween20, proteins such as BSA, bacteriostatic such as azide, protease inhibitors such as AEBSF,...). See related products below for other PBS products including additives.

*Note: Phosphate buffer should not be used in assays where competition for phosphate groups, or complex formation with a metal ion is essential for the enzyme activation. Phosphate ions will inhibit carboxypeptidase, carboxylase, urease, muscle diaminease, formase and phosphoglucomutase.*

## Directions for use

Use ultrapure water to reconstitute PBS 1x from provided powder pack, tablet and concentrated solutions.

Note: If preparing **sterile PBS** from PBS 20X #N1376A, the added water should also be sterile (or choose the ready to use PBS 1x #N13522)

- **Preparing PBS 1X from powder packs UP68723A:**

Pour 1 pack (9.88g) in a clean recipient.

Mix well with 10L of dist.water until the solution is clear and without particulates.

Note: one also may prepare 10X PBS mother solution (1pack with 1L dist.water).

- **Preparing PBS 1X from concentrate solutions N1376A(20X) and N14012(10x):**

Pipette the desired quantity into a clean recipient.

Add distilled water to 1X concentration,

Mix well until the solution is clear.

- **Preparing PBS 1X from tablets UP307157:**

Put 1 tab in a clean recipient.

Dissolve in 100ml dist.water until the solution is clear and without particulates.

- **Use PBS 1X rapidly, or store in the fridge.**

Avoid preparing larger volume that required for the same day or few days: this might affect your buffer quality.

## Other information

### Related products

- \* Other PBS buffers:
  - \* PBS pH7.4, tablets for 100-200-500-1000ml ([AWJGB0](#))
  - \* PBS pH7.4, sterile solutions (N13522, [GS3574](#))
  - \* PBS for saturation in immunodetection techniques (BSA or Tween or Milk PBS blends)

PBS with Non-Fat Powdered Milk 3%, pH 7.5 , powder packs [#GS4180](#)

PBS with BSA 1%, pH 7.5, powder packs [#GS4190](#)

PBS with Tween® 0.05%, pH 7.5 , powder packs [#GS4250](#) (BB8210) (687230)

PBS with Calcium and magnesium: Dulbecco's PBS buffer: D-PBS 20X [#RJ226A](#)

\* Other Buffers See <sup>PH</sup>[BC002b](#), including Balanced Saline Solutions , Saline media,...  
e.g. Hanks' Balanced Saline Solutions (HBSS) ([NJI661](#)) ; Dulbecco's Modified Eagle Media #63426A;...

- \* Additives:

Detergents; Tween20 ([158740](#)), TritonX100 ([521121](#)), CHAPS ([33351A](#)), ...

Bacteriostatics: Na Azide ([081125](#)), ... Chelatants: EDTA ([036290](#))...

Protease inhibitors: AEBSF ([401071](#)), PMFS,... Reducing agents: DTT ([054721](#)), TCEP ([242214](#)), ...

\*Other buffers: TBS buffer (Tris Buffer Saline) ([74004A](#)), [buffers for genomics](#)

### Ordering information

For any information, please ask : Biosciences / Interchim; Hotline : +33(0)4 70 03 88 55

[Order on-line](#) or [Contact](#) your local distributor.

**Disclaimer** : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses.  
Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.T04E-S08VB