


HECAMEG, Glucose based detergent

A very mild Glucose based detergent effective for biological applications

Product Description

Catalog number::	UP785480, 5g
Name:	Hecameg 6-O-(N-Heptylcarbamoyl)-methyl- α -D-glucopyranoside Methyl-6-O-(N-heptylcarbamoyl)- α -D-glucopyranoside
Formula:	MW: 335.40 ; C ₁₅ H ₂₉ NO ₇ CAS: [115457-83-5] CMC: 19.5 mM
Storage:	Room temperature (powder) Following reconstitution, store in the refrigerator (4°C). Stock solutions are stable for up to 3 months at 4°C.



Applications

Extraction, purification, stabilization of proteins (recombinant or natural proteins)
Surfactant for chromatography, electrophoresis and ELISA analysis
Extraction of other biomolecules (DNA and RNA) from proteous samples
Study of protein structure, crystallisation of membrane proteins, enzymes or antigens ([Sauter 1999](#))
Liposomes preparation
Sanitization of chromatography columns

Scientific and Technical Information

A high quality detergent

Hecameg is a synthetic, well defined molecule, provided as a pure material (>98%) graded for biotech applications with consistent and reliable results.

It is soluble in usual aqueous buffers at +4°C or room temperature (>100mg/ml) and stable for weeks at +4°C. Molecular Weight is 335.4 g/mol. The micellar concentration is CMC= 19.5mM (0.65%) which allow it's easy removal by dialysis. Aggregation number (H₂O₂) is ~ 92.

Stability is >2years as supplied. Stock solutions are stable for up to 3 months at +4°C.

An effective surfactant detergent that preserves proteins

Hecameg dissociates aggregated proteins, helps breaking biological membranes

Hecameg does not denature proteins, enzymes or antigens, because it is non charged

Hecameg does not interfere with their biological activity, as shown for more than hundred enzymes, antigens and receptors.

FT-785480

Extraction

Hecameg was used for extraction of proteins from Chromaffin granules of mammalian cells, at 4% (Hodel 1994), photosystem II core complex (Kouimtzoğlu 1994), and Heparan Sulfate ProteoGlycan (Kiran 1994). 50mM Hecameg with EDTA, was found to give the highest yields of active lectinic factor (involved in yeast flocculation) in comparison to other detergents (El-Behhari 1998).

Hecameg was used at several steps of the isolation of Bf6 cytochrome from tylakoids of *Chlamydomonas* alga (Pierre 1995): cells were suspended in saline buffer with 25mM Hecameg, and centrifuged. The supernatant was fractionated by centrifugation in 10-30% w/w sucrose density gradient in presence of 20mM Hecameg and 0.1g/l egg phosphatidyl choline. Lastly, affinity chromatography was performed on hydroxylapatite, eluting with 20mM Hecameg plus 0.1g/l egg phosphatidyl choline.

Purification

Addition of 0.05% w/v **Hecameg** enhances recovery of material from hydroxyapatite and Q-Sepharose columns, and decreases elution volumes (Gerngross 1994).

Analysis

Hecameg produced the best diffracting crystals of Cytochrome bc1 complex (Lee 1995), in comparison with octyl- β -D-glucopyranoside, MEGA-9, n-octanoylsucrose and octyl- β -D-maltopyranoside. This was attributed to a better stability of proteins.

Hecameg was shown effective for reconstitution procedures in which detergents must be removed by dialysis, and for the lipid solubilization and uptake of vesicle contents at concentration well below the solubilizing range (BegonaRuiz 1994)

Literature

.Begona Ruiz M., and al; An assesment of the Biochemical Applications of the Non Ionic Surfcant HECAMEG; *Biochimica et Biophysica. Acta*, 1994, 1193- 301-306

.El-Behhari and al, Comparative extraction procedure for a galactosidase-specific lectin involved in flocculation of *Kluyveromyces lactis* strains; *Appl. Microbiol. Biotechnol.* 1998; 49; 16-43

.Gerngross and al; Overexpression and Purification of the Soluble Polyhydroxyalkanoate Synthetase (PHA synthetase) from *Alcaligenes eutrophus*: evidence for a Required Posttranslational Modification for Catalitic activity; *Biochemistry* 1994, 33, 9311-9320

.Kiran B. M., and al; Purification and Characterization of Adipocyte Heparan Sulfate Proteoglycans with affinity for Lipase; *The journal of Biological Chemistry* 1994, v.269, n°38, 23838-23844

.Kouimtzoğlu ad al., Effect of Mild Trypsinization on the Properties of a (17. 23 a,d 33kDa) retaining Photosystem (PS II) core Complex; *Progress in Photosynthesis Research* 1995, in press

.Lee J.W. and al; Preliminary Cryocrystallographic study of the Mitochondrial Cytochrome bc1 Complex: Improved Crystallisation and Flash-cooling of a large Membrane Protein; *Journal of Molecular Biology* 1995, 252, 15-19

.Pierre Y. and al; Purification and Characterization of the cytochrom B6f Complex from *Chlamydomonas reinhardtii*; *J.Biol. Chem.* 1995, 270, 29342-29349

.Plusquellec, D. et al. (1989) *Anal. Biochem.* 179, 145-53 [Abstract](#)

.Sauter Claude, Joseph D. N, Bernard Lorber, Gerard Keith, Philippe Brion, Mir Wais Hosseini, Jean-Marie Lehn, Richard Giege (1996); Additives for the crystallization of proteins and nucleic acids; *Journal of Crystal Growth* 196; 365-376 [Article](#)

Other Information

For R&D use only.

Not regulated as a Hazardous Material (IATA/ADR)

Merck Index: 13, 4636

For any question, please ask Uptima or your local distributor

rev : B07E