

PRODUCT INFORMATION

AMRESCO Destaining Bags (E732)

Remove Ethidium Bromide and other dyes from solutions safely, easily.

Ethidium Bromide is an extremely useful, easy-to-use dye for visualizing nucleic acids in solution and in agarose electrophoresis. Although the recent introduction of ethidium bromide solutions into the marketplace has eliminated many of the hazards of using the powdered dye for these applications, researchers are still faced with the difficult and unpleasant task of disposing of this toxic stain. The result is that many researchers neglect to remove ethidium bromide from their staining and electrophoresis buffers before disposing of them down the drain, which contaminates municipal waste-water with this toxic, mutagenic, carcinogenic substance.

AMRESCO now offers a way for you to remove Ethidium Bromide, Coomassie Blue, and other biological dyes and stains from solution for easy, safe disposal. AMRESCO Destaining Bags effectively extract milligram quantities* of ethidium bromide from solution without exposing research personnel to corrosive agents or contaminated substrates. The special adsorbent mixture retains dye molecules in a convenient bag for removal and incineration.

AMRESCO Destaining Bags are made with materials that can be safely disposed in incinerated trash.

- A 0.05 µg/ml solution of ethidium bromide has an A_{285} of approximately 0.01 at a 1X concentration.
- A 1 µg/ml solution of ethidium bromide has an A_{285} of approximately 0.1 at a 1X concentration.
- Changing bags during destaining will increase the rate of destaining and yield faster results.
- Destaining rates are independent of salt concentration.

* Each bag will extract up to 5 mg of ethidium bromide from solution during an overnight treatment.

ORDERING INFORMATION:

CODE	SIZE
E732	25 bags

Related Products:

Ethidium Bromide Dropper Bottle, Code E406
Ethidium Bromide 10 mg/ml Solution, Code X328
Ethidium Bromide High Purity Powder, Code 0492
Coomassie Brilliant Blue G-250, Code 0615
Coomassie Brilliant Blue R-250, Code 0472

Performance Evaluation with Typical Results:

Starting Concentration (µg/ml) 1 liter of solution	Beginning A_{285}	A_{285} at 1 hour (% removed)	A_{285} at 2 hours (% removed)	A_{285} at 4 hours (% removed)	A_{285} after overnight incubation (% removed)
0.05	0.0089	0.0061 (32%)	0.0025 (72%)	N/A	0.0008 (92%)
0.5	0.0586	0.0205 (65%)	0.0142 (76%)	N/A	0.0001 (99%)
1.0	0.1186	0.0443 (63%)	0.0291 (75.5%)	0.0113 (89%)	0.0025 (98%)
8.8	0.8908	0.0692 (92%)	0.0464 (95%)	0.0413 (95%)	N/A

093 9/00



TB-900

AMRESCO, Inc. (USA)
30175 Solon Industrial Pky. • Solon, OH 44139
TEL: 800-829-2805 • INTL.: 440-349-1313 • FAX: 440-349-1182
A2000 AMRESCO, Inc. All rights reserved. All trademarks are the property of their respective owners.



amresco-inc.com
EMAIL: info@amresco-inc.com