

HiScript II 1st Strand cDNA Synthesis Kit (+gDNA Wiper)

Catalog # R212



Version 7.0

Vazyme biotech co., ltd.

Introduction

The Vazyme HiScript II Reverse Transcriptase (+gDNA Wiper) is designed for the 1st strand cDNA synthesis with genomic DNA removal treatments. The HiScript II Reverse Transcriptase is a new generation reverse transcriptase optimized from the M-MLV (RNase H-) Reverse Transcriptase. The half-life of HiScript II at 50°C is > 240 min. Even at 55°C, the HiScript II can stay stable for a long time, which significantly benefits the transcription of RNA templates with complex secondary structures. In addition, the HiScript II has a improved template affinity and cDNA synthesis efficiency. It has a good resistance to most RT-PCR inhibitors and is suitable for long-fragment cDNA amplification (as long as 20 kb).

The residual genomic DNA in RNA template can be removed rapidly and completely after a treatment (42°C for 2 min) with the 4x gDNA Wiper. The 10x RT Mix contains an optimized buffer and dNTPs. The HiScript II Enzyme Mix contains the HiScript II Reverse Transcriptase and the RNase inhibitor. The Oligo-(dT)₂₃VN has a better affinity to Ploy A⁺ RNA than Oligo-(dT)₁₈. In addition, random hexamers and gene-specific primers (GSP) are also optional.

Contents of Kits

Components	R212-01 50 rxn (20 µl/rxn)	R212-02 100 rxn (20 µl/rxn)
RNase free ddH ₂ O	1 ml	1 ml
4x gDNA Wiper Mix	200 µl	400 µl
10x RT Mix ^a	100 µl	200 µl
HiScript II Enzyme Mix ^b	100 µl	200 µl
Oligo-(dT) ₂₃ VN (50 µM)	50 µl	100 µl
Random Hexamers (50 ng/µl)	50 µl	100 µl

a. contains dNTPs.

b. contains RNase inhibitor.

Storage

All the components should be stored at -20°C.

Protocol

Note: 1. Use high quality total RNA with high integrity for reverse transcription.

- To avoid RNase contamination, please keep the experiment area clean, wear clean gloves and masks, and use RNase-free tubes and tips.
- Primer selection** (Oligo-(dT)₂₃VN, Random hexamers, or GSP)

If the cDNA product will be used for PCR

- For eukaryotic RNA templates, generally, use Oligo-(dT)₂₃VN to obtain the highest yield of full-length cDNA.
- Use gene-specific primer (GSP) to obtain the highest specificity. However, switch to Oligo-(dT)₂₃VN or random hexamers if GSP fails in the 1st-strand cDNA synthesis.
- Random hexamers with the lowest specificity can be used for RNA templates, including mRNA, rRNA, and tRNA. Use random hexamers when Oligo-(dT)₂₃VN or GSP fails in cDNA synthesis due to complex secondary structure, high GC content, or prokaryotic RNA template.

If the cDNA product will be used for qPCR

- Use the mixture of Oligo-(dT)₂₃VN or random hexamers.

1. If the cDNA product will be used for PCR

1.1. RNA Denaturation*

Mix the following components in a RNase-free PCR tube:

RNase free ddH ₂ O	to 12 µl
Oligo (dT) ₂₃ VN (50 µM)	
or Random hexamers (50 ng/µl)	1 µl
or Gene Specific Primers (2 µM)	
Total RNA	10 pg-5 µg
or Poly A ⁺ RNA	10 pg-500 ng

Incubate at 65°C for 5 min and then chill on ice immediately for 2 min.

Note: * RNA denaturation benefits the cDNA yield. However, for cDNA < 3 kb, please skip the denaturation step.

1.2. Removal of Genomic DNA

Add 4 µl of 4x gDNA Wiper to the mixture of **Step 1.1** (12 µl), mix thoroughly, and incubate at 42°C for 2 min.

1.3. Mix the following components in a RNase-free PCR tube:

Mixture of Step 1.2.	16 µl
10x RT Mix	2 µl
HiScript II Enzyme Mix	2 µl

1.4. Start the 1st-strand cDNA synthesis.

25°C*	5 min
50°C**	45 min
85°C	2 min

Note: * Only necessary when using random hexamers. Please skip this step when using Oligo-(dT)₂₃ VN or Gene Specific Primers (GSP).

** For templates with complex secondary structure or high GC-content, the temperature can be increased to 55°C, which will benefit the yield.

The products can be used for PCR immediately or be stored at -20°C for 6 months. However, it is recommended to stored at -80°C and make aliquots to avoid repeated freezing and thawing.

2. If the cDNA product will be used for qPCR

2.1. Removal of Genomic DNA

Mix the following components in a RNase-free microtube by pipetting, and incubate at 42°C for 2 min.

RNase free ddH ₂ O	to 16 µl
4x gDNA Wiper Mix	4 µl
Oligo-(dT) ₂₃ VN (50 µM)	1 µl
Random hexamers (50 ng/µl)	1 µl
Total RNA	10 pg-1 µg
or Poly A ⁺ RNA	10 pg-100 ng

2.2. Mix the following components in a RNase-free PCR tube:

Mixture of Step 1.2.	16 µl
10x RT Mix	2 µl
HiScript II Enzyme Mix	2 µl

2.3. Start the 1st-strand cDNA synthesis.

25°C	5 min
50°C*	15 min
85°C	2 min

Note: * For templates with complex secondary structure or high GC-content, the temperature can be increased to 55°C, which will benefit the yield.

The products can be used for PCR immediately or be stored at -20°C for 6 months. However, it is recommended to stored at -80°C and make aliquots to avoid repeated freezing and thawing.

