



BAC DNA Fluorescent Probes

Quick Reference – Hybridization Protocol

Manual Denaturation/Hybridization



***Optional Automated Hybrite/Thermobrite Method**



Hybrite/Thermobrite Program Guide

For Peripheral Blood Preparations:

Denature at 72-73°C for 2 minutes/Hybridize at 37°C for at least 16 hours

For Paraffin Embedded tissue sections after pretreatment:

***Denature at 83°C for 3 minutes/Hybridize at 37°C for at least 16 hours - May require troubleshooting**

Optional FISH Pretreatment:

We recommend the Abbott Molecular FISH Pretreatment Kit:

http://www.abbottmolecular.com/FISHPretreatmentKit_5792.aspx

Paraffin Pretreatment:

We recommend the Abbott Molecular Paraffin Pretreatment Kits: I, II or III:

http://www.abbottmolecular.com/PretreatmentKits_5777.aspx