



Borondipyrromethene activated dyes for fluorescent labeling of biomolecules by click chemistry BrDIPY

BorondiPyrromethene dyes (BrDIPY, BDP) have small stockes shift, environment-independent high quantum yields, brightness thanks to sharp excitation and emission peaks, and high solubility in organic solvents. These dyes lack of ionic charge resulting in minimal effect on the isoelectric point of the conjugate. The small size and relatively long excited-state lifetime are useful for fluorescence polarization-based assays, notably for studying ligand-receptor interactions. The large two-photon cross-section is perfect for multiphoton excitation. The fluorescent dyes - DBCO (azodibenzocyclooctyne) conjugates can be conjugated with azides by copper free Click chemistry reaction, yielding stable triazole conjugates.



Storage:

LIFE SCIENCES

 -20° C, protected from light. Avoid prolonged exposure to light. Desiccate. Stable for 12 months at -20°C in the dark, and +1 week at room temperature.



FT-B35T42

BrDIPY FL DBCO is a photostable borondipyrromethene for FAM (fluorescein) channel. Absorption and emission wavelengths of BrDIPY FL are compatible with light sources and filter sets for fluorescein (FAM), but BrDIPY FL is significantly more photostable.

BrDIPY R6G is a bright and photostable substitute for Rhodamine 6G (R6G). BrDIPY stands for borondipyrromethene, a versatile fluorophore scaffold that is specially tuned in this molecule to match absorption and emission of R6G.

BrDIPY 581/591 is a borondipyrromethene dye with a conjugated olefinic system. It can be used either as a regular fluorophore, or as a probe for the detection of reactive oxygen species (ROS); after oxidation, its fluorescence moves to the green part of the spectrum.

Related products

* **AF dyes** functionalized by **NHS** (<u>R08112</u>), **Azide** (<u>AXCJ91</u>), **Alkyne** (<u>AXCECA</u>), **DBCO** (<u>AXCECA</u>), **Maleimide** (820731), **Hydrazide** (846631), **Amine** (), **Carboxyl** ()

* Alternative fluorescent dyes:

- FluoProbes (Superior brightness and photostability, i.e. FP488) -NHS (BA6800), -Azide (YE4970), -Maleimide (BA6810)

- CYanine dyes -Azide (HO7250)

- Classic dyes such as FAM, R110, JOE TAMRA, and ROX.

* Other labeling/conjugation chemistries: • Click Chemistry reagents

Ordering information

Catalog size quantities and prices may be found at <u>www.interchim.com/</u> Please inquire for higher quantities (availability, shipment conditions). For any information, please ask : FluoProbes[®] / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes[®] are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes[®] is not liable for any damage resulting from handling or contact with this product.

