

HotStart *Taq* 2X PCR Master Mix

Catalog: RK20603

Size: 100 RXN / 500 RXN

Concentration: 2X

Components:

HotStart *Taq* 2X PCR Master Mix RM20351

Product Description

HotStart *Taq* DNA polymerase is an innovative antibody-modified thermostable enzyme. The activity of the enzyme is completely closed at room temperature, avoiding the non-specific amplification and primer dimer in the preparation of the sample and the first cycle reaction temperature rise, increasing the specificity of DNA amplification.

HotStart *Taq* 2X PCR Master Mix is an optimized ready-to-use solution containing HotStart *Taq* DNA Polymerase, dNTPs, MgCl₂, KCl and stabilizers. It is ideally suited to routine PCR applications from templates including pure DNA solutions, bacterial colonies, and cDNA products. It can amplify up to 4 kb from complex genomic DNA or up to 5 kb from lambda DNA. Applicable to the PCR reaction, colony PCR, primer extension, *etc.*

Storage Temperature: -20°C

Heat Inactivation: No

5' - 3' Exonuclease: Yes

3' - 5' Exonuclease: No

Strand Displacement: +

Resulting Ends: Single-base 3' Overhangs

Error Rate: ~ 285x10⁻⁶ bases

1X Master Mix Composition:

10 mM Tris-HCl, 1.5 mM MgCl₂, 50 mM KCl, 0.08% IPGAL 630, 0.05% Tween 20, pH8.6@25°C; 200 µM dNTPs, 5% Glycerol, 25 U/ml HotStart *Taq* DNA Polymerase.

Order: order@abclonal.com
Tech: support@abclonal.com

For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.
 Please visit <http://abclonal.com.cn> for a complete listing of recommended products.

Instructions

Reaction setup:

We recommend assembling all reaction components on ice and quickly transferring the reactions to a thermocycler preheated to the denaturation temperature (95°C).

Take 25 µl /50 µl system as an example.

Composition	25 µl	50 µl	Final Conc.
Nuclease-free water	to 25 µl	to 50 µl	
10 µM Forward Primer	0.5 µl	1 µl	0.2µM (0.05~1 µM)
10 µM Reverse Primer	0.5 µl	1 µl	0.2µM (0.05~1 µM)
Template DNA	variable	variable	<1 µg/50 µl
HotStart <i>Taq</i> 2X PCR Master Mix	12.5 µl	25 µl	1X

Incubated in a thermocycler as the below program:

Temperature	Time	Cycles
95°C	5 min	1
95°C	15-30s	
45-68°C	15-60s	30
68°C	1kb/min	
68°C	5min	
4-10°C	∞	1

General Guidelines:

1. **Template:**

Use of high quality, purified DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 µl reaction are as follows:

DNA	Amount
Genomic	1 ng–1 µg
Plasmid or viral	1 pg–1 ng

2. **Primers:**

Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. Computer programs such as Primer3 (<http://frodo.wi.mit.edu/primer3>)

can be used to design or analyze primers. The final concentration of each primer in a reaction may be 0.05–1 μM , typically 0.1–0.5 μM .

3. **Mg⁺⁺ and additives:**

Mg⁺⁺ concentration of 1.5–2.0 mM is optimal for most PCR products generated with *Taq* DNA Polymerase. The final Mg⁺⁺ concentration in 1X HotStart *Taq* PCR Master Mix is 1.5 mM. This supports satisfactory amplification of most amplicons. However, Mg⁺⁺ can be further optimized in 0.5 or 1.0 mM increments using MgCl₂.

Amplification of some difficult targets, like GC-rich sequences, may be improved with additives, such as DMSO or formamide.

4. **Denaturation:**

An initial denaturation of 30 seconds at 95°C is sufficient for most amplicons from pure DNA templates. For difficult templates such as GC-rich sequences, a longer denaturation of 2–4 minutes at 95°C is recommended prior to PCR cycling to fully denature the template. With colony PCR, an initial 5 minutes denaturation at 95°C is recommended.

During thermocycling a 15–30 second denaturation at 95°C is recommended.

5. **Annealing:**

The annealing step is typically 15–60 seconds. Annealing temperature is based on the T_m of the primer pair and is typically 45–68°C. Annealing temperatures can be optimized by doing a temperature gradient PCR starting 5°C below the calculated T_m.

When primers with annealing temperatures above 65°C are used, a 2-step PCR protocol is possible.

6. **Extension:**

The recommended extension temperature is 68°C. Extension times are generally 1 minute per kb. A final extension of 5 minutes at 68°C is recommended.

7. **Cycle number:**

Generally, 25–35 cycles yields sufficient product. Up to 45 cycles may be required to detect low-copy-number targets.

8. **2-step PCR:**

When primers with annealing temperatures above 65°C are used, a 2-step thermocycling protocol is possible.

Thermocycling conditions for a routine 2-step PCR:

Temperature	Time	Cycles
95°C	5 min	1
95°C 65-68°C	15-30s 1kb/min	30
65-68°C 4-10°C	5min ∞	1

9. **PCR product:**

The PCR products generated using HotStart *Taq* DNA Polymerase contain dA overhangs at the 3'-end; therefore the PCR products can be ligated to dT/dU-overhang vectors.

Order: order@abclonal.com
Tech: support@abclonal.com

For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.
Please visit <http://abclonal.com.cn> for a complete listing of recommended products.