

# T4 DNA Ligase (High Conc.)



**Catalog:** RK21500

**Size:** 80,000 U / 400,000 U

**Concentration:** 2,000,000 U/ml

**Components:**

T4 DNA Ligase (2,000,000 U/ml)	RM21500
10X T4 DNA ligase Reaction Buffer	RM20108

## Product Description

T4 DNA Ligase can catalyze the formation of a phosphodiester bond between juxtaposed 5' phosphate and 3' hydroxyl termini in duplex DNA or RNA. This enzyme joins blunt end and cohesive end termini as well as repair single-stranded nicks in duplex DNA and some DNA/RNA hybrids. T4 DNA Ligase seals nicks for these DNA substrates. T4 DNA Ligase is applicable to cloning restriction fragments and to joining linkers and adapters to blunt-ended DNA.

**Product Source:** An *E. coli* strain that carries the T4 DNA ligase gene.

**Unit Definition:** One unit is defined as the amount of enzyme required to ligate 50% of *Hind*III digestion fragments of  $\lambda$  DNA (5' DNA termini concentration of 0.12  $\mu$ M, 300  $\mu$ g/ml) in a total reaction volume of 20  $\mu$ l over 30 minutes at 16°C in 1X T4 DNA Ligase Reaction Buffer.

**Storage Conditions:** 10 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH7.4 @ 25°C

**Storage Temperature:** -20°C

**Reaction Conditions:**

1X T4 DNA Ligase Reaction Buffer.

**1X T4 DNA Ligase Reaction Buffer:**

50 mM Tris-HCl, 10 mM MgCl<sub>2</sub>, 10 mM DTT, 1 mM ATP, pH7.5 @ 25°C

**Heat Inactivation:** 65°C for 10 min.

**Order:** [service@abclonal.com](mailto:service@abclonal.com)  
**Tech:** [support@abclonal.com](mailto:support@abclonal.com)

*For research purposes only. Not for therapeutic or diagnostic purposes in humans or animals.*  
 Please visit <http://abclonal.com> for a complete listing of recommended products.

## Instructions

- ◆ Set up the following reaction in a microcentrifuge tube on ice.

Composition	Amount
10X T4 DNA Ligase Reaction Buffer*	2 $\mu$ l
Vector DNA (4 kb)	50 ng (0.02 pmol)
Insert DNA (1 kb)**	37.5 ng (0.06 pmol)
Nuclease-free dH <sub>2</sub> O	up to 19 $\mu$ l
T4 DNA Ligase ***	1 $\mu$ l
Volume	20 $\mu$ l

\*:10X T4 DNA Ligase Reaction Buffer should be thawed and resuspended at room temperature.

\*\* Insert DNA (1 kb): a ligation using a vector to insert molar ratio of 1:3 for the indicated DNA sizes.

\*\*\*:T4 DNA Ligase should be added last.

- ◆ Short centrifugation after gentle percussion.
- ◆ Gently mix the reaction by pipetting up and down and microfuge briefly.
- ◆ For cohesive (sticky) ends, incubate at 16°C overnight or room temperature for 10 minutes.
- ◆ For blunt ends or single base overhangs, incubate at 16°C overnight or room temperature for 2 hours (alternatively, high concentration T4 DNA Ligase can be used in a 10-minute ligation).
- ◆ Heat inactivate at 65°C for 10 minutes.
- ◆ Chill on ice and transform 1-5  $\mu$ l of the reaction into 50  $\mu$ l competent cells.

**QC Process:**

- Purity is above 95% detected by SDS-PAGE.
- No exonuclease, nuclease, RNase contamination.
- No residual host genomic DNA detected by PCR.