



Chromatrap

Chromatrap[®] Homogeniser Spin Column

Product Number: 500289

Quickly and efficiently homogenise cells and tissues lysates in a single step. Contained within a spin column, our unique proprietary bio-shredding Vyon[®] polymer reduces lysate viscosity and captures insoluble debris by centrifugation.

One Stop Solution

- Perfect alternative to syringe & needle homogenisation
- Quick, clean and efficient homogenisation
- Reduces viscosity of lysate
- Ideal sample preparation for nucleic acid miniprep and midiprep

Protocol

1. Position the Chromatrap[®] Homogeniser Spin Column in the collection tube provided.
2. Transfer up to 700µl cell or tissue lysate to the Chromatrap[®] Homogeniser Spin Column and close the cap.
3. Centrifuge at full speed for 2 mins at room temperature.
4. Remove and discard the homogeniser spin column from the collection tube. **Do not discard the supernatant, this contains the homogenised lysate.**
5. Carefully aspirate the suspension from the collection tube using a pipette and transfer to a fresh clean 1.5ml microcentrifuge tube.
6. Sample is now ready for downstream processing (eg. RNA extraction).

Advancements in Epigenetics

Chromatrap[®] is a brand of Porvair Sciences Ltd. Clywedog Rd. South, Wrexham Industrial Estate, Wrexham LL13 9XS, United Kingdom
Technical Support: support@chromatrap.com **Sales:** sales@chromatrap.com

Company registered in England No. 2047818 Registered office; 7 Regis Place, Bergen Way, Kings Lynn, P330 2JN, United Kingdom