


IST-Fluolid, the solid fluorescence technology




Fluolid is

Features

1. **Brighter : high quantum yield** in the solid state
2. **no photobleach** - high stability for light, heat and pH
3. **higher labeling rate** is than that of traditional dye
4. **longer stocke's shift** – lowers cross-fluorescence in multiplex



[Superior fluorescence features of IST Fluolid]






UV spectra of Fluolid labeled IgG (left: green, center: yellow, right: orange Wavelength (Ex and Em))

	IST Fluolid	Traditional Dye
Stability of Light	good (Stable over 2 years under daylight)	bad
Stability of Temperature	Good (over 200°C)	bad (degrade at r.t.)
Stability of pH	good (stable for pH1-14)	bad (usable pH5-9)
Stokes Shifts	Wide (about 100-150nm)	Narrow (about 10-50nm)
Fluorescent in the solid state	Highly Fluorescent	quenching
Storage Condition	Under room temperature*	Under -20°C

	λ_{abs} /nm	λ_{em} /nm
Commercial Dyes	Stokes shift	
Acridine Orange	460/500	526/650
Cy3	550	565
Cy5	650	↔ 670
FITC	494	518
IST Dyes	Stokes shift	
Fluolid-W Orange	440	602
Fluolid-W Yellow	410	↔ 541
Fluolid-W Green	395	522

* The fluorescence of "IST Fluolid" has high stability under various conditions. But, to keep the activity of the succinimidyl ester group, please store at 2-8 °C.

Applications - [Labeling of Oligo DNA with IST Fluolid]






Bio Rad Molecular Imager FX 488nm $\text{fmol} = 10^{-15} \text{ mol}$

- **Almost 90-100% labeling of P72 oligo DNA** (traditional dye : 20-80%)
- **Measurement at low concentration with existing micro imager**
- **Detectable of the dried sample at low concentration ($\sim 10\text{-}12\text{mol}$)**


Applications - [Labeling of BSA using IST Fluolid]

BSA (Bovine Serum Albumin) Protein




Freeze dried BSA Protein showed strong fluorescence in the solid state.

Applications - [Labeled GFP antibody]



- The labeling of GFP antibody advanced quantity.
- The labeled GFP antibody keeps the activity.
- The labeled GFP antibody was trapped by fixed EGFP.

Application - [Multi-color Immunostaining-1]




Protocol of Multi-Color Immunostaining ¹:

Tissue Section: Rat Brain(Bregma-0.26mm)
 PBS: 0.1 M Phosphate Buffer(pH7.2~7.4)
 Primary Antibody A: Anti-GFAP Mouse monoclonal IgG (Astrocyte marker)
 Primary Antibody B: Anti-Iba1 Rabbit polyclonal IgG (Microglia marker)
 Secondary Antibody C: Anti-Mouse IgG, Biotin Conjugate
 Secondary Antibody D: Anti-Rabbit IgG, Alexa Fluor® 488 Conjugate


1. Washing in PBS for 20 min is repeated 3 times
2. Blocking by 1%BSA/PBS Solution for 1 hour under room temperature
3. Reaction with primary antibodies*/PBS for 1 day under room temperature [* Mixed with A and B.]
4. Washing in PBS for 10 min is repeated 3 times
5. Reaction with secondary antibodies*/PBS for 3 hours under room temperature [* Mixed with C and D.]
6. Washing in PBS for 10 min is repeated 3 times
7. Reaction with Streptavidin, Fluolid-W Orange 600 Conjugate/PBS for 1 hours under room temperature
8. Washing in PBS for 10 min is repeated 3 times
9. Mounting in Vectashield with DAPI
10. Observation of Immunostaining

Application - [Multi-color Immunostaining-2]

(Fluolid-W Orange 600 exclusive filter)



Application - [Observation of Lectin staining]



Protocol of Lectin Immunostaining ¹:

- A / Preparation of tissue sections - Rat Kidney
Animal: Rat aged of 7 months to 8 months
1. Perfusion fixation with solution: 2.8% paraformaldehyde, 0.2% picric acid, 0.1 mol/L Phosphate buffer(pH7.2~7.4)
 2. Infiltration for 2 H to 3 H [25% sucrose/ 0.1 M Phosphate Buffer(pH7.2~7.4)]
 3. Preparation of 8μm cryosection

B / Protocol of Lectin Staining

- PBS: 0.1 M Phosphate Buffer(pH7.2~7.4)
Lectin: Peanuts Lectin, biotin Conjugate
Streptavidin: Streptavidin, Fluolid-W Conjugate
1. Washing in PBS for 20 min is repeated 3 times
 2. Reaction with Lectin/PBS for 1 day under room temperature
 3. Washing in PBS for 10 min is repeated 3 times
 4. Reaction with Streptavidin, Fluolid-W Conjugate/PBS for 3 hours at room temp.
 5. Washing in PBS for 10 min is repeated 3 times
 6. Mounting in Vectashield with DAPI
 7. Observation of Lectin staining

For ordering:

Fluolid-W Orange 600 Protein Labeling Kit	DU7720, 1 KIT
Fluolid-W Yellow 540 Protein Labeling Kit	DU7740, 1 KIT
Fluolid-W Green 520 Protein Labeling Kit	DU7780, 1 KIT
Fluolid-W Orange 600 Oligonucleotide Amine Labeling Kit	DU7710, 1 KIT
Fluolid-W Yellow 540 Oligonucleotide Amine Labeling Kit	DU7750, 1 KIT
Fluolid-W Green 520 Oligonucleotide Amine Labeling Kit	DU7770, 1 KIT
Fluolid-W Orange 600 succinimidyl ester	DU7730, 1 MG
Fluolid-W Yellow 540 succinimidyl ester	DU7760, 1 MG
Fluolid-W Green 520 succinimidyl ester	DU7790, 1 MG

[Price and technical sheet on line](#)

Related products/documents

[Products HighLights Overview](#)

Including FluoProbes labeling dyes, Desalting tools,

Information inquire

Reply by Fax : +33 (0) 4 70 03 82 60 or email at interbiotech@interchim.com

I would like to receive further information on: _____

Title : _____ First name: _____ Surname: _____ Position: _____

Company/Institute: _____ Service, Lab: _____

Adress: _____

Postcode: _____ Town: _____

Tel _____ Fax _____ Email: _____

