

FT-BB7493

NHS activated CY_{anine} agents for fluorescent labeling of biomolecules via their amine groups

Introduction

A variety of **Cyanine dyes** has been used to label proteins, nucleic acids and other biomolecules for fluorescence techniques (imaging, biochemical analysis). They replace advantageously the conventional fluorochromes such as Fluorescein(FITC) and rhodamines (TRITC, RRX).

CY_{anine}3 : An orange-to-red fluorescence-emitting dye. Excitation/emission maximum 555/567 nm

Cy_{anine}3 is one of the most broadly used fluorophores which can be detected by various fluorometers, imagers, and microscopes. Due to inherently high extinction coefficient, this dye is also easily detected by naked eye on gels, and in solution. Can replace TRITC, AF555, DL550, Chromeo546.

See also alternative superior dye: <u>FluoProbes547H</u>.

CY_{anine}**3.5** can replace SulfoRhodamine 101.

See also alternative superior dye: FluoProbes594.

CY_{anine}5 : A far-red fluorescence-emitting dye. Excitation/emission maximum 647/657 nm

 $CY_{anine}5$ fluorophore has become an incredibly popular label in life science research and diagnostics. Fluorophore emission has maximum in red region, where many CCD detectors have maximum sensitivity, and biological objects have low background. Dye color is very intense, therefore quantity as small as 1 nanomol can be detected in gel electrophores by naked eye. Can replace far red red fluorescent dyes.

See also alternative superior dye: FluoProbes647H

 CY_{anine} 5.5 : A near-infrared (IR) fluorescence-emitting dye. Excitation/emission maximum 678/694 nm can replace near infrared fluorescent dyes.

See also alternative superior dye: FluoProbes682.

LIFE SCIENCES

 CY_{anine} ?: A near-IR fluor that is invisible to the naked eye. Excitation/emission maximum 750/776 nm It is used in in vivo imaging applications. See also alternative superior dye: <u>FluoProbes752</u>.

CY_{anine}**7.5** is a near infrared red fluorophores used for *in vivo* imaging applications. See also alternative superior dye: <u>FluoProbes800</u>.

CY_{anine}**NHS esters** are reactive dyes for the labeling of amino-groups typically found in peptides, proteins, and some derivates such as aminoallyl-oligonucleotides. The reaction is carried out at physiological pH.

Sulfo NHS – CY_{anine} are water-soluble derivatives of the CY_{anine} dyes, allowing more convenient use, and to achieve higher coupling ratio and coupled ratio. DiSulfonated forms are the most classic, but some mono-, tri- and quadri-sulfonated forms are available as well, for even lower or higher hydrosolubility.

PEG-CY_{anine} derivatives are also water soluble, but they also confer hydrophylicity to the conjugate via the PEG spacer.

bis-NHS derivatives are available for higher labeling or bi-functional conjugate applications

The triethylammonium (TEA) salts are more soluble in DMSO and DMF than the corresponding potassium salts.



FT-BB7493

Products Description

The table below gives main physical and fluorescence characteristics of the activated dyes

Product name	MW	λ	mol. abs.	Comment, structure
cat.number/qty*	g∙mol ⁻¹	abs./em.		
CV 2	(+added MW)	nm	M ⁻¹ cm ⁻¹	Penlace Eluorascoin(EITC)
CY_{anine}	550.6	489 / 506	150,000	
$\begin{array}{c} \mathbf{U} \\ $	550.0	1097 500	$OY \cdot 0.12$	
FF-L V 2550, 5111101			Q 11 0112	
low solubility in water				
CY _{anine} 3				Replace Cy 3, AF555, DL550, Chromeo546
CY _{anine} 3 – NHS	590.15	555 / 570	150 000	
^(L) FP-BB7493, 1mg	(+474.2)		QY: 0.31	
FP-1H4680, 1nmol				BF4
Soluble in organic solvents (DMSO, DMF,				
dichloromethane); low solubility in water				∑ _N ^O
Sealf- CV 2 his MUS				harra 2 aming reactive groups (2 NHS) and Sulfa group
SUIIO-C Y anine 3 DISINHS				bears 2 annue reactive groups (2 1013) and Suno group
DiSulfo-CYaning3 – NHS	735.80	548 / 567	162 000	*NarO ₃ S SO ₃ -
FP-1B7980 1mg	(598.2)		QY: 0.1	
FP-1B7981 5mg				T T
Dark red nowder				
Very high solubility in water, good in polar				
organic solvents (DMF, DMSO)	820.02	555/565		(incl. 1 othyl)
DISUITO-CY anine $3 - NHS$	629.05 (TEA salt)	333/303		
FP-LQU980, Img	(1211 5410)			
				~ , î »>
CY _{anine} 3 – PEG – NHS		649/665	"	0
2000DaFP-Inquire, 1mg	~2500			CVA NH-PEG-CH2-C-O-N
3400DaFP-Inquire, 1mg	~3900			
5000DaFP-Inquire, 1mg	~5500			0 ő"
CY _{anine} 3.5	600.27	501 /604	116.000	
$CY_{anine} 3.5 - NHS$	090.27	391/004	$OY \cdot 0.35$	
Pr-PZ0/30, IIIIg			Q 1. 0.00	$\sim \sqrt{N_{+}}$
Soluble in organic solvents (DMF, DMSO,				BF4
dichloromethane), insoluble in water				ko °
				Ô-N
TetraSulfo_CV . 35_NHS	1291.66	581 / 596	"	
FP-LOV060 1mg	(TEA salt)	5017 570		
				O VILLO O O VILLO
				(incl. 1 ethyl)



FT-BB7493				
Product name	MW	λ	mol. abs.	Comment, structure
cat.number/qty*	g∙mol⁻¹	abs./em.		
	(+added MW)	nm	M ⁻¹ cm ⁻¹	
CY _{anine} 5				
CV	616.19	646 / 662	250 000	
FP-BB2070 1mg	(465.68)		OY: 0.2	
FP_BB2071_5mg			X - · · · · -	
FP_BB2074_25mg				
11-DD2074, 25mg				
MonoSulfo-CV . 5 - NHS	659 79	642 / 660	180 000	0
TD 1C2520 5	0	0.27000	OY: 0.26	
FP-1G2520, 5mg			X I . 0.20	
FP-1G2521, 25mg				
Good solubility in polar (DMSO, DMF) and				
chlorinated (DCM, chloroform) organic				L o. L
solvents, low solubility in water				
CV		645 / 660		0
minimal dve				
^(L) FP-SJH900, 5nmol				
DiSulfo CV 5 NUS	777 95	646 / 662	271.000	K0 ₃ S S0 ₃ `
$DISUII0-C I_{anine}S = IVIIS$	111.55	0407 002	OY: 0.28	
FP-100510, 1mg			Q 11 0.20	N ~ ~ ~ N+
				Ô-N
DiSulfo-CVarian 5 - NHS.	855.07	649/665	"	
TEA solt	(TEA salt)			
EP LOU000 1mg				
11-LQ0990, 111g				
				° • (incl. 1 ethyl)
TriSulfo-CY _{anine} 5 – NHS	1050.35	649 / 670	250 000	o ₃ s So ₃ H
FP_111951_5mg	(TEA sait)			
11 111991, Sing	847.97			$\langle \rangle$
	(proto-nated)			2 x Et ₃ N ^{SO₃H}
				~ N ~
D'SIG- CV 5 L:-NUIS	1038 23	640/665	"	
DISUIIO-C Y anine 5 – DISINHS	(TEA salt)	049/003		
FP-LQV120, Img				l on the second
CY _{anine} 5 – PEG – NHS		649/665	"	
2000DaFP-1I2220, 1mg	~2500			N-O-C BEG CH CH NH-C Cy5
3400DaFP-1I2230, 1mg	~3900			N=0=C=PEG=Ch ₂ Ch ₂ =NH=
5000DaFP-1I2240, 1mg	~5500			
CY _{anine} 5.5				
CY _{anine} 5.5 – NHS	716.33	674 / 707	209 000	
EP-E78750 1mg			QY: 0.2	
11 120730, 1mg				N-Correction
Soluble in organic solvents (DMSO, DMF, dichloromethane): low solubility in water				≪ · · · · · · · · · · · · · · · · · · ·
				↓ p° µ°
				° N Å
		1		0





FT-BB7493				
Product name	MW	λ	mol. abs.	Comment, structure
cat.number/qty*	g∙mol ⁻¹	abs./em.		
1 2	(+added MW)	nm	M ⁻¹ cm ⁻¹	
Tetra-Sulfo-CV	1014.14	678 / 694	190 000	SO3H \$O3H
FP_I OV160 1mg	10	0,0.1		
11-LQ V 100, 1mg			CF ₂₈₀ : 0.18	O3S SO3H
Blue powder			CF ₂₆₀ : 0.05	
Soluble in water, DMF, DMSO				$\langle \rangle$
				Cho
Tetra-Sulfo-CY_min_5.5-NHS	1114.37	675 / 694	195 000	
FP-0B8400 1mg				ко,s
11 oboroo, ing			CF ₂₈₀ : 0.11	
Blue powder			CF ₂₆₀ : 0. 09	NN+
Soluble in water, DMF, DMSO				
Totus Sulfa CV 55	1500.86	678/701	Inquire	
Tetrasuno-Cranine 5.5 –	(TEA salt)	070/701	inquire	
bisNHS	(1211 5410)			Solution of the second
FP-LQV130, 1mg				o o
				$\langle \rangle$
CY _{anine} 7				
CY _{anine} 7 – NHS	682.29			Soluble in organic solvents (DMSO, DMF, dichloromethane),
FP-JV9770, 1mg				low solubility in water
DiSulfo-CV7 – NHS	827.94	740 / 773	240 600	*NarO ₃ SSO ₃ -
				$\bigcirc +$ $+ \bigcirc$
FP-1B8001, 5mg				
Dark green powder				
Soluble in water, DMF, DMSO				
				J.N.K.
MonoSulfo-CV · 7 – NHS	881.11		Inquire	,CH3 0, Et,
$ED \perp OV140 1mc$	(TEA salt)		1	
FF-LQ V 140, 111g	· · · ·			
				o)
				o´ Ethy
			. .	l substituent and spacer C7
DiSulfo-CY _{anine} 7 – bisNHS	1064.27	749 / 776	Inquire	
FP-LQV150, 1mg	(TEA salt)			o et
_				
				Y Y
				reactive groups (2 NHS), 2 sulfo, Ethyl, and spacer C7
CY _{anine} 7.5				
CY _{anine} 7.5 – NHS	782.41	788 / 808	223 000	
FP-JO3060, 1mg	[631.9]			BF4 BF4
Soluble in organic colvents (DMSO, DME				N
dichloromethane): low solubility in water				
				¢°°
				0-N
				0

LIFE SCIENCES

Section 211 bis, avenue JF Kennedy BP 1140 - 03100 Montluçon Hotline +33 4 70 03 73 06 • interbiotech@interchim.com



FT-BB7493					
Product name	MW	λ	mol. abs.	Comment, structure	
cat.number/qty*	g∙mol ⁻¹	abs./em.			
	(+added MW)	nm	M ⁻¹ cm ⁻¹		
Tetra-Sulfo-CY _{anine} 7.5 – NHS FP-AWHGX0, 1mg Soluble in water, DMF, DMSO	1180.47 [950.2]	778 / 797	222 000		

*: search <u>online</u> or <u>ask</u> for other quantities (5mg to +100g).

The different packaging have the same catalog number as the standard size above, only the last digit is changed.

 Storage:
 -20°C, protected from light (+4°C possible for short term) (M)

 Avoid prolonged exposure to light. Desiccate.

 Stable for 24 months after receipt at -20°C in the dark. Transportation: at room temperature for up to 3 weeks.

Directions for use

Protocol for Labeling of Amino-Biomolecules

Introduction:

NHS (N-HydroxySuccinimide) esters and other activated esters (sulfo-NHS, sulfotetrafluorophenyl - STP) are highly reactive compounds suitable for the modification of amino groups. NHS is most common type of activated esters.

Usual modifications are fluorescent labels, fluorescence quenchers, and other reporter groups. Alkyne and azido group can be attached using activated esters to adapt biomolecules to Click Chemistry.

Since amino groups are nearly always contained in proteins and peptides, modification of these biopolymers is especially common. Other examples are amino-oligonucleotides, amino-modified DNA, and amino-containing sugars.

The reaction of NHS esters with amines is strongly pH-dependent: at low pH, the amino group is protonated, and no modification takes place. At higher-than-optimal pH, hydrolysis of NHS ester is quick, and modification yield diminishes. Optimal pH value for modification is 8.3-8.5.

Water is most common solvent for the labeling. If NHS ester is poorly soluble, it can be added as a solution in DMSO or DMF to a solution of protein in water, adjusted to pH 8.3-8.5. Note that DMF must not contain amines.

We recommend using the following general protocol for the labeling of biomolecules with NHS esters. See <u>related products</u> for auxiliary reagents.

Protocol:

1. Calculate required amount of NHS ester:

NHS_ester_weight [mg] =

8 × amino_compound_weight [mg] × NHS_ester_molar_weight [Da] / amino_compound_molar_weight [Da].

- 8 is molar excess of NHS ester. It is experimental value for mono-labeling, suitable for many common proteins and peptides. However, in some cases using less or more NHS ester is required. It depends on protein structure.
- For example, to label 3 mg of insulin (molar weight 69300 Dalton) with CY_{anine}5 NHS ester (molar weight 616 Dalton), and obtain maximum yield of mono-labeled product, you should calculate $10 \times 3 \text{ mg} \times 616\text{Da} / 69300\text{Da} = 0.266 \text{ mg of CY}_{anine}5$ dye NHS ester.
- The molar weights of NHS esters are displayed in the product description at first page (note that molar weights of reagents produced by other vendors may vary).
- Determine volume of reaction mixture. The labeling can be performed on any scale from nanomols to dozens of grams. When the scale is low, use minimal volume (10-20 uL). Higher concentrations (1-10 mg of aminobiomolecule per mL of mixture) are optimal.

Kenned

LIFE SCIENCES





- 3. Dissolve NHS ester in 1/10 reaction volume of DMF or DMSO. Amine-free DMF is preferred solvent. After the reaction, NHS ester can be stored in solution for 1-2 months at -20°C.
- 4. Dissolve biomolecule in 9/10 reaction volume of buffer with pH 8.3-8.5. 0.1 M Sodium bicarbonate solution has appropriate pH. Other alternatives are 0.1 M Tris buffer (although Tris has amino group, it is hindered and does not react with NHS esters), or 0.1 M phosphate buffer. Note pH is most important thing.
- When doing large-scale labeling (hundreds of milligrams of NHS ester), note that the mixture tends to acidify with time because of hydrolysis of NHS ester. Monitor pH, or use more concentrated buffer then.
- 5. Add NHS ester solution to the solution of biomolecule, and vortex well. Keep on ice overnight, or at room temperature during at least 4 hours.
- 6. Purify the conjugate using appropriate method: gel-filtration for macromolecules is most universal. Precipitation and chromatography is another alternative. Organic impurities (such as N-hydroxysuccinimide, NHS ester, acid produced by hydrolysis) are almost always easily separated.

Related / associated products and documents

* <u>CY_{anine} Dyes</u> (Cy3, Cy5, Cy5.5, Cy7...) functionalized by NHS (<u>FT-BB7493</u>), Maleimide (<u>FT-JO6660</u>), Azide (<u>FT-HO7250</u>), Alkyne (<u>FT-1A6320</u>), DBCO (<u>FT-DQP790</u>), Hydrazide (<u>FT-LQV050</u>), DBCO (<u>FT-DQP790</u>: CycloAlkynes, for strain-promoted Click reactions), Amino group (<u>FT-CY3AM0</u>), Carboxyl group (<u>FT-CY3CA0</u>). 3Dye 2D DI GE (CY2/CY3/CY5) labeling kit (FT-EV0870)

- * Superior FluoProbes fluorescent dyes
- activated by -NHS (list), i.e. FP488-NHS (FT-BA6800) activated by -Azide, i.e. FP488-Azide (FT-YE4970)
- activated by -Maleimide (list), i.e. FP488-MAL (FT-BA6810)
- * Classic dyes such as FAM, CR110, JOE TAMRA, and ROX, functionalized by –Alkyne or DBCO (<u>FT-DQP790</u>), or by -Azide (<u>FT-EV0920</u>).
- * Fluorescently labeled ligands:
- Labeled secondary antibodies
- Labeled lectins, i.e. ConA-CY_{anine}3 (#WT868/<u>FT-47496A</u>) Labeled tags, i.e. CY_{anine}3-polylysine (<u>FT-WT8550</u>)
- * Other labeling/conjugation chemistries: Click Chemistry reagents

Ordering information

Catalog size quantities and prices may be found at <u>www.interchim.com/</u> Please inquire for higher quantities (availability, shipment conditions). For any information, please ask : FluoProbes[®] / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes[®] are sold for research use only, and are not intended for food, drug, household, or cosmetic use. FluoProbes[®] is not liable for any damage resulting from handling or contact with this product.