

# VisPRO 5 minutes Protein Stain Kit

**Did you have the following problem in your experiment ?**

- **SDS-PAGE staining**
  - **Q1 A long time to staining: SDS-PAGE stained spend 2 hours to over night. All of stain can't immediately get staining result.**
  - **Q2 Low sensitivity: the Commassie blue stain, the sensitivity might be to 50ng. We hope to get more sensitivity staining?**

## Did you have the following problem in your experiment ?

- SDS-PAGE staining
  - **Q3 the complex protocol in staining preparation : Silver stain**, it's high sensitivity staining. But it has complex step, stained for a long time, and not easy to operating. It is easy to be fail and re-do it. We always hope to has high sensitivity and easy staining method.
  - **Q4 MASS not compatible** : To identify protein , we always use MALDI-TOF or LS/MS/MS, but commonly used stain can interfere with MASS result. Even mass-compatible silver staining also reduce the success rate of identification.

Did you have the following problem in your experiment ?

- **SDS-PAGE staining**
  - **Q5 Can't estimate the pattern of electrophoresis before transferring:** Before transferring, we can't check the pattern of electrophoresis. If the pattern is bad, we can't check it. And it is always waste time. We need a new staining method can make SDS-PAGE transferring.
- **Q6 Recovery activity protein in gel:** When recovery protein, we need to use CBR stain or other stain to check protein's position. But all commonly used stain cause protein loss activity. So we need a new stain can recovery protein and keep protein activity.

Did you have the following problem in your experiment ?

- SDS-PAGE staining
  - **Q7 Expensive.** To Consider the cost of reagents and the cost of time, VisPRO staining can obtain the higher C/P value.
  - **Q8** Using SYPRO Ruby stain, spotting protein must use UV to excite gel, it is a great burden on the eye. When using Cydye, It is not clear to saw the spot and inconvenient to select protein spot. Is there any suggestion in safety and convenience?

# The fastest stain

**Q1 A long time to staining: SDS-PAGE stained spend 2 hours to over night. All of stain can't immediately get staining result.**

**VisPRO Stain is the fastest stain . it just spends 15 min obtaining the result.**

Protocol for VisPRO 5 minutes Protein Stain Kit

1. add Solution 1 and shack 5 mins
2. Discard Solution 1, quickly wash it
3. add Solution 2 and shack 30 sec
4. discard Solution 2 and quickly wash twice



**VISUAL  
PROTEIN**

# Comparison staining time

## VisPRO Stain costs shorten time

Methods	VisPRO Stain	Sypro™ Ruby	Silver Stain	CBR Stain
Preparation of solutions	0 min	5 min	20 min	0 min
Fixation	0 min	1 hr	1hr~overnight	0 min
Staining procedure	5~15 min	30 min~ o/n	2.5~4 hr	30 min~ o/n
Destaining procedure	0 min	1 hour	0 min	2 hr~ o/n
Total	5 ~15 min	2.5~o/n	3.5~0/n	2.5 hr~ o/n



### ⑩ VisPRO stain


Only spends 5~15 min.

The fastest stain in the market

# Excellent Sensitivity

- Q2 **Low sensitivity:** the Commassie blue stain, the sensitivity might be to 50ng. We hope to get more sensitivity staining?


**VisPRO detects protein levels accurate to 1 ng as sensitivity as SYPRO Ruby and Sliver stain.**





# Comparison of staining sensitivity-VisPRO Stain is the best


Comparison of staining sensitivity :  
even as little as 1ng of protein can be  
detected



# Comparison of staining results :VisPRO Stain VS silver stain

VisPRO 5 minutes protein stain kit

硝酸銀染色法



# Easy to use VisPRO Stain

–Q3 the complex protocol in staining preparation : Silver stain, it's high sensitivity staining. But it has complex step, stained for a long time, and not easy to operating. It is easy to be fail and re-do it. We always hope to has high sensitivity and easy staining method.

**VisPRO Stain, Two reagents, Four steps, 5-15 mins**

Protocol for VisPRO 5 minutes Protein Stain Kit

1. add Solution 1 and shack 5 mins
2. Discard Solution 1, quickly wash it
3. add Solution 2 and shack 30 sec
4. discard Solution 2 and quickly wash twice



**VISUAL  
PROTEIN**

# The Best MASS Choice

**Q4 MASS not compatible** : To identify protein , we always use MALDI-TOF or LS/MS/MS, but commonly used stain can interfere with MASS result. Even mass-compatible silver staining also reduce the success rate of identification.


## **VisPRO Stain is a reverse staining method**

- 1.** Reagents don't combine with proteins, and the molecular weight would be change.
- 2.** Discard reagent easily, non-interference in the MASS result
- 3.** Recommended by authoritative papers, rising the success rate in identification ( ie. Journal of Proteome Research Vol. 8, No. 10, 2009 4393 )
- 4.** Recommended by many professors.


# What's Reverse Staining ?

(VisPRO 5 minutes Protein Stain Kit)

PostiveStain



Reverse Stain





Macromolecules form Space obstacles

Small insoluble white powder left Zinc-imidazole complex

**Reverse Stain: Stain non-protein part in gel and enhance the target protein by background**

# Journal of Proteome Research, Aug. 2009. Research Profile



Refer VisPRO Stain—obtaining the best result of MASS



# Papers(I)

## VisPRO 5 minutes Protein Stain Kit in Papers

- Bin Huang, Shih Chung Chen, and Danny Ling Wang. **Cardiovascular Research** 2009; 83 : 536-546. Shear flow increases S-nitrosylation of proteins in endothelial cells.
- Ching-Yu Lin, Vinchi Wang, Hao-Ai Shui, Rong-Huay Juang, Ai-Ling Hour, Pei-Sing Chen, Hui-Ming Huang, Szu-Yu Wu, Jen-Chieh Lee, Tzung-Lin Tsai and Han-Min Chen. **Proteomics** 2009; 9:696-709. A comprehensive evaluation of imidazole-zinc reverse stain for current proteomic researches.
- Bin Huang, Chung Ling Liao, Ya Ping Lin, Shih Chung Chen and Danny Ling Wang. **Journal of Proteome Research** 2009; 8: 4835–4843. S-nitrosoproteome in endothelial cells revealed by a modified biotin switch approach coupled with western blot-based two-dimensional gel electrophoresis

# Papers (II)

## VisPRO 5 minutes Protein Stain Kit in theses

- Yu-Jen Chen, Chih-Chia Yu, Shui-Tein Chen, Tai-Yuan Chen and Hui-Fen Liao. **Proteomics Clin. Appl.** 2009; 3: 563–573. Functional regulation and proteomic characterization of human natural killer cells through recombinant human granulocyte-colony stimulating factor treatment.
- Ching-Yu Lin, Hui-Ming Huang, and Han-Min Chen, **Biotechniques** 2006; 41:560-564. Use of backlit light plate to enhance visualization of imidazole-zinc reverse stained gels.
- Hui-Chung Wu, Chien-Chang Yen, Wen-Huei Tsui, Han-Min Chen. **Analytical biochemistry**, 2009; 03: 2693-2697. A red line not to cross: Evaluating the limitation and properness of gel image tuning procedures.



# Papers(I)


- **VisPRO 5 minutes Protein Stain Kit in theses**
- Yi-Jen Liao, Lisa Wen, Jei-Fu Shawc, Chi-Tsai Lin. **Journal of Biotechnology**, 2007 ; 131:84–91. A highly stable cambialistic-superoxide dismutase from *Antrodia camphorata*: *Exprcamphorata: xpression* in yeast and enzyme properties.
- Shih Chung Chen, Bin Huang, Yu Chi Liu, Kou Gi Shyu, Pen Y. Lin and Danny Ling Wang. **Biochemical and Biophysical Research Communications** 2008 ; 377:1274-1278. Acute hypoxia enhances proteins' S-nitrosylation in endothelial cells.
- Yasuhiko Kizuka, Kyoko Kobayashi, Shinako Kakuda, Yukari Nakajima, Satsuki Itoh, Nana Kawasaki, and Shogo Oka. **Glycobiology** 2008; 18(4): 331–338. Laminin-1 is a novel carrier glycoprotein for the nons

# Reference: VisPRO 5 minutes Protein Stain Kit

{The Instrument Center, Institute of Biomedical Sciences, Academia SINICA}


**July, 2007**

E.Coli-K12 lysate




12% SDS-PAGE  
Stain: Sypro Ruby 3 hours

E.Coli-K12 lysate




12% SDS-PAGE  
Stain: VisPro 5 ( < 10 mins )

# 2007.07.09 Sypro Ruby / VisPro 5 gel staining



## Sypro Ruby

Run 12% Acrylamide / bis gel  
Fixation: 50% Methanol, 7% Acetic acid. 30mins/time, 2 times  
Stain: 3 hours  
Washing: 10% Methanol, 7% Acetic acid. 30mins, once  
Rinse: ddH<sub>2</sub>O 5mins/time, 2 times  
Keep in ddH<sub>2</sub>O 2 nights, then scan the gel by Typhoon9410  
Excitation: ~532nm  
Emission: 610nm




## VisPro 5

使用染完Sypro Ruby之gel作染色

1. Balance: Solution I 5 mins
2. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
3. Stain: Solution II, shake the gel gently ( less than 30 sec )
4. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
5. Keep in ddH<sub>2</sub>O , then use scanner to store gel image.


# 2007.07.013 Sypro Ruby / VisPro 5 gel staining



## Sypro Ruby

- 1.Run 10% Acrylamide / bis gel
- 2.Sample: BSA (MW.~68KDa)
- 3.Fixation: 50% Methanol, 7% Acetic acid. 30mins/time, 2 times
- 4.Stain: 3 hours
- 5.Washing: 10% Methanol, 7% Acetic acid. 30mins, once
- 6.Rinse: ddH<sub>2</sub>O 5mns/time, 2 times
- 7.Keep in ddH<sub>2</sub>O over night, then scan the gel by Typhoon9410  
Excitation: 488nm  
Emission: 610nm

\* Use different wavelength to scan gel ( 457nm, 488nm, 532nm)  
Obtaining the similar result, but the signal(532nm) is stronger after scanning



## VisPro 5 min protein stain kit

The gel with Sypro Ruby

1. Balance: Solution I 5 mins
2. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
3. Stain: Solution II, shake the gel gently ( less than 30 sec)
4. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
5. Keep in ddH<sub>2</sub>O , then use scanner to store gel image.


**Figure a. reflect mode scan**

**Figure b. Penetration mode scan**

\* observing with naked eyes is better than scanner

# 2007.07.30 Sypro Ruby / VisPro 5 gel staining ( test BSA degradation)

0.78 1.56 3.12 6.25 12.5 25 50 100 ( ng )




## Sypro Ruby

- 1.Run 10% Acryl amide / bis gel
- 2.Sample: BSA (MW.~68KDa)
- 3.Fixation: 50% Methanol, 7% Acetic acid. 30mins/time, 2 times
- 4.Stain: over night (~20hr) in RT
- 5.Washing: 10% Methanol, 7% Acetic acid. 40mins, once
- 6.Rinse: ddH<sub>2</sub>O 5mns/time, 2 times
- 7.Keep in ddH<sub>2</sub>O, then scan the gel with Typhoon9410

Excitation: 532nm  
Emission: 610nm

0.78 1.56 3.12 6.25 12.5 25 50 100 ( ng )



## VisPro 5 min protein stain kit

The gel with Sypro Ruby(not the same gel as above)


1. Balance: Solution I 5 mins
2. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
3. Stain: Solution II, shake the gel gently ( less than 30 sec)
4. Rinse: Rinse the gel with ddH<sub>2</sub>O rapidly ( 1 time )
5. Keep in ddH<sub>2</sub>O , then use scanner to store gel image.

# In-gel digestion à PMF detected by MALDI-TOF

stain: VisPro 5

100ng BSA

Mass fit research result: 41%




# In-gel digestion à PMF detected by MALDI-TOF

stain: VisPro 5

75ng BSA

Mass fit research result: 44%






# In-gel digestion à PMF detected by MALDI-TOF

stain: VisPro 5

25ng BSA

Mass fit research result: 27%



## Recovery activity protein in gel

**Q6 Recovery activity protein in gel:** When recovery protein, we need to use CBR stain or other stain to check protein's position. But all commonly used stain cause protein loss activity. So we need a new stain can recovery protein and keep protein activity.


**Comparison of other stain, VisPRO stain don't contain any acid compounds. It can't denature protein and make protein fixing in gel. Therefore, it can recovery the active enzyme from gel.**

**VisPRO Stain's compound doesn't combine with protein; therefore, VisPRO is the compatibility different applications.**


# What's Reverse Staining ?

(VisPRO 5 minutes Protein Stain Kit)

PostiveStain




Reverse Stain



**Reverse Stain: Stain non-protein part in gel and enhance the target protein by background**

# The compatibility to most applications with VisPRO Stain



1. Application other stain
2. MS identification
3. Blotting in WB
4. Electrophoresis elution activity staining and detection
- 5.


# the method of the highest C/P stain

- **Q7 Expensive.** To Consider the cost of reagents and the cost of time, VisPRO staining can obtain the higher C/P value.

## The cost of VisPRO Stain

To Consider the cost of reagents and the cost of time, VisPRO staining can obtain the higher C/P value.

# Comparison of staining methods



Methods	VisPRO Stain	Sypro™ Ruby	Silver Stain	CBR Stain
Total time	★ 5 ~15 min	2.5~o/n	3.5~o/n	2.5 hr~ o/n
Imaging instrument	★ visible light (cheap)	fluorescent (expensive)	★ visible light (cheap)	★ visible light (cheap)
Irritant or toxic chemicals	★ no	Acetic acid, Methanol UV	Acetic acid, Silver Nitrate, Glutaraldehyde	Acetic acid, Methanol
sensitivity	★ <1 ng	1 ng	★ <1 ng	50 ng
Quantitative range	1-200 ng	★ 1-1000 ng	1-80ng	50-1000 ng
Capability to downstream applications	★ Yes Compatibility with MASS and downstream applications	★ Yes Compatibility with MASS	Limited protein cross link Low Compatibility with MASS	Limited Low compatibility with MASS

# Safe and Easy to use

**Q8 Using SYPRO Ruby, it is hurt for eyes while watching with UV.**

**It is invisible completely by cy dye.**

**It is inconvenient to select protein spot, either.**


**Is there any suggestion in safety and convenience?**

**Spending 5 minute only**

**Detect protein levels accurate to 1 ng with high sensitivity**

**None reaction with protein and MASS**

**SYPRO Ruby and Cy dye are complemented by VisPRO**



[ 2D application ]

## VisPRO Stain combine with other Fluorescent stain

the visualized result with eyes

take the point correctly and position completely


without the risk of UV

compatibility with MASS and downstream application



## The result and observation: Reverse Staining (VisPRO 5 minutes Protein Stain Kit)

- ⑩ After staining :
  - ⑩ Gel stained milk-white color
  - ⑩ Protein sample is colorless
  
- ⑩ Observing under the dark background
  - ⑩ Clearly to observe the colorless sample band changes to black.



Staining result for 2-DE gels  
by VisPRO 5 minutes  
Protein Stain Kit

# VisPRO 5 minutes Protein Stain kit

## The **three high** features

### ⑩ **High Speed**

⑩ Only 5~15 min staining process

### ⑩ **High Sensitivity**

⑩ Detect protein levels accurate to 1ng

### ⑩ **High Compatibility**

⑩ Apply to mass spectrometry analysis

⑩ Restain by other methods

⑩ Recycle proteins

⑩ Gels stained with VisPRO can be re-stained by other protein staining methods or WB.

# The features of VisPRO 5 minutes Protein Stain Kit

## ⑩ Safety

- ⑩ Non-toxic
- ⑩ non- highly toxic heavy metal
- ⑩ Visualized with naked eye ( without UV )


## ⑩ Convenience

- ⑩ Stored at room temperature for two year before opening.
- ⑩ Use it without any preparation

## ⑩ Economy

- ⑩ available price
- ⑩ Unlimited principle
- ⑩ Cheap image and accessories

# Reverse Stain's (VisPRO's) Capability to downstream applications—MASS analysis



The MALDI-TOF spectrum of spots from zinc reverse staining.

The zinc reverse staining spot of phosphorylase b was cut and subjected to mass sample preparation. The observed ionized fragments of rabbit phosphorylase b were annotated.

**⑩ Reverse stain** is perfectly suitable with mass spectrometry analysis


# Reverse Stain's (VisPRO's) Capability to downstream applications—MASS analysis

In-gel digestion à PMF detected by MALDI-TOF

stain: VisPro 5

75ng BSA


Mass fit research result: 44%




Reference : VisPRO is tested by the Instrument Center, Institute of Biomedical Sciences, Academia SINICA, in July 2007

**VISUAL  
PROTEIN**


# Comparison of linear range: Reverse Stain (VisPRO) VS others



VisPRO vs. Sypro Ruby




VisPRO vs. Silver Stain



VisPRO vs. CBR Stain

# Comparison of staining methods



Methods	Reverse Stain (VisPRO)	Sypro™ Ruby	Silver Stain	CBR Stain
Total time	★ 5 ~15 min	2.5~o/n	3.5~o/n	2.5 hr~ o/n
Imaging instrument	★ visible light (cheap)	fluorescent (expensive)	★ visible light (cheap)	★ visible light (cheap)
Irritant or toxic chemicals	★ no	Acetic acid, Methanol UV	Acetic acid, Silver Nitrate, Glutaraldehyde	Acetic acid, Methanol
sensitivity	★ <1 ng	1 ng	★ <1 ng	50 ng
Quantitative range	1-200 ng	★ 1-1000 ng	1-80ng	50-1000 ng
Compatibility to downstream application (e.g. WB & MASS)	★ Yes High compatibility with MASS and other	★ Yes Compatibility with MASS	Limited protein cross link low compatibility with MASS	Limited Low compatibility with MASS

# Why is the reverse stain not popular before 2000?

## ⑩ Image storage technology

- ⑩ Zinc-imidazole reverse stain was created in 1990
- ⑩ In 1990, adhesive save or silver photo preservation
- ⑩ After the year of 2005, Digital photo or reflective scanner
- ⑩ After the year of 2005, Visual Protein provides an see-through scanner, rising the quality and sensitivity on image.



## ⑩ Equipment

- ⑩ In 2005, Visual Protein provides **black staining box and plate which can be easy visualized with naked eye**



## ⑩ the Change of needs for the laboratory

- ⑩ Enquiry for high excellent sensitivity
- ⑩ Be **compatible with mass spectrometry**





# Protocol for VisPRO 5 minutes Protein Stain



## 1. Using Black Staining Box

2. Scan image  
See-through scanner  
EPSON V750  
PRO

### Gel in the Staining Box



Shake 5 min



Wash Gel with ddH<sub>2</sub>O



Discard Solution(2)  
Wash Gel with ddH<sub>2</sub>O



Get image



Quick Staining Protocol for  
VisPRO 5 minutes Protein Stain Kit



## VisPRO 5 minutes Protein Stain Kit



## 3. View and select spot by Gel Lighting Plate

**VISUAL**  
**PROTEIN**