



# 1-Sec Real-time Amplification





## Fast Probe Multiplex

# RealHelix™ 1-Sec qPCR Kit

- ❖ **Fast:** 30 min for 40 cycles (1-sec denature, 1-sec anneal / extension)
- ❖ **Multiplex:** Up to 5-plex probe qPCR in a reaction
- ❖ **High sensitivity:** Reliable detection as low as 10 copies of target
- ❖ **Robust hot-start enzyme:** Novel Antibody inhibited *Taq* polymerase



Comparison of amplification plots of a general qPCR enzyme (RealHelix™ qPCR kit, probe type) and 1-sec qPCR kit using either a standard cycling protocol or a fast cycling protocol. Multiplex Probe qPCR reactions were done with a purified human genomic DNA over range of 0.01 ng to 10 ng per reaction. Each 20 µl reactions were operated on CFX96 real-time PCR detection system (BioRad).





Template (ng)	Nanohelix		Competitor A		Competitor B		Δ avg Ct Nanohelix - A	Δ avg Ct Nanohelix - B
10	25.45	25.49	27.39	27.63	29.50	29.75	- 2.04	- 4.16
	25.47	27.51	29.63					
1	29.08	29.17	31.00	31.06	33.92	34.25	- 1.09	- 4.96
	29.13	31.03	34.09					
0.1	32.24	32.40	34.91	35.08	40.98	41.10	- 2.68	- 8.72
	32.32	35.00	41.04					
0.01	35.44	35.46	39.42	40.17	N/A	N/A	- 4.35	N/A
	35.45	39.80	N/A					

GloB gene fragment was amplified from 10-fold serial dilutions of human genomic DNA (10 ng – 0.1 ng/20 µl reaction) using 1-sec qPCR Kit or other competitor's fast probe qPCR kits: Competitor A (Probe Fast qPCR kit) and Competitor B (Universal Probe supermix) in fast condition. Fluorescence data were obtained in a multiplex (3-plex) qPCR reactions.

# RealHelix™ 1-Sec qRT-PCR Kit

- ❖ **Fast:** 40 min for 40 cycles (1-sec denature, 1-sec anneal / extension)
- ❖ **Multiplex:** Up to 5-plex probe qRT-PCR in a reaction
- ❖ **High sensitivity:** Reliable detection as low as 10 copies of target
- ❖ **Robust hot-start enzyme:** Novel Antibody inhibited *Taq* polymerase



Comparison of amplification plots of a general qRT-PCR enzyme (RealHelix™ qRT-PCR kit) and 1-sec qRT-PCR kit using either a standard cycling protocol or a fast cycling protocol. Multiplex Probe qRT-PCR reactions were done with the purified and mixed RNA of Avian influenza type H5 and type H7 (tenfold serial dilutions over range of 1 fg – 1 ng per reaction). Each 20 µl reactions were operated on CFX96 real-time PCR detection system (BioRad).



## **NanoHelix Co., Ltd.**

(Tel : +82-42-867-9055, Fax : +82-42-867-9057)

43-15, Techno 5-ro, Yuseong-gu, Daejeon 34014, South Korea

[www.nanohelix.net](http://www.nanohelix.net)

[info@nanohelix.net](mailto:info@nanohelix.net)



## RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)

### Kit contents

RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)		
Cat. No.	SQPU2-P200 (200 rxns)	SQPU2-P500 (500 rxns)
1-sec 2x Premix [Probe] w/ UDG (V2)	1 ml x 2ea	1 ml x 5ea
ROX Dye (25 µM)	0.2 ml	0.5 ml
Instruction for Use	1 ea	1 ea

### Description

**RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)** is designed to perform a rapid real-time quantification of DNA samples using the fluorescent probe based detection. The UDG/dUTP system prevents the carryover contamination of PCR products from previous reactions.

The kit consists of convenient 2x premix and separate ROX reference dye. The 2x premix contains antibody-inhibited hot-start *Taq* polymerase (*Ab+Taq* polymerase), thermo-labile Uracil-DNA glycosylase, dUTP, dNTPs, buffers, Mg<sup>2+</sup> and a stabilizing agent. The hot-start PCR enzyme provides high specific amplification of target DNA and minimizes the side products such as primer dimers.

### Application

Quantification of target DNA sample  
by real time PCR

### Store -20 °C

ROX Dye should be stored under dark  
conditions.

### Quality control assay data

#### Functional analysis

**RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)** was evaluated by real-time PCR using a human genomic DNA, a set of human gene-specific primers, and a dual-labeled probe with ABI 7500 real-time PCR system (Applied Biosystems, USA) and CFX96™ real-time PCR detection system (Bio-Rad, USA).

Quality authorized by Yountaek Go



### NanoHelix Co., Ltd.

45-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net < [www.nanohelix.net](http://www.nanohelix.net) [www.nanohelix.co.kr](http://www.nanohelix.co.kr) >

## Protocol

### 1. Program a real-time PCR instrument according to the recommendations below.

#### <2-step cycling protocol>

If annealing temperature (AT) of primers used in real-time PCR is between 55°C and 60°C, thermo cycling can be performed using 2-step cycling protocol as follow.

Step	Condition		Cycle(s)
[Optional] UDG reaction*		20 ~ 25°C for 5 min	1
Enzyme activation		95°C for 2 ~ 5 min	1
PCR Amplification	Denaturation	95°C for 1 ~ 10 sec**	45
	Annealing & extension	<b>55 ~ 60°C</b> for 1 ~ 30 sec** <b>Collect the fluorescence data</b>	

\* The UDG reaction step is not essential. The UDG will efficiently remove carryover contaminant DNA during sample setup and cyler ramping.

\*\* The reaction time for each steps should be optimized on the applied thermocycler.

#### <3-step cycling protocol>

If annealing temperature (AT) of primers used in real-time PCR is under 58°C or above 62°C, thermo cycling can be performed using 3-step cycling protocol as follow.

Step	Condition		Cycle(s)
[Optional] UDG reaction*		20 ~ 25°C for 5 min	1
Enzyme activation		95°C for 2 ~ 5 min	1
PCR Amplification	Denaturation	95°C for 1 ~ 10 sec**	45
	Annealing	<b>AT<sup>1</sup>°C</b> for 1 ~ 20 sec**	
	Extension	72°C for 1 ~ 30 sec** <b>Collect the fluorescence data</b>	

\* The UDG reaction step is not essential. The UDG will efficiently remove carryover contaminant DNA during sample setup and cyler ramping.

\*\* The reaction time for each steps should be optimized on the applied thermocycler.

## 1-sec qPCR Kit [P] w/ UDG (Ver.2.0)

<sup>1)</sup>AT, the annealing temperature of primers used

$$\text{Annealing Temperature} = T_m - (4 \sim 6^\circ\text{C})$$

$$\text{Where, } T_m \text{ (Melting Temp.)} = [4^\circ\text{C} \times (\text{G} + \text{C})] + [2^\circ\text{C} \times (\text{A} + \text{T})]$$

### 2. Add following components for a single 20 $\mu\text{l}$ reaction volume.\*

Components	Volumes
DNA Template	X $\mu\text{l}$
1-sec 2x Premix [Probe] w/ UDG (V2)	10 $\mu\text{l}$
Forward primer (10 $\mu\text{M}$ )	0.5 ~ 1.0 $\mu\text{l}$
Reverse primer (10 $\mu\text{M}$ )	0.5 ~ 1.0 $\mu\text{l}$
Fluorescent Probe (5 $\mu\text{M}$ )	0.5 ~ 1.0 $\mu\text{l}$
ROX Dye (25 $\mu\text{M}$ )	Optional **
DEPC-treated water	Adjust to final 20 $\mu\text{l}$

\* For multiple reactions, prepare a master mix by adding the required volumes of each above components (except the DNA template) and dispense appropriate volumes into PCR tubes or plates.

\*\* Use the recommended amount or concentration of ROX Dye (Passive Reference) depending on the instrument.

### 3. Gently mix and immediately centrifuge the reaction mix.

### 4. Perform the Real-time PCR.

## Products

Cat. No.	Products	Size
<b>SQPU2-P200</b>	RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)	200 rxns
<b>SQPU2-P500</b>	RealHelix™ 1-sec qPCR Kit [Probe] w/ UDG (Ver.2.0)	500 rxns

## NanoHelix Co., Ltd.

45-15, Techno 5-ro, Yuseong-Gu, Daejeon, 34014, South Korea. TEL : 82-42-867-9055, FAX : 82-42-867-9057

E-mail : info@nanohelix.net < [www.nanohelix.net](http://www.nanohelix.net) [www.nanohelix.co.kr](http://www.nanohelix.co.kr) >