


Test comparison for D-Malic acid between Roche and Enzytec™ Liquid

Art. No. E8270

Overview

Item	Roche	Enzytec™ Liquid																												
Art. No.	11 215 558 035	E8270																												
Kit size	3 x 11 tests	50 tests																												
Test principle	D-Malic acid is oxidized to oxaloacetate by D-malate-dehydrogenase (D-MDH). Oxaloacetate is immediately split by the same enzyme to pyruvate and CO ₂ : D-Malate + NAD ⁺ — D-MDH —> pyruvate + CO ₂ + NADH + H ⁺																													
Calculation of results	The method is an end-point reaction where results are calculated using the Lambert-Beer formula (excel table available on request)																													
Registration	This method is registered by the AIJN/IFU (fruit juices), DIN (Germany), EN (European norm), OIV (wine)																													
Reagents	<ul style="list-style-type: none">• Vial 1 = Buffer, glutamate (liquid)• Vial 2 = NAD (lyophilizate)• Vial 3 = D-MDH (lyophilizate)	<ul style="list-style-type: none">• R1 = Buffer• R2 = NAD, D-MDH																												
Stability after reconstitution	<ul style="list-style-type: none">• Vial 2 = 3 weeks• Vial 3 = 5 days	All reagents are liquid, ready-to-use and stable after opening																												
QC samples	QC included (vial 4)	QC must be purchased separately																												
Procedure	<table><tr><th>Steps</th><th>Samples (or blank)</th></tr><tr><td>Vial 1</td><td>1.000 mL</td></tr><tr><td>Vial 2</td><td>0.200 mL</td></tr><tr><td>Sample (or H₂O)</td><td>0.100 mL</td></tr><tr><td>Water</td><td>1.700 mL</td></tr><tr><td colspan="2">Incub. 6 min, read A₁</td></tr><tr><td>Vial 3</td><td>0.050 mL</td></tr><tr><td colspan="2">Incub. 20 min, read A₂</td></tr></table>	Steps	Samples (or blank)	Vial 1	1.000 mL	Vial 2	0.200 mL	Sample (or H ₂ O)	0.100 mL	Water	1.700 mL	Incub. 6 min, read A ₁		Vial 3	0.050 mL	Incub. 20 min, read A ₂		<table><tr><th>Steps</th><th>Samples (or blank)</th></tr><tr><td>R1</td><td>2.000 mL</td></tr><tr><td>Sample (or H₂O)</td><td>0.100 mL</td></tr><tr><td colspan="2">Incub. 3 min, read A₁</td></tr><tr><td>R2</td><td>0.500 mL</td></tr><tr><td colspan="2">Incub. 15 min, read A₂</td></tr></table>	Steps	Samples (or blank)	R1	2.000 mL	Sample (or H ₂ O)	0.100 mL	Incub. 3 min, read A ₁		R2	0.500 mL	Incub. 15 min, read A ₂	
Steps	Samples (or blank)																													
Vial 1	1.000 mL																													
Vial 2	0.200 mL																													
Sample (or H ₂ O)	0.100 mL																													
Water	1.700 mL																													
Incub. 6 min, read A ₁																														
Vial 3	0.050 mL																													
Incub. 20 min, read A ₂																														
Steps	Samples (or blank)																													
R1	2.000 mL																													
Sample (or H ₂ O)	0.100 mL																													
Incub. 3 min, read A ₁																														
R2	0.500 mL																													
Incub. 15 min, read A ₂																														
Linearity	→ 500 mg/L (sample = 0.100 mL)	→ 500 mg/L (sample = 0.100 mL)																												
Calculation with Lambert-Beer: <ul style="list-style-type: none">• LoD (v = 1 mL, ΔA = 0.010)• LoQ (v = 1 mL, ΔA = 0.050)	<ul style="list-style-type: none">• 0.65 mg/L• 3.15 mg/L	<ul style="list-style-type: none">• 0.7 mg/L• 3.7 mg/L																												
The real LoD and LoQ have been measured only for Enzytec™ Liquid E8270 using the method DIN 32645:2008-11 with 100 µL sample: LoD = 10 mg/L and LoQ = 30 mg/L. By increasing the sample volume to 1 mL, these values can be reduced by factor 10.																														



Test comparison for D-Malic acid between Roche and Enzytec™ *Liquid*

Art. No. E8270

Evaluation Enzytec™ *Liquid* D-Malic acid (Art. No. E8270) versus Roche

a) Evaluation program for individual labs

We suggest following evaluation program:

- Purchase one kit Enzytec™ *Liquid* D-Malic acid (E8270, 50 tests) and run in parallel with the Roche kit
- For testing the recovery, always run a QC sample (to be purchased separately)
- If available, run a QC sample specific for the laboratory and calculate the recovery
- Run the routine samples in parallel with both kits until the 50 tests are finished (method comparison)
- Keep the same sample preparation as for the Roche kit and use the same vial containing the sample solution for both tests (100 µL each)
- Perform the two assays at the same time to avoid problems with sample stability (especially for Acetaldehyde, Acetic acid, Ammonia, Ascorbic acid, Ethanol and Sulfite which are unstable)

It is possible to test more validation criteria (e.g. linearity, precision or reproducibility), but more kits will be necessary.


b) Results of internal evaluation at R-Biopharm

- Recovery of QC samples

Recovery is $100 \pm 5\%$

Note: this is part of the specifications for all Enzytec™ *Liquid* test kits and is checked for every lot produced

- Method comparison



This method comparison was performed with different control materials and different types of fruit juices (orange juice, multi-fruit juice, apple juice). All points are well aligned along the equivalence line ($y = x$) and the coefficient of correlation (R^2) is 0.99.