


**The Best Way to GENOMIC DNA**

**Never Before ! 4-28-260 !**  
**4 kinds of Products divided !**  
**28 Protocols provided !**  
**260 Samples adjusted !**  
***i-genomic series***

**PART IV**  
***i-genomic Plant***  
**DNA Extraction Mini Kit Handbook**



[www.interchim.com](http://www.interchim.com)



211 bis Avenue Kennedy - BP 1140  
03103 Montluçon - France  
33 (0) 4 70 03 88 55  
Fax 33 (0) 4 70 03 82 60  
e-mail [interchim@interchim.com](mailto:interchim@interchim.com)

Agence Paris - Normandie  
33 (0) 1 41 32 34 40  
Fax 33 (0) 1 47 91 23 90  
e-mail [interchim.paris@interchim.com](mailto:interchim.paris@interchim.com)

---

## Contents

<b>Kit Contents</b>	<b>5</b>
<b>Storage</b>	<b>6</b>
<b>Product Use Limitations</b>	<b>7</b>
<b>Precautions and Safety Information</b>	<b>7</b>
<b>Product Warranty and Satisfaction Guarantee</b>	<b>7</b>
<b>Quality Control</b>	<b>7</b>
<b>Safety Information</b>	<b>7</b>
<b>Technical Assistance</b>	<b>8</b>
<b>i-genomic series Genomic DNA Extraction Mini Kits</b>	<b>9</b>
Table 2. Four Kinds of i-genomic series Genomic DNA Extraction Mini Kits	
<b>i-genomic series Protocol Table</b>	<b>10</b>
Table 3. Protocol Table (28 kinds of protocols)	
<b>i-genomic Plant DNA Extraction Mini Kit</b>	<b>11</b>
Table 4. Seven Kinds of Protocols according to Plant Samples	
<b>Plant Samples Grouping according to Protocols</b>	<b>13</b>
Figure 2. Plant Samples Grouping (Samples tested by iNtRON)	
<b>Equipments and Reagents to be Supplied by User</b>	<b>22</b>
<b>Notice Before Use</b>	<b>23</b>
Important Points Before Starting	23
Column Information	23
<b>Important Notes</b>	<b>24</b>
Choosing the Right Protocol according to Plant Sample	24
Collection and Storage of Plant Samples	24
Disruption and Homogenization	24
<b>Notes for Sample Sizing</b>	<b>25</b>
Measuring the Sample Amount after Disruption Step	25
Sample Volume	25
Table 5. Recommended Volume of Starting Material according to Plant Sample	

---

<b>Standard Protocols</b>	<b>29</b>
<b>Leaf</b>	<b>30</b>
Lyophilized Leaf :Type A Protocol	30
Fresh, Dried, or Frozen Leaf :Type B Protocol	36
<b>Stem</b>	<b>47</b>
Type C Protocol	47
<b>Root</b>	<b>54</b>
Type D Protocol	54
<b>Fruit</b>	<b>60</b>
Type E Protocol	60
<b>Seed</b>	<b>67</b>
All seeds w/o Gramineae :Type F Protocol	67
Gramineae :Type G Protocol	74
<b>Troubleshooting Guide</b>	<b>81</b>
<b>Appendix A : Electrophoresis Results from Various Plant Samples</b>	<b>84</b>
<b>Appendix B : Determination of Yield, Purity and PCR Data of DNA</b>	<b>85</b>
<b>Leaf</b>	
Lyophilized Leaf : From Type A Protocol	
Fresh, Dried, or Frozen Leaf : From Type B Protocol & Modified Protocol	
<b>Stem</b>	
From Type C Protocol & Modified Protocol	
<b>Root</b>	
From Type D Protocols & Modified Protocol	
<b>Fruit</b>	
From Type E Protocols & Modified Protocol	
<b>Seed</b>	
All seeds w/o Gramineae : From Type F Protocol	
Gramineae : From Type G Protocol & Modified Protocol	
<b>Appendix C : Source list of genomic DNA</b>	<b>102</b>
<b>Appendix D : Summarized photo procedure</b>	<b>105</b>
<b>Ordering Information</b>	<b>111</b>
<b>Global Distributors</b>	<b>112</b>

## i-genomic Plant DNA Mini Kit Contents

**i-genomic Plant DNA Extraction Mini Kit**  
(Cat. No. 17371, 50 columns)

**Table 1. Kit Contents**

Label	Description	Contain
<b>Buffer PG<sup>1</sup></b>	Lysis Buffer	30 ml
<b>Buffer PPT</b>	Precipitation buffer	7 ml
<b>Buffer PB</b>	Binding Buffer	50 ml
<b>Buffer PWA</b>	Washing Buffer A	40 ml
<b>Buffer PWB(concentrate)<sup>2</sup></b>	Washing Buffer B	10 ml (add 40 ml of EtOH)
<b>Buffer PE<sup>3</sup></b>	Elution Buffer	20 ml
<b>Enhancer Solution</b>	Binding enhancer buffer	0.5 ml
<b>Spin Columns (Green color O-ring)</b>	Inserted into the collection tubes. (2.0 ml tubes)	50 columns
<b>Collection Tubes (2.0ml tubes)</b>	Additionally supplied.	100 tubes
<b>RNase A Solution<sup>4</sup></b>	20 mg/ml (store at -20°C)	0.3 ml
<b>Proteinase K Solution<sup>5</sup></b>	20 mg/ml (store at -20°C)	1.2 ml

<sup>1</sup> Contains a chaotropic salt. Carefully handle. See page 7 for safety information.

<sup>2</sup> Buffer PWB is supplied as concentrate. Add 40 ml of ethanol (96 ~ 100%) according to the bottle label before use.

<sup>3</sup> Buffer PE is finally 10 mM Tris-HCl (pH 8.0). You may use your lab's buffer.

<sup>4</sup> Store at -20°C. The RNase A solution is completely free of DNase activity.

<sup>5</sup> Store at -20°C. After thawing, freshly use. We recommend to aliquot to small volume of Proteinase K. Use carefully according to the instruction manual (page 23).

Figure 1. i-genomic Plant DNA Extraction Mini Kit



### Storage

We recommend that all components of i-genomic Plant DNA Extraction Mini Kit is stored dry at room temperature (15 ~ 25°C). However, two components, including RNase A, and Proteinase K as a stock solution, should be stored at -20°C, and are stable for 1 year under these conditions.

### **Product Use Limitations**

All i-genomic series Kits are developed, designed and sold for research purpose only. They are not to be used for human or animal diagnosis of diseases. Do not use internally or externally in humans or animals. Be careful in the handling of the products.

### **Precautions and Safety Information**

All chemicals should be considered as potentially hazardous. When working with chemicals, always wear a suitable lab coat and disposable glove. Some buffer contain the chaotropic salt which may be an irritant and carcinogen, so appropriate safety apparel such as gloves and eye protection should be worn. If a spill of the buffers occurs, clean with a suitable laboratory detergent and water. If the liquid spill contains potentially infectious agents, clean the affected area first with laboratory detergent and water, then with a suitable laboratory disinfectant.

Only persons trained in laboratory techniques and familiar with the principles of good laboratory practice should handle these products.

### **Product Warranty and Satisfaction Guarantee**

All products undergo extensive quality control test and are warranted to perform as described when used correctly. Immediately any problems should be reported. Satisfaction guarantee is conditional upon the customer providing full details of the problem to iNtRON within 60 days, and returning the product to iNtRON for examination.

### **Quality Control**

As iNtRON quality control program, the performance of iNtRON's products are monitored routinely on a lot-to-lot basis. The genomic DNA yield of i-genomic series Genomic DNA Mini Kit is tested by preparing various samples and assaying the genomic DNA yield spectrophotometrically. The quality of isolated genomic DNA is checked by restriction digestion, PCR, agarose gel electrophoresis, and spectrophotometry. The i-genomic Plant DNA Mini Kit is tested to ensure the absence of DNase contamination. All buffers are each incubated with 1mg pUC18 DNA for 6 hours at 37 °C. The DNA is then visualized by electrophoresis on an agarose gel and compared to a positive control to determine if any linear or nicked plasmid DNA is visible.

### **Safety Information**

When working with chemicals, always should wear a suitable lab coat, disposable gloves and protective goggles. For more information, please request the appropriate material safety data sheets (MSDS). Do not add bleach or acidic solutions directly to the waste.

Buffer PG contains a chaotropic salts, which can form highly reactive compounds when combined with bleach. If liquid containing this buffer is spilt, clean with suitable laboratory detergent and water

### **Technical Assistance**

Our Technical Assistance Team is staffed by experienced researchers with extensive practical and theoretical expertise in molecular biology and the use of iNtRON products. If you have any questions or experience any difficulties regarding the i-genomic series Genomic DNA Extraction Mini Kits or other products in general, please do not hesitate to contact us.

Your information and questions are helpful to other scientists as well as to the researchers at iNtRON. We therefore encourage you to contact us if you have any suggestions about products' performances or new applications and techniques. For technical assistance and more information, please call or send e-mail.







Hot Line e-mail : [intronbio@intronbio.com](mailto:intronbio@intronbio.com), [info@intronbio.com](mailto:info@intronbio.com)

## i-genomic series

### Genomic DNA Extraction Mini Kits

i-genomic series Genomic DNA Extraction Mini Kits provide four kinds of kits according to the type of samples as seen in Table 2. These i-genomic series Kits provide a fast and easy way to purify DNA from various samples. The kits procedures provide pure genomic DNA for reliable PCR performance Southern blotting less than 1~2 hours. Purification require no phenol or chloroform extraction or alcohol precipitation. Pure DNA extracted by i-genomic series Kit is eluted in low-salt buffer or water, ready to use in downstream applications, including PCR, RAPD analysis, AFLP analysis, RFLP analysis, Southern blotting, microsatellite analysis, SNP-genotyping and quantitative real-time PCR. Purified DNA get an  $A_{260/280}$  ratio of 1.7 ~ 1.9, indicating high purity of the DNA.

**Table 2. Four Kinds of i-genomic series Genomic DNA Extraction Mini Kits**

Product Name	Samples	Examples
<b>CAT. NO. 17341</b> <b>i-genomic CTB</b> <b>DNA Extraction Mini Kit</b>	<u>C</u> ells <u>T</u> issues Gram(-) <u>B</u> acteria	Mouse / Guinea pig / Rabbit / Chicken / Zebra fish / Shrimp / Pig / Human cultured cells / Mouse cultured cells / Insect / Animal hair / Worm / Stool / Buccal swab / Gram(-) bacteria / Others
<b>CAT. NO. 17351</b> <b>i-genomic Blood</b> <b>DNA Extraction Mini Kit</b>	Blood	Whole blood / Buffy Coat / Dried Spot / Blood Swab / Plasma / Serum / Others
	<div style="display: flex; justify-content: space-around; text-align: center;"> <div><i>Whole blood</i> </div> <div><i>Buffy Coat</i> </div> <div><i>Dried Spot</i> </div> <div><i>Blood Swab</i> </div> <div><i>Plasma</i> </div> <div><i>Serum</i> </div> </div>	
<b>CAT. NO. 17361</b> <b>i-genomic BYF</b> <b>DNA Extraction Mini Kit</b>	Gram(+) <u>B</u> acteria <u>Y</u> east <u>F</u> ungi	<i>Azotobacter sp.</i> <i>Staphylococcus sp.</i> <i>Saccharomyces sp.</i> <i>Aspergillus sp.</i> Others
<b>CAT. NO. 17371</b> <b>i-genomic Plant</b> <b>DNA Extraction Mini Kit</b>	Plant	Leaf / Root / Stem / Fruit / Seed /Others



# i-genomic series Protocol Table

**Table 3. Protocol Table (28 kinds of protocols)**

i-genomic series Kits provide 28 of different protocols for each sample, and therefore are compatible with almost all samples. Just select optimal protocol type according to your sample, and follow the procedure. Refer to each kit's instruction manual.

Protocol Name	CTB	Blood	BYF	Plant	
Protocol Type	A B C D E F G H I J K	A B C D E F	A B C D	A B C D E F G	
Sample Treatment step	Preparation step Disrupt. & Homogen. Sample Sizing step Pre-Treating step				
DNA Extraction Step	Pre-Lysis step Lysis step Precipitation step DNA Binding step Washing step A Washing step B Elution step				
<input checked="" type="checkbox"/> Perform <input type="checkbox"/> Do Not Perform					

Sample Type	Total : 11 Types	Total : 6 Types	Total : 4 Types	Total : 7 Types
Cultured Cell	Whole Blood	Gram(+)/Bacteria	Yeast	Leaf
Animal Tissue	Buffy Coat	Gram(-)/Bacteria	Fungal Tissue	Lyophilized Leaf
Animal Tail	Dried Spot	Serum	Fungi	Stem
Rodent Tail	Blood Swab	Plasma	Yeast	Root
Formalin Fixed Tissue	Formalin Embedded Tissue	Blood Swab	Gram(+)/Bacteria	Fruit
Paraffin Embedded Tissue	Animal Hairs	Buffy Coat	Gram(-)/Bacteria	Seed (Others)
Animal Hairs	Insect /Worm	Whole Blood	Serum	Seed (Gramineae)
Buccal Swab	Bone	Buffy Coat	Plasma	
Stool	Stool	Dried Spot	Blood Swab	

### Sample Treatment Step

- Preparation step** : To prepare and to pick sample
- Disrupt. & Homogen.** : To disrupt or homogenize sample
- Sample Sizing step** : To take correct quantity of sample
- Pre-Treating step** : To treat PBS solution before lysis

### DNA Extraction Step

- Pre-Lysis & Lysis** : To pre-lyse or lyse sample
- Precipitation step** : To precipitate protein, polysaccharides, etc.
- DNA Binding step** : To bind DNA to the membrane
- Washing step A & B** : To wash the columns
- Elution step** : To elute genomic DNA from the column

## i-genomic Plant DNA Extraction Mini Kit

i-genomic Plant DNA Extraction Mini Kit provides a fast and easy way to purify DNA from plant-like samples such as various leaves, stems, roots, fruits, and seeds.

Furthermore, we have tested i-genomic Plant DNA Mini Kit to get more practical data with 104 plant samples. You can see vast sample photos, vast samples and vast practical data.


i-genomic Plant DNA Mini Kit provides 7 kinds of protocols, Type A, Type B, Type C, Type D, Type E, Type F and Type G. You can also extract genomic DNA from various plant samples in addition to 75 plant samples by selecting an appropriate protocol. When you choose a protocol, please refer to Plant Sample List (Figure 2). If you need some more information for selecting a protocol, please do not hesitate to contact our Technical Assistance Team.

**Table 4. Seven Kinds of Protocols according to Plant Sample**

Plant Samples	Protocol Type
<b>Leaf</b>	
Lyophilized leaf	Type A Protocol
Fresh, Dried, or Frozen leaf	Type B Protocol
<b>Stem</b>	
Stem	Type C Protocol
<b>Root</b>	
Root	Type D Protocol
<b>Fruit</b>	
Fruit	Type E Protocol
<b>Seed</b>	
All seeds w/o <i>Gramineae</i>	Type F Protocol
Gramineae	Type G Protocol

(Ex) Seed samples





## **Plant Samples Grouping according to Protocols**

---

i-genomic Plant DNA Mini Kit provides seven kinds of protocols according to plant samples. We recommend to select an appropriate protocol for your samples. The samples show 75 samples tested with i-genomic Plant DNA Mini Kit (see Figure 2). You can extract efficiently genomic DNAs from various plant samples.

iNtRON customers are a major source of information regarding advanced or specialized use of our products. This information is helpful to other scientists as well as to the researchers at iNtRON. We therefore encourage you to contact us if you have any suggestions about product performance, new applications or techniques.

*See next pages.*

Figure 2. Plant Samples Grouping (Samples tested by iNtRON)

## **LEAF**


For Lyophilized Leaf

Lyophilized Leaf



---

**Type A Protocol**



# LEAF

For Fresh, Dried, or Frozen Leaf



## Type B Protocol

Geranium  
(*Pelargonium inquinans*)



Lettuce  
(*Lactuca sarita*)



Peanut  
(*Arachis hypogaea*)



Buckwheat  
(*Fagopyrum esculentum*)



Green Pepper  
(*Capsicum annuum*)



African Millet  
(*Sorghum bicolor*)



Spinach  
(*Spinacia oleracea*)



Rhododendron  
(*Rhododendron schlippenbachii*)




Arabidopsis  
(*Arabidopsis*)



Cabbage  
(*Brassica oleracea*)



Cabbage  
(*Brassica campestris*)



Radish (Turnip)  
(*Brassica rapa*)



Sesame  
(*Sesamum indicum*)




Black Bean  
(*Glycine max*)



Fallen Leaf  
(*Fallen leaves*)



Moss  
(*Moss*)



## **STEM**

For Stem


Stem



---

### **Type C Protocol**

Potato  
(*Solanum tuberosum*)



Onion  
(*Allium cepa*)




Elm  
(*Ulmus davidiana*)



Radish  
(*Raphanus sativus*)



Geranium  
(*Pelargonium inquinans*)



Rhododendron  
(*Rhododendron schlippenbachii*)



Sesame  
(*Sesamum indicum*)




iNtRON Waits For  
Your Data.



# ROOT

For Root



## Type D Protocol

Sweet Potato  
(*Ipomoea batatas*)



Carrot  
(*Daucus carota*  
*var sativa*)



Platycodon  
(*Platycodon grandiflorum*)



Radish (Turnip)  
(*Brassica rapa*)



Sesame  
(*Sesamum indicum*)



iNtRON Waits For  
Your Data.



iNtRON Waits For  
Your Data.



iNtRON Waits For  
Your Data.





## **FRUIT**

For Fruit

Fruit



---

**Type E Protocol**

Melon  
(*Cucumis melo*)



Peach  
(*Prunus persica*)



Apple  
(*Malus pumila*)



Muscat  
(*Vitis spp*)



Tomato  
(*Lycopersicon esculentum*)



Pine Apple  
(*Ananas comosus*)



Grape  
(*Vitis vinifera* L.)



*iNtRON* Waits For  
Your Data.



# SEED

For All Seeds w/o Gramineae

## Seed (Others)



## Type F Protocol

### Leguminosae

Black Bean  
(*Glycine max*)



Mung Bean  
(*Phaseolus radiatus*)



Yellow Bean  
(*Glycine max*)



Pea  
(*Pisum sativum*)



### Cruciferae

Radish Type 1  
(*Raphanus sativus*)



Radish Type 2  
(*Raphanus sativus*)



Radish Type 3  
(*Brassica rapa*)



iNtRON Waits For Your Data.



### Solanaceae

Green Pepper  
(*Capsicum annuum*)



Paprika  
(*Capsicum annuum*)



Petunia  
(*Petunia hybrida*)



Tomato Type 1  
(*Lycopersicon  
esculentum Type 1*)



Tomato Type 2  
(*Lycopersicon  
esculentum Type 2*)



Tomato Type JPN  
(*Lycopersicon  
esculentum Type JPN*)



iNtRON Waits For Your Data.



iNtRON Waits For Your Data.



## SEED

For All Seeds w/o Gramineae

### Seed (Others)



---

### Type F Protocol

#### Cucurbitaceae

Pumpkin  
(*Cucurbita moschata*  
*duchesne*)



Melon  
(*Cucumis melo*  
*var. makuwa*)



Cucumber  
(*Cucumis sativus*)



Water Melon  
(*Citrullus vulgaris*)



Gourd  
(*Lagenaria leucantha*)



Melon  
(*Cucumis melo*)



iNtRON Waits For  
Your Data.



iNtRON Waits For  
Your Data.




#### Others


Carrot  
(*Daucus carota* *var. sativa*)



Spinach  
(*Spinach oleracea*)



Sesame  
(*Sesamum indicum*)



Sunflower  
(*Helianthus annuus*)



# SEED

For Gramineae

Seed (Gramineae)



---

## Type G Protocol

Grass  
(*Zoysia japonica*)



Rice  
(*Oryza sativa*)



African Millet  
(*Sorghum bicolor*)



Wheat  
(*Triticum aestivum vulgare*)



Barley  
(*Hordeum vulgare var. hexastichon*)



iNtRON Waits For  
Your Data.



iNtRON Waits For  
Your Data.



iNtRON Waits For  
Your Data.



## Equipments and Reagents to Be Supplied by User

---

i-genomic Plant DNA Mini Kit provides almost all reagents for extracting DNA, including RNase A and Proteinase K stock solutions. However, you should prepare some equipments and reagents as follows for a fast and easy extraction. When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

- ✓ Equipment for disruption and homogenization, including grinding jar set (mortar)
- ✓ Pipettes and pipette tips
- ✓ Water bath or heating block
- ✓ Vortex mixer
- ✓ Microcentrifuge with rotor for 2.0 ml tubes.
- ✓ Microcentrifuge tubes (1.5 ml)
- ✓ Liquid nitrogen
- ✓ Absolute ethanol (EtOH, 96 ~ 100%)
- ✓ Ice
- ✓ Other general lab equipments

## Notice Before Use

### Important Points Before Starting

- **Buffer PWB (Washing Buffer B)**  
: Buffer PWB is supplied as concentrate. Before using for the first time, be sure to add 40 ml of absolute ethanol (96 ~ 100% EtOH) to obtain a working solution
- **Enhancer Solution**  
: We recommend keeping the Enhancer Solution at 2 ~ 8°C upon arrival, and then it is stable for 1 year. The Enhancer Solution is a yellow color solution with precipitates upon storage. Before use, warm it to 65°C to redissolve. After use, store at 2 ~ 8°C.
- **Proteinase K Solution (20 mg/ml)**  
: Proteinase K possesses a high specific activity which remains stable over a wide range of temperature and pH values with substantially increased activity at higher temperatures. Proteinase K Solution shows a milk-white color, since it is supplied as concentrate. After thawing, freshly use. DO NOT heat to redissolve. We recommend to aliquot to small volume of Proteinase K Solution.
- **Pre-heat a water bath or a heating block to 65°C.**
- **Centrifugation**  
: All centrifugation steps are carried out at RT (15 ~ 25°C ) in a microcentrifuge.

### Column Information

- **i-genomic series Spin Column**

<b>Column Membrane</b> <sup>1</sup>	Silica-based membrane
<b>Spin Column</b> <sup>1</sup>	Individually, is inserted in a 2.0 ml collection tube <sup>2</sup> .
<b>Loading Capacity</b>	Maximum 800 $\mu\text{l}$
<b>DNA Binding Capacity</b>	Maximum 45 $\mu\text{g}$
<b>Recovery</b>	85 ~ 95% depending on the elution volume
<b>Elution Volume</b>	Generally, eluted with 30 ~ 200 $\mu\text{l}$ of elution buffer

<sup>1</sup> After use, seal the pack containing spin columns tightly without getting dry. Then, the spin columns are stable for over 1 years under these conditions. It's not good for DNA binding to be dried completely.

<sup>2</sup> Additional collection tubes (100 ea) are also supplied for your convenient handling.

## Important Notes

---

### **Choosing the Right Protocol according to Plant Sample**

Seven kinds of different protocols in this handbook provide detailed instructions to use i-genomic Plant DNA Mini Kit for purifying genomic DNA from various plant samples (see Figure 2). These protocols are optimized for use. Especially, with iNtRON's 75 samples more than five categories, we have verified practically by several experiments to ensure the quality and the application of i-genomic Plant DNA Mini Kit. We recommend you will choose the right protocol according to your plant sample. For more information, please contact iNtRON Technical Assistance Team.

### **Collection and Storage of Plant Samples**

Generally, for higher quality of genomic DNA, it is preferable to collect from young plant materials. The fresh plant tissues contain more viable cells per weight and therefore result in higher yields and purity of genomic DNA. If plant tissue will not be used freshly, after harvesting, we recommend it should be quickly frozen in liquid nitrogen, and then stored at -80°C. When grind plant tissues, use liquid nitrogen as in detailed instruction manual. Ground plant tissue powder can also be stored at -80°C.

Alternatively, plant tissues can be dried or lyophilized after harvesting to allow storage at room temperature. To ensure high quality of genomic DNA, plant tissues should be completely dried within 24 hours of collection.

### **Disruption and Homogenization**

Almost all samples can be disrupted on a mortar after freezing in liquid nitrogen without Buffer PG. Namely, disruption can be performed without lysis buffer by keeping the sample submerged in liquid nitrogen before and during disruption on a mortar. Especially hard tissues, such as roots or seeds, are relatively difficult to be disrupted, and therefore be careful to use a mortar in liquid nitrogen. Alternatively, fresh leaf can be directly disrupted in lysis buffer (Buffer PG) without using liquid nitrogen, but this may cause shearing of high molecular-weight DNA. We do not recommend frozen materials to disrupt in lysis buffer as this can result in low yields and degraded DNA.

In case of dried or lyophilized leaf, this disruption step is omitted, and therefore do not add liquid nitrogen. After Sample Sizing step, directly treat lysis buffer (Buffer PG) to the sample for lysis step. For optimal results, we recommend to keep the disruption time as short as possible.

When disruption for more than 1 minute it may lead to shearing of genomic DNA.

## Notes for Sample Sizing

### Measuring the Sample Amount after Disruption Step

We recommend to measure the amount of starting material after disrupting plant tissues in i-genomic Plant DNA Mini Kit, since it brings a loss of starting material. If the plant samples are disrupted on a mortar submerged in liquid nitrogen, the samples will turn into powder-like form. It makes them to measure conveniently the amount of starting material.

Table 5 shows a recommended amount of starting material after disruption step according to plant samples. Please follow the manual instruction not to be over the required amounts.

### Sample Volume

i-genomic Plant DNA Mini Kit procedures are optimized for 5 mg ~ 100 mg of wet-weight starting material. Table 5 provides guidelines according to the plant tissues.

Exceeding the recommended amount of starting material will result in inefficient lysis, resulting in low yield and purity. In the large, DNA yields and purity will be varied by depending on genome size, sample viscosity and age of sample.

**Table 5. Recommended Volume of Starting Material according to Plant Sample**

#### ▪ Type A Protocol : Lyophilized Leaf

Plant Tissue	Amount	Plant Tissue	Amount
Persimmon ( <i>Diopyros kaki</i> )	5 mg	Amaranth ( <i>Amarantus</i> )	5 mg
Millet, Hog Millet ( <i>Panicum millacerm</i> )	5 mg	Corn ( <i>Zea mays</i> )	5 mg
Mung Bean ( <i>Phasedus radiutus</i> )	5 mg	Leek ( <i>Allium tuberosum</i> )	5 mg
Perilla ( <i>Perilla fru tescens</i> )	5 mg	Italian Millet ( <i>Setania italica</i> )	5 mg
Paddy ( <i>Oryza sativa</i> )	5 mg	Adlay ( <i>Coix lachrymajobi</i> )	5 mg
Ginger ( <i>Zingiber officinale</i> )	5 mg	Ginseng ( <i>Panax ginseng</i> )	5 mg



▪ **Type B Protocol : Leaf**

*Fresh, Dried, or Frozen Leaf*

Plant Tissue	Amount
Geranium ( <i>Pelargonium inquinans</i> )	50 mg
Lettuce ( <i>Lactuca sarita</i> )	50 mg
Peanut ( <i>Arachis hypogaea</i> )	50 mg
Buckwheat ( <i>Fagopyrum esculentum</i> )	50 mg
Green Pepper ( <i>Capsicum annum</i> )	50 mg
African Millet ( <i>Sorghum bicolor</i> )	50 mg
Spinach ( <i>Spinacia oleracea</i> )	50 mg
Rhododendron ( <i>Rhododendron schlippenbanchii</i> )	50 mg

Plant Tissue	Amount
<i>Arabidopsis (Arabidopsis)</i>	50 mg
Cabbage ( <i>Brassica oleracea</i> )	50 mg
Cabbage ( <i>Brassica campestris</i> )	50 mg
Radish, Turnip ( <i>Brassica rapa</i> )	50 mg
Sesame ( <i>Sesamum indicum</i> )	50 mg
Black Bean Leaves ( <i>Glycine max</i> )	50 mg
Fallen Leaf ( <i>Fallen leaves</i> )	50 mg
Moss ( <i>Moss</i> )	50 mg

▪ **Type C Protocol : Stem**

Plant Tissue	Amount
Potato ( <i>Solanum tuberosum</i> )	50 mg
Onion ( <i>Allium cepa</i> )	50 mg
Elm ( <i>Ulmus davidiana</i> )	50 mg
Radish ( <i>Raphanus sativus</i> )	50 mg

Plant Tissue	Amount
Geranium ( <i>Pelargonium inquinans</i> )	50 mg
Rhododendron ( <i>Rhododendron schlippenbanchii</i> )	50 mg
Sesame ( <i>Sesamum indicum</i> )	50 mg

▪ **Type D Protocol : Root**

Plant Tissue	Amount
Sweet Potato ( <i>Ipomoea batatas</i> )	50 mg
Carrot ( <i>Daucus carota var.sativa</i> )	50 mg
Platycodon ( <i>Platycodon grandiflorum</i> )	50 mg

Plant Tissue	Amount
Radish, Turnip ( <i>Brassica rapa</i> )	50 mg
Sesame ( <i>Sesamum indicum</i> )	50 mg

▪ **Type E Protocol : Fruit**

Plant Tissue	Amount
Melon ( <i>Cucumis melo</i> )	100 mg
Peach ( <i>Prunus persica</i> )	100 mg
Apple ( <i>Malus pumila</i> )	100 mg
Muscat ( <i>Vitis spp</i> )	100 mg

Plant Tissue	Amount
Tomato ( <i>Lycopersicon esculentum</i> )	100 mg
Pine Apple ( <i>Ananas comosus</i> )	100 mg
Grape ( <i>Vitis vinifera L.</i> )	100 mg

▪ **Type F Protocol : Seed**

**Leguminosae**

Plant Tissue	Amount
Black Bean ( <i>Glycine max</i> )	50 mg
Mung Bean ( <i>Phaseolus radiatus</i> )	50 mg
Yellow Bean ( <i>Glycine max</i> )	50 mg
Pea ( <i>Pisum sativum</i> )	50 mg

**Cruciferae**

Plant Tissue	Amount
Radish Type 1 ( <i>Raphanus sativus</i> )	50 mg
Radish Type 2 ( <i>Raphanus sativus</i> )	50 mg
Radish Type 3 ( <i>Brassica rapa</i> )	50 mg

**Solanaceae**

Plant Tissue	Amount
Green Pepper ( <i>Capsicum annuum</i> )	50 mg
Paprika ( <i>Capsicum annuum</i> )	50 mg
Petunia ( <i>Petunia hybrida</i> )	50 mg

Plant Tissue	Amount
Tomato Type 1 ( <i>Lycopersicon esculentum Type 1</i> )	50 mg
Tomato Type 2 ( <i>Lycopersicon esculentum Type 2</i> )	50 mg
Tomato JPN ( <i>Lycopersicon esculentum JPN</i> )	50 mg

**Cucurbitaceae**

Plant Tissue	Amount
Pumpkin ( <i>Cucurbita moschata duchesne</i> )	50 mg
Melon ( <i>Cucumis melo var. makuwa</i> )	50 mg
Cucumber ( <i>Cucumis sativus</i> )	50 mg

Plant Tissue	Amount
Water Melon ( <i>Citrullus vulgaris</i> )	50 mg
Gourd ( <i>Lagenaria leucantha</i> )	50 mg
Melon ( <i>Cucumis melo</i> )	50 mg

**Others**

Plant Tissue	Amount
<b>Umbelliferae</b>	
Carrot ( <i>Daucus carota var. sativa</i> )	50 mg
<b>Chenopodiaceae</b>	
Spinach ( <i>Spinacia oleracea</i> )	50 mg

Plant Tissue	Amount
<b>Pedaliaceae</b>	
Sesame ( <i>Sesamum indicum</i> )	50 mg
<b>Compositae</b>	
Sunflower ( <i>Helianthus annuus</i> )	50 mg

**• Type G Protocol : Seed****Gramineae**

Plant Tissue	Amount
Grass ( <i>Zoysia japonica</i> )	100 mg
Rice ( <i>Oryza sativa</i> )	10 mg
African Millet ( <i>Sorghum bicolor</i> )	10 mg

Plant Tissue	Amount
Wheat ( <i>Triticum aestivum vulgare</i> )	10 mg
Barley ( <i>Hordeum vulgare var. hexastichon</i> )	10 mg

# ***Standard Protocols***



**Plant**  
**Type A** *Type A Protocol*  
*Bench Protocol* **For Lyophilized Leaf**


A-1. Sample Treatment Step

▪ **Lyophilized Leaf**

- I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step



I. Preparation step

**1**  Prepare lyophilized leaf sample

III. Sample Sizing step

**2**  Measure 5 mg of lyophilized leaf




 Transfer to 1.5 ml tube

A-2. DNA Extraction Step


## A-2. DNA Extraction Step

▪ **Lyophilized Leaf**


- I. Pre-Lysis step   
  II. Lysis step   
  III. Precipitation step   
  IV. DNA Binding step  
 V. Washing step A   
  VI. Washing step B   
  VII. Elution step




V. Washing setp A



VI. Washing step B



VII. Elution step



## A-1. Sample Treatment Step

**Lyophilized Leaf**

- I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step

**I. Preparation step****1. Prepare lyophilized leaf sample.**

Plant tissue can be lyophilized after harvesting to allow storage at room temperature (15 ~ 20°C). To ensure DNA quality, we recommend that samples should be completely lyophilized within 24 hours of collection. Generally, lyophilized leaves are fine powder form, and therefore don't need to have any special disruption & homogenization steps. Furthermore, when using lyophilized tissue, the samples do not need to be frozen in liquid nitrogen. Directly, you can perform the Sample Sizing Step.

**III. Sample Sizing step****2. Measure 5 mg of lyophilized leaf, and then transfer into 1.5 ml tube using a spatula.**

It's difficult to handle to measure the fine powder sample due to its static electricity. It can be inhibited by previously chill the spatula and 1.5 ml tube in liquid nitrogen. Exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.

## A-2. DNA Extraction Step

**Lyophilized Leaf**

- I. Pre-Lysis step    II. Lysis step    III. Precipitation step    IV. DNA Binding step  
 V. Washing step A    VI. Washing step B    VII. Elution step

Equilibrate samples at room temperature (15 ~ 25°C).  
Heat a water bath or heating block to 65°C for use in step 2.  
All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu$ l Buffer PG, 7  $\mu$ l Enhancer Solution, 20  $\mu$ l Proteinase K, and 5  $\mu$ l RNase A Solution into sample tube, and vortex vigorously.**

With lyophilized leaf it absorbs lysis buffer, and becomes swollen. It may be difficult to handle plant tissue due to its viscosity. Always keep the recommended amount of starting material. Furthermore, vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.



**2. Incubate the lysate at 65°C for 30 min.**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate. Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended sample amount.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min.**

Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step.

**IV. DNA Binding step****5. Transfer carefully 200  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

Although the supernatant is typically 350 ~ 400  $\mu\text{l}$ , we recommend to recover only 200  $\mu\text{l}$  of lysate. More lysate can result in shearing of the DNA and contaminating the next step with impurities. When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane. A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA Mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge at 13,000 rpm (RT) for 1 min, and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

**V. Washing step A****9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through and reuse the collection tube in step 10.**

## VI. Washing step B

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge at 13,000 rpm for 1 min .**

**Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane.**

**Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

## VII. Elution Step

**11. Place the spin column into a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate at room temperature for 1 min, and then centrifuge at 13,000 rpm for 1 min to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.


Plant  
**Type B**  
Bench Protocol

**Type B Protocol**  
**For Fresh, Dried, or Frozen Leaf**

B-1. Sample Treatment Step

▪ **Fresh Leaf**

- I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step




• **Dried Leaf**

■ I. Preparation step ■ II. Disrupt.& Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step


Type B Protocol  
Leaf





I. Preparation step

1  Prepare dried leaf sample


II. Disruption & Homogenization step


2  Place the sample into mortar


 Add liquid nitrogen for sample freezing

 Pouring liquid nitrogen slowly, Disrupt and homogenize completely

III. Sample Sizing step

3  Measure 50 mg of sample powder

 Transfer to 1.5 ml tube



B-2. DNA Extraction Step


Type B Protocol  
Leaf

▪ **Frozen Leaf**


■ I. Preparation step ■ II. Disrupt.& Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step





I. Preparation step

1  Prepare from frozen leaf sample


II. Disruption & Homogenization step


2  Place the sliced or whole sample into mortar


 Add liquid nitrogen for sample freezing

 Pouring liquid nitrogen slowly. Disrupt and homogenize completely

III. Sample Sizing step

3  Measure 50 mg of sample powder

 Transfer to 1.5 ml tube




 B-2. DNA Extraction Step

B-2. DNA Extraction Step




▪ **Fresh, Dried, or Frozen Leaf**

- I. Pre-Lysis step   
  II. Lysis step   
  III. Precipitation step   
  IV. DNA Binding step  
 V. Washing step A   
  VI. Washing step B   
  VII. Elution step

Type B Protocol  
 Leaf



Type B Protocol  
Leaf



## B-1. Sample Treatment Step

▪ **Fresh Leaf**

- I. Preparation step   ■ II. Disrupt.& Homogen.   ■ III. Sample Sizing step   □ IV. Pre-Treating step

Type B Protocol  
Leaf**I. Preparation step****1. Prepare fresh leaf sample.**

If possible, it is preferable to collect young fresh leaf since they contain more viable cells per weight and therefore result in higher yields. In addition, young fresh leaves contain smaller amounts of polysaccharides and polyphenolics. Therefore, they are easier to handle than other plants. For storage of harvested fresh leaf, in general, when genomic DNA is to be isolated, plant leaves from most species should be frozen and kept at -80°C after harvesting. It is good for disruption and homogenization if the sample is sliced off when it stores at -20°C or -80°C.

**II. Disruption & Homogenization step****2. Slice off the prepared sample to small pieces by the scalpel or scissor.**

To reduce disruption and homogenization time, we recommend to slice it off.

**3. Place the sliced sample material into grinding mortar and add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately step 4.**

Disruption and homogenization time depends on the leaf samples. We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of fresh leaf sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen.

Generally, it is a fine powder form after disruption and homogenization.

**III. Sample Sizing step****4. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5 ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation.

And more, exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.



▪ **Dried Leaf**

■ I. Preparation step ■ II. Disrupt. &amp; Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step

**I. Preparation step****1. Prepare dried leaf sample.**

If possible, it is preferable to collect young leaves since they contain more viable cells per weight and therefore result in higher yields. Alternatively, leaf can be also dried after harvesting to allow storage at room temperature (15 ~ 20°C).

**II. Disruption & Homogenization step****2. Place the sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately step 3.**

Disruption and homogenization time depend on the leaf samples. We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of dried sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Generally, it is a fine powder form after disruption and homogenization.

**III. Sample Sizing step****3. Measure 20 ~ 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

If too much dry, it absorbs all of Buffer PG (lysis buffer) to difficult to handle to lysis. And in that time, some dried sample becomes swollen. To prevent it, reduce the amount of starting sample material if it is. If it is, we recommended below 30 mg of amount starting material. Especially, the ocean dried grass organ. In case of organ, normally use 10 mg for amount of starting material. But in case of fallen leaves, use 50 mg for extracting DNA. Exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity. To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5 ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA

**▪ Frozen Leaf**

I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step

Type B Protocol  
Leaf

**I. Preparation step****1. Prepare frozen leaf sample.**

It is important to keep the frozen leaf sample frozen in liquid nitrogen during all of sample treatment step to inhibit low DNA yields and degraded DNA. More long-term storage of the sample, it should be frozen and kept at -80°C after harvesting. It is good for disruption and homogenization if the sample is sliced off when store at -20°C or -80°C.

Young fresh leaf is good for extraction of DNA, if possible, we recommend to collect young fresh leaf since they contain more viable cells per weight and therefore result in higher yields.

**II. Disruption & Homogenization step****2. Place the sliced or whole sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately to step 3.**

Disruption and homogenization time depends on the leaf samples. We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of fresh leaf sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Generally, it is a fine powder form after disruption and homogenization.

**III. Sample Sizing step****3. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transfer it, previously pre-chill the spatula and 1.5 ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation.

And more, exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.

## B-2. DNA Extraction Step

▪ **Fresh, Dried or Frozen Leaf**

- I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step
- V. Washing step A     VI. Washing step B     VII. Elution step

- Equilibrate samples to room temperature (15 ~ 25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu$ l Buffer PG, 7  $\mu$ l Enhancer Solution, 20  $\mu$ l Proteinase K, and 5  $\mu$ l RNaseA Solution into sample tube, and vortex vigorously.**

Like lyophilized or dried leaf, ground fresh and frozen leaf absorbs Buffer PG (lysis buffer), and becomes swollen. It may be difficult to handle plant tissue due to its viscosity. Always keep the recommended amount of starting material. Furthermore, vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

**2. Incubate the lysate at 65°C for 30 min .**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu$ l Buffer PPT to the lysate, mix well, and incubate on ice for 5 min .**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate. Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min .**

Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step. Plant leaf contains chloroplast, including much have chlorophyll, therefore shows green color. So, after centrifugation, the color of supernatant shows clear green color. But disappear after pass through washing process two times.

#### IV. DNA Binding step

**5. Transfer carefully 200  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

Although the supernatant is typically 400 ~ 450  $\mu\text{l}$ , we recommend to recover only 200  $\mu\text{l}$  of lysate. More lysate can result in shearing of the DNA and contaminating the next step with impurities. When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane. A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA Mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge at 13,000 rpm (RT) for 1 min, and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

#### V. Washing step A

**9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through and reuse the collection tube in step 10.**

#### VI. Washing step B

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

### **VII. Elution Step**

**11. Place the spin column into the a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate at room temperature for 1 min, and then centrifuge at 13,000 rpm for 1 min to elute.**

With 50  $\mu\text{l}$  of Buffer PE increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

Plant  
**Type C**  
Bench Protocol


**Type C Protocol**  
For Stem

C-1. Sample Treatment Step

• Stem

- I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step

Type C Protocol  
Stem









C-2. DNA Extraction Step

▪ **Stem**

- I. Pre-Lysis step   
  II. Lysis step   
  III. Precipitation step   
  IV. DNA Binding step  
 V. Washing step A   
  VI. Washing step B   
  VII. Elution step

Type C Protocol  
Stem





Type C Protocol  
Stem



## C-1. Sample Treatment Step

▪ **Stem**

■ I. Preparation step   ■ II. Disrupt.& Homogen.   ■ III. Sample Sizing step   □ IV. Pre-Treating step

Type C Protocol  
Stem

**I. Preparation step****1. Prepare the stem sample.**

We recommend to collect the fresh stem, but if it is impossible, it is preferable to collect non dried stem sample. The more stem sample is dried, the more it is difficult to lysis because it absorbs all of lysis buffer. Also, plant stem sample is very tough and thick, therefore it is difficult to disrupt and homogenize. Same as other plant samples, although its organization is similar to leaf's organization, it is so hard plant tissue. For storage of harvested stem sample, in general, when genomic DNA is to be isolated, to be not dried, we recommend to freeze or keep in plastic bag containing a wet paper towel after harvesting.

It is good for disruption and homogenization if the sample is sliced off when it store.

**II. Disruption & Homogenization step****2. Cut off the prepared sample to small pieces by the blade or scissor.**

To disruption and homogenization perfectly in trunk of tree, peel the shell of trunk. If not, the shell powder is floated above the supernatant after the precipitation step. It has an effect on low DNA yield and purity.

To reduce disruption and homogenization time, we recommend to cut it off. Some stem sample is very tough and thick, to cut off it by the saw.

**3. Place the sliced sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately step 4.**

In stem sample, it is difficult to disrupt and homogenize to the fine powder using liquid nitrogen with mortar and pestle. But, we recommend to be disrupted completely until no tissue clumps are not visible. Clumps of stem sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen.

**III. Sample Sizing step****4. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5 ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation. And more, exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.

## C-2. DNA Extraction Step

▪ **Stem**

- I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step
- V. Washing step A     VI. Washing step B     VII. Elution step

- Equilibrate samples to room temperature (15 ~ 25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu\text{l}$  Buffer PG, 7  $\mu\text{l}$  Enhancer Solution, 20  $\mu\text{l}$  Proteinase K, and 5  $\mu\text{l}$  RNase A Solution into sample tube, and vortex vigorously.**

Although ground stem sample is not fine powder, it absorbs Buffer PG (lysis buffer). Especially, dried stem sample absorbs all of Buffer PG. To prevent this case, always keep the recommended amount of starting material. Be careful. This step is cell lysis step. After adding them, immediately vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

**2. Incubate the lysate at 65°C for 30 min .**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate. Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min .**

Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step.

Due to the nature color of stem, sometimes supernatant shows green or dark brown color. Nevertheless, it will disappear after washing step.

#### **IV. DNA Binding step**

**5. Transfer carefully 200  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

In trunk of tree, if not peel the shell when the disruption and homogenization step, you can see the shell in above supernatant. In this situation, the shell is attached the external of yellow tip, while taking the supernatant with yellow tip by pipette. Be careful of taking the supernatant. Although the supernatant is typically 400 ~ 450  $\mu\text{l}$ , too much dried stem sample supernatant is 100 ~ 200  $\mu\text{l}$ . If it is, reduce the amount of starting material.

Always keep the recommended amount of starting material.

Normally, we recommend to recover only 200  $\mu\text{l}$  of lysate in stem sample. More lysate can result in shearing of the DNA and contaminate the next step with impurities.

When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane. A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge at 13,000 rpm (RT) for 1 min, and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

After centrifugation, sometimes spin column membrane becomes light brown or light green color. Do not worry about that. It will be disappeared after washing step.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

#### **V. Washing step A**

**9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through and reuse the collection tube in step 10.**

#### **VI. Washing step B**

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

### **VII. Elution Step**

- 1. Place the spin column into a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate for 1 min at room temperature, and then centrifuge for 1 min at 13,000 rpm to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

Plant  
**Type D**  
Bench Protocol


**Type D Protocol**  
For Root

D-1. Sample Treatment Step

• Root

■ I. Preparation step ■ II. Disrupt. & Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step




Type D Protocol  
Root



D-2. DNA Extraction Step




• **Root**

- I. Pre-Lysis step   
  II. Lysis step   
  III. Precipitation step   
  IV. DNA Binding step  
 V. Washing step A   
  VI. Washing step B   
  VII. Elution step



Type D Protocol  
**Root**

Type D Protocol  
Root



## D-1. Sample Treatment Step

▪ **Root**

■ I. Preparation step ■ II. Disrupt. & Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step

**I. Preparation step****1. Prepare root sample.**

In root sample, we recommend to collect same condition of the fresh stem sample. Although root is very tough and thick as stem, it is not difficult to disrupt and homogenize. Same as stem, root sample does not keep wet condition, therefore we recommend to frozen condition (-20°C or -80°C) for storage of harvested of root sample. Before freezing, remove impurities on surface of root sample that are inhibited DNA extraction by washing. And then, store at the frozen condition after removing the wetness. As ever the other plant sample, it is good for disruption and homogenization if the sample is sliced off when it stores.

In frozen root, it is important to keep frozen state in liquid nitrogen during all of sample treatments step to inhibit low DNA yields and degraded DNA. Do not repeat freezing-thawing.

Type D Protocol  
Root

**II. Disruption & Homogenization step****2. Cut off the prepared sample to small pieces by the blade or scissor.**

To reduce disruption and homogenization time, we recommend to cut it off.  
Some root sample is very tough and thick, to cut off it by the saw.

**3. Place the sliced sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely.****Allow the liquid nitrogen to evaporate, and proceed immediately to step 4.**

We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of sample will not lyse properly and will be resulted in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Generally, it is a fine powder form after disruption and homogenization. Be careful to handle liquid nitrogen.

**III. Sample Sizing step****4. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation.

And more, exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.



## D-2. DNA Extraction Step

▪ **Root**

I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step

V. Washing step A     VI. Washing step B     VII. Elution step

- Equilibrate samples to room temperature (15 ~ 25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu\text{l}$  Buffer PG, 7  $\mu\text{l}$  Enhancer Solution, 20  $\mu\text{l}$  Proteinase K, and 5  $\mu\text{l}$  RNase A Solution into sample tube, and vortex vigorously.**

In ground root sample, it is properly fine powder form. Like ground fresh leaf, it absorbs Buffer PG (lysis buffer) after adding them. It is not difficult to mix until any root clumps are not visible if only keep the recommended amount of starting material and immediately vortex or pipette them vigorously. If not, it may be difficult to handle it due to its viscosity. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA.

A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

During lysis step, the lysate color express to mix the cloudy and the nature color of root sample.

**2. Incubate the lysate at 65°C for 30 min.**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate.

Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min .**

Plant materials can create very viscous lysates and large amounts of precipitates during this step.

If you keep our recommended amounts of starting material, optimal results are obtained.

If not, you should perform one more centrifugation step.

**IV. DNA Binding step****5. Transfer carefully 200  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

Although the supernatant is typically 350 ~ 400  $\mu\text{l}$ , we recommend to recover only 200  $\mu\text{l}$  of lysate.

More lysate can results in shearing of the DNA and contaminating the next step with impurities.

When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane.

A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge for 1 min at 13,000 rpm (RT), and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

#### V. Washing step A

**9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through and reuse the collection tube in step 10.**

#### VI. Washing step B

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge at 13,000 rpm for 1 min.**

**Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

#### VII. Elution Step

**11. Place the spin column into a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate at room temperature for 1 min, and then centrifuge at 13,000 rpm for 1 min to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

Plant

**Type E**

**Bench Protocol**

**Type E Protocol**

**For Fruit**

E-1. Sample Treatment Step


• **Fruit**

■ I. Preparation step ■ II. Disrupt. & Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step


Type E Protocol  
**Fruit**





I. Preparation step


1  Prepare fruit sample

II. Disruption & Homogenization step


2  Slice off the prepared sample to small pieces


3  Place the sliced sample into mortar

 Add liquid nitrogen for sample freezing

 Pouring liquid nitrogen slowly, Disrupt and homogenize completely

III. Sample Sizing step

4  Measure 100 mg of sample powder




 Transfer to 1.5ml tube

E-2. DNA Extraction Step

E-2. DNA Extraction Step




• **Fruit**

- I. Pre-Lysis step   
  II. Lysis step   
  III. Precipitation step   
  IV. DNA Binding step  
 V. Washing step A   
  VI. Washing step B   
  VII. Elution step



Type E Protocol  
**Fruit**

Type E Protocol  
Fruit



## E-1. Sample Treatment Step

▪ **Fruit**

■ I. Preparation step   ■ II. Disrupt.& Homogen.   ■ III. Sample Sizing step   □ IV. Pre-Treating step

**I. Preparation step****1. Prepare fresh fruit sample.**

While other plant's organ contains a little water, fruit organ contains water about 80 ~ 95%, therefore increase the amount of starting material during the experiment. When keep fruit sample long at RT, It should denaturalize easily and will therefore result in a lower yield of DNA. For storage of fruit, we recommend to frozen and keep at -80°C or lyophilize it and store at room temperature (15°C ~ 25°C) after harvesting. Before freezing, peel the shell the fruit sample and slice off it to suitable size. When extracting DNA from fruit, generally use the flesh of fruit.

In frozen fruit, it is important to keep frozen state in liquid nitrogen during all of sample treatment step to inhibit low DNA yields and degraded DNA. Do not repeat freezing-thawing.

**II. Disruption & Homogenization step****2. Slice off the prepared sample to small pieces by the blade or scissor.**

To reduce disruption and homogenization time, we recommend to slice off.

**3. Place the sliced sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately to step 4.**

We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Generally, it is a fine powder form after disruption and homogenization.

**III. Sample Sizing step****4. Measure 100 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation. Because fruit organ contains water about 80 ~ 95%, increase the amount of starting material. But exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.

## E-2. DNA Extraction Step

▪ **Fruit**

I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step

V. Washing step A     VI. Washing step B     VII. Elution step

- Equilibrate samples to room temperature (15 ~ 25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu\text{l}$  Buffer PG, 7  $\mu\text{l}$  Enhancer Solution, 20  $\mu\text{l}$  Proteinase K, and 5  $\mu\text{l}$  RNase A Solution into sample tube, and vortex vigorously.**

With fruit sample absorbs lysis buffer, and becomes swollen. When apply exceeding the recommended amount of starting material, it may be difficult to handle plant tissue due to its viscosity. Therefore, always keep the recommended amount of starting material. Furthermore, vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

**2. Incubate the lysate at 65°C for 30 min .**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate.

Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min.**

Because plant fruit sample contain a few amount of fat, protein and other material, fruit sample create a small size and few amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step.

**5. Transfer carefully 300  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

Although the fruit supernatant is typically over 450 ~ 500  $\mu\text{l}$ , we recommend to recover only 300  $\mu\text{l}$  of lysate. Generally other plant sample take 200  $\mu\text{l}$  of supernatant. But, take more quantitative supernatant than other plant sample, because fruit sample contains an amount of water. When pipette the supernatant, please be careful without disturbing the cell-debris pellet. If much pellet is transferred with the supernatant, try one more centrifugation with the transferred supernatant. More lysate can result in shearing of the DNA and contaminating the next step with impurities. When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 975  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane. A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA Mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge at 13,000 rpm (RT) for 1 min and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm. After centrifugation, sometimes spin column membrane becomes light brown or light green color. Do not worry about that. Disappear after washing step.

**8. Repeat step 7 with remaining sample (maximum 625  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again at 13,000 rpm for 1 min.

**V. Washing step A****9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge for 1 min at 13,000 rpm. Discard the flow-through and reuse the collection tube in step 10.****VI. Washing step B****10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge for 1 min at 13,000 rpm.**

**Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE:** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.



### **VII. Elution Step**

- 11. Place the spin column into a new 1.5 ml tube (not supplied), and 50  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate at room temperature for 1 min , and then centrifuge for 1 min at 13,000 rpm to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but it will be reduced overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

Plant  
**Type F**  
*Bench Protocol*

**Type F Protocol**  
**For Seed(Others)**


F-1. Sample Treatment Step

▪ **Seed (Others)**


- I. Preparation step    II. Disrupt.& Homogen.    III. Sample Sizing step    IV. Pre-Treating step





I. Preparation step

**1**  Prepare dried or fresh seed sample


II. Disruption & Homogenization step


**2**  Place the sample into mortar

 Add liquid nitrogen for sample freezing

 Pouring liquid nitrogen slowly, Disrupt and homogenize completely

III. Sample Sizing step

**3**  Measure 50 mg of sample powder

 Transfer to 1.5 ml tube

F-2. DNA Extraction Step




Type F Protocol  
Seed (Other)




F-2. DNA Extraction Step

▪ **Seed (Others)**

- I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step
- V. Washing step A     VI. Washing step B     VII. Elution step

Type F Protocol  
Seed (Other)





Type F Protocol  
Seed (Other)

F-1. Sample Treatment Step

▪ **Seed (Others)**

■ I. Preparation step ■ II. Disrupt.& Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step

**I. Preparation step**

**1. Prepare dried or fresh seed sample.**

Plant seed sample corresponds to dormancy state at cycle of plant and consists embryo and albumen. Seed is kept in dry well in dry state, for albumen consist starch or fat. Even if seed keep at room temperature (15°C ~ 25°C) after dry, there is no problem when using later.

If not, store at -20°C or -80°C for long term storage. Before storage, to isolate pure DNA, wash the surface of the sample with distilled water.

If frozen seed sample, it is important to keep frozen state in liquid nitrogen during all of sample treatment step to inhibit low DNA yields and degraded DNA. Do not repeat freeze-thaw.

**II. Disruption & Homogenization step**

**2. Place the sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately to step 3.**

In seed sample, break outer shell of seed due to the hard outer shell after putting one, two or more dried or fresh seed sample in the mortar. It is difficult to disrupt and homogenize. Slowly pouring the liquid nitrogen and pressing the power, repeat again repeat again. We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Finally, It is a fine powder form after disruption and homogenization.

**III. Sample Sizing step**

**4. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

To prevent thawing the frozen sample during transferring it, previously pre-chill the spatula and 1.5 ml tube in liquid nitrogen. The freezing-thawing repetition of frozen sample will result in the DNA degradation. And more, exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity.

## F-2. DNA Extraction Step

**Seed (Others)**

I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step

- Equilibrate samples to room temperature (15 ~ 25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu\text{l}$  Buffer PG, 7  $\mu\text{l}$  Enhancer Solution, 20  $\mu\text{l}$  Proteinase K, and 5  $\mu\text{l}$  RNase A Solution into sample tube, and vortex vigorously.**

With seed sample absorbs lysis buffer, and becomes swollen. It may be difficult to handle plant tissue due to its viscosity. Always keep the recommended amount of starting material.

Furthermore, vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and therefore it will result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

**2. Incubate the lysate at 65°C for 30 min.**

For complete lysis, mix 5 ~ 6 times during incubation by inverting tube. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate.

Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased.

Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min at room temperature.**

Seed sample is contain much fat, so let you see thick fat layer on supernatant. Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step.

Type F Protocol  
Seed (Other)

#### **IV. DNA Binding step**

**5. Transfer carefully 200  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

As seed sample is containing much fats and it forms fatty layer, so it is difficult to take the supernatant in fatty layer. When pipetting the supernatant, must observe so that much fat should not stick to tip and do not take pellet together. If much pellet and fat are transferred with the supernatant together carelessly, debris prevent binding the column membrane.

Although the supernatant is typically over 250 ~ 300  $\mu\text{l}$ , we recommend to recover only 200  $\mu\text{l}$  of lysate. More lysate can results in shearing of the DNA and contaminating the next step with impurities. When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 to 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane.

A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA Mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge for 1 min at 13,000 rpm (RT), and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

After centrifugation, sometimes spin column membrane becomes light brown or light green color. Do not worry about that. It will be disappeared after washing step.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

#### **V. Washing step A**

**9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge for 1 min at 13,000 rpm. Discard the flow-through and reuse the collection tube in step 10.**

#### **VI. Washing step B**

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge at 13,000 rpm for 1 min .**

**Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

### **VII. Elution Step**

- 1. Place the spin column into a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate for 1 min at room temperature, and then centrifuge for 1 min at 13,000 rpm to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

Type F Protocol  
Seed (Other)



Plant

**Type G**

**Bench Protocol**


**Type G Protocol**

**For Seed (Gramineae)**

G-1. Sample Treatment Step

▪ **Seed (Gramineae)**

■ I. Preparation step   ■ II. Disrupt.& Homogen.   ■ III. Sample Sizing step   □ IV. Pre-Treating step






Type G Protocol  
Seed (gramineae)

## G-2. DNA Extraction Step




▪ **Seed (Gramineae)**

I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step

V. Washing step A     VI. Washing step B     VII. Elution step



Type G Protocol  
Seed (gramineae)



## G-1. Sample Treatment Step

▪ **Seed (Gramineae)**

I. Preparation step
  II. Disrupt.& Homogen.
  III. Sample Sizing step
  IV. Pre-Treating step

**I. Preparation step****1. Prepare dried or fresh gramineae seed sample.**

It is hardly difficult to get good result more than other seed sample as gramineae seed sample, it has too much secondary product and polyphenolic compounds. Even if seed keep at room temperature(15°C ~ 25°C) after drying, there is no problem when use later. If not, store at -20°C or -80°C for long-term storage.

If frozen seed sample, it is important to keep frozen state in liquid nitrogen during all of sample treatment step to inhibit low DNA yields and degraded DNA. Do not repeat freeze-thaw.

**II. Disruption & Homogenization step****2. Place the sample material into a grinding mortar. Add liquid nitrogen to the mortar and freeze. Keep the sample submerged in liquid nitrogen, and disrupt carefully until the sample is homogenized completely. Allow the liquid nitrogen to evaporate, and proceed immediately to step 3.**

Because seed sample of gramineae is thinner than other seed sample, easily grind sample. It's very important to keep the sample frozen in liquid nitrogen during disruption and then homogenization step to inhibit low DNA yields and degraded DNA. We recommend to be disrupted completely until no tissue clumps are not visible. Clumps of sample will not lyse properly and will therefore result in a lower yield of DNA. It's very important to keep the sample frozen in liquid nitrogen during disruption and homogenization step to inhibit low DNA yields and degraded DNA. Be careful to handle liquid nitrogen. Generally, It is a fine powder form after disruption and homogenization.

Type G Protocol  
Seed (gramineae)

**III. Sample Sizing step****3. Measure 10 mg of sample powder, and then transfer into 1.5 ml tube using a spatula.**

Although the thin seeds are easily ground, but rice, wheat etc are so thick to difficult to grind. Rice and wheat like this sample, when lysis incubation time, It absorbs all of Buffer PG (lysis buffer) to difficult to handle to lysis. And in that time, becomes swollen. To prevent it, reduce the amount of starting sample material. In gramineae seed sample, previously, use below 10 mg for amount of starting material. But, in case of grass sample, use 100 mg for amount of starting material. Exceeding the recommended optimal amount of starting material will result in inefficient lysis, resulting in low DNA yield and purity. To prevent thawing the frozen sample during transferring it, previously pre-chilling the spatula and 1.5ml tube in liquid nitrogen.

The freezing-thawing repetition of frozen sample will result in the DNA

## G-2. DNA Extraction Step

▪ **Seed (*Gramineae*)**

I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step

V. Washing step A     VI. Washing step B     VII. Elution step

- Equilibrate samples to room temperature (15~25°C).
- Heat a water bath or heating block to 65°C for use in step 2.
- All centrifugation steps should be carried out at room temperature.

**II. Lysis step****1. Add 390  $\mu\text{l}$  Buffer PG, 7  $\mu\text{l}$  Enhancer Solution, 20  $\mu\text{l}$  Proteinase K, and 5  $\mu\text{l}$  RNase A Solution into sample tube, and vortex vigorously.**

With seed sample absorbs lysis buffer, and becomes swollen. It may be difficult to handle plant tissue due to its viscosity. Always keep the recommended amount of starting material. Furthermore, vortex or pipette vigorously to remove any clumps until any plant tissue clumps are not visible. Clumps of plant tissue will not lyse adequately and will therefore result in a lower yield of DNA. A disposable micropestle may be used if clumps are not be removed by pipetting and vortexing.

**2. Incubate the lysate at 65°C for 30 min.**

Like rice, it becomes boiled rice when incubation for lysis. You must mix by tapping every 3 minutes for complete lysis. Be careful this. In case of grass sample, mix 5 ~ 6 times during incubation by inverting tube for complete lysis. The incubation time can be prolonged for more yields of DNA. The complete lysis let you see clear lysate.

**III. Precipitation step****3. Add 100  $\mu\text{l}$  Buffer PPT to the lysate, mix well, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides. During incubation on ice, please mix 5 ~ 6 times by inverting tube. The reaction makes clear lysate into opaque slushy lysate.

Generally, plant tissues contain large amounts of polysaccharides and polyphenolics, and are therefore relatively not easier to handle. These impurities may be present in the purified DNA if the amount of starting material is increased. Always keep the recommended amounts of samples.

**4. Centrifuge the lysate at 13,000 rpm (RT) for 5 min .**

In case of rice, although rice becomes the boiled rice, you can see the clear supernatant after centrifugation if you use 10 mg for amount of starting material. Plant materials can create very viscous lysates and large amounts of precipitates during this step. If you keep our recommended amounts of starting material, optimal results are obtained. If not, you should perform one more centrifugation step.

#### **IV. DNA Binding step**

**5. Transfer carefully 300  $\mu\text{l}$  of supernatant from step 4 into a new 1.5 ml tube.**

Although the supernatant is typically over 300 ~ 350  $\mu\text{l}$ , we recommend to recover only 200  $\mu\text{l}$  of lysate. More lysate can result in shearing of the DNA and contaminating the next step with impurities. When pipetting, please be careful without disturbing the cell-debris pellet.

**6. Add 650  $\mu\text{l}$  Buffer PB to the lysate, and mix well by gently inverting 5 ~ 6 times or by pipetting. DO NOT vortex.**

This step is an equilibration step for binding genomic DNA to column membrane.

A precipitate may form after the addition of Buffer PB, but this will not affect the i-genomic Plant DNA Mini Kit procedure.

**7. Pipette 650  $\mu\text{l}$  of the mixture from step 6, including any precipitate that may have formed, into the spin column inserted in a 2.0 ml collection tube. Centrifuge for 1 min at 13,000 rpm (RT), and discard the flow-through. Reuse the collection tube in step 8.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

After centrifugation, sometimes spin column membrane becomes light brown or light green color. Do not worry about that. Disappear after washing step.

**8. Repeat step 7 with remaining sample (maximum 200  $\mu\text{l}$ ). Discard flow-through and collection tube altogether.**

If a small amount will not pass through, please centrifuge again for 1 min at 13,000 rpm.

#### **V. Washing step A**

**9. Place the spin column into a new 2.0 ml collection tube (additionally supplied), add 700  $\mu\text{l}$  Buffer PWA, and centrifuge at 13,000 rpm for 1 min. Discard the flow-through and reuse the collection tube in step 10.**

#### **VI. Washing step B**

**10. Add 700  $\mu\text{l}$  Buffer PWB to the spin column, and centrifuge for 1 min at 13,000 rpm.**

**Discard the flow-through, and again centrifuge for additional 1 min to dry the membrane. Discard the flow-through and collection tube altogether.**

It is very important to dry the membrane of the spin column since residual ethanol may inhibit subsequent reactions. Following the centrifugation, remove carefully the spin column from the collection tube without contacting with the flow-through, since this will result in carryover of ethanol.

**NOTE :** Ensure that 40 ml of ethanol (EtOH) has been added to Buffer PWB.

### **VII. Elution Step**

- 11. Place the spin column into a new 1.5 ml tube (not supplied), and 100  $\mu\text{l}$  Buffer PE directly onto the membrane. Incubate for 1 min at room temperature, and then centrifuge at 13,000 rpm for 1 min to elute.**

Elution with 50  $\mu\text{l}$  (instead of 100  $\mu\text{l}$ ) increases the final DNA concentration, but reduces overall DNA yield conventionally. Alternatively, if you need larger amounts of DNA, eluting with 200  $\mu\text{l}$  increases generally overall DNA yield.

**NOTE :** A new 1.5 ml tube can be used for the second elution step to prevent dilution of the first eluate.

Alternatively, the tube can be reused for the second elution step to combine the eluates.

---

## Troubleshooting Guide

When working with i-genomic Plant DNA Extraction Mini Kit, always follow the description of each protocols. Nevertheless, if it causes problems upon extracting DNA, please refer to the following Troubleshooting Guide. This Troubleshooting Guide may be helpful in solving any problems that may arise. For more information, please contact our Technical Assistance Team. Our Technical Assistance Team is staffed by experienced researchers with extensive practical and theoretical expertise in molecular biology and the use of iNtRON products.

---

### Comments and Suggestions

---

#### ***Low flow rate in column***

---

- ✓ Clogged spin column by particulate material
  - (1) Completely perform the Disruption & Homogenization step.
  - (2) Increase the incubation time at 65°C in Lysis step.
  - (3) Ensure that no particulate material is transferred following precipitation step or when supernatants transferred to new 1.5 ml tube prior to addition of Buffer PB.
- ✓ High viscosity of Lysate
  - (1) Reduce the amounts of starting material.
  - (2) Increase the incubation time at 65°C in Lysis step.
- ✓ Problem in centrifugation
  - (1) Check your centrifuge, and then speed up or increase the centrifugation time.

---

#### ***Low DNA yield***

---

- ✓ Not enough disruption & homogenization
    - (1) Check again your Disruption & Homogenization step, and always follow the protocol. The complete disruption and homogenization of the starting material is very important to get high DNA yield.
  - ✓ Inadequate lysis
    - (1) Reduce the amounts of starting material.
    - (2) Increase the incubation time at 65°C in Lysis step.
  - ✓ Error in DNA binding
    - (1) After adding Buffer PB in DNA Binding step, please mix well by gently inverting.
    - (2) Check that the amount of Buffer PB is added correctly to the supernatant.
  - ✓ Incorrect Washing step
    - (1) Check again that the amount of ethanol (EtOH) is added correctly to Washing buffer.
    - (2) When storing Washing Buffer, always keep a lid shut tightly without evaporation.
-



continued from **Low DNA yield.**

- ✓ Insufficient DNA elution
  - (1) Increase the volume of Buffer PE or water to 200  $\mu\text{l}$ .
  - (2) Increase the incubation time on the column to 5 ~ 10 min at room temperature prior to centrifugation.

---

**Low DNA concentration**

---

- ✓ Excess addition of elution buffer
  - (1) Reduce the amount of Buffer PE.
  - (2) Increase the incubation time on the column to 5 ~ 10 min at room temperature prior to centrifugation.

---

**DNA sheared**

---

- ✓ Incorrect storage of plant tissue
  - (1) When store frozen tissues, always keep the samples frozen below  $-80^{\circ}\text{C}$ .
  - (2) If possible, it is preferable to use fresh plant tissues.
- ✓ Vigorously vortex
  - (1) Do not vortex the mixture after adding Buffer PB as described in protocol.
  - (2) Increase the incubation time on the column to 5 ~ 10 min at room temperature prior to centrifugation.
- ✓ Debris of precipitates in lysate
  - (1) Perform the optional centrifugation step before loading a large amount of the lysate onto the spin column.
  - (2) Always use the recommended amounts of starting material.

---

**Problems in downstream experiments**

---

- ✓ Ethanol contamination
  - (1) Ensure that during Washing Step B, the column membrane should be dried completely. Please centrifuge at full speed for 5 ~ 10 min to dry the membrane.
  - (2) During Washing Step B, after centrifugation, remove carefully the spin column from the collection tubes without contacting with the flow-through. This careless contact will result in contamination of ethanol.
- ✓ Salt contamination
  - (1) Check again to add EtOH previously into Buffer PWB.
  - (2) Store Buffer PWB at room temperature (15 ~  $20^{\circ}\text{C}$ ).
- ✓ Amount of DNA used in experiments.
  - (1) Optimize the amount of DNA used in your downstream experiments.

---

**For more questions, please contact us without hesitation.**

---

# ***Data Information***

## ***Determination of Yield and Purity***




## Appendix A

### Electrophoresis Results from Various Plant Samples

i-genomic Plant Mini Kit provides a reliable and practical method to purify efficiently genomic DNA from all kinds of plant tissues. DNA purified by i-genomic Plant Mini Kit is up to over 40 Kb, and has an  $A_{260/280}$  ratio of 1.7 ~ 1.9, indicating high purity of the DNA. The following Figure 3 shows overall electrophoresis data from representative various samples in each plant groups.

**Figure 3. Agarose Gel Electrophoresis of Eluted Genomic DNA (1.0 %)**



## Appendix B


### Determination of Yield and Purity Data of DNA(Standard Protocol Applied)

- Type A Protocol  
: Lyophilized Leaf

#### (A) DNA Yield and Purity

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Persimmon	5	3 ~ 5	1.73
2	Adray	5	2 ~ 4	1.83
3	Perilla	5	3 ~ 6	1.72
4	Millet (Hog millet)	5	3 ~ 6	1.79
5	Ginseng	5	1.5 ~ 4	1.88
6	Millet	5	5 ~ 9	1.82
7	Mung bean	5	2 ~ 5	1.81
8	Leak	5	2 ~ 5	1.85
9	Amaranth	5	2 ~ 5	1.77
10	Ginger	5	2 ~ 5	1.89
11	Corn	5	2 ~ 5	1.90
12	Paddy	5	2 ~ 5	1.88


#### (B) DNA Purification and Enzyme Digestion (RE)



#### Fig 4-1. Results of DNA purification and enzyme digestion with *EcoRI*

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

#### (C) Genomic DNA PCR



#### Fig 4-2. PCR Amplification


The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)

**• Type B Protocol****: Fresh, Dried, or Frozen Leaf****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Spinach	50	8 ~ 9	1.85
2	Arabidopsis	50	8 ~ 11	1.81
3	Rhododendron	50	8 ~ 10	1.83
4	Lettuce	50	10 ~ 12	1.79
5	Cabbage ( <i>Brassica oleracea</i> )	50	8 ~ 10	1.78
6	Green pepper	50	8 ~ 12	1.83
7	African millet	50	8 ~ 10	1.81
8	Peanut	50	8 ~ 12	1.81
9	Buckwheat	50	10 ~ 12	1.83
10	Sesame	50	8 ~ 10	1.77
11	Cabbage ( <i>Brassica campestris</i> )	50	8 ~ 10	1.85
12	Radish	50	10 ~ 12	1.79
13	Black Bean leaves	50	8 ~ 10	1.75
14	Geranium	50	8 ~ 10	1.77
15	Fallen leaves	50	2 ~ 5	1.82
16	Moss	50	2 ~ 4	1.85

**(B) DNA Purification and Enzyme Digestion (RE)**


M	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
RE	-	+	+	-	+	-	+	-	+	-	+	-	+	-	+	-


**Fig 5-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR**

M	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16


**Fig 5-2. PCR Amplification**


The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (*i-StarTaq*, Cat. 25165) for PCR amplification reaction.

**Additional Data****(Type B protocol-based modification)**

When we have been developing i-genomic Plant DNA Mini Kit, iNtRON customers requested to be tested previously in their various samples for DNA yield, purity and PCR amplification.

The additional data are included in the following section. We have performed to extract DNA from various plant leaves supplied by customers with Type B protocol (occasionally, we slightly have modified Type B protocol). You can show good results from various plant leaves.

**Fig 5-3. Additional Data (sample morphology)**






**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Paddy	50	4 ~ 6	1.79
2	Cucumber	50	10 ~ 11	1.87
3	Chinese matrimony vine (type 1)	50	8 ~ 9	1.82
4	Chinese matrimony vine (type 2)	50	10 ~ 11	1.79
5	Japanese lady bell	42.2	9 ~ 12	1.82
6	Phyteuma	10	3 ~ 4.2	1.79
7	Adenophora	10	5 ~ 6	1.75
8	Bell flower	5	2 ~ 4	1.75
9	Platycodon	80	7 ~ 9	1.77
10	Sea weed	50	8 ~ 10	1.78
11	Pine needles	100	12 ~ 15	1.70
12	Sweet potato	50	12 ~ 13	1.87
13	Ginko	100	7 ~ 8	1.75
14	Korea box tree	100	9 ~ 11	1.77
15	Rhododendron	300	8 ~ 9	1.76
16	Cabbage ( <i>Brassica oleracea</i> )	100	10 ~ 12	1.74
17	Forsythia	50	11 ~ 13	1.88

**(B) Genomic DNA PCR**



**Fig 5-4. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction.

(Lane M : 100 bp DNA Ladder, Cat. 24012)


• **Type C Protocol**  
**: Stem**

**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Potato	50	6 ~ 8	1.80
2	Onion	50	4 ~ 6	1.81
3	Elm	50	4 ~ 6	1.74
4	Geranium	50	4 ~ 5	1.78
5	Radish	50	5 ~ 6	1.78
6	Rhododendron	50	4 ~ 6	1.85
7	Sesame	50	4 ~ 7	1.80

**(B) DNA Purification and Enzyme Digestion (RE)**

M	1	2	3	4	5	6	7	
RE	-	+	-	+	-	+	-	+




**Fig 6-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR**

M	1	2	3	4	5	6	7



**Fig 6-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng). We used Maxime PCR Premix Kit (*i-StarTaq*, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)




**Additional Data****(Type C protocol-based modification)**

When we have been developing i-genomic Plant DNA Mini Kit, iNtRON customers requested to be tested previously in their various samples for DNA yield, purity and PCR amplification.

The additional data are included in the following section. We have performed to extract DNA from various plant stems supplied by customers with Type C protocol (occasionally, we slightly have modified Type C protocol). You can show good results from various plant stems.

**Fig 6-3. Additional Data (sample morphology)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Agelica	10	5 ~ 7	1.82
2	Potato	100	8 ~ 10	1.75
3	Cactus	100	4 ~ 6	1.75
4	Sea weed fusiforme	30	3 ~ 5	1.72

**(B) Genomic DNA PCR****Fig 6-4. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)


▪ **Type D Protocol**

: **Root**

**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Platycodon	50	5 ~ 7	1.83
2	Radish	50	6 ~ 10	1.79
3	Sweet potato	50	5 ~ 9	1.82
4	Carrot	50	8 ~ 10	1.85
5	Sesame	50	5 ~ 8	1.81


**(B) DNA Purification and Enzyme Digestion (RE)**



**Fig 7-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR**




**Fig 7-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)


**Additional Data****(Type D protocol-based modification)**

When we have been developing i-genomic Plant DNA Mini Kit, iNtRON customers requested to be tested previously in their various samples for DNA yield, purity, and PCR amplification.

The additional data are included in the following section. We have performed to extract DNA from various plant roots supplied by customers with Type D protocol (occasionally, we slightly have modified Type D protocol). You can show good results from various plant roots.

**Fig 7-3. Additional Data (sample morphology)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Ginseng (Type 1)	50	5 ~ 7	1.78
2	Ginseng (Type 2)	50	3 ~ 4	1.80
3	Ginseng (Type 3)	50	4 ~ 6	1.81

**(B) Genomic DNA PCR****Fig 7-4. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction.

(Lane M : 100 bp DNA Ladder, Cat. 24012)


• **Type E Protocol**

: **Fruit**

**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Muscat	100	3 ~ 4	1.87
2	Tomato	100	3 ~ 5	1.82
3	Pineapple	100	3 ~ 5	1.79
4	Apple	100	2 ~ 4	1.83
5	Melon	100	2 ~ 4	1.75
6	Peach	100	2 ~ 5	1.83
7	Grape	100	2 ~ 5	1.85


**(B) DNA Purification and Enzyme Digestion (RE)**



**Fig 8-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 50  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR**



**Fig 8-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)

**Additional Data**

**(Type E protocol-based modification)**

When we have been developing i-genomic Plant DNA Mini Kit, iNtRON customers requested to be tested previously in their various samples for DNA yield, purity and PCR amplification. The additional data are included in the following section. We have performed to extract DNA from various plant fruits supplied by customers with Type E protocol (occasionally, we slightly have modified Type E protocol). You can show good results from various fruits.


**Fig 8-3. Additional Data (sample morphology)**



**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Tomato	50	2 ~ 3	1.83
2	Green pepper	50	2 ~ 4	1.85

**(B) Genomic DNA PCR**



**Fig 8-4. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)


• **Type F Protocol**

[1] **Seed (*Leguminosae*)**

**(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Pea	50	12 ~ 14	1.87
2	Black bean	50	10 ~ 14	1.88
3	Yellow bean	50	12 ~ 15	1.84
4	Mung bean	50	10 ~ 12	1.86


**(B) DNA Purification and Enzyme Digestion (RE)**

M	1		2		3		4	
RE	-	+	-	+	-	+	-	+
								

**Fig 9-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR**


M	1	2	3	4	5	6	7
							

**Fig 9-2. PCR Amplification**


The housekeeping gene (18s, 222 bp) was amplified with the purified DNAs as template DNA (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)

**[2] Seed (*Solanaceae*)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Green pepper	50	8 ~ 10	1.85
2	Paprika	50	10 ~ 12	1.85
3	Tomato type(normal)	50	9 ~ 12	1.78
4	Tomato JPN	50	8 ~ 12	1.77
5	Petunia	50	8 ~ 10	1.83

**(B) DNA Purification and Enzyme Digestion (RE)****Fig 10-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR****Fig 10-2. PCR Amplification**


The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng).

We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction.


(Lane M : 100 bp DNA Ladder, Cat. 24012)

**[3] Seed (Cruciferae)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Radish type 1 ( <i>Raphanus sativus</i> )	50	8 ~ 10	1.78
2	Radish type 2 ( <i>Raphanus sativus</i> )	50	10 ~ 12	1.84
3	Radish type3 ( <i>Brassica rapa</i> )	50	10 ~ 12	1.81

**(B) DNA Purification and Enzyme Digestion (RE)****Fig 11-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)


**(C) Genomic DNA PCR****Fig 11-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)




**[4] Seed (*Cucurbitaceae*)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Green pepper	50	8 ~ 10	1.85
2	Paprika	50	10 ~ 12	1.85
3	Tomato type(normal)	50	9 ~ 12	1.78
4	Tomato JPN	50	8 ~ 12	1.77
5	Petunia	50	8 ~ 10	1.83

**(B) DNA Purification and Enzyme Digestion (RE)****Fig 12-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA marker, Cat. 24022)

**(C) Genomic DNA PCR****Fig 12-2. PCR Amplification**


The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng).

We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction.


(Lane M : 100 bp DNA Ladder, Cat. 24012)

**[5] Other Seeds****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Spinach	50	8 ~ 10	1.74
2	Sunflower	50	8 ~ 10	1.75
3	Carrot	50	5 ~ 8	1.82
4	Sesame	50	10 ~ 12	1.77

**(B) DNA Purification and Enzyme Digestion (RE)****Fig 13-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

**(C) Genomic DNA PCR****Fig 13-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)


### •Type G Protocol

:Seed (*Gramineae*)

#### (A) DNA Yield and Purity

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Wheat	10	4 ~ 5	1.75
2	Grass	100	4 ~ 6	1.82
3	Barley	10	5 ~ 8	1.79
4	Rice	10	6 ~ 8	1.85
5	African millet	10	4 ~ 8	1.81


#### (B) DNA Purification and Enzyme Digestion (RE)



**Fig 14-1. Results of DNA purification and enzyme digestion with *EcoRI***

After eluting genomic DNA with 100  $\mu\text{l}$  Buffer PE, each 100 ng of DNA were used in DNA electrophoresis and enzyme digestion, respectively. Restriction enzyme digestion was performed with 5 units of *EcoRI* at 37°C for 1 hour. (Lane M : 1 Kb DNA Ladder, Cat. 24022)

#### (C) Genomic DNA PCR



**Fig 14-2. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNA as template DNA (10 ng).


We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction.

(Lane M : 100 bp DNA Ladder, Cat. 24012)


**Additional Data****(Type G protocol-based modification)**

When we have been developing i-genomic Plant DNA Mini Kit, iNtRON customers requested to be tested previously in their various samples for DNA yield, purity and PCR amplification.

The additional data are included in the following section. We have performed to extract DNA from various plant seeds supplied by customers with Type G protocol (occasionally, we slightly have modified Type G protocol). You can show good results from various plant seeds.

**Fig 10-3. Additional Data (sample morphology)****(A) DNA Yield and Purity**

Lane	Samples	Amounts (mg)	DNA Yield ( $\mu\text{g}$ )	$A_{260/280}$
1	Unpolished rice	50	3 ~ 5	1.73
2	Rice (1 of grain)	17.5	1 ~ 2.5	1.81
	Rice	50	2 ~ 3	1.75
3	Rice seed (1 of grain)	25.6	2 ~ 4	1.82
	Rice seed	250	2 ~ 4	1.71

**(B) Genomic DNA PCR****Fig 10-4. PCR Amplification**

The housekeeping gene (18s, 222 bp) was amplified with the purified DNAs as templates (10 ng). We used Maxime PCR Premix Kit (i-StarTaq, Cat. 25165) for PCR amplification reaction. (Lane M : 100 bp DNA Ladder, Cat. 24012)

## Appendix C : Source List of Genomic DNA

Scientific Name	Common Name			
	Korean	English	Japanese	Chinese
<i>Gallus domesticus</i>	닭	fowl	ニワトリ	雞
<i>Bos taurus</i>	소	Cow	うし	牛
<i>Canis familiaris</i>	개	Dog	いぬ	犬
<i>Drosophilidae</i>	초파리	Drosophila	ショウジョウバエ	醃鷄
<i>Lumbricidae</i>	지렁이	Earth worm	ミミズ	蚓
<i>Cavia porcellus</i>	기니피그	Guinea pig	モルモット	天竺鼠
<i>Homo sapiens</i>	인간	Human	にんげん	人間
<i>Gryllotalpa orientalis</i>	땅강아지	Mole cricket	ケラ	东方蝼蛄
<i>Culicidae</i>	모기	Mosquito	か	蚊
<i>Mus musculus</i>	생쥐	Mouse	マウス	小家鼠
<i>Sus domesticus</i>	돼지	Pig	ブタ	豚
<i>Oryctolagus cuniculus</i>	토끼	Rabbit	ウサギ	兔
<i>Rattus</i>	쥐	Rat	ねずみ	鼠
<i>Pandalus hypsinotus</i>	도화새우	Humpback shrimp	ボタンエビ	牡丹蝦
<i>Zebra fish</i>	제브라피시	Zebra fish	ゼブラフィッシュ	狂放类型线
<i>Bladder</i>	방광	Bladder	ぼうこう	膀胱
<i>Brain</i>	뇌	Brain	のう	腦
<i>Bronchus</i>	기관지	Bronchus	きかんし	氣管支
<i>Hair</i>	털	Hair	ヘア	毛
<i>Head</i>	머리	Head	あたま	頭
<i>Heart</i>	심장	Heart	しんぞう	心臟
<i>Kidney</i>	신장	Kidney	じんぞう	腎臟
<i>Liver</i>	간	Liver	レバ	肝
<i>Lung</i>	폐	Lung	はい	肺
<i>Muscle</i>	근육	Muscle	きんにく	筋肉
<i>Ovary</i>	난소	Ovary	らんそう	卵巢
<i>Pancreas</i>	췌장	Pancreas	すいぞう	胰臟
<i>Placenta</i>	태반	Placenta	たいばん	胎盤
<i>Skin</i>	피부	Skin	ひふ	皮膚
<i>Spleen</i>	비장	Spleen	ひぞう	脾臟
<i>Stomach</i>	위	Stomach	い	胃
<i>Subcutaneous fat</i>	피하지방	Subcutaneous fat	ひかしぼう	皮下脂肪
<i>Tail</i>	꼬리	Tail	テ-ル	尾
<i>Thymus</i>	흉선	Thymus	きょうせん	胸腺

<i>Adenophora remotiflora</i>	모시대	Adenophora	ソバナ	莽莨
<i>Adenophora triphylla</i>	잔대	Japanese lady bell	つりがねこんじん	山羊百科
<i>Agelica gigas</i>	당귀	Agelica	オニノダケ	當歸
<i>Allium cepa</i>	양파	Onion	タマネギ	洋葱
<i>Allium tuberosum</i>	부추	Korean leek	ニラ	韭
<i>Amarantus</i>	아마란스	Amaranth	ガンライコウ	反枝苋
<i>Ananas comosus</i>	파인애플	Pine apple	パイナップル	菠萝
<i>Arabidopsis</i>	애기장대	Arabidopsis	シロイヌナズナ	拟南芥
<i>Arachis hypogaea</i>	땅콩	Peanut	らっかせい	花生
<i>Brassica campestris</i>	배추	Cabbage	ハウサイ	菘
<i>Brassica oleracea</i>	양배추	Cabbage	タマナ	洋菘
<i>Brassica rapa</i>	순무	Radish (Turnip)	カブ	菁
<i>Campanula punctata</i>	초롱꽃	Bell flower	ホタルブクロ	山小菜
<i>Capsicum annuum</i>	고추	Green pepper	トウガラシ	苦椒
<i>Capsicum annuum</i>	파프리카	Paprika	ピマン	甘唐芥子
<i>Citrullus vulgaris</i>	수박	Water melon	スイカ	水瓜
<i>Coix Lachrymajobi</i>	율무	Adlay	ハトムギ	薏
<i>Cucumis melo</i>	멜론	Melon	メロン	甜瓜
<i>Cucumis melo var. makuwa</i>	참외	Melon	マクワウリ	香瓜
<i>Cucumis sativus</i>	오이	Cucumber	キウリ	黄瓜
<i>Cucurbita moschata</i>	호박	Pumpkin	カボチャ	南瓜
<i>Daucus carota var. sativa</i>	홍당무	Carrot	ニンジン	紅蘿蔔
<i>Diospyros kaki</i>	감	Persimmon	かき	柿
<i>Ecklonia cava</i>	감태	Sea wood	カジメ	槁布
<i>Fagopyrum esculentum</i>	메밀	Buckwheat	ソバ	蕎
<i>Fallen leves</i>	낙엽	Fallen leves	らくよう	落葉
<i>Forsythia korean</i>	개나리	Forsythia	チョウセンレン	ギョウ
<i>Ginkgo biloba</i>	은행나무	Ginkgo	いちょう	銀杏
<i>Glycine max</i>	검은콩	Black bean	くろまめ	黒大豆
<i>Glycine max</i>	메주콩	Yellow bean	ダイズ	大豆
<i>Helianthus annuus</i>	해바라기	Sunflower	ヒマワリ	向日葵
<i>Hizikia fusiforme</i>	툇	Sea weed fusiforme	かいそう	海藻
<i>Hordeum vulgare</i>	보리	Barley	ムギ	麥
<i>Ipomoea batatas</i>	고구마	Sweet potato	サツマイモ	甘薯
<i>Lactuca sativa var. capitata</i>	양상추	Lettuce	ガーデンレタス	莴苣

<i>Lagenaria leucantha</i>	박	Gourd	ウリ科	扁蒲
<i>Lycopersicon esculentum</i>	토마토	Tomato	トマト	番茄
<i>Lycopersicon esculentum JPN</i>	토마토 JPN	Tomato JPN	ソバナ	番茄
<i>Lyticum chinense</i>	구기자	Chinese matrimony vine	クコ	枸杞子
<i>Malus pumila</i>	사과	Apple	リンゴ	沙果
Moss	이끼	Moss	こけ	苔
<i>Opuntia ficus-indica</i>	선인장	Cactus	サボテン	仙人掌
<i>Oryza sativa</i>	벼	Paddy	いね	稻
<i>Oryza sativa</i>	벼	Rice	いね	米
<i>Oryza sativa</i>	현미	Unpolished rice	げんまい	玄米
<i>Panax ginseng</i>	인삼	Ginseng	こうらいこんじん	人蔘
<i>Panicum millacerm</i>	기장	Millet (Hog millet)	キビ	粱
<i>Pelargonium inquinans</i>	제라늄	Geranium	ペラルゴニウム	小花天竺葵
<i>Perilla frutescens</i>	들깨	Green perilla	エゴマ	紫苏
<i>Petunia hybrida</i>	피튜니아	Petunia	サントリー「サフィニア	矮牽牛
<i>Phasedus radiatus</i>	녹두	Mung bean	リョクトウ	綠豆
<i>Phyteuma japonicum</i>	영아자	Phyteuma	キキョウ科 シデシヤジン	牧根草属
<i>Pinus densiflora</i>	솔잎	Pine needles	まつば	松葉
<i>Pisum sativum</i>	완두콩	Pea	エンドウ	豌豆
<i>Platycodon grandiflorum</i>	도라지	Platycodon	キキョウ	桔梗
<i>Prunus persica</i>	복숭아	Peach	もも	桃
<i>Raphanus sativus</i>	무	Radish	だいこん	萝卜
<i>Rhododendron schlippenbachii</i>	철쭉	Rhododendron	クロフネツツジ	灌木
<i>Sesamum indicum</i>	참깨	Sesame	ごま	胡麻
<i>Setania italica</i>	조	Italian millet	アワ	粟
<i>Solanum tuberosum</i>	감자	Potato	ジャガタライモ	马铃薯
<i>Sorghum bicolor</i>	수수	African millet	モロコシ	高粱
<i>Spinacia oleracea</i>	시금치	Spinach	ホウレンソウ	菠菜
<i>Triticum aestivum (vulgare)</i>	밀	Wheat	コムギ	小麥
<i>Ulmus davidiana</i>	느릅나무	Elm	ハレニレ	榆
<i>Vitis vinifera L.</i>	포도	Grape	ぶどう	葡萄
<i>Vitis spp</i>	청포도	Muscat	あおぶどう	青葡萄
<i>Zea mays</i>	옥수수	Corn	とうもろこし	玉米
<i>Zingiber officinale</i>	생강	Ginger	ショウガ	生薑
<i>Zoysia japonica</i>	잔디	Grass	しば	结缕草

## Appendix D

### Summarized Photo Procedure

This Summarized Photo Procedure is a representative figure for Type B Protocol in i-genomic Plant DNA Extraction Mini Kit Handbook v1.0. This figure is only reference for beginners, and therefore please follow your set-up method. For more detailed information, please contact us.

#### B-1. Sample Treatment Step

##### ▪ Fresh Leaf

■ I. Preparation step ■ II. Disrupt.& Homogen. ■ III. Sample Sizing step □ IV. Pre-Treating step

#### I. Preparation step

1. Prepare fresh leaf sample.



3. Add liquid nitrogen to the mortar and freeze. Disrupt and homogenize, pouring liquid nitrogen slowly.




#### II. Disruption & Homogenization step

2. Slice off the prepared sample to small pieces by the scalpel or scissor and place them the mortar.



#### III. Sample Sizing step

4. Measure 50 mg of sample powder, and then transfer into 1.5 ml tube.



See next page for DNA Extraction step.



B-2. DNA Extraction Step

- I. Pre-Lysis step     II. Lysis step     III. Precipitation step     IV. DNA Binding step  
 V. Washing step A     VI. Washing step B     VII. Elution step

**II. Lysis step**

1. Add 390  $\mu\text{l}$  of Buffer PG,  
7  $\mu\text{l}$  of Enhancer Solution,  
20  $\mu\text{l}$  of Proteinase K,  
5  $\mu\text{l}$  of RNase A Solution.



2. Incubate the lysate at 65°C for 30 min.  
(invert 5 ~ 6 times for more lysis)



Vortex vigorously.



**III. Precipitation step**

3. Add 100  $\mu\text{l}$  of Buffer PPT.



Then, mix well.



Incubation for 5 min on ice.



4. Centrifuge at 13,000 rpm at RT for 5 min.



**IV. DNA Binding step**

5. Transfer 200  $\mu\text{l}$  of supernatant to a new 1.5 ml tube.



After centrifugation



*See next page.*

6. Add 650  $\mu\ell$  of Buffer PB and mix by inverting or pipetting. DO NOT vortex !



Discard flow-through.



8. Transfer the rest of mixture to spin column (maximum 200  $\mu\ell$ ).



7. Transfer 650  $\mu\ell$  of the mixture to spin column. Do not remove the rest of mixture.



Centrifuge at 13,000 rpm (RT) for 1 min.



Centrifuge at 13,000 rpm (RT) for 1 min .



Discard flow-through.



**V. Washing step A**

9. Place spin column into a new 2.0 ml collection tube.



Add 700  $\mu\text{l}$  of Buffer PWA to the spin column.



Centrifuge at 13,000 rpm (RT) for 1 min .



Discard flow-through.



**VI. Washing step B**

10. Add 700  $\mu\text{l}$  of Buffer PWB to column.

**NOTE:** Ensure that 40 ml of ethanol has been added to Buffer PWB.



Centrifuge at 13,000 rpm (RT) for 1 min .



Discard flow-through.



Centrifuge at 13,000 rpm (RT) for 1 min .



*See next page.*

### **VII. Elution step**

11. Transfer the spin column to a new 1.5 ml tube.



Centrifuge at 13,000 rpm (RT) for 1 min .



Add 100  $\mu$ l of Buffer PE directly onto the membrane.



Discard the spin column.



Incubate at RT for 1 min .



Finished DNA Extraction.



## Ordering Information

Product Name	Samples	Examples
CAT. NO. 17341 (50 Columns) i-genomic CTB DNA Extraction Mini Kit	<u>C</u> ells <u>T</u> issues Gram(-) <u>B</u> acteria	Human cultured cells / Mouse cultured cells / Mouse / Guinea pig / Rabbit / Chicken / Zebra fish / Shrimp / Pig / Insect / Animal hair / Worm / Stool / Buccal swab / Gram(-) bacteria / Others
CAT. NO. 17351 (50 Columns) i-genomic Blood DNA Extraction Mini Kit	Blood	Whole blood / Buffy Coat / Dried Spot / Blood Swab / Plasma / Serum / Others
CAT. NO. 17361 (50 Columns) i-genomic BYF DNA Extraction Mini Kit	Gram(+) <u>B</u> acteria <u>Y</u> east <u>F</u> ungi	<i>Azotobacter sp.</i> <i>Staphylococcus sp.</i> <i>Saccharomyces sp.</i> <i>Aspergillus sp.</i> Others
CAT. NO. 17371 (50 Columns) i-genomic Plant DNA Extraction Mini Kit	Plant	Leaf / Root / Stem / Fruit / Seed / Others

**Larger kit sizes available; please inquire.**

## Related Products

Product Name	Cat. No.	Size
<i>Maxime</i> <sup>TM</sup> PCR PreMix (i-Taq)	25025 / 25026	96 Tubes / 480 Tubes
<i>Maxime</i> <sup>TM</sup> PCR PreMix (i-StarTaq)	25165 / 25167	96 Tubes / 480 Tubes
<i>Maxime</i> <sup>TM</sup> PCR PreMix (i-pfu)	25185	96 Tubes
<i>Maxime</i> <sup>TM</sup> PCR PreMix (i-MAX II)	25265	96 Tubes
i-Taq <sup>TM</sup> DNA Polymerase	25021 / 25022	250 Units / 500 Units
i-StarTaq <sup>TM</sup> DNA Polymerase	25161 / 25162	250 Units / 500 Units
i-MAX <sup>TM</sup> II DNA Polymerase	25261	250 Units
i-pfu DNA Polymerase	25181	250 Units
easy-Labeler <sup>TM</sup> Random Primed DNA Labeling PreMix for dCTP	16043	30 Rxns
PROBER <sup>TM</sup> Probe DNA Purifying System	17072	100 Columns
MicroPROBER <sup>TM</sup> Probe DNA Purifying System	17073	50 Columns

## Global Distributors

---

Australia  
Scientifix  
Unit C3 Adamco Business Park 2A Westall Rd, Clayton, VIC  
Phone : +61-3-9562-7488  
Fax : +61-3-9548-7177  
URL : [www.scientifix.com.au](http://www.scientifix.com.au)

Belgium  
Goffin Meyvis Analytical & Medical Systems  
Park Rozendal Terhulpensesteenweg 6a 1560 Hoeilaart  
Phone : +32-2-769-7474  
Fax : +32-2-769-7499  
URL : [www.goffinmeyvis.com](http://www.goffinmeyvis.com)

China  
Innovating Santaibio INC.  
Room 322 Building A, Shengdi Mansion,  
No.29 Xueyuan Road,  
Haidian District, Beijing 100083  
Phone : +86-10-8232-3705, +86-10-8232-3706  
Fax : +86-10-8232-3707  
URL : [www.santaibio.com](http://www.santaibio.com)

China  
Chinagen INC.  
FL1, Building one, 2nd Industrial District, Kefa Road,  
Science&Technology Park Central,  
Shenzhen 518057  
Phone : +852-2646-5368  
Fax : +852-2646-5037  
URL : [www.chinagen.com.cn](http://www.chinagen.com.cn)

China  
Tech Dragon Limited  
Rm 808, Topsail Plaza, 11 On Sum Street,  
Shatin, N.T., Hong Kong  
Phone : +86-755-2601-4525  
Fax : +86-632-418-5877  
URL : [www.techdragon.co.hk](http://www.techdragon.co.hk)

Hungary  
Kasztel-Med Ltd.  
1118 Budapest Kelenhegyi 22  
Phone : +36-1-385-3887  
Fax : +36-1-381-0695  
URL : [www.kasztel.com](http://www.kasztel.com)

India  
Biogene  
B-41, 1st Floor, Moti Nagar, New Delhi-110015  
Phone : +91-11-2592-0047/48  
Fax : +91-11-2592-0047  
URL : [www.biogene-india.com](http://www.biogene-india.com)

Israel  
Talron Biotech Ltd.  
17 Hazait street Rehovot 76349  
Phone : +97-2-8-947-2563  
Fax : +97-2-8-947-1156  
URL : [www.talron.co.il](http://www.talron.co.il)

Indonesia  
PT.Blue Sky Biotech  
Blok J 3 No.46 Taman Tekno BSD  
Serpong-Tangerang 15314  
Phone : +62-21-7028-3685  
Fax : +62-21-756-5165

Italy  
Li StarFISH S.r.l  
Via Silvo Pellico, 3. 20061 Carugate (MI)  
Phone : +39-02-9215-0794  
Fax : +39-02-9215-7285  
URL : [www.listarfish.it](http://www.listarfish.it)

Japan  
Cosmo Bio Co. Ltd.  
Toyo-Ekimae Bldg, 2-20, Toyo 2-Chome, Koto-ku,  
Tokyo, 135-0016  
Phone : +81-3-5632-9610/9620  
Fax : +81-3-5632-9619  
URL : [www.cosmobio.co.jp](http://www.cosmobio.co.jp)

Japan  
Funakoshi Co. Ltd.  
9-7, Hongo 2-Chome, Bunkyo-ku, Tokyo, 113-0033  
Phone : +81-3-5684-1620  
Fax : +81-3-5684-1775  
URL : [www.funakoshi.co.jp](http://www.funakoshi.co.jp)

Japan  
Nacalai Tesque, Inc.  
272 higashitsuchikawa-cho, Kuze Higashi, Minami-ku,  
Kyoto, 601-8204  
Phone : +81-75-251-1723  
Fax : +81-75-251-1762  
URL : [www.nacalai.co.jp](http://www.nacalai.co.jp)

Malaysia  
NHK BIOSCIENCE SOLUTIONS SDN. BHD.  
4 Lorong Checkor, 3rd mile off Jalan Kelang  
Lama, 58000, Kuala Lumpur  
Phone : +603-7987-8218  
Fax : +603-7987-8213

Netherlands  
Goffin Meyvis Analytical & Medical Systems  
Ecustraet 11 4879 NP Etten-Leur  
Phone : +31-76-508-6000  
Fax : +31-76-508-6086  
URL : [www.goffinmeyvis.com](http://www.goffinmeyvis.com)

Pakistan  
Lab Line  
98, National Medicine Market, Arambagh Road,  
Karachi 74200  
Phone : +92-21-262-6097  
Fax : +92-21-262-2194

Philippines  
Hebborn Analytics INC.  
A-7 Lester building, Quirino Highway, Lagro Quezon city  
Phone : +632-461-7173  
Fax : +632-418-5877

Singapore  
ProGen Scientific  
32 Maxwell Road #03-13, Whitehouse Singapore 069115  
Phone : +65-9557-0665  
Fax : +65-6363-4381

Spain  
LABOTAQ, S.C  
C1 Aviacion, 5 Poligono Industrial Calonge  
41007 Seville  
Phone : +34-954-31-7216  
Fax : +34-954-31-7360  
URL : [www.labotaq.com](http://www.labotaq.com)

Thailand  
Pacific Science Co. Ltd.  
90 Soi Charansanitwong 49/1 Charansanitwong rd,  
Bandbunru Bangplad, Bangkok 10700  
Phone : +66-2-433-0068/9  
Fax : +66-2-434-2609  
URL : [www.Pacificscience.co.th](http://www.Pacificscience.co.th)

Turkey  
MedSanTek Ltd. Co.  
Millet Cad. Findikzade sok. Emre Apt. No : 19/5  
Findikzade, Istanbul  
Phone : +90-212-635-8546  
Fax : +90-212-635-8350  
URL : [www.Medsantek.com.tr](http://www.Medsantek.com.tr)

Taiwan  
Hong-jing Co., Ltd.  
5F, No. 174, CHUNG-SHAN Rd., Sec 1,  
YOUNG-HO CITY, Taipei  
Phone : +886-2-3233-8585  
Fax : +886-2-3233-8686

U.A.E.  
Emphor FZCO.  
Shed No. MB-6, P.O. Box 61232  
Jebel Ali Free Zone, Dubai  
Phone : +971-4-883-0233  
Fax : +971-4-883-0133

United Kingdom  
Molecular Solutions Europe Ltd.  
2nd floor LPL, 145-157 St. John Street,  
London, EC1V 4PY?Phone : +44-0870-199-5067  
Fax : +44-0207-253-9040  
URL : <http://www.mseu.co.uk>

U. S. A.  
Boca Scientific  
1140 Holland Drive, Suite 11, Boca Raton, FL 33487  
Phone : +1-561-995-5017  
Fax : +1-561-995-5018  
URL : [www.bocascientific.com](http://www.bocascientific.com)

U. S. A.  
Lambda Biotech  
3830 Washington Avenue, St. Louis, Mo 63108  
Phone : +1-314-533-8650  
Fax : +1-314-533-4007  
URL : [www.lambdabiotech.com](http://www.lambdabiotech.com)

U. S. A.  
Nacalai USA, Inc.  
6640 Lusk Blvd, Suite A 200,  
San Diego, CA 92121  
Phone : +1-858-404-0403  
Fax : +1-858-404-0408



## **Note**

---

iNtRON Biotechnology, Inc.

***Plant***

i-genomic BYF DNA Extraction Mini Kit

