

FavorPrepTM Plasmid DNA Extraction Mini Kit

User Manual

**Cat. No.: FAPDE 001 (100 Preps)
FAPDE 001-1 (300 Preps)**

For Research Use Only

v.1103

Introduction

FavorPrep Plasmid Extraction Mini Kit is an excellent tool offering a speed and economic method to purify plasmid DNA from bacteria cultures. This technology is based on binding DNA to silica-based membranes in chaotropic salts, washing DNA with ethanol-contained Wash Buffer. Compare with other harmful and time-consuming procedure, such as phenol/ chloroform extraction and ethanol precipitation, FavorPrep Plasmid extraction kit shortens the handling time to about 25 minutes. The high quality plasmid DNA can be used directly for the downstream application.

Specification

Sampling: 1~5 ml overnight culture

Plasmid Size: < 12Kb

Yield: 20~30 µg of high-copy plasmid

Handling time: about 25 min

Kit Contents

	FAPDE 001	FAPDE 001-1
FAPD1 Buffer	30 ml	90 ml
FAPD2 Buffer	30 ml	90 ml
FAPD3 Buffer	40 ml	120 ml
W1 Buffer (concentration)*	35 ml	98 ml
Wash Buffer (concentration)**	20 ml	50 ml
Elution Buffer	15 ml	35 ml
RNase A (50mg/ml)	60 µl	180 µl
FAPD Column	100 pcs	300 pcs
Collection Tube	100 pcs	300 pcs
User Manual	1	1


* Add 13 ml/ 36 ml ethanol (96 ~ 100%) to W1 Buffer when first open.

** Add 80 ml/ 200 ml ethanol (96 ~ 100%) to Wash Buffer when first open.

Important Notes

1. Buffer provided in this kit contain irritants. Wear gloves and lab coat when handling these buffer.
2. Brief spin RNase A tube to remove drops from the inside of the lid. Add 1 ml of FAPD1 Buffer into RNase A tube and mix well. Transfer the mixture into FAPD1 Buffer bottle and store at 4 ° C.
3. Check FAPD2 Buffer before use. Warm FAPD2 Buffer at 55 ° C for 10 minutes if any precipitation formed. Don't shake FAPD2 Buffer vigorously.
4. To avoid acidification of FAPD2 Buffer from CO₂ in the air, close the bottle immediately after use.
5. For FAPDE 001, add 13 ml ethanol (96~100%) to W1 Buffer when first open. For FAPDE 001-1, add 36 ml ethanol (96~100%) to W1 Buffer when first open.
6. For FAPDE 001, add 80 ml ethanol (96~100%) to Wash Buffer when first open. For FAPDE 001-1, add 200 ml ethanol (96~100%) to Wash Buffer when first open.
7. All centrifuge steps are done at full speed (14,000 rpm or 10,000 x g) in a microcentrifuge.

Brief Procedure



General Protocol

- 1. Transfer 1-5 ml of well-grown bacteria culture to a microcentrifuge tube (not provided).**
- 2. Descend the bacteria by centrifuging for 1-2 min and discard the supernatant completely.**
- 3. Add 250 μ l of FAPD1 Buffer to the pellet and resuspend the cells completely by pipetting.**
 - Make sure that RNase A has been added into FAPD1 Buffer when first open.
 - No cell pellet should be visible after resuspension of the cells.
- 4. Add 250 μ l of FAPD2 Buffer and gently invert the tube 5 times to lyse the cells and incubate at room temperature for 2 min.**
 - Do not vortex, vortex may shear genomic DNA. If necessary, continue inverting the tube until the lysate become clear.
 - Do not proceed this step over 5 min.
- 5. Add 350 μ l of FAPD3 Buffer and invert the tube 5 times immediately but gently.**
 - Invert immediately after adding FAPD3 Buffer will avoid asymmetric precipitation.
- 6. Centrifuge for 10 min. During centrifuging, place a FAPD Column in a Collection Tube.**
- 7. Transfer the supernatant carefully to FAPD Column. Centrifuge for 1 min then discard the flow-through.**
 - Do not transfer any white pellet into the column.
- 8. Add 400 μ l of W1 Buffer to FAPD Column. Centrifuge for 1 min then discard the flow-through.**
 - Make sure that ethanol (96-100 %) has been added into W1 Buffer when first open.

9. Add 750 µl of Wash Buffer to FAPD Column. Centrifuge for 1 min then discard the flow-through.

- Make sure that ethanol (96-100 %) has been added into Wash Buffer when first open.

10. Centrifuge for an additional 5 min to dry the column.

- Important step ! This step will remove the residual liquid completely that will inhibit subsequent enzymatic reaction.

11. Place FAPD Column to a new 1.5 ml microcentrifuge tube (not provided).

12. Add 50 µl ~ 100 µl of Elution Buffer or ddH₂O to the membrane center of FAPD Column. Stand the column for 1 min.

- Important step ! For effective elution, make sure that the elution solution is dispensed on the membrane center and is absorbed completely.

- Important : Do not Elute the DNA using less than suggested volume (50ul). It will lower the final yield.

13. Centrifuge for 1 min to elute plasmid DNA.

14. Store plasmid DNA at 4 °C or -20 °C.

Troubleshooting

Low yield

Bacterial cells were not lysed completely

- Too many bacterial cells were used ($OD_{600} > 10$). Separate the bacterial culture into multiple tubes.
- After FAPD3 Buffer addition, break up the precipitate by inverting to ensure higher yield.

Overgrown of bacterial cells

- Incubation time should not longer than 16 hours.

Bacterial cells were insufficient

- Ensure that bacterial cells have grown to an expected amount ($OD_{600} > 1$) after incubation under suitable shaking modes.

Incorrect DNA Elution Step

- Ensure that Elution Buffer was added and absorbed to the center of FAPD Column Martix.

Incomplete DNA Elution

- If size of DNA fragments is larger than 10 kb, use preheated Elution Buffer (60~70°C) on Elution Step to improve the elution efficiency.

Incorrect Wash Buffer

- Ensure that Ethanol was added to Wash Buffer prior to use.

Eluted DNA does not perform well

Residual ethanol contamination

- After Wash Step, dry FAPD Column with additional centrifugation at top speed for 5 minutes or incubation at 60°C for 5 minutes.

Genomic DNA Contaminates

Lysate prepared improperly.

- Gently invert the tube after adding FAPD2 Buffer. And the incubation time should not longer than 5 minutes.
- Do Not use overgrown bacterial culture.

Troubleshooting

RNA Contaminates Plasmid DNA

Insufficiency of RNase A activity in FAPD1 Buffer because of long-term storage

- Prior to using FAPD1 Buffer, ensure that RNase A was added. If RNase A added FAPD1 Buffer is out of date, add additional RNase A into FAPD1 Buffer to a concentration of 50µg/ml then store 4°C.
- Too many bacterial cells were used, reduce sample volume.

Smearing or degrading of Plasmid DNA

Nuclease contamination

- If used host cells have high nuclease activity (e.g., *enA*⁺ strains), perform this Optional Wash Step to remove residuary nuclease.
- After DNA Binding Step, add 400µl of W1 Buffer into FAPD column and column and incubate for 2 minutes at room temperature.
- Centrifuge at full speed (14,000 rpm or 10,000 xg) for 30 seconds.
- Followed using standard Wash Step.

Plasmid DNA is not adequate for enzymatic digestions

Eluted plasmid DNA contains residual ethanol

- Make sure you have discarded the flow-through after washing with Wash Buffer (Step 9) and centrifuged for an additional 3 minutes (Step 10).

Denatured Plasmid DNA migrate faster than supercoiled form during electrophoresis

Incubation in FAPD2 Buffer is too long

- Do not incubate longer than 5 minute in FAPD2 Buffer