User Manual



FavorPrep[™] Blood/ Cultured Cells Genomic DNA Extraction Mini Kit

 For extraction of genomic DNA from fresh blood, forzen blood, cultured cells and fungus

Kit Contents:

For Research Use Only

Cat. No:	FABGK 004 (4 preps_sample)	FABGK 100 (100 preps)	FABGK 300 (300 preps)
RBC Lysis Buffer	7 ml	135 ml	405 ml
FATG Buffer	1.5 ml	30 ml	75 ml
FABG Buffer	1.5 ml x 2	40 ml	100 ml
W 1 Buffer	1.3 ml	45 ml	130 ml
Wash Buffer * (concentrate)	1 ml	25 ml	50 ml
Elution Buffer	1 ml	30 ml	75 ml
FABG Mini Column	4 pcs	100 pcs	300 pcs
Collection Tube	8 pcs	200 pcs	600 pcs
User Manual	1	1	1

Preparation of Wash Buffer by adding ethanol (96 ~ 100%)					
* Ethanol volume for Wash Buffer	4 ml	100 ml	200 ml		

Specification:

Principle: mini spin column (silica matrix) Operation time: 30 ~ 60 minutes Binding capacity: up to 50 µg DNA/ column Typical yield: 15 ~35 µg/ prep Column applicability: centrifugation and vaccum Minimum elution volume: 50 µl Sample size: up to 300 µl of Whole blood up to 200 µl of frozen blood up to 200 µl of frozen blood up to 200 µl of buffy coat up to 1 x 10^7 of Cultured animal cells up to 1 x 10^9 of Cultured bacterial cells up to 5 x 10^7 of Fungus cells

Important Notes:

- Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
- 2. Add ethanol (96-100%) to Wash Buffer when first open.
- 3. Preheat the Elution Buffer to 70 °C for step 13.
- 4. All centrifuge steps are done at full speed (~18,000 x g) in a microcentrifuge.

Brief procedure:



General Protocol: Isolation of DNA from Fresh Human Blood Please Read Important Notes Before Starting Following Steps.

RBC Lysis

- 1. Collect fresh human blood in an anticoagulant-treat collection tube.
- 2. Transfer up to 300 µl of fresh blood to a 1.5 ml microcentrifuge tube (not provided). If the sample is more than 300 µl (up to 1 ml), add the sample to a sterile 15 ml centrifuge tube.
- 3. Add 3x the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex. **For example:** add 900 µl of RBC Lysisi Buffer to the 300 µl of the blood sample.
- 4. Incubate the sample mixture at room temperature for 10 min.
- Note: Make sure that the sample mixture become deep-red and transparant after incubation.
- 5. Centrifuge at 3,000 x g for 5 min. And completely remove the supernatant.
- 6. Add 100 µl of RBC Lysis Buffer to the pellet and resuspend the cells by pipetting.

Cell Lysis

- 7. Add 200 µl of FABG Buffer to the sample mixture. And mix well by vortexing.
- 8. Incubate the sample mixture at room temperature for 10 min or until the sample mixture is clear. During incubation, invert the tube every 3 min.
- 9. Preheat required Elution Buffer (for Step 5 DNA Elution) in a 70 °C water bath.
- 10. (Optional Step): If RNA-free genomic DNA is required, add 5 µl of 10 mg/ml RNase A and mix by vortexing. Incubate for 5 min at room temperature.

DNA Binding

- 11. Add 200 µl of ethanol (96~100%) to the sample and vortex for 10 sec. Pipette the sample to mix well if there is any precipitate formed.
- 12. Place a FABG Column to a Collection Tube. Transfer the sample mixture carefully to FABG Column. Centrifuge at speed 14,000 rpm or 18,000 x g for 1 min. Discard the Collection Tube and place the FABG Column to a new Collection Tube.

Column Washing

- 13. Add 400 µl of W1 Buffer to the FABG Column and centrifuge for 30 sec at speed 14,000 rpm or 18,000 x g. Discard the flow-through and place the FABG Column back to the Collection Tube.
- 14. Add 600 µl of Wash Buffer to the FABG Column and centrifuge for 30 sec at speed 14,000 rpm or 18,000 x g. Discard the flow-through and place the FABG Column back to the Collection Tube.
 - Make sure that ethanol has been added to Wash Buffer when first open.
- Make sole that enabled in a speed 14,000 rpm or 18,000 x g to dry the column.
- - Important Step! This step will avoid the subsequent enzymatic reactions from being inhibited by residual liquid.

Elution

- 16. Place the dry FABG Column to a new 1.5 ml microcentrifuge tube.
- 17. Add 100 µl of Preheated Elution Buffer or TE to the membrane center of FABG Column.
- -- Important Step! For effective elution, make sure that the elution solution is dispensed onto the membrane center and absorbed completely.
- 18. Incubate the FAGB Column at 37 °C for 10 min in an incubator.
- 19. Centrifuge for 1 minute at full speed 14,000 rpm or 18,000 x g to elute the DNA .
 - - Standard volume for elution is 100 μl. If higher DNA yield is required, repeat the DNA Elution step to increase DNA recovery and the total volume could be 200 μl.
- 20. Store the DNA fragment at 4°C or -20°C.

Protocol: Isolation of DNA from Fresh Non-Human Blood

Please Read Important Notes Before Starting Following Steps.

- i. The sample volume of mammalian blood (non-nucleated) can be up to 50 µl; the sample volume of nucleated erythrocytes (eg. bird or fish) can be up to 10 µl.
- ii. Add 150 µl of FATG Buffer and the blood sample into a 1.5 ml microcentrifuge tube (not provided). Mix by vortexing.

Cell Lysis

- 1. Add 200 μl of FABG Buffer to the sample and vortex for 5 seconds.
- 2. Incubate the sample mixture 70 °C for 10 min or until the sample mixture is clarified. During incubation, invert the tube every 3 min.
- 3. Preheat required Elution Buffer in a 70 $^{\circ}\mathrm{C}$ water bath for DNA Elution step.
- 4. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and mix by vortexing. Then incubate for 5 min at room temperature.
- 5. Follow the General Protocol starting from Step 11 (DNA Binding).



Protocol: Isolation of DNA from Frozen Blood

Please Read Important Notes Before Starting Following Steps.

Sample Preparation

- 1. Transfer up to 200 µl blood to a 1.5ml microcentrifuge tube (not provided). If the sample volume is less than 200 µl, add the appropriate volume of PBS.
- 2. Add 30 µl Proteinase K (10 mg/ml, not provided) to the sample and briefly mix. Then incubate for 15 min at 60 °C.

Cell Lysis

- 3. Add 200 µl FABG Buffer to the sample and mix by vortexing.
- 4. Incubate in a 70 °C water bath for 15 min to lyse the sample. During incubation, invert the sample every 3 min.
- 5. Preheat required Elution Buffer in a 70 °C water bath for DNA Elution.
- (Optional Step): If RNA-free genomic DNA is required, add 5 µl of 10 mg/ml RNase A to the sample and mix by vortexing. Incubate for 5 min at room temperature.
- 7. Follow the General Protocol starting from Step 11 (DNA Binding).

Protocol: Isolation of DNA from Buffy Coat

Please Read Important Notes Before Starting Following Steps.

Sample Preparation

Centrifuge whole blood at 3,300 x g for 10 min at room temperature and you will get three different fractions: the upper clear layer is plasma; the intermediate layer is buffy coat, containing concentrated leukocytes; the bottom layer contains concentrated erythrocytes. Extraction total DNA from buffy coat will yield 5-10 times more DNA than an equivalent volume of whole blood.

RBC Lysis

- 1. Transfer up to 200 µl buffy coat to a 1.5 ml microcentrifuge tube (not provided).
- 2. Add 3 x the sample volume of RBC Lysis Buffer and mix by inversion. For example: add 600 µl of RBC Lysisi Buffer to the 200 µl of the buffy coat.
- 3. Incubate at room temperature for 10 min. During incubation, invert the tube every 3 min.
- 4. Centrifuge for 1 minutes at speed 14,000 rpm or 18,000 x g and completely remove the supernatant.
- 5. Add 500 µl of RBC Lysis Buffer and resuspend the pellet by pipetting. Centrifuge for 1 min at 14,000 rpm or 18,000 x g and completely remove the supernatant.
- Add 200 µl of RBC Lysis Buffer and resuspend the pellet by vortexing. Note! Make sure the pellet is completely resuspended.

Cell Lysis

- 7. Add 250 µl of FABG Buffer to the sample and mix by vortex.
- 8. Incubate for 30 min at room temperature or until the sample lysate is clear. During incubation, invert the tube every 3 min.
- 9. Preheat required Elution Buffer in a 70°C water bath for DNA Elution.
- 10. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and mix by vortexing. Then incubate for 5 min at room temperature.

DNA Binding

- 11. Add 250 µl of ethanol (96~100%) to the sample and vortex for 10 sec. Pipette the sample to mix well if there is any precipitate formed.
- 12. Place a FABG Column to a Collection Tube. Transfer the sample mixture carefully to FABG Column. Centrifuge at speed 14,000 rpm or 18,000 x g for 1 min. Discard the Collection Tube and place the FABG Column to a new Collection Tube.
- 13. Follow the General Protocol starting from Step 13 (Column Washing).

Protocol: Isolation of DNA from Cultured Cells

Please Read Important Notes Before Starting Following Steps.

Sample Preparation

- i. Trypsinize the adherent cells before harvesting.
- ii. Transfer the appropriate number of cells (up to 1×10^{7}) to a 1.5 ml microcentrifuge tube (not provided) and centrifuge at 6,000 x g for 20 seconds.
- iii. Remove the supernatant and resuspend the cells with 150 μl of RBC Lysis Buffer.

Cell Lysis

- 1. Add 200 µl of FABG Buffer to the sample and vortex for 5 seconds.
- 2. Incubate for 10 min at 70 °C or until the sample lysate is clear. During incubation, invert the tube every 3 min.
- 3. Preheat required Elution Buffer (for Step 5 DNA Elution) in a 70 °C water bath.
- 4. (Optional Step): If RNA-free genomic DNA is required, add 5µl of 10 mg/ml RNase A to the sample and mix by vortexing. Incubate for 5 min at room temperature.
- 5. Follow the General Protocol starting from Step 11 (DNA Binding).

Protocol: Isolation of DNA from Bacterial cells

Please Read Important Notes Before Starting Following Steps.

Sample Preparation

A. For Gram-negative bacteria:

- i. Transfer the appropriate number of bacterial cells (up to 1 x 10^9) to a 1.5 ml microcentrifuge tube (not provided) and centrifuge at 14,000 rpm or 18,000 x g for 1 min. Discard the supernatant.
- ii. Add 200 µl of FATG Buffer and resuspend the pellet by vortexing or pipetting. Incubate for 5 min at room temperature.
- iii. Follow the Cultured Cell Protocol starting from Step 1 (Cell Lysis).

B. For Gram-positive bacteria:

- i. Transfer the appropriate number of bacterial cell (up to 1 x 10^9) to a 1.5ml microcentrifuge tube (not provided) and centrifuge at speed 14,000 rpm or 18,000 x g for 1 min. Discard the supernatant.
- ii. Add 200 µl of lysozyme buffer (20 mg/ml lysozyme; 20 mM Tris-HCl; 2 mM EDTA; 1% Triton X-100, pH 8.0; prepare fresh lysozyme buffer immediately prior to use) and resuspend the pellet by vortex or pipetting.
- iii. Incubate for 10 minutes at room temperature. During incubation, invert the tube every 2-3 min.
- iv. Follow the Cultured Cell Protocol starting from Step 1 (Cell Lysis).

Protocol: Isolation of DNA from Fungus

Please Read Important Notes Before Starting Following Steps.

Sample Preparation

- i. Harvest appropriate number of fungus cell (up to 5 x 10^7) to a 1.5ml microcentrifuge tube (not provided) and centrifuge at 5,000 x g for 10 min. Discard the supernatant.
- ii. Add 600 µl of sorbitol buffer (1.2 M sorbitol; 10 mM CaCl; 0.1 M Tris-HCl pH 7.5; 35mM β-mercaptoethanol) and resuspend the pellet.
- iii. Add 200 U of lyticase or zymolase. Incubate for 30 min at 30°C.
- iv. Centrifuge the mixture at 2,000 x g for 10 min to harvest the spheroplast, and then remove the supernatant.
- v. Add 200 µl of FATG Buffer to the tube and resuspend the cell pellet by vortex or pipetting.
- vi. Incubate at room temperature for 5 min, and then follow the Cultured Cell Protocol starting from Step 1 (Cell Lysis).

