

FavorPrep™
Plant Genomic
DNA Extraction Mini Kit

User Manual

Cat. No.: FAPGK 001 (50 Preps)
FAPGK 001-1 (100 Preps)

For Research Use Only

v.0905

Introduction

Genomic DNA Mini Kit provides a fast and simple method to isolate total DNA (genomic DNA, mitochondrial and chloroplast) from plant tissue and cells. In the process, sample is disrupted by grinding in liquid nitrogen and lysis buffer incubation. The Lysate is treated with RNase A to degrade RNA and filtrated by filter column to remove cell debris and salt precipitations. In the presence of binding buffer with chaotropic salt, the genomic DNA in the lysate binds to glass fiber matrix in the spin column. The contaminants are washed with an ethanol contained wash buffer and finally, the purified genomic DNA is eluted by low salt elution buffer or water. The protocol does not require DNA phenol extraction and alcohol precipitation. The entire procedure can be completed in 60 minutes. The purified genomic DNA is ready for PCR, real-time PCR, Southern blotting and RFLP.

Quality Control

The quality of Plant Genomic DNA Mini Kit is tested on a lot-to-lot basis. The Kits are tested by isolation of genomic DNA from 50 mg young leave. More than 10µg of genomic DNA could be quantified with spectrophotometer and checked by agarose gel.

Sample: 100 mg of plant tissue

Yield: 5~40µg

Operation time: <60 min

Kit Contents

	FAPGK 001 (50 preps)	FAPGK 001-1 (100 preps)
FAPG1 Buffer	24ml	50ml
FAPG2 Buffer	8ml	15ml
FAPG3 Buffer*	15ml	30ml
W1 Buffer**	22ml	44ml
Wash Buffer***	10ml	20ml
Elution Buffer	15ml	30ml
RNase A (50mg/ml)	500µl	840µl
Filter Column	50 pcs	100 pcs
FAPG Column	50 pcs	100 pcs
1.5 ml Elution tube	50 pcs	100 pcs
2ml Collection tube	100 pcs	200 pcs

*Add 30 ml/60 ml ethanol to FAPG3 Buffer before first use.


**Add 8 ml/16ml ethanol to W1 Buffer before first use.

***Add 40 ml/80ml ethanol to Wash Buffer before first use.

Caution

The component contains irritant agent. During operation, always wear a lab coat, disposable gloves, and protective goggles.

Procedure



Protocol

<p>Step 1 Tissue Dissociation</p>	<ul style="list-style-type: none"> • Cut off 50mg (up to 100mg) of fresh or frozen plant tissue or 5 mg (up to 100 mg) of dried sample. • Grind the sample under liquid nitrogen to a fine powder. Transfer it into a microcentrifuge tube (not provided). For some plant sample, we can destruct it without liquid nitrogen.
<p>Step 2 Lysis</p>	<ul style="list-style-type: none"> • Add 400 µl FAPG1 Buffer and 8µl RNase A (50 mg/ml) into the sample tube and mix by vortexing. Do not mix FAPG1 Buffer and RNase A before use. • Incubate at 65°C for 10 minutes. During incubation, invert the tube every 5 minutes. At the same time, preheat required Elution Buffer (200µl per sample) at 65°C. • Add 130µl FAPG2 Buffer and mix by vortexing. • Incubate at ice for 5 minutes. • Place a Filter Column in a 2 ml Collection Tube. • Apply the mixture from previous step to the Filter Column. Centrifuge for 3 minutes at full speed (13,000 rpm). • Discard the Filter Column and carefully transfer clarified supernatant in Collection Tube to a new microcentrifuge tube (not provided).
<p>Step 3 DNA Binding</p>	<ul style="list-style-type: none"> • Add 1.5 volumes of FAPG3 Buffer (ethanol added) to the cleared lysate and mix immediately by vortexing for 5 seconds. For example, add 750µl FAPG3 Buffer to 500µl lysate. • Place a FAPG Column in a 2 ml Collection Tube. • Apply 750µl the mixture (including any precipitate) from previous step to the FAPG Column. • Centrifuge at full speed (13,000 rpm) for 2 minute. • Discard flow-through in Collection Tube and apply remaining mixture to FAPG Column. • Centrifuge at full speed (13,000 rpm) for 2 minute. • Discard flow-through in Collection Tube.

<p>Step 4 Wash</p>	<ul style="list-style-type: none"> • Add 500µl of W1 Buffer (ethanol added) into the column. • Add 750µl of Wash Buffer (ethanol added) into the column. • Centrifuge at full speed (13,000 rpm) for 30 seconds. • Discard the flow-through and place the FAPG Column back in the Collection Tube. • Centrifuge at full speed for 3 minutes to dry the column matrix. ---Important Step! The residual liquid can affect the quality of DNA and inhibit subsequent enzymatic reactions.
<p>Step 5 DNA Elution</p>	<ul style="list-style-type: none"> • Transfer dried FAPG Column into a clean 1.5 ml microcentrifuge tube (not provided). • Add 50-200µl of preheated Elution Buffer into the center of the column matrix. • Stand for 3 minutes until Elution Buffer absorbed by the matrix. • Centrifuge full speed (13,000 rpm) for 2 minutes to elute purified DNA. <p>---Important Step! For effective elution, make sure that the elution solution is dispensed on the membrane center and is absorbed completely.</p> <p>---Standard elution volume is 200µl. If less sample to be used, reduce the elution volume (50-150µl) to increase DNA concentration.</p>

Troubleshooting

Problem	Possible Reasons/ Solution
Low yield	Insufficient Lysis <ul style="list-style-type: none"> • Prolong the incubation time in lysis buffer to obtain higher yields of DNA.
	Insufficient disruption <ul style="list-style-type: none"> • For most of species we recommend grinding with liquid nitrogen. Homogenization should be done thoroughly until the plant material is ground to a fine powder.
	DNA still bound to the membrane <ul style="list-style-type: none"> • The DNA can be either eluted in higher volumes or by repeating the elution step up to three times. Elution buffer should be preheated to 60°C prior to elution. • To ensure correct pH, use supplied elution buffer.
DNA is degraded	Sample was contaminated with DNase <ul style="list-style-type: none"> • Preheat elution buffer to 60°C for 5 minutes to eliminate any possible Dnase
	Centrifugation speed was too high <ul style="list-style-type: none"> • Higher velocities may cause shearing of the DNA. The centrifugation maximum speed is at 11,000xg.
Column clogged	Too much tissue was used. <ul style="list-style-type: none"> • Too much tissue was used. Reduce the amount of sample material or separate it into multiple tubes.
	Insufficient centrifugation <ul style="list-style-type: none"> • Centrifuge again and extend centrifugation time.
	Precipitate was formed at DNA Binding Step <ul style="list-style-type: none"> • Reduce the sample material. • Before loading the column, break up the precipitate in ethanol-added lysate by pipetting.