## Water Soluble Tetrazolium Salts (WSTs) WST-1,-3,-4,-5,-8,-9

Color changes when reduced. Usefull for metabolism enzymatic assays, staining, cell viability monitoring and immunoassays in conjunction with an electron mediator.

Water-soluble tetrazolium salts (WSTs) were developed by introducing positive or negative charges and hydroxy groups to the phenyl ring of the tetrazolium salt. Positive charges, such as trialkylammonio groups, improve the water-solubility of the formazan dye. However, a large cation is easily precipitated out with organic anions such as carboxylate or phosphate. Though a hydroxy group also improved the water-solubility of the tetrazolium salt, it corresponding formazan dye was not sufficiently water-soluble.

Following WSTs have sulfonate groups added directly or indirectly to the phenyl ring to improve watersolubility.

We also offer several newly developed phenlyazo-type tetrazolium salts, which are easily reduced with NADH or other reducing agents to give orange or purple formazan dyes. Due to the phenylazo group, the color changes with heavy metal ion. Since the water solubility of WSTs are high, 10 mM to 100 mM solutions can be prepared.



Absorption spectra of WST formazan dyes:

Storage: store at 0-5°C (L) or  $-20^\circ C$  for long term.

	PMS	1-Methoxy PMS	1-Methoxy- PES
Water Solubility	0	0	0
Light Stability	х	0	0
pH Stabiity	$\nabla$	$\nabla$	0
ph Stability	V	V	0



Catalog #:	Name / Features	
F98880, 100mg F98881, 500mg F98884, 1g F98885, 10g	WST-1 Syn.: 2-(4-Iodophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt M.W.: 651.35 g/mol CAS: 150849-52-8	$\lambda_{max}$ (formazan): 438nm EC: 21 600(244nm) EC: 37 000(formazan) Solubility: 10 mg/ml H <sub>2</sub> O, 0.65 mg/ml 50mM Tris buffer, pH 8.0

The tetrazolium salt WST-1 is cleaved to formazan by a complex cellular mechanism that occurs primarily at the surface. This bioreduction is mostly dependent on the glycolytic NAD(P)H production of viable cells. Therefore the amount of formazan dye formed directly correlates to the munber of metabolically active cells in the culture. The stability of WST-1 allows for its formulation as a ready-to-use solution. Suitable for cell proliferation and cytotoxcity assays.



P.2/4





the stability of PMS is very poor. However, 1-Methoxy PMS solution can be stored at room temperature for over 3 months without protection from light. Therefore, it is a useful regent for NAD(P)H-tetrazolium-based assay systems.

B2XHY0, 100mg	Methoxy-PES 1-Methoxy-5-ethylphenazinium ethylsulfate M.W.: 364.42 g/mol CAS: 133395-54-7	O N
	-Higher solution stability than 1-Methoxy PMS -Stable in wide range of pH -Higher solution stability than diaphorase	N <sup>+</sup> 0 0,5,0,0,0,0,0,0,0,0,0,0,0,0,0,0,0,0,0,

211bis Avenue J.F. Kennedy - BP 1140 - 03103 Montluçon Cedex - France

The 1-Methoxy PES solution can be stored long term.

erchin

LIFE SCIENCES

P.3/4

## **Related products:**

*WST based cell assays
Cell Counting Kit (WST-8 based, not toxic)
Cell Counting Kit –SK (high sensitivity)
SOD assay kit (WST-1 based)
Total Glutathione Quantification Kit

# <u>899650</u>
# <u>LL0851</u>
# <u>899650</u>
#T33220

## **Other information**

## For R&D use in vitro only.

Please contact InterBioTech – Interchim for any other information

Rev.V01E-J08E

