

# Taq DNA Polymerase

Cat No. MB101-0500

Size: 500 units

Conc. 5 units/ $\mu$ l

Store at -20°C

## Description

*Taq* DNA Polymerase is purified from *E. coli* expressing a *Thermus aquaticus* DNA polymerase gene. This enzyme has a 5'  $\rightarrow$  3' DNA polymerase and a 5'  $\rightarrow$  3' exonuclease activity but lacks a 3'  $\rightarrow$  5' exonuclease activity. The enzyme consists of a single polypeptide with a molecular weight of approximately 94 kDa. *Taq* DNA polymerase is heat-stable and will synthesize DNA at elevated temperatures from single-stranded templates in the presence of a primer.

*Taq* Polymerase is recommended for use in routine PCR reactions. The buffer system is optimized for high specificity and guarantees minimal by-product formation. We supplied *Taq* Polymerase with appropriate buffers. Usually 1-1.5 u of *Taq* DNA Polymerase are used in 50  $\mu$ l of reaction mix. Higher *Taq* DNA Polymerase concentrations may cause synthesis of nonspecific products. However, if inhibitors are present in the reaction mix (e.g., if the template DNA used is not highly purified), higher amounts of *Taq* DNA Polymerase (2-3 u) may be necessary to obtain a better yield of amplification products.

## Storage Buffer

The enzyme is supplied in a storage buffer consisting of 20 mM Tris-HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 50% (v/v) glycerol, and 1% Triton X-100.

## 10X PCR Buffer

Taq DNA Polymerase(5 units/ $\mu$ l)	500 units
10X PCR buffer	1.5ml

## Unit Definition

One unit is defined as the amount of enzyme that will catalyze the incorporation of 10 nmol of dNTP into acid-insoluble form in 30 min at 74°C in a reaction containing 25 mM TAPS (tris-[hydroxymethyl]-methyl-amino-propane-sulfonic acid, sodium salt), pH 9.3 at 25°C, 50 mM KCl, 2 mM MgCl<sub>2</sub>, 1 mM  $\beta$ -mercaptoethanol, 0.2 mM dATP, dGTP, and dTTP, 0.1  $\mu$ M [ $\alpha$ -<sup>32</sup>P] dCTP, and activated salmon sperm DNA.

## Standard PCR with Taq DNA Polymerase

1. For each 50 $\mu$ l reaction, assemble the following in a 0.5 ml PCR tube on ice just prior to use:

Volume	Final Conc.	
add to 50 $\mu$ l	-	PCR Grade Water
1 $\mu$ l	200uM	dNTP Mix (10 mM each dATP, dCTP, dGTP, dTTP)
1 $\mu$ l	0.1-0.2uM	Forward primer, 5-10 uM
1 $\mu$ l	0.1-0.2uM	Reverse primer, 5-10 uM
5 $\mu$ l	2mM MgCl <sub>2</sub>	10X PCR Buffer
0.25 $\mu$ l	1.25 units	<i>Taq</i> DNA Polymerase (5 units/ $\mu$ l)
xul	10ng	DNA template
50 $\mu$ l		Total volume

2. Mix gently. If necessary, centrifuge briefly. Cap tubes and place in thermal cycler.

3. Process in thermal cycler for 25-35 cycles as follows:

Initial Denaturation	2-5 minutes at 94°C	} 30 cycles
Denaturation	20-40 seconds	
Annealing	1 min at the proper annealing temperature	
Extention	2 min at 72°C	
Final extension	5 min at 72°C	

Note: Optimal conditions for amplification will vary depending on the primers and thermal cycler used. It may be necessary to optimize the system for individual primers, template, and thermal cycler.

[www.interchim.com](http://www.interchim.com)

