

## FluoProbes® 488


This technical sheet presents comparative datas and applications of **FluoProbes488 dye**, that is the most popularly used green label, matching conventional fluorophores and competitors dyes. These highlight its superior features (brightness, photostability,...)

For more information, and information about other FluoProbes labels, please ask FluoProbes at [info@fluoprobes.com](mailto:info@fluoprobes.com).

### Fluorescence characteristics

Product name cat.number	$\lambda_{exc}/\lambda_{em}$ max. (nm)	mol. abs. (M <sup>-1</sup> cm <sup>-1</sup> )	Comments
<b>FluoProbes® 488</b>	493 / 519	85 000	<ul style="list-style-type: none"> <li>Bright green fluorescence, compatible with standard filters for FITC/Cy<sup>TM</sup>2</li> <li>pH-independent fluorescence between pH 5 and 9</li> <li>Ultimate photostability, hence minimal fading</li> <li>Superior alternative to FITC, A488</li> <li>Ideal for confocal microscopy, but suits also any other microscopy or technique including microplate readers and FCM.</li> </ul>
<b>FITC</b>	494 / 515	70 000	<ul style="list-style-type: none"> <li>photostability is insufficient with pH, and upon continuous illumination</li> </ul>
<b>Cy<sup>TM</sup>2</b>	489 / 506	150 000	<ul style="list-style-type: none"> <li>brighter than FITC, not photostable</li> </ul>
<b>A488</b>	494 / 517	71 000	<ul style="list-style-type: none"> <li>very bright, photostable, expensive</li> </ul>

### Abs. & Fluorescence Spectra



### FluoProbes488 use by FCM & Microscopy

compared **by CFSM** (confocal microscopy) with competitor labels:

**FP488 demonstrates superior photostability:** 15 min photostability whereas competitors have only 1-2min photostability (see below). FP488 produces also brighter images (see below), and has been used in multiplex combined to blue, orange, red and IR emitting dyes (see below 3-color experiment with Hoechst and FP505).

compared **by flow cytometry (FCM):**

FP488 has a **superior signal and lower background**. It can be used alternatively to conventional dyes (FITC, DTAF) and competitor dyes (A488, Cy2). I.e. +30% higher signal than FITC was shown analysing RBCs with labeled streptavidins using a BD FACS scan. See also comparisons of AnnexinV-FP488 to AnnexinV-FITC (NT-BH9390).

: Streptavidin labeled to detect red blood cells with BD FACSscan. FP488 has a signal superior to FITC (+30%) and lower signal (not shown).

[Info@fluoprobes.com](mailto:Info@fluoprobes.com)

[Technical-support@fluoprobes.com](mailto:Technical-support@fluoprobes.com)

[Order-online@fluoprobes.com](mailto:Order-online@fluoprobes.com)

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140  
03103 Montluçon Cedex - France  
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60


P.1

## FluoProbes488 superior photostability


**FluoProbes®488 demonstrates superior photostability** compared by CFSM (confocal microscopy) with competitor labels: it shows incredible 15 min lasting signal upon continuous illumination whereas competitors fade within 1-2min photostability. A superior signal is also shown.

Myc tagged CHO cells were labeled by a biotinylated anti myc Ab followed by green fluorescent streptavidins, then imaged with a Nikon TE2000E equipped with FITC filter Chroma 41001 at times 0, 30sec, 5min, 10min and 15min durations of light exposure. After quantitation of the signal, normalized fluorescence is plotted against light illumination time.

### High and lasting signal to noise of FluoProbes488



### Ultimate photostability of FluoProbes488 (FP488)




#### Material & methods – CFSM Microscopy:

Nikon TE200E  
Exfo W-cite light source-no attenuation  
Camera: cooled CCD digital camera Roper CoolSnap HQ  
Exposure time: 170 msec  
Digitization: 12 bits (levels of gray 0-4095, full scale)

CHO-KiSS-1 receptor, Myc tagged – fixed on polylysine coated slides, mounted in 50% PBS 50% glycerol.  
Anti myc clone 9<sup>E</sup>10 at 10µg/ml  
Streptavidin-FluoProbes488 (or labeled by competitor dyes AF488, V2, FITC) 2µg/ml  
Filter set: Chroma 41001 (for FITC... Ex 480/40, DM 505, Em 535/50)

## FluoProbes488 great brightness in Microscopy



HEK-293-Myc-GPR54 cells stained with anti myc tag mouse antibody then with anti mouse IgG labeled by FP547H, Phalloidin-FP505 and Hoechst33342. DIC: observed by differential interference contrast; FITC: observed with a FITC filter; Overlay= merge of images with FITC, Hoechst and TexasRed filters

[Info@fluoprobes.com](mailto:Info@fluoprobes.com)  
[Technical-support@fluoprobes.com](mailto:Technical-support@fluoprobes.com)  
[Order-online@fluoprobes.com](mailto:Order-online@fluoprobes.com)

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140  
93103 Montluçon Cedex - France  
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

FT-FP488c(FluoProbes488-Comparison)

*Technical tip: Photostability*

**Fluorescence fading** raises important limitations of use, i.e. during direct observation of anatomy and pathology by fluorescent microscopy, the declining fluorescence strains eyes and demands to occlude the beam and wait fluorescence restoration (switch off pas ok). More decidedly, fading ruins the benefits of advanced microscope imaging systems, as digital camera and confocal microscopy: the fluorescence is recorded during important times to acquire more photons and thus increase the signal (provided background stays low), or to allow autofluorescence to fade (especially with visible light), or because the sample is illuminated over a long lasting scanning.

## FluoProbes Green dyes Applications

### [NT-FP647c](#)      **Using Annexin-FP488 in Confocal Microscopy application.**

AnnexinV-FP488 greatly detects apoptotic Jurkat cells, multiplexed to GPR54 detection with FP647H.

### [NT-FP547c](#)      **Using Phalloiding-FP505 in Confocal Microscopy application in 3-color experiment.**

The cytoskeleton of HEK293 cells are visualized, including stress fiber due to apoptosis, multiplexed to GPR54 detection with FP547H and the nucleus detected using Hoechst33342.

#### Other applications

Use of FP488 dye for **microArrays** – FP488 excellent photostability with excellent brightness are taken to good account in production of diagnostic arrays - InoMust serology multiplexed diagnostic. [Inquire](#)

Use of FP488 dye for **polarisation anisotropy fluorescence** – FP488 is one of the rare dyes with a negative value of anisotropy. [Inquire](#)  
Comparison of FP488 in Confocal Microscopy. FP488 labeled Secondary antibody is compared to A488 in double staining with Cy5 – No crosstalk. [Inquire](#)

Using secondary Abs labeled by FP488 in Confocal Microscopy - insulin secretory granules. American Society for Biochemistry and Molecular Biology ; Brunner ; [Inquire](#)

Using Annexin-FP488 in Confocal Microscopy application - ischemia/reperfusion injured mouse heart ; Shows superior quantum yield and improved stability compared with FITC, A488; yields a 1:1 stoichiometric complex with the reporter group [Inquire](#)

## Other information

Other FluoProbes dyes: see [FT-FPlist](#)

Other products: please consult the [BioSciences catalogue](#)

Catalog size quantities and prices may be found at <http://www.interchim.com/> .

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

FluoProbes® is a registered trademark from Interchim

AlexaFluor® is a trademark from molecularProbes (Invitrogen) ; Cy™ is a registered trademark from Amersham Biosciences Corp.

Rev.109E

[Info@fluoprobes.com](mailto:Info@fluoprobes.com)  
[Technical-support@fluoprobes.com](mailto:Technical-support@fluoprobes.com)  
[Order-online@fluoprobes.com](mailto:Order-online@fluoprobes.com)

Contact your local distributor

FluoProbes®, powered by



213 Avenue J.F. Kennedy - BP 1140  
03103 Montluçon Cedex - France  
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60