

FEATURES

1. Effective:

Favorgen EtBr Destroyer is a laboratory reagent intended for the removal and destruction of Ethidium Bromide contamination. This effect can be monitored and confirmed by UV light exposure whereby once the EtBr has been destroyed the fluorescence will disappear.

2. Fast:

For general protection of uncontaminated area, spray the EtBr Destroyer on the entire working area, leave for about 5 minutes, then wipe it dry with paper towel. When dealing with more seriously contaminated liquid samples the EtBr Destroyer bag is an ideal solution. Three liters of double distilled H₂O containing 30mg of Ethidium Bromide can be destroyed in a matter of days using Favorgen easy-to-follow up protocol with the bag format.

3. Specific:

Two specific designs are for different purposes : Ethidium Bromide Destroyer bags and Ethidium Bromide Destroyer Spray

4. Economical:

Favorgen EtBr Destroyer Sprayer bottles can be refilled after use and the bags are available in both a 10 and 20 pack. In comparison to the costs associated with properly disposing of hazardous materials the EtBr Destroyer products save your money in the long run. You can even re-use your buffer after treatment! One of the big economical benefits of using this product is that after using the Favorgen EtBr Destroyer, the once-contaminated buffer can be filtrated and use again without worry.

5. Safe:

Favorgen EtBr Destroyer has been demonstrated by Ames test to safely and effectively destroy EtBr without any mutagenic effect. The AMES test result data supports that the genotoxicity of EtBr contaminated solution, treated by both the Favorgen EtBr Destroyer Band and Sprayer was reduced significantly and a negative result was determined.



FAVORGEN EtBr DESTROYER — BAG Degrades EtBr Degrades SYBR Green Degrades SYBR Gold Degrades SYBR Orange

Description	Cat. No.
Favorgen EtBr Destroyer Bag: 20 pack	FAEBD 001
Favorgen EtBr Destroyer Bag: 20 pack	FAEBD 001-1
Favorgen EtBr Destroyer Sprayer Pack: 2 x 200ml bottles and 1 sprayer head	FAEBD002
Favorgen EtBr Destroyer Supplement: 2 X 200ml refill bottles, no sprayer head	FAEBD002-1
Favorgen EtBr Destroyer Supplement: 1 x 1L refill bottle, no sprayer head	FAEBD002-2
Favorgen EtBr Destroyer Bag & Sprayer Starter Kit: 5 bags, 1 x 200ml bottle, & 1 sprayer Head	FAEBD 003



FAVORGEN EtBr DESTROYER—BAG PROTOCOL

Favorgen EtBr Destroyer is a specifically designed reagent effectively degrade and destroy Ethidium Bromide and result in both non-fluorescence and non-mutagenic remain. The product also has a proven effectiveness at destroying SYBR Dye from materials after use.

Favorgen EtBr Destroyer is provided in two different formats for the treatment of both solid and liquid Ethidium Bromide contaminant. The ready pack Favorgen EtBr Destroyer Bag is idea for the treatment of liquid Ethidium Bromide contaminant while Sprayer can be used for the treatment of solid contaminant waste including electrophoresis gels, glassware, paper towels, gloves, laboratory equipment, bench surface etc.



211 bis Avenue Kennedy - BP 1140
03103 Montluçon - France
33 (0) 4 70 03 88 55
Fax 33 (0) 4 70 03 82 60
e-mail interchim@interchim.com

Agence Paris - Normandie
33 (0) 1 41 32 34 40
Fax 33 (0) 1 47 91 23 90
e-mail interchim.paris@interchim.com

Ethidium Bromide

Ethidium Bromide containing buffer :

Step 1

Tear open the package and take out the Favorgen EtBr Destroyer Bag. Use the product immediately or seal it into a zip-lock bag to avoid prolonged exposure to the air.

Step 2

Put the EtBr Destroyer into the EtBr-containing buffer. Completely submerge the bag so that the entire Destroyer bag is in the solution.

Step 3

Destroyer Working Timetable

EtBr-containing Quantity	2mg/3L	5mg/3L	10mg/3L	20mg/3L	30mg/3L
Working Time	16hrs	30hrs	3 days	5 days	6 days

Step 4

When the treated buffer becomes clear, please use a UV light (254nm) to observe the EtBr-content remaining. If exposure to the UV light displays no fluorescence and the O.D. value is equal to that of d.d.H₂O it means that the EtBr containing buffer has been totally destroyed.

Step 5

The Favorgen EtBr Destroyer Bag has absorbed and destroyed all of the EtBr content in the liquid. Following treatment the bag can be thrown away as a non-hazardous item in compliance with the local waste removal regulations.

Ethidium Bromide containing gels :

Step 1

Tear open the package and take out the Favorgen EtBr Destroyer Bag onto the container which fills with 3L d.d.H₂O

Step 2

30 minutes later, you can drop the EtBr-contaminated gels into this container filling with the solution of Favorgen EtBr Destroyer Bag.

Step 3

You can drop as many as gels into the container during the 6 months period, but the total amount of EtBr in all gels can't be over 30mg and the solution of Favorgen EtBr Destroyer Bag has to cover the whole gels. The last gel dropped into the container only can be thrown away after 24 hours treatment.

Step 4

When the treated gel becomes clear, please use the UV light box (UV Wavelength 254nm) to ensure the Ethidium Bromide has been totally destroyed.

Step 5

After 6 months, no matter how much EtBr has been destroyed, the solution of Favorgen EtBr Destroyer Bag can't be used any more. The users need to clean up the container for another new solution.

Step 6

After gels containing EtBr which have been destroyed can be collected as the non-hazardous item in compliance with the local waste removal regulations.

SYBR Green I

Buffer containing SYBR Green I :

Use the same procedure as for the Ethidium Bromide containing buffer. However please note that 1.) The decontamination time is different to that of Ethidium Bromide containing buffer, and 2.) You must use UV light box (wavelength 488nm) to observe SYBR Green I.

Gels containing SYBR Green I:

The procedure is the same as that of the Ethidium Bromide containing gels. However please note that 1.) The decontamination

time is different to that of Ethidium Bromide containing gels, and 2.) You must use UV light box (wavelength 488nm) to observe SYBR Green I.

SYBR Gold

Buffer containing SYBR Gold :

Use the same procedure as for the Ethidium Bromide containing buffer. However please note that 1.) The decontamination time is different to that of Ethidium Bromide containing buffer, and 2.) You must use UV light box (wavelength 488nm) to observe SYBR Gold.

Gels containing SYBR Gold :

The procedures are the same as the usage for Ethidium Bromide containing gels. However please note that 1.) The decontamination time is different to that of Ethidium Bromide containing gels, and 2.) Please use the UV light box (wavelength 488nm) to observe the SYBR Gold.

SYPRO Orange

Buffer containing SYPRO Orange :

Use the same procedure as for the Ethidium Bromide containing buffer. However please note that 1.) The decontamination time is different to that of Ethidium Bromide containing buffer, and 2.) You must use UV light box (wavelength 488nm) to observe SYPRO Orange.

Gels containing SYPRO Orange :

The procedure is the same as that of the Ethidium Bromide containing gels. However please note that 1.) The decontamination time is different to that of Ethidium Bromide containing gels, 2.) Please use the UV light box (wavelength 488nm) to observe the SYPRO Orange.