

FT-JQ6940



## Caspase Substrates

### Product Description

Product name cat. number	MW (g·mol <sup>-1</sup> )	$\lambda_{exc} / \lambda_{em}$ max. (nm)	Function (for Caspases)	Soluble in
<b>Ac-DEVD-AFC</b> FP-R1147A, 5 mg	729,61	380 / 500	Caspase 3/7	DMSO
<b>Ac-DEVD-AMC</b> FP-46764A, 100 mg	675,64	351 / 430	Caspase 3/7	DMSO
<b>Ac-DEVD-CHO</b> FP-AN1842, 1 mg	502,47	N/A	Caspase 3/7 inhibitor	DMSO
<b>Ac-DEVD-pNA</b> FP-BI0930, 5 mg	638,58	408 / N/A	Caspase 3/7	DMSO
<b>FITC-C6-DEVD-FMK</b> FP-JQ6930, 1 mg	994,99	492 / 516	Caspase 3/7 (cell permeable)	DMSO
<b>FITC-C6-LEHD-FMK</b> FP-JQ6940, 1 mg	1031,07	492 / 516	Caspase 9 (cell permeable)	DMSO
<b>Ac-IETD-AFC</b> FP-HG804A, 5 mg	729,65	380 / 500	Caspase 8	DMSO
<b>Ac-IETD-AMC</b> FP-HH2900, 5 mg	675,68	351 / 430	Caspase 8	DMSO
<b>Ac-IETD-CHO</b> FP-JQ6970, 5 mg	502,52	N/A	Caspase 8 inhibitor	DMSO
<b>Z-IETD-pNA</b> FP-JQ6980, 5 mg	730,72	380 / 500	Caspase 8	DMSO
<b>Z-IETD-AFC</b> FP-R1147A, 5 mg	821,71	380 / 500	Caspase 3/7	DMSO
<b>Z-DEVD-AMC</b> FP-AM355A, 5 mg	767,74	351 / 430	Caspase 3/7	DMSO
<b>Z-DEVD-pNA</b> FP-HG269A, 5 mg	730,68	408 / N/A	Caspase 3/7	DMSO
<b>Z-IETD-AFC</b> FP-HG830A, 5 mg	821,75	380 / 500	Caspase 8	DMSO
<b>Ac-LEDH-AMC</b> FP-F98962, 5 mg	711,72	351 / 430	Caspase 9	DMSO
<b>(Ac-LEDH)2-R110</b> FP-CJF520, 5 mg	1403,41	498 / 520	Caspase 9	DMSO
<b>(Z-DEVD)2-R110</b> FP-HG9031, 1 mg	1515,44	498 / 520	Caspase 3/7	DMSO

FT-JQ6940

Product name cat. number	MW (g·mol <sup>-1</sup> )	$\lambda_{exc}/\lambda_{em}$ max. (nm)	Function (for Caspases)	Soluble in
<b>(Ac-IETD)2-R110</b> FP-HG9281, 1 mg	1331,34	498 / 520	Caspase 8	DMSO
<b>Z-DEVD-ProRed 620</b> FP-GCZ830, 1 mg	1565,50	534 / 619	Caspase 3/7	DMSO
<b>Z-IETD-ProRed 620</b> FP-GCZ840, 1 mg	1565,59	534 / 619	Caspase 8	DMSO
<b>Z-LEHD-ProRed 620</b> FP-GCZ850, 1 mg	1453,47	534 / 619	Caspase 9	DMSO

**Storage:** -20°C Protect from light and moisture

## Introduction

Caspases are essential in cells for apoptosis, or programmed cell death, in development and most other stages of adult life, and have been termed "executioner" proteins for their roles in the cell. Some caspases are also required in the immune system for the maturation of lymphocytes. Failure of apoptosis is one of the main contributions to tumour development and autoimmune diseases.

Caspases are also involved in ischemia and Alzheimer's disease. There are two types of apoptotic caspases: initiator and effector caspases. Initiator caspases (e.g., CASP2, CASP8, CASP9, and CASP10) cleave inactive pro-forms of effector caspases, thereby activating them. Effector caspases (e.g., CASP3, CASP6, CASP7) in turn cleave other protein substrates within the cell, to trigger the apoptotic process. The initiation of this cascade reaction is regulated by caspase inhibitors.

FluoProbes offers both color and a group of blue, green, and red fluorescent substrates for monitoring caspase activities. The AMC-, AFC-, R110- and ProRed™-derived protease substrates are colorless and non-fluorescent. Cleavage of blocking protease-cleavable peptide residue by caspases generates the strongly blue, green, or red fluorescent respectively that can be monitored by fluorescence instruments. Our proprietary ProRed™-derived caspase substrates are the most sensitive red indicators for the fluorimetric detection of various caspase activities. In general, R110 and ProRed™ substrates are much more sensitive than the AMC-, AFC- or 4-nitroaniline-based substrates.

## Directions for use

### Guidelines for use

Following protocol only provides a guideline, and should be modified according to your specific needs.

1. General Solution Caspase Assays Using AMC, AFC, pNA, R110 and ProRed Substrates
  - 1.1. Prepare a 10 mM stock solution in DMSO.
  - 1.2. Prepare a 2X caspase substrate (50  $\mu$ M) assay solution as the following:
    - 50  $\mu$ L substrate stock solution
    - 100  $\mu$ L DTT (1M)
    - 400  $\mu$ L EDTA (100 mM)
    - 10 mL Tris Buffer (20 mM), pH =7.4
  - 1.3. Mix equal volume of the caspase standards or samples with 2X caspase substrate assay solution (from Step 1.2), and incubate the solutions at room temperature for at least 1 hour.
  - 1.4. Monitor the fluorescence using a fluorescence microplate reader, or absorbance using an absorbance microplate reader.
2. Cell Caspase Assays Using Cell-Permeable FMK Caspase Probes
  - 2.1. Prepare a 2-5 mM stock solution in DMSO.
  - 2.2. Treat cells as desired.
  - 2.3. Prepare a 2X permeable caspase substrate (20  $\mu$ M) assay solution by diluting the DMSO stock solution (from Step 2.1) in Hanks with 20 mM Hepes buffer (HHBS).
  - 2.4. Mix equal volume of the treated cells with 2X caspase substrate assay solution (from Step 2.3), and incubate the cells in a 37°C, 5% CO<sub>2</sub> incubator for at least 1 hour.
  - 2.5. Wash the cells with HHBS for at least once.

FT-JQ6940

- 2.6. Monitor the fluorescence intensity by a flow cytometer, a fluorescence microscope or a fluorescence microplate reader.
  
3. Cell Caspase Assays Using Cell-Permeable FMK Caspase Probes (For #13470-13476 only)
  - 3.1. Prepare a 150X stock solution by adding 50  $\mu$ L DMSO into the vial.
  - 3.2. Treat cells as desired.
  - 3.3. Add 150 X DMSO stock solution (from Step 3.1) into the cell solution at a 1:150 ratio, and incubate the cells in a 37°C, 5% CO<sub>2</sub> incubator for 1 hour.
  - 3.4. Wash the cells with HHBS for at least once.
  - 3.5. Monitor the fluorescence intensity by flow cytometer, fluorescence microscopy or fluorescent microplate reader.

## References

1. Imre G, Heering J, Takeda AN, Husmann M, Thiede B, Zu Heringdorf DM, Green DR, van der Goot FG, Sinha B, Dotsch V, Rajalingam K. (2012) Caspase-2 is an initiator caspase responsible for pore-forming toxin-mediated apoptosis. *EMBO J*.
2. Zhang F, Lau SS, Monks TJ. (2012) A Dual Role for Poly(ADP-ribose) Polymerase-1 (PARP-1) During Caspase-dependent Apoptosis. *Toxicol Sci*.
3. Xiao XL, Peng J, Su Q, Xiang SL, Tang GH, Huang YS, Zhou XT. (2006) [Diallyl Trisulfide Induces Apoptosis of Human Gastric Cancer Cell Line MGC803 Through Caspase-3 Pathway.]. *Ai Zhong*, 25, 1247.
4. Sakaue M, Motoyama Y, Yamamoto K, Shiba T, Teshima T, Chiba K. (2006) Quantitative measurement of caspase-3 activity in a living starfish egg. *Biochem Biophys Res Commun*, 350, 878.
5. Kume T, Taguchi R, Katsuki H, Akao M, Sugimoto H, Kaneko S, Akaike A. (2006) Serofendic acid, a neuroprotective substance derived from fetal calf serum, inhibits mitochondrial membrane depolarization and caspase-3 activation. *Eur J Pharmacol*, 542, 69.
6. Fennell M, Chan H, Wood A. (2006) Multiparameter measurement of caspase 3 activation and apoptotic cell death in NT2 neuronal precursor cells using high-content analysis. *J Biomol Screen*, 11, 296.
7. Wu X, Simone J, Hewgill D, Siegel R, Lipsky PE, He L. (2006) Measurement of two caspase activities simultaneously in living cells by a novel dual FRET fluorescent indicator probe. *Cytometry A*, 69, 477.
8. Jiang DJ, Jia SJ, Dai Z, Li YJ. (2006) Asymmetric dimethylarginine induces apoptosis via p38 MAPK/caspase-3-dependent signaling pathway in endothelial cells. *J Mol Cell Cardiol*, 40, 529.

## Technical and scientific information

### Related products

- Dulbecco's PBS 20X (Sterile), RJ226A

## Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.