

## Amplite™ Fluorimetric NAD/NADH Assay Kit

*\*Red Fluorescence\**

**Ordering Information:**

Product Number: #15257 (400 assays)

**Instrument Platform:**

Fluorescence microplate readers


**Storage Conditions:**

Keep in freezer and avoid light

### Introduction

Nicotinamide adenine dinucleotide (NAD<sup>+</sup>) and nicotinamide adenine dinucleotide phosphate (NADP<sup>+</sup>) are two important cofactors found in cells. NADH is the reduced form of NAD<sup>+</sup>, and NAD<sup>+</sup> is the oxidized form of NADH. It forms NADP with the addition of a phosphate group to the 2' position of the adenyl nucleotide through an ester linkage. NADP is used in anabolic biological reactions, such as fatty acid and nucleic acid synthesis, which require NADPH as a reducing agent. In chloroplasts, NADP is an oxidizing agent important in the preliminary reactions of photosynthesis. The NADPH produced by photosynthesis is used as reducing power for the biosynthetic reactions in the Calvin cycle of photosynthesis. The traditional NAD/NADH and NADP/NADPH assays are done by monitoring of NADH or NADPH absorption at 340 nm. This method suffers low sensitivity and high interference since the assay is done in the UV range that requires expensive quartz microplate.

This Amplite™ NAD/NADH Assay Kit provides a convenient method for sensitive detection of NAD, NADH and their ratio. The enzymes in the system specifically recognize NAD/NADH in an enzyme cycling reaction. There is no need to purify NAD/NADH from sample mix. The enzyme cycling reaction significantly increases detection sensitivity. In addition, this assay has very low background since it is run in the red visible range that significantly reduces the interference from biological samples. The assay has demonstrated high sensitivity and low interference with 570 nm excitation 590 nm emission.



The Amplite™ Fluorimetric NAD/NADH Assay Kit provides a sensitive, one-step fluorimetric assay to detect as little as 10 picomoles of NAD (H) in a 100  $\mu\text{L}$  assay volume (100nM; Figure 1). The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation with no separation steps required. Its signal can be easily read by either fluorescence microplate reader with Ex/Em

### Kit Key Features

**Broad Application:** Can be used for quantifying NAD/NADH in solutions, in cell extracts.

**Sensitive:** The kit detect as low as 10 picomoles of NAD/NADH in solution.

**Continuous:** Easily adapted to automation with no separation required.

**Convenient:** Formulated to have minimal hands-on time. No wash is required.

**Non-Radioactive:** No special requirements for waste treatment.

---

Component	Amount
Component A: NAD/NADH cycling enzyme mixture	2 bottles (lyophilized powder)
Component B: NADH sensor buffer	1 bottle (20 mL)
Component C: NADH standard (FW: 709)	1 vial (142 µg)
Component D: DMSO	1 vial (0.5 mL)

### Assay Protocol for one 96-well plate

#### Brief Summary

**Prepare NAD/NADH reaction mixture (50 µL) → Add NADH standards or test samples (50 µL)  
→ Incubate at room temperature for 15 min-2hr → Read fluorescence at Ex 570 nm/Em 590 nm**

*Note: Thaw 1 vial (or bottle) each of all the kit components to room temperature before starting the experiment.*

#### 1. Prepare NADH stock solution:

Prepare NADH standard stock solution: Add 200 µL of DMSO (Component D) into the NADH standard vial (Component C) to have 1 mM (1nmol/µL) stock solution.

*Note: The unused NADH solution should be divided as single use aliquots and stored at -20°C.*

#### 2. Prepare NAD/NADH reaction mixture:

Prepare the NAD/NADH reaction mixture: Add 10 mL of NAD/NADH sensor buffer (Component B) to the bottle of NAD/NADH cycling enzyme mixture (Component A), mixed well.

*Note: This solution is enough for two 96-well plates. The unused NAD/NADH mixture should be divided as single use aliquots and stored at -20°C.*

#### 3. Prepare Serial NADH (0 to 10 µM) solutions

3.1 Add 10 µL of NADH standard stock solution (from step 1) to 990 µl PBS (pH 7.4) buffer to generate 10 µM (10 pmol/µL) standard.

*Note: Diluted NADH standard solution is unstable, should be used within 4 hours.*

3.2 Take 200 µL of 10 µM solution to perform 1:3 serial dilutions to get 3, 1, 0.3, 0.1, 0.03, 0.01 and 0 standard NADH solutions.

3.3 Add NADH standards and NAD/NADH-containing test samples into a 96-well solid black microplate as described in Tables 1 and 2

*Note: Prepare your cell or tissue samples as desired.*

---

NS1	NS1	....	....	....	....							
NS2	NS2											
NS3	NS3											
NS4	NS4											
NS5	NS5											
NS6	NS6											
NS7	NS7											

Note: NS= NADH standards, BL=Blank control, TS=test samples.

**Table 2.** Reagent composition for each well:

NADH Standard	Blank Control	Test Sample
Serial dilutions* (50 µL)	PBS: 50 µL	50 µL

\*Note: Add the serially diluted NADH standards from 0.01 µM to 10 µM into wells from NS1 to NS7 in duplicate.

High concentration of NADH (e.g., >100 µM, final concentration) may cause reduced fluorescence signal due to the overoxidation of NADH sensor (to a non-fluorescent product).

#### 4. Run NAD/NADH Assay in supernatants Reaction

4.1 Add 50 µL of NADH reaction mixture (from step 2) to each well of the NADH standard, blank control, and test samples (see step 3.3) so that the total NADH assay volume is 100 µL/well

Note: For a 384-well plate, add 25 µL sample, 25 µL of NADH reaction mixture per well.

4.2 Incubate the reaction for 15 minutes to 2 hours at room temperature, protected from light.

4.3 Monitor the fluorescence increase with 530-570 nm (optimal at 570) excitation and 590-600 nm emission by using a fluorescence plate reader.

Note1: The contents of the plate can also be transferred to a white clear bottom plate and read by absorbance microplate reader at the wavelength of 576±5 nm. The absorption detection has lower sensitivity compared to fluorescence reading.

Note2: To detect NADH only, aliquot 200 µL samples into Eppendorf tubes. Heat samples to 60°C for 30 min in a heating block or a water bath. All NAD will be decomposed while NADH will be still intact under the conditions. Cool samples on ice, and quick spin samples if precipitates occur. Transfer 50 µL of NADH samples into the wells as indicated in Table 1 and 2.

#### 5. Run Data Analysis

The fluorescence in blank wells (with the PBS buffer only) is used as a control, and is subtracted from the values for those wells with the NADH reactions. The typical data are shown in Figure 1 (NADH standard curve).

Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point.

www.interchim.com



211 bis Avenue Kennedy - BP 1140  
03103 Montluçon - France  
33 (0) 4 70 03 88 55  
Fax 33 (0) 4 70 03 82 60  
e-mail interchim@interchim.com

Agence Paris - Normandie  
33 (0) 1 41 32 34 40  
Fax 33 (0) 1 47 91 23 90  
e-mail interchim.paris@interchim.com




Figure 1. NADH dose response on 96-well black plate was measured with Amplite™ NAD/NADH Assay Kit using a BMG LabTech NOVOSTar microplate reader. As low as 100 nM (10 pmol/well) of NADH can be detected with 1 hour incubation time (n=3) while there is no response from NADPH.

### Related Products

15257	Amplite™ Fluorimetric NAD/NADH Assay Kit *Red Fluorescence*	1 kit
15258	Amplite™ Colorimetric NAD/NADH Assay Kit	1 kit
15259	Amplite™ Fluorimetric NADP/NADPH Assay Kit *Red Fluorescence*	1 kit
15260	Amplite™ Colorimetric NADP/NADPH Assay Kit	1 kit
15261	Amplite™ Fluorimetric NADH Assay Kit *Red Fluorescence*	1 kit
15262	Amplite™ Fluorimetric NADPH Assay Kit *Red Fluorescence*	1 kit

### References:

1. Ziegenhorn J, Senn M, Bucher T. (1976) Molar absorptivities of beta-NADH and beta-NADPH. Clin Chem, 22, 151.
2. Ikegami T, Kameyama E, Yamamoto SY, Minami Y, Yubisui T. (2007) Structure and Properties of the Recombinant NADH-Cytochrome b(5) Reductase of Physarum polycephalum. Biosci Biotechnol Biochem.
3. Kimura N, Fukuwatari T, Sasaki R, Shibata K. (2006) Comparison of metabolic fates of nicotinamide, NAD<sup>+</sup> and NADH administered orally and intraperitoneally; characterization of oral NADH. J Nutr Sci Vitaminol (Tokyo), 52, 142.
4. O'Donnell JM, Kudej RK, LaNoue KF, Vatner SF, Lewandowski ED. (2004) Limited transfer of cytosolic NADH into mitochondria at high cardiac workload. Am J Physiol Heart Circ Physiol, 286, H2237.