






PEGylation silane reagents (for surface modification)

Products Description

Silane (Si) functionalized PEG can be used to modify glass, silica, silicone and other surfaces via the reaction between hydroxyl groups and triethoxyl silane. They can thus graft these support with functional groups that will confer specific physico-chemical properties (such as with nPEG-Silane, or electronic charges of amine, COOH), or binding properties (e.g. with Biotin groups), or allow further chemical modifications by chemistry (e.g. with amine, COOH) or direct reaction in aqueous samples (such as with SH, Maleimide, succinimidyl (NHS)).

PolyEthylene Glycol (**PEG**, also known as **PEO**: PolyEthylOxy) is a convenient structure used to create a link biomolecules and supports. It is available in a variety of lengths, providing a flexible spacer, and it is highly hydrophilic. PEG have found great application to create drug delivery systems, hydrophilic gels, special coated surfaces, biocaptors, microscopy slides, microarrays, beads, vaccines, coatings, biomaterials or biocompatible materials such as optical lenses, purification supports,...

Product name cat.number	Cat.number Qty: 100mg	MW (g·mol ⁻¹)	
mPEG-silane (M,L)			 IO5230
PEG6-0102 PG1-SL-1K	OO7020	1 000	
PEG6-0101 PG1-SL-2K	EV4980	2 000	
PEG6-0103 PG1-SL-5K	DY8030	5 000	
PEG6-350 PG1-SL-10K	GV7350	10 000	
PEG6-350 PG1-SL-350	GV7360	350	
PEG6-550 PG1-SL-550	GV7370	550	
PEG6-750 PG1-SL-750	GV7380	750	
Silane-PEG-amine (NH2) (M)			
PEG6-0013 PG2-AMSL-2K	WT9700	2 000	
PEG6-0011 PG2-AMSL-3K	DZ5160	3 400	
PEG6-0012 PG2-AMSL-5K	DZ5380	5 000	
PEG6-0010 PG2-AMSL-10K	DZ5390	10 000	
Silane-PEG-Maleimide (M)			
PEG6-0033 PG2-MLSL-2K	GV7310	2 000	
PEG6-0032 PG2-MLSL-3K	WT9720	3 400	
PEG6-0031 PG2-MLSL-5K	FO9490	5 000	
PEG6-0031 PG2-MLSL-10K	LO4660	10 000	
Silane-PEG-NHS (M)			
PEG6-0043 PG2-NSSL-2K	GV7260	2 000	
PEG6-0042 PG2-NSSL-3K	WT9730	3 400	
PEG6-0041 PG2-NSSL-5K	GV7250	5 000	
PEG6-0044 PG2-NSSL-10K	LO8210	10 000	
Silane-PEG-COOH (M)			
PEG6-0053 PG2-CASL-5K	GV7320	2 000	
PEG6-0052 PG2-CASL-3K	WT9740	3 400	
PEG6-0051 PG2-CASL-5k	GV7330	5 000	
PEG6-0054 PG2-CASL-10k	LO5310	10 000	
Silane-PEG-SH (M)			
PEG6-0062 PG2-SLTH-3K	WT9750	3 400	

Contact your local distributor

Uptima, powered by



interbiotech@interchim.com

P.1

FT-LO5310

Silane-PEG-biotin (M)				
	PG2-BNSL-1K	IV2270	1 000	
PEG6-0023	PG2-BNSL-2K	GV7270	2 000	
PEG6-0022	PG2-BNSL-3K	WT9710	3 400	
PEG6-0021	PG2-BNSL-5K	GV7280	5 000	
	PG2-BNSL-10K	LO5320	10 000	
DSPE-PEG-Silane				
	PG2-DSSL-5K	1J9920	5 000	

- Physical Properties:

Appearance: Off white solid or viscous liquid depends on molecule weight;

Solubility: Soluble in regular aqueous solution (H₂O) as well as most organic solvents (DMSO).

- Storage Conditions:

Store: -20°C, desiccated. Silane PEG tends to hydrolyze from moisture. Avoid frequent thaw and frozen. (M)

- Reaction Procedures:

Generally, silane PEG is first dissolved in a mixture of ethanol/water solution, then allow the solution to contact with oxidized glass or silica surface for 30 min to 2 hours. Under this condition, silane PEG should bind to the material covalently.

- Materials Required:

Pegylation solution: Ethanol/water (w/w, 95%/5%);

Silane PEG stock solution, 10~50 mg/mL in pegylation solution, prepare in fresh.

Washing solution: Distilled water.

- Reaction Steps:

1. Dissolve targeted materials in Pegylation buffer.

2. Add silane PEG stock solution to the targeted conjugation materials with the final concentration keep at least 10 mg/mL. 10~50 molar excess of silane PEG needed for optimal conjugation;

3. Allow mixture agitates at room temperature for 30 min at room temperature or 2 hours at 4 -20 OC.

4. Wash out unreacted materials either by distilled water.

Related / associated products and documents

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

Other PEO/PEG reagents, including purified PEGs, synthetic PEG products (PEO3 to PEO36):

- Heterobifunctional crosslinkers: [NHS-PEO-MAL \(AL6581\)](#) (altern.to SMCC 17412A)

- Homobifunctional crosslinkers: [NHS-PEO-NHS \(BH8811\)](#) (altern.to DSS 54940A)

- [MAL-PEO-MAL \(L7736A\)](#) (altern.to BMOEL7730A)

- PEO Linkers & modifiers:

- [Maleimide-PEG-COOH \(AZ4170\)](#) (altern.to BMPA 43064A)

- [NHS-PEG-COOH \(AN1280\)](#)

- [Maleimide-PEG-Amine \(FK3520\)](#)

- [Azide-PEG-COOH \(WU0930\)](#)

- PEG modifiers:

- [mPEG reagents \(DZ3531\)](#): mPEG-NHS and others mono-fonctionnal (MAL) (-SH, -OH,...)

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use.

Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.R10E-P01E

Contact your local distributor

Uptima, powered by



interbiotech@interchim.com