




PRODUCT INFORMATION SHEET



MARKERGENE™ OPTI-KLEAR™ LIVE CELL IMAGING BUFFER Product M1898

**MARKER GENE TECHNOLOGIES, INC.
UNIVERSITY OF OREGON RIVERFRONT RESEARCH PARK
1850 MILLRACE DRIVE
EUGENE, OREGON 97403-1992
1-888-218-4062
EMAIL: TECHSERVICE@MARKERGENE.COM**



I. OVERVIEW

Physiology, the study of living tissues, has gained wide spread application in microscopy with the advent of fluorescent dyes and inverted microscopes dedicated to their viewing. A critical component is maintaining physiological conditions in media formulations, including pH and osmolarity, without contributing to background fluorescence. Most cells are grown in media formulations that use bicarbonate buffers, requiring the continuous presence of CO₂ to maintain physiological pH. Just 30 minutes outside a CO₂ incubator can change the pH from 7.3 to 7.8, and to 8.5 after an hour.

Many of these formulations include phenol red which can have estrogenic effects, and quench fluorescence signals. Most media formulations also contain autofluorescent components like riboflavin and folic acid, that while quenched when phenol red is present, will autofluoresce when phenol red-free solutions are used. And while phenol red itself is not fluorescent, there may be contaminants present in the dye that are.

Marker Gene Technologies, Inc. now offers HEPES based **MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer** as an ideal substitute for growth media for all extended live cell imaging sessions.

HEPES has a pKa of 7.66 making it a much stronger buffer in the pH 7.2-7.6 range than media formulations based on bicarbonate buffer. Unlike the bicarbonate buffering system that requires the use of a CO₂ incubator, the HEPES buffering system may be used with or without a CO₂ atmosphere.

Shipped in amber Nalgene bottles to prevent ambient light-induced toxicity, the sterile filtered **MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer** is formulated to:

- i) Maintain proper pH and osmolarity for hours
- ii) Lower background fluorescence
- iii) Provide energy source and preserve fluorescent signals



II. MATERIALS

MarkerGene™Opti-Klear™ Live Cell Imaging Buffer:
30 mL of 5X stock. pH 7.5

Storage and Handling.

Store HEPES based **MarkerGene™Opti-Klear™ Live Cell Imaging Buffer** in the amber Nalgene bottle provided to avoid creation of toxic hydrogen peroxide which can occur when HEPES based buffers are exposed to ambient light for long periods. Solution is sterile filtered – and should be stored at 4°C or, for extended periods, stored frozen (-20°C). In case of contact with skin or eyes, wash thoroughly with soap and cold-water. Solution should be stable for at least 6 months following purchase.

III. PROTOCOL:

This bottle provides sufficient reagent to make 150 mL of imaging solution.

Standard protocol:

1. To prepare 10 mL of 1X imaging solution, combine 2mL of 5X **MarkerGene™Opti-Klear™ Live Cell Imaging Buffer** with 8 mL sterile water and briefly mix.
2. Remove media, wash cells 1-2 times with sterile PBS (if desired), and replace with 1X **MarkerGene™Opti-Klear™ Live Cell Imaging Buffer**.
3. Proceed to imaging as normal.



FIGURES:

Figure 1: Reduction of background fluorescence by switching from DMEM with 10% fetal calf serum (Fisher) to **MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer**. HeLa cells were treated with JC-1 (Marker Gene Technologies, Inc. Product M1842) for 60 minutes to detect J-aggregate formation in healthy, polarized mitochondria and imaged with a 40X oil objective and photographed with identical exposure times and light intensity.




Figure 2: Reduction of background fluorescence by switching from DMEM with 10% fetal calf serum (Fisher) to **MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer**. HeLa cells were treated with JC-1 (Product M1842) for 60 minutes to detect total mitochondria in healthy, polarized mitochondria and imaged with a 40X oil objective and photographed with identical exposure times and light intensity.






Figure 3: AG06173 Normal Human Fibroblasts were stained with LysoTracker™ (Life Technologies) for 1 hour then washed 3X with PBS. Media was replaced with A) complete media B) serum-free media C) **MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer** and imaged with 40X dry long working distance objective and identical exposure times and light intensity.



PRODUCT M1898			
DESCRIPTION	QUANTITY	PART NO.	STORAGE
SAMPLE PREPARATION			
MarkerGene™ Opti-Klear™ Live Cell Imaging Buffer - sterile filtered	30 mL	M1898	C, L
DOCUMENTATION			
MSDS SHEETS	1		
PRODUCT INFORMATION SHEET	1		

Notes: C=store cold (4°C); L=light sensitive



REFERENCES:

1. Brown PS, Wang E, Aroeti B, Chapin SJ, Mostov KE, Dunn KW. (2000) "Definition of distinct compartments in polarized Madin-Darby canine kidney (MDCK) cells for membrane-volume sorting, polarized sorting and apical recycling." *Traffic*. 2000 1(2):124-40.
2. Gandhi SP, Stevens CF. (2003) "Three modes of synaptic vesicular recycling revealed by single-vesicle imaging." *Nature*. 423(6940):607-13.
3. Diaz G, Isola R, Falchi A, Diana A. (1999) "CO₂-Enriched Atmosphere on the Microscope Stage." *BioTechniques* 27:292-294.
4. Spierenburg GT, Oerlemans FT, van Laarhoven JP, de Bruyn CH. (1984) "Phototoxicity of N-2-hydroxyethylpiperazine-N'-2-ethanesulfonic acid-buffered culture media for human leukemic cell lines." *Cancer Res*. 44(5):2253-4.
5. Zigler JS Jr, Lepe-Zuniga JL, Vistica B, Gery I. (1985) "Analysis of the cytotoxic effects of light-exposed HEPES-containing culture medium." *In Vitro Cell Dev Biol*. 21(5):282-7.
6. Lepe-Zuniga JL, Zigler JS Jr, Gery I. (1987) "Toxicity of light-exposed Hepes media." *J Immunol Methods*. 103(1):145.
7. Berthois Y, Katzenellenbogen JA, Katzenellenbogen BS. (1986) "Phenol red in tissue culture media is a weak estrogen: implications concerning the study of estrogen-responsive cells in culture." *Proc Natl Acad Sci U S A*. 83(8):2496-500.
8. Lelong IH, Rebel G. (1998) "pH drift of "physiological buffers" and culture media used for cell incubation during in vitro studies." *J Pharmacol Toxicol Methods* 39(4):203-10.
9. Giandomenico AR, Cerniglia GE, Biaglow JE, Stevens CW, Koch CJ. (1997) "The importance of sodium pyruvate in assessing damage produced by hydrogen peroxide." *Free Radic Biol Med*. 23(3):426-34.
10. Babich H, Liebling EJ, Burger RF, Zuckerbraun HL, Schuck AG. (2009) "Choice of DMEM, formulated with or without pyruvate, plays an important role in assessing the in vitro cytotoxicity of oxidants and prooxidant nutraceuticals." *In Vitro Cell Dev Biol Anim*. 45:226-33.
11. Pawley, JB. (1990) *Handbook of biological confocal microscopy* / edited by James B. Pawley New York : Plenum Press 232 pages.



NOTES:



Contact and Support

For questions or comments on this or any product from Marker Gene Technologies, Inc., you may contact us by phone or via our website. We welcome customer feedback and we make every effort to improve our products based on input from our clients.

To ask a question or make a comment or suggestion, you can call us at **1-888-218-4062** or fax to **541-342-1960**.

For more information on our products and services, please visit our website at www.markergene.com, where you can find:

- Secure online ordering
- Product Information
- MGT Scientific Newsletters
- Corporate Information
- Custom Synthesis Info

We want to thank you for your purchase and hope that you will continue to order from Marker Gene Technologies, Inc.

MARKER GENE TECHNOLOGIES, INC.
UNIVERSITY OF OREGON RIVERFRONT RESEARCH PARK
1850 MILLRACE DRIVE
EUGENE, OREGON 97403-1992
1-888-218-4062
EMAIL: TECHSERVICE@MARKERGENE.COM