



S200R

Keep cards refrigerated;
Keep other kit components at
room temp

Safety information

Slightly hazardous (irritant, sensitizer) in case
of skin and/or eye contact, always wear gloves and
safety glasses.

Description:

RNA^{Sound}™ RNA Sampling Card
unprecedentedly stabilizes RNA on
filter paper by the proprietary
impregnating lysis buffer. It features:

- Room temperature RNA sample
collection, storage and
transportation;
- Immediate inactivation and thus
safe handling of infectious
agents;
- Easy RNA recovery;
- Integrated RNA sample
collection and RNA recovery

ReadyPunched™ format (Patent
pending) eliminates the tedious card
punching and the risk of cross
contamination.

Kit contents

item	description	quantity
RNA Sampling Cards	Individually packaged in zip bag with desiccant	25
Polyester swabs	For sample application	25
Card Elution Solution	For DNA elution	2.5 mL

Protocol

1. Sample preparation

- 1) Serum, saliva, nasal fluids,
environmental water samples
 - Applied directly

- 2) Cells or bacteria cultures:
 - (For adherent cells) Detach cells and
inactivate trypsin;
 - Cells pelleted down;
 - Cells washed with 1XPBS;
 - Cells resuspended in 1XPBS

2. Sample application on card

- 1) Directly drop sample on the perforated
discs on the card;
 - The two perforated discs on each
card are for more volume of the
same sample.
- 2) Or, collect sample on a cotton swab,
and press and roll the contents of the
swab onto the perforated discs on the
card;
 - The stability of sample RNA is not
guaranteed outside the two
perforated discs.
- 3) Dry the card on a portable Card Drying
station (Cat. # U100) for about 10
minutes or at room temperature for
about an hour.


- 4) Return the card to its original zip bag
with desiccant.
 - RNA are stable at room temp for
at least one week;
 - If accessible, store cards at 4 °C or
lower for longer storage.

2. RNA elution

- 1) Take the card out of the zip bag;
- 2) Further dry the card if necessary on
portable Card Drying station (Cat. #
U100) for 10 min, or in air for an hour;
 - Cards need to be thoroughly dry
to avoid inhibitors to be eluted
off with RNA;
- 3) Push out the two perforated discs into
a 1.5 mL eppendorf tube using a sterile
pipette tip;
 - If the disc hangs on the card, push
the disc against the tube wall, and
pull the card to detach the disc.

- 4) Add 100 µL of provided Card Elution
Solution to each well;
 - Vortex the Card Elution Solution tube
to suspend the beads in solution;
 - Pipette up and down twice before
each aspiration
- 5) Use the pipette tip to submerge the discs
in water;
- 6) Vortex for 1 min;
- 7) Spin the tubes at top speed for 1 min;
- 8) Transfer the supernatant to a new tube;
- 9) Use elute for immediate RT-PCR or store
at -20° C.

Brief workflow



- Squeeze and roll the
polyester swab on
perforated discs;
- Dry the card on portable
Card Drying station (Cat.
U100) for 10 minutes;
or in air for 1 hour

- RNA are stable for
one week;
- Push out discs to a
1.5 mL eppendorf
tube

Add 100 µL of
provided Card
Elution
Solution, vortex
for 1 min, top
speed spin for 1
min to elute
RNA