

UptiTherm DNA Polymerase

A superior polymerase with lower error rate (1/10 Kbp) to use for routine applications


Product Description

UptiTherm Product format (buffer)	Qty	Cat. No. 1 U/μl	Cat. No. 5 U/μl
Standard buffer (including MgCl ₂)	1000 U	UPS53663	UPS53881
Mg free buffer + MgCl ₂ 50 mM	1000 U	UPS53703	UPS53921

Notes:

1- UptiTherm DNA polymerase shows optimal specificity at **2 mM MgCl₂** final concentration (rather than 1.5 mM) in reaction buffer.

2- UptiTherm is also available as 500 U and 250 U package , and as a [gel form](#) (stable at room temperature!).

Store at - 20°C 

Technical and Scientific Information

Highly thermostable DNA polymerase. It is a recombinant, modified form of the enzyme from the thermophilic bacterium *Thermus thermophilus* expressed in *E. coli* (see Note 1).

UptiTherm DNA polymerase is suitable for applications which require a highly thermostable and processive enzyme capable of synthesising DNA strands at elevated temperatures in amplification reactions or similar (e.g. primer extension), thus resolving the most complex secondary structures.

The enzyme is free of unspecific endo- or exonucleases activities, as well as nicking activities. It does not either exhibit significant reverse-transcriptase activity. Terminal transferase activity inherent to the enzyme renders A-tailed amplification products suitable to be further used in T/A cloning approaches.

The enzyme is supplied at a concentration of 5 U/μl and 1 U/μl in a storage buffer. 1 U/μl concentration allows accurate pipetting of small amounts of the DNA polymerase, so that it is not necessary to perform ulterior dilutions.

Unit Definition

One unit is defined as the amount of enzyme which incorporates 10 nanomoles of dNTPs into acid-insoluble DNA within 30 minutes at 72 °C.

Storage buffer

10 mM Tris-HCl (pH 8.0), 50 mM KCl, 1 mM EDTA, 0.1% Triton X-100, 50% glycerol (v/v). Other storage buffers available (contact our Technical Dpt.).

Reaction buffer

Recommended reaction buffer is: 75 mM Tris HCl (pH 9.0), 2 mM MgCl₂ (see Note 2), 50 mM KCl, 20 mM (NH₄)₂SO₄. This reaction buffer (the so-called Standard Buffer, Ordering Information at the end) is supplied at 10X concentration together with the enzyme.

Reaction buffer can be supplied MgCl₂ free (the so-called Free Buffer, see Ordering Information): Mg 2+ ion, being the enzyme cofactor, plays a key role on polymerase activity, this is why its concentration must be optimized in certain amplification-based experiments. In this case, the MgCl₂ is supplied as a separate vial at 50 mM concentration. This solution must be completely thawed, vigorously vortexed and spun down in a bench-top centrifuge before use.

Storage Conditions

Store at -20°C in a constant temperature freezer (i.e. do not use frost-free freezers). Under these conditions the activity of the enzyme remains unaltered over 18 months of storage. The glycerol in the storage buffer prevents freezing at -20°C.

Reaction Conditions

After thawing the reaction buffer (and MgCl₂ solution, in case the “free buffer” choice is adopted), shake all vials (buffer, enzyme, 50 mM MgCl₂ solution) by gentle vortexing, later spin them down in a bench-top centrifuge, and eventually pipette desired volumes.

Keep all reagents on ice while they remain out of the -20°C storage freezer, otherwise enzyme activity will decrease over the time. Wear disposable gloves and make use of sterile, DNase- and RNase-free pipette tips and tubes in order to avoid contaminations and false negative results.

Recommended enzyme volumes to be added to the reaction mix

Final reaction volumes	Recommended enzyme volumes	Volume using 1U/μl items (or a 1U/μl master mix)
100 μl	up to 2.5U	Up to 2.5 μl
50 μl	1-1.25U	1-1.25 μl
25 μl	0.5-0.75U	0.5-0.75 μl

These volumes may vary depending on the reaction conditions and on the particular nucleotide composition of the product to be amplified. Further information may be obtained from our Technical Department.

It is recommended to increase the enzyme units (up to three times) in order to perform certain applications such as PRINS (Primed In Situ Synthesis) or when working on long DNA fragment amplifications (longer than 2 Kb from genomic DNA).

The dNTP final concentration recommended is 200 μM (50 μM each), but this figure may be decreased (e.g. when unspecific amplimers occur), enlarged (e.g. long amplifications) or even unbalanced in favour of any dNTP in particular (e.g. in vitro mutagenesis experiments) depending on the intended approach. This enzyme also accepts modified dNTPs (e.g. radioactively or fluorescein labelled) as substrate.

List of UptiTherm products

Quantity	Standard Buffer 2 mM		Free MgCl ₂ + 50 Mm MgCl ₂	
	1 U/μl	5 U/μl	1 U/μl	5 U/μl
5000 units	S53665	S53885	S53705	S53925
1000 units	UPS53663	UPS53881	UPS53703	UPS53921
500 units	UPS53662	UPS53882	UPS53702	UPS53922

Notes

Note 1: this enzyme is not recommended for certain experiments dealing with amplification of sequences homologous to those found in E. coli or very low-annealing temperature amplification approaches (e.g. RAPDs, Random Amplified Polymorphic DNAs). This enzyme is not recommended for sequencing.

Note 2: at difference with the vast majority of the thermostable DNA polymerases existing in the market, UptiTherm DNA polymerase shows optimal specificity at 2 mM MgCl₂ final concentration (rather than 1.5 mM) in reaction buffer.

References

- Bailly X. et al.**, Globin gene family evolution and functional diversification in annelids, *FEBS J.*, 274: 2641 - 2652 (2007)
- Boutet I. et al.**, A transcriptomic approach of salinity response in the euryhaline teleost, *Dicentrarchus labrax*, *Archimer*, 379:40-50 (2006)
- Boutet I. et al.**, Molecular characterisation of prolactin and analysis of extrapituitary expression in the European sea bass *Dicentrarchus labrax* under various salinity conditions, *Comparative Biochemistry and Physiology Part D: Genomics and Proteomics*, Volume 2, Issue 1, Pages 74-83 (2007) [Article](#)
- Bourin M. et al.**, Sex- and tissue-specific expression of "similar to nothepsin" and cathepsin D in relation to egg yolk formation in *Gallus gallus*, *Poult. Sci.*, 91: 2288 - 2293 (2012) [Abstract](#)
- Charrier G. et al.**, An investigation of the population genetic structure of pollack (*Pollachius pollachius*) based on microsatellite markers, *ICES Journal of Marine Science: Journal du Conseil* 63(9):1705-1709 (2006) [Article](#)
- Guardia S. et al.**, Effects of stocking density on the growth performance and digestive microbiota of broiler chickens, *Poult Sci* 90:1878-1889 (2011) [Article](#)
- Gilbert F. et al.**, Evaluation of tandem repeats for MLVA typing of *Streptococcus uberis* isolated from bovine mastitis, *BMC Veterinary Research*, 2:33 (2006)
- Malki S. et al.**, Prostaglandin D₂ induces nuclear import of the sex-determining factor SOX9 via its cAMP-PKA phosphorylation, *EMBO J.* 24:1798-1809. (2005)
- Marchand J. et al.** : Molecular Identification and Expression of Differentially Regulated Genes of the European Flounder, *Platichthys flesus*, Submitted to Pesticide Exposure, *Marine Biotechnology*, Vol.8, No 3 :275-294 (2006)
- Moussard H. et al.** : Uncultured Archaea in a hydrothermal microbial assemblage: phylogenetic diversity and characterization of a genome fragment from a euryarchaeote, *FEMS Microbiology Ecology*, 57 :452 (2006)
- Roux M. et al.** Characterization of the bovine PRKAG3 gene: structure, polymorphism, and alternative transcripts, *Mammalian Genome*, Vol. 17, No 1 :83-92 (2006)
- Roux M. et al.** Cotranscription and intergenic splicing of the PPARG and TSEN2 genes in cattle, *BMC Genomics*, 7:71 (2006)
- Tesseraud S. et al.**, Insulin and amino acid availability regulate atrogen-1 in avian QT6 cells, *Biochemical and Biophysical Research Communications*, Volume 357, Issue 1, Pages 181-186 (2007) [Abstract](#)
- Uzbekova S. et al.**, Glycogen synthase kinase 3B in bovine oocytes and granulosa cells: possible involvement in meiosis during in vitro maturation, *Reproduction*, 138: 235 - 246 (2009) [Abstract](#)

Related Products

- Mineral oil, #UPS5422A
- dNTP Mix (dATP, dGTP, dCTP, dTTP), #127823
- dNTP Set 1 (dATP, dGTP, dCTP, dTTP), #UP968640
- Mix Sequencing Term. (dNTP+ddNTP), #UP996780
- dNTP Set 2 (dATP, dGTP, dCTP, dUTP), #UP968660
- Agarose, regular uses, #[31272L](#)

Notice to buyers/users:

Some of the applications which may be performed with this product are covered by applicable patents in certain countries. The purchase of this product does not include or provide a license to perform patented applications. Users may be required to obtain a license depending on the country and/or application.

Rev.D0825E