



FT-WT9771

Thiol – PEG_x – Carboxy reagents

Heterobifunctional crosslinkers

HeteroBifunctional SH- PEG – COOH reagents

Product name synonymes	Cat.number Qty 1-100mg, 2-1g 3-5g	MW (g·mol ⁻¹)	Structure
Thiol-PEG_x-NH₂	GV9882	200	 <p>Ask also Thiol-PEO_n-Carboxyl (monodisperse) :</p>  <p>(B4CFS1 : PEG7)</p> <p>#8X331B: Thiol-PEO₃-Carboxyl MW: 194.25 CAS : 1379649-73-6 #8X553B: Thiol-PEO₃-Carboxyl MW: 238.30 CAS : 1347750-82-6 #RPX67B: Thiol-PEO₄-Carboxyl MW: 282.35 CAS : 749247-06-1 #8X482B: Thiol-PEO₅-Carboxyl MW: 326.4 CAS : 1449390-67-3 #AWWUTB:Thiol-PEO₆-Carboxyl MW: 370.46 CAS : 1347750-77-9 #B4CFS1: Thiol-PEO₇-Carboxyl MW: 414.51 . #B4CFT1: Thiol-PEO₈-Carboxyl MW: 458.56 CAS :866889-02-3 #B4CFU1: Thiol-PEO₉-Carboxyl MW: 502.63 . #B4CFV1: Thiol-PEO₁₀-Carboxyl MW: 546.68 . #B4CFW1: Thiol-PEO₁₂-Carboxyl MW:634.77 CAS : 1032347-93-5 #B4CSE0 Thiol-PEO₁₆-Carboxyl MW:810.99 #B4CSF0 Thiol-PEO₂₀-Carboxyl MW:987.2 #B4CFX0: Thiol-PEO₂₄-Carboxyl MW:1163.40 . Colorless liquid (PG2) to White Solid or Colorless Liquid (PG7-9) or Viscous Liquid or White Solid (PEG10-24)</p>
Thiol PEG Carboxyl, HS-PEG-COOH ⁽¹⁾	AWKL32	400	
	Inquire	500	
PG2-CATH	7A4052	600	
HE003017	GV9872	800	
RPW570 – Pack of 2KDa, 3KDa, 5KDa, 10KDa	1N9042	1 000	
	GV9892	2 000	
	WU0602	3 400	
	WT9772	5 000	
	LV7912	7 500-	
	WU0612	10 000	
	WU0622	20 000	
		30 000	
		40 000	

Description:

Heterobifunctional PEG derivative that can be used with Amine or Sulfhydryl reactive chemical groups. PEGylation can modify peptides and proteins and other materials, to create conjugates or to increase solubility and stability and reduce immunogenicity. It can also suppress the non-specific binding of charged molecules to the modified surfaces.

Physical Properties:

Off-white/white solid or viscous liquid depending on molecule weight;
Soluble in regular aqueous solution as well as most organic solvents: water, ethanol, chloroform, DMSO, etc

Storage Conditions:

Store at -20°C (J,M). Keep desiccated. Protect from light. Stable for +12months at -20°C.

Handling and Use:

For best use, material should always be kept in low temperature in dry conditions and under inert gaz for best stability. Prepare fresh solution right before use. Avoid frequent thaw and freezing.

FT-WT9771

Related products

See or [ask](#) for other PEG and PEO reagents
4Arm-PEG-Thiol and 4Arm-PEG-Amine

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Uptima is not liable for any damage resulting from handling or contact with this product. Rev. T09<E-T04E-R04E