

## Fluorimetric NAD/NADH ratio Assay Kit


### Product Description

<b>Name :</b>	<b>Fluorimetric NAD/NADH ratio Assay Kit, red fluorescence</b>	
<b>Catalog Number :</b>	FP-BZ8120	250 assays
<b>Components :</b>	Component A: NAD/NADH Recycling Enzyme Mixture	2 bottles (lyophilized powder)
	Component B: NADH Sensor Buffer	1 bottle (20 mL)
	Component C: NADH Standard	1 vial (142 µg)
	Component D: NADH Extraction Solution	1 bottle (10 mL)
	Component E: NAD Extraction Solution	1 bottle (10 mL)
	Component F: NAD/NADH Control Solution	1 bottle (10 mL)
	Component G: NAD/NADH Lysis Buffer	1 bottle (10 mL)

**Storage:** -20°C      Protect from light and moisture

### Introduction

Nicotinamide adenine dinucleotide (NAD<sup>+</sup>) and nicotinamide adenine dinucleotide phosphate (NADP<sup>+</sup>) are two important cofactors found in cells. NADH is the reduced form of NAD<sup>+</sup>, and NAD<sup>+</sup> is the oxidized form of NADH. It forms NADP with the addition of a phosphate group to the 2' position of the adenyl nucleotide through an ester linkage. NADP is used in anabolic biological reactions, such as fatty acid and nucleic acid synthesis, which require NADPH as a reducing agent. In chloroplasts, NADP is an oxidizing agent important in the preliminary reactions of photosynthesis. The NADPH produced by photosynthesis is used as reducing power for the biosynthetic reactions in the Calvin cycle of photosynthesis. The traditional NAD/NADH and NADP/NADPH assays are done by monitoring of NADH or NADPH absorption at 340 nm. This method suffers low sensitivity and high interference since the assay is done in the UV range that requires expensive quartz microplate.



This NAD/NADH Ratio Assay Kit provides a convenient method for sensitive detection of NAD, NADH and their ratio. The enzymes in the system specifically recognize NAD/NADH in an enzyme cycling reaction. There is no need to purify NAD/NADH from sample mix. The enzyme cycling reaction significantly increases detection sensitivity. In addition, this assay has very low background since it is run in the red visible range that significantly reduces the interference from biological samples. The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation with no separation steps required. Its signal can be easily read by either fluorescence microplate reader with Ex/Em = 530 to 570/590 to 600 nm (maximum Ex/Em = 540/590 nm) or absorbance microplate reader at 576±5 nm. This kit also provides cell lysis buffer for convenience.

## Kit Key Features

<b>Broad Application:</b>	Can be used for quantifying NAD/NADH in solutions, and in cell extracts.
<b>Sensitive:</b>	The kit detects as low as 10 picomoles of NAD/NADH in solution.
<b>Continuous:</b>	Easily be adapted to automation with no separation required.
<b>Convenient:</b>	Formulated to have minimal hands-on time.
<b>Non-Radioactive:</b>	No special requirements for waste treatment.

## Directions for use

### Assay Protocol (for one 96-well plate) :

#### Brief Summary

Prepare NADH standards or test samples (25  $\mu$ L) → Add 25  $\mu$ L of NADH or NAD Extraction Solution  
 → Incubate at room temperature for 15 min → Add 25  $\mu$ L of NAD or NADH Extraction Solution  
 (25  $\mu$ L) → Add NAD/NADH reaction mixture (75  $\mu$ L) → Incubate at RT for 15 min – 2 hr  
 → Read fluorescence at Ex/Em = 540/590 nm

*Note: Thaw 1 vial (or bottle) each of all the kit components to room temperature before starting the experiment.*

#### 1. Prepare NADH stock solution:

- 1.1 Prepare NADH standard stock solution: Add 200  $\mu$ L of PBS buffer into the NADH standard vial (Component C) to have 1 mM (1 nmol/ $\mu$ L) stock solution.

*Note: The unused NADH solution should be divided as single use aliquots and stored at -20 °C.*

#### 2. Prepare NAD/NADH reaction mixture:

- 2.1 Prepare the NAD/NADH reaction mixture: Add 10 mL of NADH Sensor Buffer (Component B) to the bottle of NAD/NADH Recycling Enzyme Mixture (Component A), mixed well.

*Note: This solution is enough for two 96-well plates. The unused NAD/NADH mixture should be divided as single use aliquots and stored at -20 °C.*

#### 3. Prepare serial NADH (0 to 10 $\mu$ M) solutions

- 3.1 Add 30  $\mu$ L of 1 mM NADH standard stock solution (from step 1) to 970  $\mu$ L PBS (pH 7.4) buffer to generate 30  $\mu$ M (30 pmols/ $\mu$ L) standard.

*Note: Diluted NADH standard solution is unstable, should be used within 4 hours.*

- 3.2 Take 200  $\mu$ L of 30  $\mu$ M NADH standard solution (from step 3.1) to perform 1:3 serial dilutions to get 10, 3, 1, 0.3, 0.1, 0.03 and 0 standard NADH solutions.

- 3.3 Add NADH standards and NAD/NADH containing test samples into a 96-well solid black microplate as described in Tables 1 and 2.

*Note: Prepare your cell or tissue samples as desired. Component G (NAD/NADH Lysis Buffer) can be used for lysing the cells for convenience.*

**Table 1.** Layout of NADH standards and test samples in a solid black 96-well microplate:

BL	BL	TS	TS	TS (NADH)	TS (NADH)	TS (NAD)	TS (NAD)						
NS1	NS1	....	....	....	....	....	....						
NS2	NS2												
NS3	NS3												
NS4	NS4												
NS5	NS5												
NS6	NS6												
NS7	NS7												

FT-ZE8120

Note: NS= NAD/NADH Standards; BL=Blank Control; TS=Test Samples; TS (NADH) = Test Samples treated with **NADH Extraction Solution** for 10 to 15 min, then neutralized by **NAD Extraction Solution**; TS (NAD) = Test Samples treated with **NAD Extraction Solution** for 10 to 15 min, then neutralized by **NADH Extraction Solution**.

**Table 2.** Reagent composition for each well:

NADH Standard	Blank Control	Test Sample (NAD/NADH)	Test Sample (NADH Extract)	Test Sample (NAD Extract)
Serial dilutions*: 25 µL	PBS: 25 µL	Test Sample: 25 µL	Test Sample: 25 µL	Test Sample: 25 µL
Component F: 25 µL	Component F: 25 µL	Component F: 25 µL	Component D: 25 µL	Component E: 25 µL
Incubate at room temperature for 10 to 15 min				
Component F: 25 µL	Component F: 25 µL	Component F: 25 µL	Component E: 25 µL	Component D: 25 µL
Total: 75 µL	Total: 75 µL	Total: 75 µL	Total: 75 µL	Total: 75 µL

\*Note: Add the serially diluted NADH standards from 0.03 µM to 30 µM into wells from NS1 to NS7 in

duplicate. High concentration of NADH (e.g., >300 µM, final concentration) may cause reduced fluorescence signal due to the over oxidation of NADH sensor (to a non-fluorescent product).

3.4 **For NADH Extraction (NADH)**, add 25 µL NADH Extraction Solution (Component D) into the wells of NAD/NADH containing test samples. Incubate at room temperature for 10 to 15 minutes, then add 25 µL of NAD Extraction Solution (Component E) to neutralize the NADH extracts as described in Tables 1 & 2.

**For NAD Extraction (NAD)**, add 25 µL NAD Extraction Solution (Component E) into the wells of NAD/NADH containing test samples. Incubate at room temperature for 10 to 15 minutes, then add 25 µL of NADH Extraction Solution (Component D) to neutralize the NAD extracts as described in Tables 1 & 2.

**For Total NAD and NADH**, add 25 µL NAD/NADH Control Solution (Component F) into the wells of NADH standards and NAD/NADH containing test samples. Incubate at room temperature for 10 to 15 minutes, and then add 25 µL of Control Solution (Component F) as described in Tables 1 and 2.

Note: Prepare your cell or tissue samples as desired. Component G (NAD/NADH Lysis Buffer) can be used for lysing the cells.

#### 4. Run NAD/NADH assay in supernatants reaction

- 4.1 Add 75 µL of NADH reaction mixture (from step 2) to each well of the NADH standard, blank control, and test samples (from step 3.4) so that the total NADH assay volume is 150 µL/well.
- 4.2 Incubate the reaction for 15 minutes to 2 hours at room temperature, protected from light.
- 4.3 Monitor the fluorescence increase with Ex/Em = 540/590 nm using a fluorescence plate reader.  
Note1: The contents of the plate can also be transferred to a white clear bottom plate and read by absorbance microplate reader at the wavelength of 576±5 nm. The absorption detection has lower sensitivity compared to fluorescence reading.

#### Data Analysis

The fluorescence in blank wells (with the PBS buffer only) is used as a control, and is subtracted from the values for those wells with the NADH reactions. The typical data are shown in Figure 1 (Total NAD and NADH vs. NAD or NADH Extract).

FT-ZE8120




Figure 1. Total NADH and NAD, and their extract dose response on 96-well black plate were measured with NAD/NADH Ratio Assay Kit using a Gemini microplate reader (Molecular Devices). 25  $\mu$ L of equal amount of NAD and NADH was treated with or without NADH or NAD extraction solution for 15 minutes, and then neutralized with extraction solutions at room temperature. The signal was acquired at Ex/Em = 540/590 nm (cut off 570 nm) after 30 minutes of the addition of 75  $\mu$ L of NADH reaction mixture. The blank signal was subtracted from the values for those wells with the NADH reactions (*Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point*).

## References

- 1- **Ziegenhorn J**, *et al.* Molar absorptivities of beta-NADH and beta-NADPH. *Clin Chem*, 22, 151 (1976)
- 2- **Ikegami T**, *et al.* Structure and Properties of the Recombinant NADH-Cytochrome b(5) Reductase of *Physarum polycephalum*. *Biosci Biotechnol Biochem.* (2007)
- 3- **Kimura N**, *et al.* Comparison of metabolic fates of nicotinamide, NAD<sup>+</sup> and NADH administered orally and intraperitoneally; characterization of oral NADH. *J Nutr Sci Vitaminol (Tokyo)*, 52, 142 (2006)
- 4- **O'Donnell JM**, *et al.* Limited transfer of cytosolic NADH into mitochondria at high cardiac workload. *Am J Physiol Heart Circ Physiol*, 286, H2237 (2004)

## Technical and scientific information

### Related products

See [BioSciences Innovations catalogue](#) and [e-search tool](#).

- Fluorimetric NADP/NADPH Assay Kit, [JQ7300](#)
- Fluorimetric NADH Assay Kit, [JQ7320](#)
- Fluorimetric NADPH Assay Kit, [JQ7330](#)
- Fluorimetric NADP/NADPH Ratio Assay Kit, [ZE8130](#)

## Ordering information

[Catalog size quantities and prices may be found at www.interchim.com/](http://www.interchim.com/)

Please inquire for higher quantities (availability, shipment conditions).

[interbiotech@interchim.com](mailto:interbiotech@interchim.com)

FluoProbes®, powered by  
 213 Avenue J.F. Kennedy - BP 1140  
 03103 Montluçon Cedex - France  
 Tél. 04 70 03 88 55 - Fax 04 70 03 82 60

P.4

FT-ZE8120

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.

**Warning: This kit is only sold for the end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of kit components is strictly prohibited.**