

Amplite™ Fluorimetric NADP/NADPH Ratio Assay Kit

Red Fluorescence

Ordering Information:

Product Number: 15264 (250 assays)

Storage Conditions:


Keep in freezer and avoid light.

Instrument Platform:

Fluorescence microplate readers

Introduction

Nicotinamide adenine dinucleotide (NAD⁺) and nicotinamide adenine dinucleotide phosphate (NADP⁺) are two important cofactors found in cells. NADH is the reduced form of NAD⁺, and NAD⁺ is the oxidized form of NADH. It forms NADP with the addition of a phosphate group to the 2' position of the adenyl nucleotide through an ester linkage. NADP is used in anabolic biological reactions, such as fatty acid and nucleic acid synthesis, which require NADPH as a reducing agent. In chloroplasts, NADP is an oxidizing agent important in the preliminary reactions of photosynthesis. The NADPH produced by photosynthesis is used as reducing power for the biosynthetic reactions in the Calvin cycle of photosynthesis. The traditional NAD/NADH and NADP/NADPH assays are done by monitoring of NADH or NADPH absorption at 340 nm. This method suffers low sensitivity and high interference since the assay is done in the UV range that requires expensive quartz microplate.



This Amplite™ NADP/NADPH Ratio Assay Kit provides a convenient method for sensitive detection of NADP, NADPH and their ratio. The enzymes in the system specifically recognize NADP/NADPH in an enzyme cycling reaction. There is no need to purify NADP/NADPH from sample mix. The enzyme cycling reaction significantly increases detection sensitivity. In addition, this assay has very low background since it is run in the red visible range that significantly reduces the interference from biological samples. The assay can be performed in a convenient 96-well or 384-well microtiter-plate format and easily adapted to automation with no separation steps required. Its signal can be easily read by either fluorescence microplate reader with Ex/Em = 530 to 570/590 to 600 nm (maximum Ex/Em = 540/590 nm) or absorbance microplate reader at 576±5 nm. This kit also provides cell lysis buffer for convenience.

Kit Key Features

- | | |
|---------------------------|--|
| Broad Application: | Can be used for quantifying NADP/NADPH in solutions, and in cell extracts. |
| Sensitive: | The kit detects as low as 10 picomoles of NADP/NADPH in solution. |
| Continuous: | Easily be adapted to automation with no separation required. |
| Convenient: | Formulated to have minimal hands-on time. |
| Non-Radioactive: | No special requirements for waste treatment. |

Note: NS= NADP/NADPH Standards; BL=Blank Control; TS=Test Samples; TS (NADPH) = Test Samples treated with NADPH Extraction Solution for 10 to 15 min, then neutralized by NADP Extraction Solution; TS (NADP) = Test Samples treated with NADP Extraction Solution for 10 to 15 min, then neutralized by NADPH Extraction Solution.

Table 2. Reagent composition for each well:

| NADPH Standard | Blank Control | Test Sample (NADP/NADPH) | Test Sample (NADPH Extract) | Test Sample (NADP Extract) |
|---|--------------------|--------------------------|-----------------------------|----------------------------|
| Serial dilutions*: 25 µL | PBS: 25 µL | Test Sample: 25 µL | Test Sample: 25 µL | Test Sample: 25 µL |
| Component F: 25 µL | Component F: 25 µL | Component F: 25 µL | Component D: 25 µL | Component E: 25 µL |
| Incubate at room temperature for 10 to 15 min | | | | |
| Component F: 25 µL | Component F: 25 µL | Component F: 25 µL | Component E: 25 µL | Component D: 25 µL |
| Total: 75 µL | Total: 75 µL | Total: 75 µL | Total: 75 µL | Total: 75 µL |

*Note: Add the serially diluted NADPH standards from 0.01 µM to 3 µM into wells from NS1 to NS7 in duplicate. High concentration of NADPH (e.g., >100 µM, final concentration) may cause reduced fluorescence signal due to the over oxidation of NADPH sensor (to a non-fluorescent product).

3.4 **For NADPH Extraction (NADPH)**, add 25 µL NADPH Extraction Solution (Component D) into the wells of NADP/NADPH containing test samples. Incubate at room temperature for 10 to 15 minutes, then add 25 µL of NADP Extraction Solution (Component E) to neutralize the NADPH extracts as described in Tables 1 & 2.

For NADP Extraction (NADP), add 25 µL NADP Extraction Solution (Component E) into the wells of NADP/NADPH containing test samples. Incubate at room temperature for 10 to 15 minutes, then add 25 µL of NADPH Extraction Solution (Component D) to neutralize the NADP extracts as described in Tables 1 & 2.

For Total NADP and NADPH, add 25 µL NADP/NADPH Control Solution (Component F) into the wells of NADPH standards and NADP/NADPH containing test samples. Incubate at room temperature for 10 to 15 minutes, and then add 25 µL of Control Solution (Component F) as described in Tables 1 and 2.

Note: Prepare your cell or tissue samples as desired. Component G (NADP/NADPH Lysis Buffer) can be used for lysing the cells.

4. Run NADP/NADPH assay in supernatants reaction

4.1 Add 75 µL of NADPH reaction mixture (from step 2) to each well of the NADPH standard, blank control, and test samples (from step 3.4) so that the total NADPH assay volume is 150 µL/well.

4.2 Incubate the reaction for 15 minutes to 2 hours at room temperature, protected from light.

4.3 Monitor the fluorescence increase with Ex/Em = 540/590 nm using a fluorescence plate reader.

Note1: The contents of the plate can also be transferred to a white clear bottom plate and read by absorbance microplate reader at the wavelength of 576±5 nm. The absorption detection has lower sensitivity compared to fluorescence reading.

Data Analysis

The fluorescence in blank wells (with the PBS buffer only) is used as a control, and is subtracted from the values for those wells with the NADPH reactions. The typical data are shown in Figure 1 (Total NADP and NADPH vs. NADP or NADPH Extract).




Figure 1. Total NADPH and NADP, and their extract dose response on 96-well black plate were measured with Amplite™ NADP/NADPH Ratio Assay Kit using a Gemini microplate reader (Molecular Devices). 25 μ L of equal amount of NADP and NADPH was treated with or without NADPH or NADP extraction solution for 15 minutes, and then neutralized with extraction solutions at room temperature. The signal was acquired at Ex/Em = 540/590 nm (cut off 570 nm) after 30 minutes of the addition of 75 μ L of NADPH reaction mixture. The blank signal was subtracted from the values for those wells with the NADPH reactions (*Note: The fluorescence background increases with time, thus it is important to subtract the fluorescence intensity value of the blank wells for each data point*).

References:

1. Ziegenhorn J, Senn M, Bucher T. (1976) Molar absorptivities of beta-NADH and beta-NADPH. Clin Chem, 22, 151.
2. Ikegami T, Kameyama E, Yamamoto SY, Minami Y, Yubisui T. (2007) Structure and Properties of the Recombinant NADH-Cytochrome b(5) Reductase of Physarum polycephalum. Biosci Biotechnol Biochem.
3. Kimura N, Fukuwatari T, Sasaki R, Shibata K. (2006) Comparison of metabolic fates of nicotinamide, NAD⁺ and NADH administered orally and intraperitoneally; characterization of oral NADH. J Nutr Sci Vitaminol (Tokyo), 52, 142.
4. O'Donnell JM, Kudej RK, LaNoue KF, Vatner SF, Lewandowski ED. (2004) Limited transfer of cytosolic NADH into mitochondria at high cardiac workload. Am J Physiol Heart Circ Physiol, 286, H2237.

Warning: This kit is only sold for the end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at info@aatbio.com if you have any questions.



211 bis Avenue Kennedy - BP 1140
03103 Montluçon - France
33 (0) 4 70 03 88 55
Fax 33 (0) 4 70 03 82 60
e-mail interchim@interchim.com

Agence Paris - Normandie
33 (0) 1 41 32 34 40
Fax 33 (0) 1 47 91 23 90
e-mail interchim.paris@interchim.com