



interchim®

Technical Information

Contents

• Column Handling	218-220
• Mobile Phases for RP-Columns	221
• HPLC Column Performance	222
• Inspection Reports	223
• FAQ	224-227
• Troubleshooting	228-231
• Basic Data	232
• Solvent Miscibility Table	233
• Column Selection Guide	234-235
• Linear Scale-Up	236
• Preparative Column Selection Guide	237
• Substance Index	238-246

Introduction

Technical Information

YMC produces chromatography packing materials and HPLC columns under very strict quality control procedures and delivers to customers only those products that pass the strict Quality Assurance tests prior to shipment. In order to ensure the best performance and long column life, the following instructions should be followed for all packed columns.

Column Handling

- shipping solvent
- mobile phase considerations
- mobile phase replacement and column cleaning
- guard columns
- column back pressures
- temperature



1. Shipping solvent

The solvent used for shipping the column is described on the COLUMN INSPECTION REPORT or in the COLUMN CARE AND USE INSTRUCTIONS leaflet which is included with each column. Please determine the miscibility of this solvent with the mobile phase to be used in

your analysis to prevent immiscibility problems. If you intend to store the column for any length of time, you should replace the mobile phase in the column with the shipping solvent or solvent specified in the column inspection report.

2. Mobile phase considerations

Reversed phase columns can be used with both aqueous and nonaqueous solvents. However, repeated alternating between solvents with extremely different polarities can result in loss of column performance. Typical general organic solvents include acetonitrile, methanol and THF.

Cyano columns can be used in both normal and reversed phase modes. However, a column should be dedicated for use in only one separation mode and not switched between normal and reversed phase modes as this can result in loss of column performance. When using the column in a normal phase mode, replace the solvent in the column with isopropanol. (Make sure that the flow rate is set so that the pressure does not exceed 15 MPa during solvent exchange.)

Silica columns are usually used with nonaqueous solvents such as n-hexane, chloroform or other weak solvents with the addition of isopropanol,

ethyl acetate or similar as appropriate to allow elution of high polarity components.

All Amino columns (ie both Polyamine II and YMC-Pack Amino) can be used with both aqueous and nonaqueous solvents. However, repeated alternating of solvents with extremely different polarities can result in loss of column performance

Solvent should flow in the direction of the arrow (as indicated on the column label) for normal use, although reversed flow for washing will not affect column stability.

The pH ranges for stability of every type of column varies by product. For specific information please refer to the instruction manual included with each column. Should this instruction manual be available, please contact YMC or your local distributor. A general overview of the characteristics of each phase is shown on the inside of the envelope.

3. Mobile phase replacement and column cleaning (general methods)

a) Reversed phase columns

When a mobile phase which contains no buffers or salts is used, wash the column with an eluent consisting of the same solvents as that of the mobile phase, but with a higher organic solvent concentration.

When a mobile phase containing buffers or salts is used, this should first be replaced with an eluent containing the same ratio of water and organic solvent as the mobile phase but which has no buffer or salt components. If the concentra-

tion of buffer or salts used is less than 100 mM, it can be replaced directly with approximately 60% acetonitrile in water.

After using a column near the usable pH limit, washing the column with water alone may cause column deterioration. Instead wash the column with a mixture of water and organic solvent containing no buffer or salt components or alternatively 60% acetonitrile in water to remove the aggressive pH eluents.

Column Handling

Should the column back pressure increase, wash the column in the reverse direction (the opposite direction of the arrow shown on the column label) making sure that the detector is not in line with the solvent stream. A solution having the same composition as that of the mobile phase, but with a higher organic solvent concentration and no added salts or buffers is usually used as the cleaning solution. However consideration should be given to the characteristics of sample so that a solvent which easily dissolves the sample is chosen.

When macromolecules, including proteins and sugars, adsorb to the column, it is usually difficult to wash them off with organic solvents. When columns are used to analyse samples containing

such macromolecules, it is preferable to pretreat the sample and/or use a guard column

b) Normal phase columns

Wash the column with a solution having the same solvents as that of the mobile phase, but with an increased content of high polar component concentration. If polar compounds absorb on the column, flush with isopropanol or similar solvent.

Before storing a column used with a mobile phase containing acid or alkali is used, replace the eluent with a simple solvent or solvent:water mixture. (for example replace n-hexane/isopropanol/acetic acid (90/10/0.1) with n-hexane/isopropanol (90/10) for storage).

RECOMMENDED COLUMN CLEANING AND REGENERATING PROCEDURES

Use the cleaning routine that matches the properties of the column and what you believe is contaminating it. Flush columns with 20 column volumes (80 ml total for 4.6 x 250 mm column) of HPLC-grade solvents. Run columns in reverse flow direction, with the outlet disconnected from the detector. Cleaning efficiency is increased by increasing mobile phase temperature to 35–55°C. If the column performance is poor after regenerating and cleaning, call us.

Silica-based particles

Non-polar-bonded phases (Carotenoid, C18, Octyl, YMCbasic™, J'sphere™, Phenyl, Butyl, TMS):

Polar Samples	Non-polar Samples	Proteinaceous Samples
1. Water	1. Isopropanol	Option 1: Inject repeated allquots of DMSO
2. Methanol	2. THF	Option 2: Gradient of 10 to 90% B where:
3. THF	3. Dichloromethane	A = 0.1% TFA in water
4. Methanol	4. Hexane	B = 0.1% TFA in CH ₂ CN
5. Water	5. Isopropanol	Option 3: Flush column with 7M guanidine
6. Mobile phase	6. Mobile phase	HCl, or 7M urea

Polar-bonded phases (Cyano, Diol, Amino, PVA-sil™, Silica):

Polar Samples	Non-polar Samples
1. Water	1. Chloroform
2. Methanol	2. Methanol
3. THF	3. Dichloromethane
4. Methanol	4. Heptane or Isocyanate
5. Water	5. Isopropanol
6. Mobile phase	6. Mobile phase

Polymer-based particles: Polymer C18™

1. Flush column with mobile phase but omit buffers or salts (i.e. just organic and water, acetonitrile is preferable)
2. Run a gradient to 100% organic
3. Flush with twenty column volumes of THF
4. Flush with twenty column volumes of acetonitrile
5. Run a gradient back to starting mobile phase conditions, omitting buffers and salts
6. Re-equilibrate in mobile phase

4. Guard columns

YMC recommends that you always use a guard column with the same packing material and of the recommended inner diameter for your column (see table).

A YMC guard column is normally composed of a cartridge holder and a guard cartridge. The cartridge holder can be used repeatedly.

Where different cartridge lengths are available, only chose the longer cartridge when samples containing high levels of contaminants are present to increase the time between cartridge changes. Guard cartridges should be changed frequently in order to maximise their protection of the main column. Cartridge holders should be connected to the main column using the shortest length of tubing possible. This tubing should be of an appropriate inner diameter for the flow rate and pressure to be used.

Samples containing particulate matter MUST always be pre-filtered (at least 0.45 µm but 0.2 µm is preferred) before being injecting onto a column.

Column ID (mm)	Recommended Guard Cartridge ID (mm)
1.0	1.0
2.1	2.1
3.0	3.0
4.0	4.0
4.6	4.0* (4.6*)
10	10
20	20
30	30

* Formerly, external guard columns of 4.6 mm ID exhibit dead volume resulting in poor performance of the column system. Additionally, these guards were rather expensive.

In the course of intense testing of the compatibility of different hardware concepts no negative influence of a 4.0 mm ID guard cartridge combined with a 4.6 mm ID main column was observed. Thus, we recommend to use 4.0 mm ID guards with a 4.6 mm ID main column.

Column Handling

5. Column Back Pressures

Column back pressure is a function of several parameters, including :-

Particle size and distribution

Packing porosity and bonded phase coating levels

Column length and inner diameter

Solvent flow rate, viscosity and temperature

Typically for a column packed with 12 nm, 5 µm ODS phase and pumped at ambient temperature with methanol:water (70:30) at 1 ml/min the back pressure should be less than 25 MPa (250 bar, 3750 psi) for 250 x 4.6 mm ID.

For wide pore (20 or 30 nm) 5 µm ODS phase and pumped at ambient temperature with methanol:water (70:30) at 1 ml/min the back pressure should be less than 17 MPa (170 bar, 2550 psi) for 250 x 4.6 mm ID.

We recommend using a column at below the maximum operating pressure to ensure maximum column life.

	Maximum Operating Pressure		
	(MPa)	(bar)	(psi)
YMC-Triart 1.9 µm	100	1000	15000
YMC UltraHT 2.0 µm	50	500	7500

Column ID (mm)	Maximum Operating Pressure		
	(MPa)	(bar)	(psi)
1.0	20/25 ¹	200/250 ¹	3000/3750 ¹
2.1	20/25	200/250	3000/3750
3.0	20/25	200/250	3000/3750
4.0	20/25	200/250	3000/3750
4.6	20/25	200/250	3000/3750
6.0	20/25	200/250	3000/3750
8.0	20/25	200/250	3000/3750
10	10	100	1500
20	10	100	1500
30	10	100	1500
50 and above	7.0/10	70/100 ²	1000/1500 ²
YMC-Actus 20/30	30	300	4500
YMC-Actus 50	10	100	1500

¹ The first figure is for up to 150 mm column length, the second figure is for 250 mm column length

² depending on certification of the column; >10 µm packings should be kept below 3 MPa (30 bar, 450 psi)

Note: PolymerC18 has a pressure limit of 14.5 MPa (145 bar, 2000 psi)

6. Temperature

The upper temperature limit for silica and bonded phases is 50°C (70°C for YMC-Triart at pH=7 or lower). However YMC recommends using columns between 20 and 40°C because certain conditions of pH or mobile phase composition

may affect column lifetime. For recommended column temperatures for other column types, please refer to the instruction manual included with each column.

Mobile Phases for RP-Columns

Mobile phases for reversed phase columns

The composition of mobile phase greatly affects the separation in HPLC. To optimise a separation, it is necessary to consider the interaction between the solutes, stationary (or solid) phase, as well as the mobile phase.

For reversed phase columns, the most commonly used in HPLC, various mobile phases are available. Attention needs to be paid to a number of points when deciding on the mobile phase composition. The variables factors to be considered include:

- miscibility of solvents
- effects on detection methods (eg., UV or MS)
- effects on the column
(column deterioration due to pressure or pH)
- separation reproducibility
- stability of solutes

Typical solvents for ODS columns and some helpful tips for establishing optimum separation conditions are described below.

General solvents

Water, acetonitrile, methanol and tetrahydrofuran (THF) are the important solvents for use with reversed phase columns.

It is important to use high purity water purified by ion-exchange, distillation, reverse osmosis, etc. The presence of organic substances or ionic impurities may cause problems, including ghost peaks during short wavelength UV detection or high sensitivity gradient elution systems.

Acetonitrile is frequently used as an HPLC solvent, due to its low UV absorption and low viscosity. Methanol has a higher viscosity and often shows different separation selectivity to

that obtained using acetonitrile. THF is used occasionally to influence selectivity in conjunction with acetonitrile and methanol, due to the cyclic ether structure of THF. THF has several adverse properties for a solvent for HPLC; it has:-

- significant UV absorption
- high viscosity
- a tendency to form peroxides, especially as the use of antioxidants can give rise to ghost peaks

Appropriate separating conditions can be obtained by using these three solvents plus water individually or in combination.

Buffers and reagents

Acetic acid, formic acid, phosphoric acid and trifluoroacetic acid (TFA) are generally used as acidic modifiers. The buffers normally used include phosphate and acetate buffers (sodium, potassium, ammonium). Monobasic phosphates provide a pH of 4.6 and are used as convenient pH adjusters rather than buffers.

In order to separate ionic compounds, such as amines and carboxylic acids, with good repeatability, the pH of mobile phase must be adjusted so that it is 1 (or preferably 2) units away from the pKa of the solute. At or near the pKa, peak broadening or splitting may be observed as the free acid/base and its salt coexist.

Most buffers are used at a concentration of about 10 mM. However, depending on dissociation of solutes and interactions with the stationary phase, this can be raised to 50 to 100 mM.

When acids or alkalis which degrade reversed phases are used, caution must be taken regarding their concentrations and pH. TFA and phosphoric acid are usually used at concentrations of 0.1% or less. Acetonitrile/water (approx 60/40) solution is a convenient storage solvent after use of acids or buffers (salts).

Tetrabutylammonium salts and sodium perchloric acid may be used as ion pair reagents for retention of highly polar compounds on reversed phase columns or for improvement of separation and peak shape. When these additives are used, it is necessary to use a reagent with the shortest alkyl chains available. If sodium dodecylsulfate, (SDS; which contains long alkyl chains) is used, it may be retained on the column phase and can cause problems with repeatability.

Other solvents for HPLC

Ethanol, 2-propanol, ethyl acetate, or chloroform may be used in the mobile phase (particularly in normal phase separations) in order to improve retention or separation of solutes. In some cases, hexane is used as a mobile phase. When a hy-

drophobic solvent is added to a mobile phase, care must be taken with regards the miscibility with the mobile phase existing in the column and a separate wash stage included before changing the eluent.

t_0	Void volume, Column dead-time	
t_R	Retention time	
h	Peak height	
$W_{0.5h}$	Peak width at half-height	
N	Theoretical plate number	$N = 5.54 \times (t_R / W_{0.5h})^2$
k'	Capacity factor	$k' = (t_R - t_0) / t_0$
α	Separation factor	$\alpha = k'_2 / k'_1$
R_s	Resolution	$R_s = 1.18 \times (t_{R2} - t_{R1}) / (W_{0.5h1} + W_{0.5h2})$
A_s	Asymmetry factor	$A_s = b/a$
T_f	Tailing factor	$T_f = W_{0.05h} / 2f$

Inspection Reports

Formation

YMC employs strict quality control of packing materials to ensure lot-to-lot and column-to-column reproducibility. All packed columns are subject to performance tests and only those columns which meet strict specifications are shipped to customers.

A test report (see below) is shipped with each column. The test method shown on this inspec-

tion report is not only a method used for column performance evaluation but can also be used as a test method for determination of column life. We provide full details of all the analytical conditions used test method, including compounds analysed, sample concentration, eluent composition, etc. to allow the end user to reproduce these tests.

YMC HPLC COLUMN INSPECTION REPORT

Serial No.: 0315002014 (W)

Product code: TA12S05-1503WT

Name, Particle: YMC-Triart C18/ S-5 μm / 12 nm

Gel lot: 9233

Column size: 150 \times 3.0 mm I.D.

Eluent: acetonitrile/water (60/40)

Flow rate: 0.4 mL/min


Temperature: ambient

Detection: UV at 270 nm

Injection volume: 2 μL

Pressure: 3.8 MPa

Shipping solvent: acetonitrile



Sample Components	Retention time [min]	Capacity factor $[k' = (t_R - t_R1)/t_R1]$	Theoretical plates $[5.54 \times (t_R/W0.5)^2]$	USP Tailing factor $[TF = W0.05/2f]$
1. Uracil (0.01 mg/ml)	1.60	0.00	2648	1.36
2. Methyl benzoate (0.1 μL /ml)	4.93	2.09	11325	1.05
3. Naphthalene (0.04 mg/ml)	9.75	5.10	12835	0.99
4. Butyl benzoate (0.3 μL /ml)	12.99	7.13	13252	1.01

[System No.104] [Inspected by T.OKADA]

YMC Co., Ltd.

Frequently Asked Questions



What is “Endcapping”?

Conventional ODS (C18) packing materials are silica gel bonded with octadecyl groups. This is the result of reaction between silanol groups on the silica surface and octadecyl groups. However some active silanol groups remain after the reaction. It is impossible for all the silanol groups to react because of steric hindrance of octadecyl groups. Such residual silanol groups create a secondary interaction in chromatography, which, in many cases, affects on chromatograms by causing peak tailing of basic compounds or irreversible absorption to the column. Therefore, a secondary silanisation reaction with residual silanol groups using a small reagent (typically trimethylsilane) should be performed. This process is called “endcapping”.



Are there ODS columns used with 100% aqueous mobile phase?

Hydrosphere C18 and YMC-Pack ODS-AQ columns can be used with 100% aqueous mobile phase. With conventional ODS columns, retention time becomes shortened due to the incompatibility between water molecules and the silica bonded surface with high hydrophobicity. Water tends to be expelled from the pores on material and the C18 chains “collapse” onto themselves. The retention time is hardly affected for Hydrosphere C18 and YMC-Pack ODS-AQ columns because the silica surface is capable of solvation between mobile phase and hydrophilic silica surface as a result of the reduced C18 functional group density and the proprietary derivatisation process.



What is the upper limit of column pressure?

Column length of 150 mm or less and diameters less than 10 mm:
20 MPa, (200 bar, 3000 psi)

Column length of 250 mm or greater and diameters less than 10 mm:
25 MPa, (250 bar, 3750 psi)



How should we store the columns?

When columns are not used for a long time, they should be stored in a cool place after replacing the eluent with the shipping solvent as described in the Inspection Report. Do not store the column in the mobile phase with salt or acid, even for very short times. Close the airtight stopper tightly to prevent the solvent from evaporating.



How can we evaluate the performance of columns?

Repeat the performance test using exactly the same conditions as the Inspection Report which accompanies the column at the time of purchase. Columns which show no change in retention time, theoretical plate number, peak asymmetry, etc are acceptable for further use.

Columns which show no change in these parameters after several years of use may, however, have changes in separation characteristics for certain types of compounds such as ionic species. It is advisable to avoid using such columns for method development as reproducibility compared to new columns may not be possible.

Frequently Asked Questions

1. To remove strongly adsorbed hydrophobic material; pump the column in the reverse direction with eluent with a greater elution ability than mobile phase. For example, for cleaning reversed phase columns, use an eluent with increased ratio of organic modifier and flush the column with at least 10 column volumes.

2. To recondition of gel surface condition caused as a damage resulting in generation of active silanol groups, and observed as irregularities in peak asymmetry and retention time. Washing with acidic solvents can be effective in such cases. Typically a mixed solvent of 0.1% aqueous phosphoric acid solution and organic solvent (between 10 and 60% organic content) can return the silanol groups to the dissociation state.

**How do I
clean the
columns?**



YMC recommend the use of guard columns, particularly if the samples being analysed contain a high level of contaminants. This will extend the useful lifetime of a column, particularly if replaced at frequent intervals. We recommend that guard columns are packed with the same packing materials as the analytical column. Guard columns with different material may cause abnormalities in peak asymmetries and reproducibility. YMC guard cartridges are particularly economic when frequent replacement is required.

**Do we need
guard
columns?**



Recommended flow rates for semi-micro column (1.0 to 3.0 mm inner diameters) are:-

1.0 mm ID	0.05 ml/min
2.1 mm ID	0.2 ml/min
3.0 mm ID	0.4 ml/min

This can be increased if the column length is short and the system back pressure is low.

Such columns can be used in conventional HPLC systems but it is advisable to use short lengths of smaller diameter connection tubing and detector flow cells which are optimised for low flow rates.

**What is
required in
system and
flow rate for
using
semi-micro
columns?**



Step 1: Determine separation conditions by using analytical columns.

Step 2: Study the preparative scale. Select the particle size of the packing material and the inner diameter of column appropriate for the sample volume.

Step 3: Optimise the separation conditions using analytical columns with inner diameter of 4.6 mm or 6.0 mm packed with the packing material selected for the preparative separation (scout column). If the particle size of the packing material is the same as in the Step 1, this process can be omitted. If the preparative column is more than 100 mm ID, it is advisable to insert another step with a scout column of 20 mm ID in order to accurately predict loadability and calculate the running costs.

Step 4: Proceed with the preparative separation.

**How do I
carry out a
scale-up of a
method?**



Frequently Asked Questions



What should I do when the column pressure rises up?

Depending on the reasons for increased pressure, the following procedures are recommended:-

Blocked frits: Flush the column in reverse flow as described on page 342-343. Reduce the flow rate in order to keep the column pressure within recommended limits whilst flushing the column.

Contamination of the packing material: Wash the column in reverse flow as described on page 218-219.

If pressure increases occur frequently despite treatment as above, it is recommended that sample pretreatment or the use of guard columns is employed to prevent the problem occurring in the first place.



What are the solutions for poor peak shapes?

The following solutions are recommended, depending on the cause.

Inappropriate Mobile phase: If pKa of the analyte and pH of mobile phase are close for ionic analytes, it will result in poor peak shape. Set the pH of mobile phase at least 1 (or better 2) units from pKa.

Effect of solvent used to dissolve sample: If the dissolving solvent of sample and mobile phase are not the same, it causes defects in the peak shape. Dilute the sample solution with mobile phase or reduce the injection volume.

Overloading sample injection: Overloading the column will cause defects in the peak shape. Reduce injection volume and/or the sample concentration.

Insufficient equilibration time: When the difference in pH between the current and a previous mobile phase is wide or the buffer concentration of mobile phase is low, column equilibration may take some time

Column contamination and degradation: If the column is contaminated, wash the column as described on page 225. If the column is degraded, it is not possible to regenerate it and it should be replaced.

System problems: Dispersion of the sample may occur within the tubing between the injector and the column or within the flow cell of detector which can result in peak tailing and/or broadening. Optimisation of the system for use with semi-micro use should be performed.



What are the solutions for ghost peaks?

The following solutions are recommended, depending on the cause.

Injector fouling: If the ghost peak(s) appears when injecting only mobile phase (no sample), wash the injector.

Gradient Analysis: When hydrophobic impurities are eluted by a stronger solvent, they appear as ghost peaks. Clean the column as described in the Instruction Manual. If this does not eliminate them, they are probably due to impurities of solvent. Use a higher grade solvent, purified specifically for HPLC or alternatively install a guard column between the solvent delivery pump and the mixing chamber or injector.



What should I do if columns dry out?

Flush the column with a solvent such as methanol for all bonded phase silica or hexane for non bonded silica and remove trapped air using a flow rate such that the column pressure is about half that normally used for analysis. When all the air has been removed, check the column performance by running a test chromatogram under the conditions stated on the original Column Inspection Report.

Frequently Asked Questions

This can arise for a number of reasons:-

Inappropriate mobile phase conditions:

It may become difficult to obtain reproducibility when analysing ionic compounds if the pH of mobile phase is not controlled or the buffer concentration is low. Increase the buffer concentration.

Retention time can fluctuate widely due to a slight variance of pH when the pH of a mobile phase is set too close to the pKa of analyte. Set the pH of the mobile phase to be at least 1 (or preferably 2) units away from the pKa.

System variance: It may be difficult to obtain reproducibility in chromatograms when using different HPLC systems. Where possible the manufacture of pumps, detectors and injectors should be the same, otherwise differences in extra column volume from mixing chamber, detector cell and plumbing will result in poor reproducibility between systems. Also, with column heaters from different manufacturers, there may be an effect on the retention time due to the set temperature being different between the 2 systems. Use of the same system throughout a sequence of analysis is to be recommendable.

Column histories: Reproducibility between chromatograms may not be obtained when using different columns of the same type. This is due to differences in the columns' prior histories. For example, changes in the chemistry of the surface of the packing material can arise by use of mobile phases containing ion pair reagents or when strongly hydrophobic material (especially proteins) becomes adsorbed on the column. Dedicating a column to a specific application is recommended.

Using 100% aqueous mobile phase: Reproducibility of chromatograms obtained on conventional ODS columns will not be obtained when using 100% aqueous mobile phase due to the short retention times obtained. Columns which can be used in 100% aqueous mobile phase are recommended. YMC recommends the use of either Hydrosphere C18 or ODS-AQ which are designed to be used in 100% aqueous mobile phase.

Grade difference in mobile phase: Reproducibility between chromatograms may not be obtained when using different grade of solvent in a mobile phase. Impurities contained in a solvent can act like salts in mobile phase and affect the separation. Solvent in HPLC grade is recommendable.

What should I do if the column fails to provide reproducibility?



This is caused by excess of ion pair reagent. In general, the higher the concentration of ion pair reagent, the greater the retention. However if the concentration of ion pair reagent is above a certain level, the retention may become poor because of micelle formation. Good separations are achieved when the concentration of ion pair reagent is between 5 mM to 20 mM. YMC recommend that the lowest possible concentration is used to avoid short column life.

I still have poor retention after adding ion pair reagent to mobile phase. Why?



Troubleshooting

1. Consideration of solvent grade for reversed phase LC

Reversed phase liquid chromatography frequently employs organic solvents such as methanol, acetonitrile or tetrahydrofuran. Although HPLC grade products of these types of solvents are available, it seems some users have trouble when using a reagent grade solvent instead of HPLC grade. This results in them wasting considerable amounts of time. How do the two solvent grades differ?

Methanol and acetonitrile

Reagent grade solvents contain larger quantity of impurities UV absorbing than HPLC grade solvents do, which makes it difficult to use them for gradient elution or trace analysis, especially when the detection requires short wavelength. This gives rise to significant increases in baseline noise

or detection sensitivity. In some cases (or at some wavelengths) it might be possible to use a reagent grade solvent, but we recommend the use of HPLC grade solvents whenever possible.


Tetrahydrofuran

Tetrahydrofuran easily generates peroxides. To prevent this, the solvent generally contains antioxidants which can cause ghost peaks. As a result solvent containing no antioxidants should be used in HPLC. The peroxides in tetrahydrofuran also have a marked effect on the baseline stability (with differences between grades and between different suppliers being greater than for other organic solvents), which leads to the recommendation that HPLC grade solvent with little or no impurities should be used.

2. Eluent conditions

Although a column is frequently thought to be the cause of HPLC analysis not providing the correct trace, many failures are attributed to causes other than the column, including improper maintenance operations. This discussion illustrates the case in which the grade of a solvent affects the peak shapes. In the chromatogram for basic compound analysis, using an eluent of acetonitrile/water, Peak 2 represents the basic compound.

Figures right show chromatograms from two identical operations except that the acetonitrile used was of different grades. One was HPLC grade (Figure 1); the other was reagent grade (Figure 2). While the peak shape was broadened with HPLC grade acetonitrile, it was much improved when using reagent grade. The differences in peak shapes which were observed were also found to be dependant on the different makers even though they were of the same specific grade. This may be the effect of traces of impurities contained in acetonitrile behaving in the same way as modifiers added to an eluent. Replacing eluent with acetonitrile / 5 mM ammonium acetate produced



the chromatogram shown in Figure 2 irrespective of the grade of solvent. To avoid the influence of different grades, solvents specifically made for HPLC must be used. Even compounds which have groups which can dissociate can be analysed with eluent containing no acid or salt, although eluents with additives such as salt must be used when reproducibility is important.



Troubleshooting

3. Peak shape anomaly

A common problem encountered during HPLC operations is peak shape anomalies such as peak tailing and double peaks. In order to remove these problems, the cause must be precisely determined. The majority of cases are the result of inappropriate conditions for the separation; including inappropriate selection of column or solvent, or use of an old column which has a void to the packing at the top of the packing. Here we discuss the method of determining the cause of the problem with peak shapes.

The simplest way is to test the column performance using the "shipping inspection criteria" as

described in the column inspection report which is included with every column. If the examination reveals no peak shape anomaly, then the cause will be the result of inappropriate selection of separation condition. The separation condition such as eluent selection must be reconsidered. If, on the contrary, the same examination reveals any anomaly, the column may be the problem. Flushing (to remove the impurity could have accumulated on the column) or replacement of the column is necessary. We recommend examining column performances on a regular basis and always under the identical conditions.



YMC provide sufficient analytical information, including sample concentration, in the column inspection reports to allow customers to evaluate the performance of the column using standard compounds.

4. Column Pressure Increases

Pressure increase is a common problem in HPLC. Some of the reasons for pressure increases in reversed phase chromatography are discussed below.

If the system pressure increases, you should first disconnect the column and run the system to determine the line pressure. If the line pressure is high, the tubing may be clogged or damaged. If there is no excessive line pressure, then the column pressure may be high and the column needs cleaning. Cleaning by pumping in the reversed direction can be very effective. Generally the relative proportion of the organic solvent in a mobile phase should be increased when washing, to speed up removal of bound material. However the key consideration is to choose, in accordance with the characteristics of the sample, an appropriate solvent that will easily dissolve the adsorbed

material and not cause precipitation. Reversed phase separations often cause protein to be adsorbed by the packing material which results in high pressure. This problem can be overcome effectively by gradient washing with acetonitrile/water containing 0.1% TFA, rather than washing with an organic solvent. If the cause of high back pressure is believed to be the result of insoluble material in samples or precipitation of a sample during separation, washing or replacing the inlet frit may be successful.

However, once high back pressures occur, it frequently becomes difficult to restore performance despite washing, etc. It is far better to prevent increased column pressure from occurring by simple sample preparation such as protein removal or filtration and using a guard column to protect the analytical column.

Troubleshooting

5. The Cause of the Ghost Peaks

As part of a test of a gradient method a chromatogram was run without a sample being injected. A number of peaks were obtained, as in trace (A). When a similar test was performed, but with the column disconnected, the ghost peaks disappeared, as in trace (B). This led to the idea that the column was at fault.


However, despite flushing and replacing the column, the baseline could not be improved. Several other factors were then examined; the cause was found to be water used to prepare the mobile phase. Standard distilled water (which is inadequate for HPLC) had been mistakenly used. When HPLC grade distilled water was used, an excellent baseline as in trace (B).

Water purity can have a great impact on gradient elution. Even HPLC grade distilled water will become contaminated with time, causing ghost peaks. This will have no significant influence on isocratic elution methods but it will cause problems in gradient elution methods.

In gradient elution methods, a column is equilibrated with an eluent with low organic content. This allows impurities in the eluent to be adsorbed and concentrated in the column. After starting the analysis, the amount of organic solvent increases and impurities begin to elute from the column,

resulting in ghost peaks. The heights of the ghost peaks are dependent on the duration of equilibration (the amount of contaminant adsorbed during equilibration).

Such ghost peaks do not appear when the column is disconnected because there is nothing to adsorb and concentrate the impurities during the equilibration stage. Therefore in gradient analysis, the grade and storage conditions of all solvents requires great care.



6. pH Adjustment of Eluents

Analysis of ionic compounds by reversed phase HPLC has to be performed with the pH of eluent controlled using acid or buffering agent. However, separation at a pH which is not the optimum for the compound of interest can cause problems such as double peak or peak broadening. Even if the peak shape is satisfactory, retention time reproducibility may not be obtained in some cases. The relation between retention of benzoic acid and pH value is shown in the figure below. Although the retention time (measured as k') varies little when the pH is in the range 2 - 3.5, it varies widely when the pH ranges is in the range 3.5 - 4.5. The pKa of benzoic acid is 4.2 and it is noticeable that the region where the retention time varies most widely is near the

pKa. If the eluent pH is adjusted to a value near the pKa, the results may not be reproducible due to very small variations of the pH adjustment having a large impact on the retention time. In fact, the eluent pH variation of just 0.1 will affect the separation significantly. Therefore, it is recommended that the eluent pH should be more than 1 unit away from the pKa.

If the pKa of the analyte is unknown, the eluent pH should be adjusted to within the value where the impact on the separation seems minimal, after having evaluated the relation between the eluent pH and the retention time by using several eluents with their pH values adjusted to be slightly different from each other.

Troubleshooting


7. Regenerating Columns

In reversed phase HPLC, column deterioration can cause poor peak shapes and/or reduced retention times. The column deterioration is the result of changes in the packing material's structure, such as the loss of bonded phase (eg C18 chains) or dissolution of the silica gel base material. Should this occur, columns are difficult to restore and reuse.

If 100% aqueous mobile phase is used in an ODS column, a sharp reduction in retention times of compounds can arise (as in the figure below). Whilst some may think this reduction in retention time is due to column deterioration, this is not the case. In this case, the cause is due to the

decrease of apparent hydrophobicity of the packing material due to polarity difference between the water and the C18 functional groups, leading the C18 chains to collapsing onto themselves. In some cases where this occurs, the initial retention times can be restored by flushing the column with 10 times its volume of mobile phase containing 50% organic solvent.

This decreases the repulsion between the eluent and the C18 chains and allows them to return to their normal pendant state. However YMC recommend that columns specifically intended for 100% aqueous eluents should be used to prevent this problem arising.



Basic Data

Conversion factors

Pressure

MPa	bar	psi	kgf/cm ²	atm
1	10	145.04	10.20	9.87
0.1	1	14.504	1.020	0.987
6.90×10^{-3}	0.069	1	0.070	0.068
0.0981	0.981	14.223	1	0.968
0.101	1.013	14.696	1.033	1

Length

m	in	ft	yd	mile
1	39.37	3.28	1.094	6.21×10^{-4}
0.025	1	0.083	0.028	0.15×10^{-4}
0.305	12	1	0.33	1.89×10^{-4}
0.91	36	3	1	5.68×10^{-4}
1609.3	63360	5280	1760	1

Weight

kg	oz	lb
1	35.274	2.204
0.0283	1	0.0625
0.454	16	1

Volume

l	gal(UK)	gal(US)
1	0.22	0.26
4.55	1	1.201
3.79	0.83	1

Temperature

K	°F	°C
0	-459.67	-273.15
255.37	0	-17.8
273.15	32	0
298.15	77	25
310.93	100	37.8
373.15	212	100

formula: $^{\circ}\text{C} = (^{\circ}\text{F} - 32) \times 5/9$ $^{\circ}\text{F} = ^{\circ}\text{C} \times 9/5 + 32$

Ratio Scale

ppb	ppm	%
1	10^{-3}	10^{-7}
10^3	1	10^{-4}
10^7	10^4	1

SI Prefixes

da (deca)	h (hecto)	k (kilo)	M (mega)	G (giga)	T (tera)	P (peta)	E (exa)	Z (zetta)	Y (yotta)
10^1	10^2	10^3	10^6	10^9	10^{12}	10^{15}	10^{18}	10^{21}	10^{24}

d (deci)	c (centi)	m (milli)	μ (micro)	n (nano)	p (pico)	f (femto)	a (atto)	z (zepto)	y (yocto)
10^{-1}	10^{-2}	10^{-3}	10^{-6}	10^{-9}	10^{-12}	10^{-15}	10^{-18}	10^{-21}	10^{-24}

1 Å (ångström) = 0.1 nm = 10^{-10} m


Column Area Ratio

Inner Diameter	1.0	2.0	3.0	4.6	10.0	20.0	30.0	50.0
Ratio	0.0473	0.189	0.425	1	4.73	18.90	42.53	118.15


Solvent Miscibility Table

[illegible]

Column Selection Guide



Column Selection Guide



Linear Scale-Up

In order to simplify your Scale-Up the three most important scale-up factors are summarised.

Scalable factor SF	ID "Linear Scale-Up"	Column length	Column length and ID "Volume"
	$SF = r_{ID, prep}^2 / r_{ID, anal.}^2$	$SF = l_{ID, prep} / l_{ID, anal.}$	$SF = (r_{ID, prep}^2 / r_{ID, anal.}^2) / (l_{ID, prep} / l_{ID, anal.})$
Impact	Flow rate Eluent composition	Retention time Cycle time Plate number	Amount of adsorbent

Linear Scale-Up

In most cases it is beneficial to develop a semi-preparative method on an analytical scale column.

The analytical separation carried out on a 150 x 4.6 mm ID column has to be scaled up to 150 x 20 mm ID. Therefore the chromatographic parameters such as flow rate and column load have to be adjusted according to the following equation:


Linear Scale-Up

$$SF = \frac{\Pi_{ID, prep}^2}{\Pi_{ID, anal.}^2} = \frac{m_{prep}}{m_{anal.}}$$

$$SF = r_{ID, prep}^2 / r_{ID, analytical}^2$$

$$SF = 20^2_{ID, prep} / 4.6^2_{ID, analytical}$$

$$SF = 18.9$$



Inner diameter → loadability / compound

Guideline for Sample Load according to column ID

Column ID (mm)	Scale-Up factor	Loadability (mg)
4.6	1	1-4
10	4.7	5-20
20	18.9	20-80
50	118	80-350
75	266	270-980
100	472	470-1900
150	1060	1000-4200

Preparative Column Selection Guide

Optimisation of preparative chromatography!

The main task for a preparative chromatographer is to find the suitable system. In order to simplify the considerations YMC developed a "Preparative Column Selection Guide".

		High ← Column efficiency ^{※1} Pressure ^{※1} Cost ^{※1} → Low					
Standard sample load	Particle size (μm) Inner diameter (mm ID)	spherical					Irregular
		5 N = 90000 ^{※2}	10 N = 40000 ^{※2}	15 N = 20000 ^{※2}	20 N = 10000 ^{※2}	50 N = 5000 ^{※2}	230/70 mesh N = 2500 ^{※2}
For investigation	4.6, 6.0						
tens-of-mg	10, 20						
hundreds-of-mg	50						
g	100, 150						
tens-of-g	200						
hundreds-of-g	300 or more						

Most appropriate
 Appropriate
 According to the purpose

※1 Value per unit length
 ※2 Standard theoretical plate number per m

The "Preparative Column Selection Guide" will help to select:-

1. the column ID for the required sample loading
2. the particle size for optimum efficiency
3. the column length for the necessary resolution

Scale-Up

The YMC Scale-Up is defined by 4 steps:

1. Analytical Scale: Method Development

Determine separation conditions by using analytical columns packed with different stationary phases and various conditions.

2. Study the preparative scale. Select the particle size of the packing material and the inner diameter of column appropriate for the sample volume.

3. Optimise the separation conditions and perform loadability studies using analytical columns with inner diameter of 4.6 mm or 6.0 mm packed with the packing material selected for the preparative separation (scout column). If the particle size of the packing material is the same as in the Step 1, this process can be omitted. If the preparative column is more than 100 mm ID, it is advisable to insert another step with a scout column of 20 mm ID in order to accurately predict loadability and calculate the running costs.

4. Proceed with the preparative separation with scale-up of chromatographic parameters such as flow rate/ column ID/ sample load as necessary.

From all the given steps above the most demanding step will be the scale-up of the chromatographic parameters in order to meet the preparative demands.

There are a number of scalable parameters: flow rate, column ID, sample load, tubing ID, sample injection concentration, volume of sample loop, consumption of solvent, dead volume, fraction mass, size of the detector cell.

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
A					
Acacetin	56	74	Angiotensin I	20,77,145,163	42,43,44,45
Acenaphthene	69,78,192,193	98,128	[Asn ¹ ,Val ⁸] Angiotensin I		42
Acenaphthylene	192,193	128	[Val ⁸] Angiotensin I	163,77,20	42
Acesulfam-K	19		Angiotensin II	20,77,145,153	42,43,45,47
Acetaldehyde 2,4-DNPH	103	127,128	Angiotensin III	20,145	42,43,46
Acetaminophen		90,91,92	Anilofos		121
Acetic acid	38,86,95	58	<i>o</i> -Anisic acid		130
Acetoacetyl-CoA		18	<i>m</i> -Anisic acid		130
Acetone 2,4-DNPH	103	127,128	<i>p</i> -Anisic acid		130
Acetylacetone	197	135	Annatto pigment		66
5-Acetylamino-6-amino-3-methyluracil		107	Anthracene	192,193	
N-Acetyl-chitobiose		28	Anthrathene		128
N-Acetyl-chitopentaose		28	α -Antitrypsin	151	53
N-Acetyl-chitotetraose		28	Apigenin	56	73,74
N-Acetyl-chitotriose		28	Apronamide		92
Acetyl-CoA		18	Arabitol		21
6''- <i>o</i> -Acetyldaidzin	30,57,58,59	60	Arachidonic acid		59
6''- <i>o</i> -Acetylgenistin	30,57,58,59	60	Arbutin	87	70,71,77,78
N-Acetylglucosamine		28	Arginine		41
6''- <i>o</i> -Acetylglycitin	30,57,58,59	60	L-Arginine	91	41,42
N-Acetylprocainamide	71,97,195	105	L-Ascorbic acid	30,33,105	13,15,16,17,19,90
Aciclovir	28		L-Ascorbic acid 2,4-DNPH		19
Acid Red		66	L-ascorbic acid 2-glycoside	30	
Acrolein 2,4-DNPH	103	127,128	Asparagine		41
Acylovir		105	L-Asparagine		42
Acylsterylglucosides	117		Aspartame	107	67
Adenine	72,79,87	34	Aspartic acid		41
Adenosine	32,87	34,35	L-Aspartic acid	91	41,42
Adenylate kinase	147,152,154	56,57	Astaxanthin	191	
5'-ADP	91	35	Asulam	17	121,122
ADP-D-glucose	89,101	36,37	Atenolol	21	
Adrenaline	86		Atorvastatin calcium hydrate	29	
Aflatoxin B1	198	78,79	5'-ATP	91	35
Aflatoxin B2	198	78,79	Atraton	197	122,123
Aflatoxin G1	198	78,79	Atrazine	197	122,123
Aflatoxin G2	198	78,79	Atrazine-desethyl		122,123
Alachlor		124	Atrazine-desisopropyl		122,123
Alanine		41	Atrazine-2-hydroxy		123
L-Alanine	91	41,42	Atropine sulfate	25,58	69,70
Albiflorin		71	Azamethiphos		121
Albumin	137,150,151	52,53	Azoxystrobin	17	
Aldicarb		124	AZT	103	104
Allantoic acid		108			
Allantoin		17,108	B		
Allopurinol		106,107	Bacampicillin hydrochloride		118
Aloin A		76	Baicalein	27,56,87	73,74
Aloin B		76	Baicalin		70,71,75
Amaranth		66	BAM-12P		45
Ametryn		122	Barbaloin	87	70,71,76
<i>p</i> -Aminobenzoic acid (PABA)	101,103,105	14,15,17	Barbital	25	
<i>o</i> -Aminophenol		136	Beclomethasone	89	80,82
<i>m</i> -Aminophenol		136	Beclomethasone dipropionate		82
<i>p</i> -Aminophenol		136	Bendiocarb		124
Aminopterin		103	Benfotiamine	197	16,17,91
Amitriptyline HCl	69,78	97,98	Bensulide		121,122
Amlodipine besilate	28,29		Benzaldehyde 2,4-DNPH	103	127,128
Amoxicillin		118,119	Benzalkonium chloride		94,95
5'-AMP	91	35,36	Benzene		80,128
Ampicillin		111,118	Benzethonium chloride		93,94,95
<i>n</i> -Amylbenzene	65		Benzo [a] anthracene	192,193	
Amyloid β -protein	160	47	Benzo [a] pyrene	192,193	128
Androsterone	199	84	Benzo [b] fluoranthene	192,193	128
<i>trans</i> -Androsterone	199	84	Benzo [e] pyrene		128

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
Benzo [g,h,i] perylene.....	192,193	128	CDP-D-glucose	89,101	36,37
Benzo [k] fluoranthene.....	192,193	128	Cefaclor.....	198	112,113
Benzocaine.....		102	Cefadroxil.....		113,114
Benzofenap.....		121	Cefamandole Na.....		115,116
Benzoic acid.....	71,77,93,107	62,63,67,105,129	Cefatrizine propylene glycol.....		113,114
Benzylpenicillin potassium.....		118	Cefazolin Na.....	198	113
Bepiridil hydrochloride.....		101	Cefoperazone Na.....		115,116
Betamethasone.....	89,199	80,82	Cefotaxime Na.....	198	113,116,117
Betamethasone valerate.....		82	Cefoxitin Na.....	198	113,117,118
Bezafibrate.....		99	Cefsulodin Na.....	198	113
1,1-Bi-2-naphthol.....	179,180	137	Ceftriaxone Na.....		116,117
BHA [3(2)-tert-Butyl-4-hydroxyanisol].....		65	Cefuroxime Na.....	198	113
BHT (2,6-Di-tert-butyl-4-methylphenol).....		65	Cephalexin.....	198	111,112,113,114,119
Biotin.....	101,103,105	14,15	Cephaloglycin.....	198	112,113,117
Bisbentiamine.....	197	17	Cephaloridine.....	198	112,113,114
Bisdemethoxycurcumin.....		62	Cephalosporin C potassium.....		113,114,117
Bombesin.....		44	Cephalothin Na.....		115,116
Bradykinin.....	77,145,163	44,45,46	Cephapirin Na.....	198	113,117
Brassicasterol.....	93	87	Cephradine.....	198	113
Bright Blue FCF.....		66	Cerebroside.....	117	
Bromovalerylurea.....		92,93	Chenodeoxycholic acid (CDAC).....		86
BSA.....	95,102,136,137,150,161,163	48,49,51	Chloramphenicol.....		111
Bucetine.....		91	Chlorhexidine hydrochloride.....		95
Buckminsterfullerene (C60).....	101		Chloridazon.....		121
Bufalin.....	93	89	o-Chloroacetophenone.....		134
Butacaine.....		102	m-Chloroacetophenone.....		134
Butachlor.....		124	p-Chloroacetophenone.....		134
Butafenacil.....		121	o-Chlorobenzoic acid.....		130,131
p-Buthoxybenzoic acid.....		129,130	m-Chlorobenzoic acid.....		130,131
n-Butylbenzene.....	65,78		p-Chlorobenzoic acid.....		92,130,131
Butyl benzoate.....	54		Chlorogenic acid.....		61,62
p-tert-Butylbenzoic acid.....		129,130	o-Chlorophenol.....	77	132,133
tert-Butylhydroquinone (TBHQ).....		64,65	m-Chlorophenol.....	77	132,133
Butyl p-hydroxybenzoate.....		63,78,96	p-Chlorophenol.....	77	132,133
C			Chlorophyll a.....	191	79
trans-Caffeic acid.....		61,62	Chlorophyll b.....	191	79
Caffeine.....	107	13,14,67,70,90,91,92, 93,95,107,117,118,119	Chlorophyll c ₁ +c ₂	191	
Campesterol.....	93	87	Chloropropham.....		124
Candesartan cilexetil.....	29		Chlorotetracycline.....	21,23	111
Canthaxanthin.....	191		Chlorpheniramine.....	16,18,35,37,38,79,97	
Capronaldehyde 2,4-DNPH.....	103	127,128	Chlorpheniramine (maleate).....		90,91,93,94,95,96
Capsaicin.....	47	79	Chlorpromazine hydrochloride.....		98
Capsanthin.....	191		Chlortetracycline.....	23	
Carbadox.....	89,199	109	Chlorthalidone.....		102,103,139
Carbamazepine.....	69	99	Cholecalciferol.....	93	8,9,10
Carbendazim.....		121	Cholesterol.....	93	87
Carbonic anhydrase.....	145,151	53	Cholic acid (CA).....		86
Carnitine chloride.....		106	Chondroitin sulfate A, Na salt.....		31
trans-α-Carotene.....	191	8	Chondroitin sulfate B, Na salt.....		31
trans-β-Carotene.....	191	8	Chondroitin sulfate C, Na salt.....		32
δ-Carotene.....	191		Chondroitin, Na salt.....		31
9-cis β-carotene.....	191		Chromafenozide.....		121
13-cis β-carotene.....	191		Chrysene.....	192,193	128
15-cis β-carotene.....	191		Chrysin.....	56	73,74
Carvedilol.....	29		Chrysoeriol.....		74
Casein.....		57	α-Chymotrypsinogen A.....	161	51
(+)-Catechin.....	198	75,76	Cinchonidine.....		69
Catechol.....	101,121	131	Cinchonine.....	25,58	69,70,102
CBZ-Phenylalanine (Z-Phe-OH) / CBZ-Phenylalanine.....	180	138,140,142	trans-Cinnamic acid.....	27	61,62
CKK-Octapeptide.....	145	45	Cinobufagin.....	93	89
5'-CDP.....	91	35	Cinoxacin.....		111,112
			Citric acid.....	86	58
			Clemastine.....	39	

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
Clemastine fumarate.....	46		D		
Clindamycin	28		D4T	103	104
Clofibrate.....		99	Daidzein	27,30,57,58,59	60
Cloquintocet-mexyl		121	Daidzin	30,57,58,59	60
Clothianidin		121	5'-dAMP	91	35
Cloxacillin sodium.....		118	5'-dCMP	91	35
5'-CMP	91	35,36	DDC.....	103	104
Codeine.....		93	DDI	103	104
Coenzyme A.....		18	Dehydroacetic acid.....	77	62,63
Coenzyme A glutathione disulfide.....		18	Dehydrocerebroside.....	117	
Coenzyme A oxidized.....		18	Dehydroisoandrosterone.....	199	84
Color additives (Acid red).....		66	Delphinidin	31	
Color additives (Amaranth).....		66	Delphinidin-3-O-arabinoside.....	31	
Color additives (Brilliant Blue FCF).....		66	Delphinidin-3-O-galactoside.....	31	
Color additives (Erythrosine).....		66	Delphinidin-3-O-glucoside.....	31	
Color additives (Fast Green FCF).....		66	Demethoxycurcumin.....		62
Color additives (Indigo Carmine).....		66	Deoxyadenosine	87	34,35
Color additives (New Coccine).....		66	Deoxycorticosterone	107,119	80,81,82,83
Color additives (Phloxine).....		66	Deoxycorticosterone acetate.....		82
Color additives (Rose Bengale).....		66	Deoxycytidine	87	34,35
Color additives (Sunset Yellow FCF).....		66	Deoxyguanosine	87	34,35
Color additives (Tartrazine)		66	Deoxyinosine	87	34,35
Conalbumin.....	102,161	51	Deoxyuridine		103
Copper 8-quinolinolate.....		121	Desipramine hydrochloride		97,98
Coronene		128	Dexamethasone.....	89,199	80,82
Corticosterone.....	89,107,119,199	80,82,83	Dextromethorphan HBr	16,18,35,37,38	90
Corticosterone acetate.....		80,82	5'-dGMP	91	35
Cortisone	89,199	80,82	Δ Di-OS.....		21
Cortisone acetate		80,82	Δ Di-4S.....		21
p-Coumaric acid.....		61,62	Δ Di-6S.....		21
Creatinine.....	151	53,127	Δ Di-diSD.....		21
o-Cresol	184	133	Δ Di-diSE.....		21
m-Cresol	184	133	Δ Di-triS.....		21
p-Cresol	184	133	Diacylglycerol.....	117	
Crimidine.....		122	Diadinoxanthin	191	
Crotonaldehyde 2,4-DNPH.....	103	127,128	Diatoxanthin	191	
β-Cryptoxanthin	191		Diaveridine		109
5'-CTP	91	35	Dibenzo [a,h] anthracene.....	192,193	128
Cuminic acid		129,13	Dicethamine hydrochloride.....	197	17
Curcumin 1 (Curcumin).....		62	2,4-Dichlorophenol		132
Curcumin 2 (Demethoxycurcumin).....		62	2,6-Dichlorophenol		132
Curcumin 3 (Bisdemethoxycurcumin).....		62	Dicloxacillin sodium		118
Cyanazine	197	122,123	Diethanolamine.....	95	135
Cyanidin	31		Diethylaniline.....	97	
Cyanidin-3-O-arabinoside.....	31		1,1-Diethylurea.....		136
Cyanidin-3-O-galactoside.....	31		1,3-Diethylurea.....		136
Cyanidin-3-O-glucoside.....	31		Digalactosyldiacylglycerol.....	117	
Cyanocobalamin.....	30,101,103,105	13,14,15,16	Digitoxigenin		88
Cyazofamid		121	Digitoxin	87	70,71,88
α-Cyclodextrin	89,153	27,28	Digoxigenin		88
β-Cyclodextrin	89,153	27,28	Digoxigenin-bis-digitoxoside		88
γ-Cyclodextrin.....	89,153	27,28	Digoxin		88
Cyclamat-Na ⁺	19		Dihydrocapsaicin	47	79
Cyflufenamid.....		121	Dihydrocodeine phosphate		90,91
L-Cystine.....	91	41,42	Dihydroquinidine		69
Cytidine.....	32,87	34,35	Dihydroquinine.....	23,58	69
Cytochrome c	103,135,139,145,147,152,154,160,161,162,163	47,48,50,51,56,57	2,3-Dihydroxybenzoic acid.....		129
Cytochrome c (dimer).....		50	2,4-Dihydroxybenzoic acid.....		129
Cytochrome c (hexamer).....		50	2,5-Dihydroxybenzoic acid.....		129
Cytochrome c (trimer).....		50	3,4-Dihydroxybenzoic acid.....		129
Cytochrome c (trimer).....		50	3,5-Dihydroxybenzoic acid.....		129
Cytosine.....	31,32,72,79,87	34	3,4-Dihydroxyphenylacetic acid.....	79	106
			3,4-Dihydroxyphenylalanine (DOPA).....	79,86	106

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
DL-3,4-Dihydroxymandelic acid.....	86		β-Estradiol.....	99,185	85,86
Diltiazem hydrochloride.....	72,73,77	100,101,102	Estradiol benzoate.....		85,86
Dimethirimol.....		121	Estradiol dipropionate.....		85,86
7,12-Dimethylbenzo [a] anthracene.....		128	Estriol.....	99	85
α-Dimethyl-β-cyclodextrin.....	89	28	Estrone.....	99	85
Dimethylaniline.....	97		Ethanolamine.....	95	135
N,N-Dimethylethanolamine.....	95	135	Ethenzamide.....	46	76,90,91,92,93
2,3-Dimethylphenol.....		132	Ethinylestradiol.....	99	85,86
2,4-Dimethylphenol.....		132	Ethiofencarb.....		124
2,5-Dimethylphenol.....		132	Ethopropazine hydrochloride.....		98
2,6-Dimethylphenol.....		132	Ethosuximide.....	69	99
3,4-Dimethylphenol.....		132	p-Ethoxybenzoic acid.....		129,130
3,5-Dimethylphenol.....		132	Ethylbenzene.....		82
1,7-Dimethyluric acid.....		107	Ethyleneglycol.....		126
1,9-Dimethyluric acid.....		107	Ethyl p-hydroxybenzoate.....		63
2,4-Dinitrophenol.....		132	Ethylurea.....		136
2,6-Dinitrophenol.....		132	Etizolam.....	22	
Diphenhydramine hydrochloride.....		96			
Diphenyl.....		86,124,125,128	F		
Diquat.....	21		Fast Green FCF.....		66
Diuron.....	197	122,123	Fenofibrate.....		99
1kb DNA ladder (75-12,126 bp).....	140	39	Fenoprofen.....		97,139
Docosahexaenoic acid.....		59	Fenoxycarb.....		121
n-Dodecyl gallate (DG).....		64,65	(E)-Ferimzone.....		121
DOPA (3,4-Dihydroxyphenylalanine).....	86	106	(Z)-Ferimzone.....		121
Dopamine hydrochloride.....	79,86	106	Ferritin.....	163	
Doxepin hydrochloride.....		97	trans-Ferulic acid.....		61,62
Doxorubicin hydrochloride.....		111	Fibrinogen.....	145,151	53
Doxycycline.....	23		Flavin adenine dinucleotide disodium.....		18
Dynorphin A (Human,1-13).....		44	Flavin mononucleotide sodium salt (Flavin sodium phosphate).....		13,14,18
E			Flazasulfuron.....	17	
α-Ecdysone.....		87	Fluocinolone acetonide.....		80,81,82
β-Ecdysone.....		87	Fluocinonide.....		81,82
Cu (II)-EDTA.....		63	Fluoranthene.....	192,193	128
Fe (III)-EDTA.....		63	Fluorobiprofen.....	183	
Echinone.....	191		Fluorene.....	192,193	128
Eicosapentaenoic acid.....		59	5-Fluoro-2'-deoxyuridine.....		103
Eicosatrienoic acid.....		59	5-Fluorocytosine.....		34
Eledoisin related peptide.....		44	5-Fluoroorotic acid.....		35
α-Endorphin.....		45	5-Fluorouracil.....		34,103
β-Endorphin (Human).....	31	44,45	Fluoxymesterone.....	115,199	83,84
γ-Endorphin.....	31	44,45	Fluphenazine hydrochloride.....		98
[D-Ala ² , Met ⁵]-Enkephalin.....		45	Flurbiprofen.....	77	97,100,139,141
Leu-Enkephalin.....	31,77,95,145,163	44,45,46,48	Folic acid.....		13,103
Met-Enkephalin.....	31,95,145,163	44,45,46,48	Formaldehyde 2,4-DNPH.....	103	127,128
[D-Ala ² , Met ⁵]-Enkephalinamide.....		45	Formic acid.....	38	
Enolase.....	147,152,154	56,57	1-Fructofuranosyl-D-nystose.....		23
Enoxacin.....		111	Fructose.....		20,22
(-)-Epicatechin.....	198	75,76	Fucose.....		20
(-)-Epicatechin gallate.....	198	75,76	Fucosterol.....	93	87
(-)-Epigallocatechin.....	198	75,76	Fucoxanthin.....	191	
(-)-Epigallocatechin gallate.....	198	75,76	(5,6)-Fullerene.....	101	
Epinephrine.....	79	106	Fullerene (C ₆₀).....		134
Ergocalciferol.....	93	8,9,10	Fullerene (C ₇₀).....		134
Erythorbic acid.....	30,33	16,19	Fumaric acid.....	86,95	58,91
Erythritol.....		21,22	Furathiocarb.....		121
Erythromycin.....	20,26,76	119	Furazolidone.....	89,199	109
Erythromycin estolate.....	26	119	Furosemide.....		102,103
Erythromycin ethylsuccinate.....	26	119			
Erythrosine.....		66	G		
Esculine.....	87	70,71	Galactose.....		20
α-Estradiol.....	185		Gallic acid.....		61,62,77,78

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
5'-GDP.....	91	35	Ibuprofen.....	46	90,93,97,140
GDP-D-mannose.....	89,101	36,37	IgA.....	145,151	
Gemfibrozil.....		99	IgG.....	137,145,150,151	
Genistein.....	30,57,58,59	60	IgM.....	145,151	
Genistin.....	30,57,58,59	60	IgY.....	138	
Ginsenoside-Rb ₁		72,73	Imazalil.....		124,125
Ginsenoside-Rg ₁		72,73	Imipramine hydrochloride.....		97,98
γ-Globulin.....	145,151	53	5'-IMP.....	91	35,36
Glucose (G1).....	89,146	20,22,24,25,26	Indapamide.....		140
Glutamate dehydrogenase.....	147,152,154	56,57	Indeno(1,2,3-c,d)pyrene.....	192,193	
Glutamic acid.....		41	Indigo Carmine.....		66
L-Glutamic acid.....	91	41,42	Indometacin (Indomethacin).....	93	89,92
Glutamine.....		41	Indoxacarb.....		121
Glutathione (GSH).....		42	Inosine.....	87	34,35
Glutathione disulfide (GSSG).....		42	Inositol.....		21
Glycerol.....		21	Insulin.....	95	44,48
Glycine.....	91,136,145,153	41,42,47	Insulin (Bovine).....	20,136,145,153,160,163	43,44,47
Glycitein.....	30,57,58,59	60	Insulin (Equine).....		43
Glycitin.....	30,57,58,59	60	Insulin (Human).....	20	43,44
Glycolic acid.....	86	58	Insulin (Porcine).....	20	43,44
Glycyrrhizic acid / Glycyrrhizin.....	46,79,87	68,70,71,95	Insulin (Sheep).....		43,44
5'-GMP.....	91	35,36	Insulin B Chain.....	145	45
5'-GTP.....	91	35	Iprodion.....		121
Guanine.....	72,79,87	34	Iprovalicarb.....		121
Guanosine.....	32,87	34,35	Isobarbaloin.....		76
H			Isobutyl p-hydroxybenzoate.....	28	63
Halosulfuron methyl.....	17		Isobutyraldehyde 2,4-DNPH.....	103	127,128
Heparan sulfate, Na salt.....		32	Isoleucine.....		41
Heparin, Na salt.....		32	L-Isoleucine.....	91	41,42
Hesperidin.....		75	Isomalt.....		21
Hexazinone.....		122	Isonicotinic acid.....	99,198	104,105
Hexobarbital.....	25,179,180,183	138,139	Isonicotinic acid hydrazide.....	99,198	104,105
Hinokitiol.....	35,38	135	Isophthalic acid.....		131
Hippuric acid.....		127	p-Isopropylbenzoic acid.....		129,130
Histidine.....		41	Isopropyl p-hydroxybenzoate.....		63
L-Histidine.....	91	41,42,106	Isoproturon.....	197	122,123
Homovanillic acid.....	79,86	106	Isoquercitrin.....		74
Hordenine sulfate.....		69,70	Isoxaflutole.....		121
Human IgA.....	151	53,54,55	K		
Human IgG.....	151,152	52,53,54,55	Kaempferol.....	27,56	73
Human IgM.....	151,152	53	Keratan polysulfate, Na salt.....		33
Human IgM (κ), myeloma protein.....	152	54,55	Keratan sulfate, Na salt.....		33
Human IgM (λ), myeloma protein.....		55	1-Kestose.....		22,23
Human serum albumin (HSA).....	145,151	53	Ketamine HCl.....	181	138
Hyaluronic acid, Na salt.....		32	Ketoprofen.....	183	96,97,140
Hydrochlorothiazide.....	29	102,103	L		
Hydrocortisone.....	89,107,119,199	80,82,83	α-Lactalbumin.....	145,160	47
Hydrocortisone acetate.....		80,82	Lactate dehydrogenase.....	147,152,154	56,57
Hydrocortisone succinate.....		82	Lactic acid.....	86,95,177	58
Hydroquinone.....	101,121	77,78,131	Lactitol.....		21,22
4'-Hydroxyacetophenone (p').....	97,195	71,105,120	β-Lactoglobulin A.....	161	51
m-Hydroxybenzoic acid.....		129	β-Lactoglobulin B.....	161,163	48,51
p-Hydroxybenzoic acid.....		61,62,129	Lactose.....		20
p-Hydroxybenzoic acid n-butyl ester.....		92	Lactosyl fructoside.....		24
5-Hydroxyindol-3-acetic acid.....	86		Lanatoside B.....		88
4-Hydroxy-3-methoxymandelic acid.....	79	106	Lanatoside C.....		88
4-Hydroxy-3-methoxyphenylglycol.....	79,86	106	Leucine.....		41
4-Hydroxymethyl-2,6-di-tert-butylphenol (HMBP).....		64,65	L-Leucine.....	91	41,42
3-Hydroxy-3-methylglutaryl-CoA (HMG-CoA).....		18	Lidocaine.....	72,73,77	93,100,101,102
Hyodeoxycholic acid (HDCA).....		86	Linoleic acid.....		59
Hypoxanthine.....	87	34,106,107			

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
γ -Linolenic acid.....		59	Methylethylketone 2,4-DNPH.....	103	127,128
Linomycin HCl.....	28		<i>o</i> -Methylhippuric acid.....		127
Linuron.....	197	122,123	<i>m</i> -Methylhippuric acid.....		127
Lipoxidase.....	102,161	51	<i>p</i> -Methylhippuric acid.....		127
Lovastatin.....		100	Methyl <i>p</i> -hydroxybenzoate.....	27	63,74,93
Lutein.....	191		Methyl linoleate.....		59
Luteolin.....	27	74	Methyl linolenate.....		59
Lycopene.....	191		Methyl oleate.....		59
Lysine.....		41	Methyl γ -linolenate.....		59
L-Lysine.....	91	41,42	Methylmalonyl-CoA.....		18
Lysophosphatidylcholine.....	117		<i>N</i> -Methylnicotinamide.....	99,198	104,105
Lysozyme.....	103,135,139,151,160,161,163	47,50,51,53	Methylprednisolone.....	89,115,199	80,82,83
M			Methyltestosterone.....	199	84
Maleate.....	97		1-Methyluric acid.....		107
Maleic acid.....	79,86,95	58,90,93,95	1-Methylxanthine.....		107
L-Malic acid.....	86,95	58	Metobromuron.....	197	122,123
Malonic acid.....	86	58	Metolachlor.....		122,124
Malonyl-CoA.....		18	Metoxuron.....	197	122,123
6''- <i>o</i> -Malonyldaidzin.....	30,57,58,59	60	Metoprolol.....	21	
6''- <i>o</i> -Malonylgénistin.....	30,57,58,59	60	Monoclonal antibodies.....	139,141	
6''- <i>o</i> -Malonylglycitin.....	30,57,58,59	60	Monogalactosyldiacylglycerols.....	117	
Maltitol.....		21,22	Mouse IgG Fab fragment.....	152	54
Maltododecaose (G ₁₂).....		25,26	Mouse IgG Fc fragment.....	152	53,54
Maltoheptaose (G ₇).....	146	24,25,26	Mouse monoclonal IgG1.....	152	53
Maltohexaose (G ₆).....		24,25	Myoglobin.....	103,136,145,160,161	47,50
Maltopentadecaose (G ₁₅).....	146	25,26	Myricetin.....	56	73
Maltopentaose (G ₅).....	146	24,25,26	N		
Maltose (G ₂).....	146	20,22,24,25	NAD.....		18
Maltotetraecaose (G ₁₄).....		25,26	NADH.....		18
Maltotetraose (G ₄).....		24,25	Nadolol.....	21	
Maltotridecaose (G ₁₃).....		25,26	NADP.....		18
Maltotriose (G ₃).....	146	24,25,26	NADPH.....		18
Maltoundecaose (G ₁₁).....	146	25,26	Nalidixic acid.....	89,199	109,111,112
Malvidin.....	31		Naphazoline hydrochloride.....	79	93,94,95
Malvidin-3- <i>O</i> -arabinoside.....	31		Naphthacene.....		128
Malvidin-3- <i>O</i> -galactoside.....	31		Naphthalene.....	54,76,192,193	128
Malvidin-3- <i>O</i> -glucoside.....	31		1-Naphthol.....		135
Mandelic acid.....	86	127	1-(1-Naphthyl)-ethylalcohol.....	184	141
Mannitol.....		21	Naproanilide.....		121
Mannose.....	89	20	Naproxen.....		96
α -Mating factor.....	95,145,163	44,48	Narcotine.....		69,70,91
Mecoprop.....	17	121	Naringin.....		75
Medroxyprogesterone acetate.....	99	85	Narirutin.....		75
Menadione.....		8	Neohesperidin.....		75
Menaquinone.....		8,12	Neostigmine methylsulfate.....	79	95
Mestranol.....	99	85	Neoxanthin.....	191	
Metazachlor.....		122,124	Neurotensin.....	31,136,145,153	44,47
Metformin HCl.....	32		New Coccine.....		66
Methicillin sodium.....		118	Nicardipine hydrochloride.....	72,73,77	101,102
Methionine.....		41	Nicotinamide.....	30,99,105,198	13,14,15,16,17,18,104,105
L-Methionine.....	91	41,42	Nicotinamide adenine dinucleotide phosphate oxidized form.....		18
Methotrexate.....		103	Nicotinamide adenine dinucleotide phosphate reduced form.....		18
DL-3-Methoxy-4-Hydroxymandelic acid.....	86		β -Nicotinamide adenine dinucleotide reduced form.....		18
<i>o</i> -Methoxybenzoic acid.....		130	β -Nicotinamide adenine dinucleotide.....		18
<i>m</i> -Methoxybenzoic acid.....		130	Nicotinic acid.....	28,30,99,198	13,16,17,18,104,105
<i>p</i> -Methoxybenzoic acid.....		130	Nifedipine hydrochloride.....		101
Methoxyfenozide.....		121	Niguldipine hydrochloride.....		101
3-Methoxytyramine hydrochloride.....	79	106	<i>o</i> -Nitrophenol.....	185	132,133
Methoxyverapamil hydrochloride.....		101	<i>m</i> -Nitrophenol.....	185	132,133
Methyl benzoate.....	35,38,54,77	95,100	<i>p</i> -Nitrophenol.....	185	132,133
dl-Methylephedrine hydrochloride.....		90,91	Noradrenaline.....	86	
<i>N</i> -Methylethanolamine.....	95	135	Nordihydrocapsaicin.....	47	79

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
Nordihydroguaiaretic acid (NDGA).....		64,65	Petunidin-3-O-glucoside.....	31	
Norepinephrine hydrochloride.....	79	106	Phaeophorbid b.....	191	
Norethisterone.....	99	85	Phenanthrene.....	192,193	128
Norfloxacin.....		111	Phenethicillin potassium.....		118
Norgestrel.....		86	Phenformin HCl.....	32	
Noscapine.....		91	Phenobarbital.....	25	
Nystose.....		22,23	Phenobarbital sodium salt.....	69	99
			Phenol.....	71,93,97,101,121	70,105,117,131
O			Phenoxyethylpenicillin potassium.....		118
<i>n</i> -Octyl gallate (OG).....		64,65	Phenoxyethylpenicillin potassium.....		118
Ofloxacin.....		111	Phenylalanine.....	86	41
Olaquinox.....	89,199	109	L-Phenylalanine.....	91	41,42
Oligonucleotides.....		37,38,39	<i>o</i> -Phenylenediamine.....		136
Orotic acid.....	105	15,35	<i>m</i> -Phenylenediamine.....		136
Ouabagenin.....		88	<i>p</i> -Phenylenediamine.....		136
Ovalbumin.....	95,136,145,161,163	48,51,52	Phenylephrine hydrochloride.....		95
Oxamyl.....		124	1-Phenylethanol.....	179	
Oxine-copper.....	17	121	Phenylethylalcohol.....	184	141
Oxipurinol.....		106,107	2-Phenylethyl β-cellobioside.....		76
Oxolinic acid.....	89,199	109,111,112	2-Phenylethyl β-gentiobioside.....		76
Oxytetracycline hydrochloride.....	21,23	111,119	2-Phenylethyl β-lactoside.....		76
Oxytocin.....	31,77,145,163	44,45	2-Phenylethyl β-maltoside.....		76
			2-Phenylethyl β-melibioside.....		76
P			Phenylglyoxylic acid.....		127
PA-diheptadecanoyl.....	117		<i>o</i> -Phenylphenol.....		124,125
Paeoniflorin.....		71	Phenylpropanolamine hydrochloride.....	97	95
Paeonol.....		78	Phenytol.....	69,72,73	99,101,102
Palatinol.....		21	Phloroglucinol.....	121	
Palatinose.....		21,22	Phloxine.....		66
Pantothenic acid, Na salt.....		13	Phosphatidylcholine.....	117	
Pantothenyl alcohol (Panthenol).....		13,17,18	Phosphatidylethanolamine.....	117	
Papaverine chloride.....		69,7	Phosphatidyl glycerols.....	117	
Paraquat.....	21		Phosphatidylinositol.....	117	
Paraxanthine.....		107	Phosphoenolpyruvic acid monopotassium salt.....		58
PA-Sugar Chain 001.....		28,29	Phthalic acid.....		131
PA-Sugar Chain 002.....		28,29	Phylloquinone.....		8,12
PA-Sugar Chain 014.....		28,29	PI, PC, PS-diheptyldecanoyl.....	117	
PCB 28.....	193		Pilocarpine chloride.....		69,7
PCB 52.....	193		Pindolol.....	21	
PCB 101.....	193		Pipemidic acid.....		111
PCB 138.....	193		Pirimicarb.....		124
PCB 153.....	193		Piromidic acid.....	89,199	109,111,112
PCB 180.....	193		Piroxicam.....		96
PE-dipalmitoyl.....	117		Plasmid pBR322.....	140,149	39,4
PEG #1,000 (MW 1,000).....		126	Poly (A).....		37
PEG #20,000 (MW 20,000).....		126	Polydextrose II.....		30
PEG #300 (MW 300).....		126	Potassium guaiacol sulfonate.....		91
PEG #4,000 (MW 3,350).....		126	Prednisolone.....	115	80,81,82,83
PEG #6,000 (MW 8,400).....		126	Prednisolone acetate.....		81,82
Penicillin G potassium.....		118	Prednisolone succinate.....		82
Penicillin V potassium.....		118	Prednisone.....		80,82,83
Pentaglycine.....	153	42	Prenylamine lactate.....		101
Pentobarbital.....	25		Primidone.....	69	99
Peonidin.....	31		Procainamide.....	97,195	
Peonidin-3-O-arabinoside.....	31		Procainamide hydrochloride.....	71	105
Peonidin-3-O-galactoside.....	31		Procaine hydrochloride.....		102
Peonidin-3-O-glucoside.....	31		Prochlorperazine maleate.....		98
Peptide mapping.....		48,49	Progesterone.....	69,99,107,119	83,85
Perphenazine.....		98	Proline.....		41
Perylene.....	193	128	L-Proline.....	91	41,42
Petunidin.....	31		Promazine hydrochloride.....		98
Petunidin-3-O-arabinoside.....	31		Prometon.....	197	122,123
Petunidin-3-O-galactoside.....	31		Prometryn.....	197	122,123

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
Propazine	197	122,123	Saccharin Na	19	62,90
Propham		124	Saikosaponin a	199	72
Propicillin potassium		118	Saikosaponin b2	199	72
Propionaldehyde 2,4-DNPH	103	127,128	Saikosaponin c	199	72
Propionamide		136	Saikosaponin d	199	72
Propionic acid	38,86,95	58	Salicin	87	70,71
Propionyl-CoA		18	Salicylic acid	27,71,93	61,62,105,129
p-Propoxybenzoic acid		129,13	Saponin		72
Propranolol HCl	21,36,72,73,180,185	101,102,138,141	Sccinyl-CoA		18
n-Propyl gallate (PG)		64,65	Scopolamine	25,58,97	
Propyl p-hydroxybenzoate	35,37,38	63,68	Scopolamine Br		69,70
Propylbenzene	69,78		Sebuthylazine	197	122,123
Propylparabene	16,18		Secobarbital	25	
Proscillaridin A		88	Sennoside A		76,77
Protocatechuic acid		61,62	Sennoside B		76,77
Pseudoephedrine	97		Serine		41
Pullulan P-5 (MW 5,800)	146	30,31	L-Serine	91	41,42
Pullulan P-10 (MW 12,200)	146		Siduron	17	121
Pullulan P-20 (MW 23,700)	146	30,31	Simazine	197	122,123
Pullulan P-50 (MW 48,000)	146	30,31	Simeconazole		121
Pullulan P-100 (MW 100,000)	146		Sinapic acid		61,62
Pullulan P-200 (MW 186,000)	146		β-Sitosterol	93	87
Pullulan P-400 (MW 380,000)	146		Sodium decylbenzenesulfonate		125,126
Pullulan P-800 (MW 853,000)	146	30,31	Sodium dodecylbenzenesulfonate		125,126
Pyrazolynate		121	Sodium tetradecylbenzenesulfonate		125,126
Pyrene	192,193	128	Sodium tridecylbenzenesulfonate		125,126
Pyridoxal hydrochloride		16,17	Sodium undecylbenzenesulfonate		125,126
Pyridoxine hydrochloride	30,79,105	13,14,15,16,17,95	Sorbic acid	77	62,63
Pyrititalid		121	Sorbitol		21,22
Co (II)-1-(2-Pyridylazo)-2-naphthol		135	Sparteine sulfate		69,7
Fe (III)-1-(2-Pyridylazo)-2-naphthol		135	Sphingomyelin	117	
Ni (II)-1-(2-Pyridylazo)-2-naphthol		135	Spiarmycin	20	
Zn (II)-1-(2-Pyridylazo)-2-naphthol		135	Stachyose		23
Pyrimethamine		109	Sterol Esters	117	
Pyrogallol	101,121	131	Sterols	117	
Pyroloquinoline quinone (PQQ)	105	15	Steryl glycosides	117	
Pyruvic acid		58	Steryl oleate	117	
			Stevioside		68
Q			Stigmasterol	93	87
Quercetin	56	73,74	Stigmasterol / sitosterol	117	
Quercitrin		74	G-Strophanthin	87	70,71,88
Quinidine	72,73	69,101,102	Substance P		44
Quinine	25,58	69,70	Succinic acid	86,95	58
8-Quinolol	197	135	Sucralose		22,23
			Sucrose		20,22,23
R			Sulfachloropyridazine		110
Rabbit IgG	152	54	Sulfadiazine		109,110,111
Raffinose		23	Sulfadimethoxine	89,199	109,11
Rebaudioside A		68	Sulfadimidine	89,199	109,110,111
Resibufogenin	93	89	Sulfadoxine		110
Resorcinol	99,121	18,117,129	Sulfaguandine		109,11
α-Resorcylic acid	99	129	Sulfamerazine	89,199	109,110,111
β-Resorcylic acid	99	129	Sulfamethazine		109,11
γ-Resorcylic acid	99	129	Sulfamethizole		109,11
Retinol		8	Sulfamethoxazole		110,111
Retinol acetate		8	Sulfamethoxypyridazine		109,11
Retinol palmitate	115	8	Sulfanomomethoxine	89,199	109,110,111
Rhamnose		20	Sulfanilamide		109,11
Ribitol		21	Sulfapyridine		109,11
Riboflavin	30,101,103,105	13,14,15,16,17,18	Sulfaquinoxaline	89,199	109,11
Ribonuclease A	103,135,139,145,161,163	50,51	Sulfathiazole		109,110,111
Rose Bengale		66	Sulfisomidine		109,11
Rutin		74	Sulfisoxazole		110

Substance Index

	Refer to this catalogue on page	Please refer to Application Data Collection on page		Refer to this catalogue on page	Please refer to Application Data Collection on page
Sulfoquinovosyldiacylglycerol	117		2,2,2-Trifluoro-1-(9-anthryl) ethanol	179,181	137
Sulindac		96	Triglycine	153	42
Sulphamerazin	20		Trigonerline		70
Sulphamethaxole	20		2,4,5-Trihydroxybutyrophenone (THBP)		64,65
Sulphathiazole	20		Trimethoprim		109,110
Sunset Yellow FCF		66	o-Trimethyl- β -cyclodextrin	89	28
Swep		124	2,4,6-Trinitrophenol		132
T			Triphenylene	65	128
TAG / Tripentadecanoin	117		Triton X-100		126
Tartaric acid		58	Trogar's base		142
Tartrazine		66	Trypsin inhibitor (soybean)	136,145	52
Taurine		41,42	Tryptophan		41
5'-TDP	91	35	L-Tryptophan		41
Terbumeton	197	122,123	5'-TTP	91	35
Terbutylazine	197	122,123	Tyrosine	86	41
Terbutryn	197	122,123	L-Tyrosine	91	41,42
Terephthalic acid		131	U		
o-Terphenyl	65	134	UDP-D-glucose	89,101	36,37
m-Terphenyl		134	UDP-N-acetyl-D-glucosamine	89,101	36,37
p-Terphenyl		134	5'-UMP	91	35,36
Testosterone	78,199	84,85	Uracil	20,32,54,65,71,72,79,87,97,195,197,231	34,35,103,105,135
Testosterone enanthate		84,85	Urea		17,136
Testosterone propionate		84,85	Uric acid	151	53,106,107,108
Tetracaine hydrochloride		102	Uridine	32,87	34,35
Tetracycline hydrochloride	21,23,136,145	111,119	5'-UTP	91	35
Tetraglyzine	136,145		V		
Tetrahydrozoline		93	n-Valeraldehyde 2,4-DNPH	103	127,128
Theobromine		70,107	Valine		41
Theophylline	71,93	70,105,107	L-Valine	91	41,42
Thiabendazole		121,124,125	Valsartan	29	
Thiacloprid		121	Vancomycin hydrochloride		120
Thiamethoxam		121	Vanillic acid	86	61,62
Thiamine hydrochloride	30	13,14,17,90	[Arg ⁸]-Vasopressin	77	45
Thiamphenicol	89,199	109	Verapamil hydrochloride	36,72,73,77	100,101,102
Thiophanate-methyl		121	Violaxanthin	191	
Thioridazine hydrochloride		98	Vitamin A		8
Thiram	17	121,122	Vitamin A palmitate		8
Threonine		41	Vitamin B ₁		13,14,17,90
L-Threonine	91	41,42	Vitamin B ₁₂		13,14,15,16
Δ -Threonine		41	Vitamin B ₁₃		15,35
Thymidine	87	34,35	Vitamin B ₂ (Flavin mononucleotide)		13,14,18
Thymine	72,79,87	34	Vitamin B ₂ (Riboflavin)		13,14,15,16,17,18
Thyroglobulin	145		Vitamin B ₆ (Pyridoxal hydrochloride)		16,17
5'-TMP	91	35,36	Vitamin B ₆ (Pyridoxinel hydrochloride)		13,14,15,16,17,95
α -Tocopherol	115,123	8,10,11,12	Vitamin C		13,15,16,17,19,90
β -Tocopherol	115,123	10,11,12	Vitamin D ₂		8,9,10
γ -Tocopherol	115,123	10,11,12	Vitamin D ₃		8,9,10
δ -Tocopherol	115,123	8,10,11,12	Vitamin E		8
α -Tocopherol acetate		8,12	Vitamin H		14,15
Tolmetin sodium		96	Vitamin K ₁		8,12
p-Tolualdehyde 2,4-DNPH	103	127,128	Vitamin K ₂		8,12
Toluene	69,78,97		Vitamin K ₃		8
Transferrin	137,145,150,151	52,53	X		
Trehalose		22	Xanthine	87	16,34,106,107
Triacylglycerol		59	Xanthosine	87	34,35
Triamcinolone		81	Xylitol		21,22
Triamcinolone acetoneide	89,199	80,81,82	Xylose		20
Triazolam	22		Z		
2,4,6-Trichlorophenol		132	Zeaxanthin	191	
Triclopyr	17				
Triethanolamine	95	135			
Trifluoperazine hydrochloride		98			

Ordering Information

The previous product listing represents commonly used standard column dimension. In order to identify any specific product version and order number, please see the example and the table below.

Full listing of all chemistries and dimensions

Gel Code								Hardware Code						
Chemistry Code		Pore size [nm]		Particle shape		Particle size [μm]		Length [mm]		Inner diameter [mm]		Column Type		
YMC30	CT	6	06	spherical	S	3	03	10	01	0.05	E5	Quick Seal Cartridge Waters type	QT	
Triart C18	TA	8	08			4	04	20	02	0.075	E8		QC	
Pro C18	AS	12	12			5	05	33	03	0.1	F0		WT	
Pro C18 RS	RS	20	20			6	06	50	05	0.2	G0			
Hydrosphere C18	HS	30	30					75	L5	0.3	H0			
ODS-A	AA	100	A0			10	11	100	10	0.5	J0			
ODS-AM	AM	proprietary	99			15	16	125	R5	0.8	M0			
ODS-AQ	AQ	non-porous	00			20	21	150	15	1.0	01			
J'sphere ODS-H80	JH					50	50	250	25	2.1	02			
J'sphere ODS-M80	JM					75	75			3.0	03			
J'sphere ODS-L80	JL							300	30	4.0	04			
ODS-AL	AL					63/210	A4	500	50					
PAH	YP					150	A5	1000	A0	4.6	46			
PolymerC18	PC									6.0	06			
Triart C8	TO									8.0	08			
Pro C8	OS									10	10			
C8 (Octyl)	OC									20	20			
YMCbasic	BA									30	30			
Ph (Phenyl)	PH													
Pro C4	BS									50 (2000 psi)	52			
C4 (Butyl)	BU									70 (2000 psi)	72			
Protein-RP	PR									100 (2000 psi)	A2			
TMS (C1)	TM									150 (2000 psi)	B2			
PVA-Sil	PV									200 (2000 psi)	C2			
Polyamine II	PB													
NH ₂ (Amino)	NH													
CN (Cyano)	CN													
Triart Diol-HILIC	TDH													
Diol	DL													
SIL (Silica)	SL													
BioPro-QA	QA													
BioPro-SP	SP													
BioPro-QA-F	QF													
BioPro-SP-F	SF													
Chiral NEA (R)	NR													
Chiral NEA (S)	NS													
Chiral CD BR α	DA													
Chiral CD BR β	DB													
Chiral CD BR γ	DG													
Chiral Prep CD ST	ST													
Chiral Prep CD PM	PM													

Example

Choose your column and fill in the “Gel and Hardware Code” or detailed description (The part number consists of the “Gel Code” and the “Hardware Code”).

YMC-Pack ODS-A

12 nm

spherical

3 μm

250 mm

1.0 mm

Quick Seal

AA

12

S

03

25

01

QT

Example

Choose your column and fill in the "Gel and Hardware Code" or detailed description (The part number consists of the "Gel Code" and the "Hardware Code").

YMC-Pack ODS-A	12 nm	spherical	3 µm	250 mm	1.0 mm	Quick Seal
AA	12	S	03	25	01	QT

Your column part number: **AA12S03-2501QT (Example)**

Please note that combinations of features cannot be selected at random, but only from the possible specifications for a chosen stationary phase. These can be determined from the individual product sections in this catalogue or from our homepage www.ymc.de.

For more details



contact your local distributor or

YMC Europe GmbH, Schöttmannshof 19, D-46539 Dinslaken, Phone: +49 (0) 2064 / 427-0,

Fax +49 (0) 2064 / 427-222, e-mail: info@ymc.de, homepage: www.ymc.de

The following brands, trademarks, or service marks are the property of the listed company and/or its subsidiaries. Every effort has been taken to ensure this list is accurate at the time of printing this brochure.

SureFit is trademark of IDEX Corporation.

Capcell Pak MG-II, Capcell Pak C18 AQ are trademarks of Shiseido Co., Ltd.

Intertsil ODS-3V, Intertsil ODS-3, Intertsil ODS-4, Inertsil C8-3 are trademarks of GL Sciences, inc..

Acquity UPLC BEH C18, Acquity UPLC CSH C18, Atlantis dC18, Atlantis T3, SunFire C18, Symmetry C18, Symmetry C8, X-Bridge C18, XTerra MS C18, OBD are trademarks of Waters Corp.

Gemini, Gemini NX, Kinetex C18 are trademarks of Phenomenex, Inc.

Hypersil GOLD and ProPac WCX-10 are trademarks of Thermo Fisher Scientific.

Kromasil Eternity is a trademark of EKA Chemicals AB.

TSK-gel ODS-100S, TSKgel DEAE-NPR, TSKgel SP-NPR, TSKgel BioAssist Q and TSKgel BioAssist S are trademarks of Tosoh Bioscience LLC.

Zorbax Extend-C18, Zorbax StableBond, Zorbax Eclipse Plus, Zorbax Eclipse XDB-C18, Zorbax SB-C18 and Zorbax XDB-C18 are trademarks of Agilent Technologies, Inc..

L-column is a trademark of Chemicals Evaluation and Research Institute.

Mini Q, Mono Q and Mono S are trademarks of General Electric Company as GE Healthcare.

ReproSil-Pur C18-AQ is a trademark of Dr. Maisch GmbH.

ECO, ECOplus and Kronlab are trademarks of YMC Europe GmbH.

Bio-Rad is a trademark of Bio-Rad Laboratories Inc..

GE is a trademark of General Electric Company.

Pall is a trademark of Pall Corporation.

Tosoh is a trademark of Tosoh Corporation.

Analytical stationary phases routinely available from YMC

		PRODUCT	PAGE	PHASE (silica-based unless stated)	END-CAPPED	USP CLASS NO.	PARTICLE SIZE (µm spherical)	PORE SIZE (nm)	CARBON LOAD (%C)	pH	TYPICAL APPLICATIONS	PAGE	PRODUCT
Reversed Phase	C30	C30	192	proprietary polymeric bonding chemistry	—	—	3, 5	proprietary	—	2.0-7.5	isomeric carotenes, retinols, steroids, fat-soluble vitamins	192	C30
	C18	Triart C18	13-40 216	multi-layer hybrid particle	yes	L1	1.9, 3, 5	12	20	1.0-12.0	acid, neutral, basic and chelating compounds, metabolites, "versatile" stationary phase	13 216	Triart C18
		Pro C18	70	very low residual non-specific interactions	yes	L1	3, 5	12	17	2.0-8.0	fat-soluble vitamins, antioxidants, metabolites, acidic, neutral, basic and chelating compounds	70	Pro C18
		UltraHT	54	2 µm Pro C18 for fast and ultra fast separations			2					54	UltraHT
		Pro C18 RS	76	high carbon load with polymeric bonding C18	yes	L1	3, 5	8	22	1.0-10.0	acidic and basic compounds	76	Pro C18 RS
		Hydrosphere C18	80	can be used in 100% aqueous eluent	yes	L1	3, 5	12	12	2.0-8.0	strong polar compounds, antibiotics, nucleic acids, water-soluble vitamins, acidic, neutral, basic and chelating compounds	80	Hydrosphere C18
		UltraHT	54	2 µm Hydrosphere C18 for fast and ultra fast separations			2					54	UltraHT
		ODS-A	90	one of the YMC's international bestsellers	yes	L1	3, 5	12, 20, 30	17, 12, 7	2.0-7.5	general purpose phase	90	ODS-A
		ODS-AM	92	high performance C18 column for validated methods operation	yes	L1	3, 5	12	17	2.0-7.5	purines, phenols, PTC-amino acids, angiotensins, alkaloids	92	ODS-AM
		ODS-AQ	86	"hydrophilic" endcapping, for 100% aqueous eluent systems	yes	L1	3, 5	12, 20	14, 10	2.0-7.5	strong polar compounds	86	ODS-AQ
		J'sphere	196	C18-family with differently controlled hydrophobicity for method development	yes	L1	4	8	22, 14, 9 (JH, JM, JL)	1.0-9.0 (JH) 2.0-7.5 (JM+JL)	positional isomers, complexing agents, pharmaceuticals	196	J'sphere
		ODS-AL	94	traditional C18 for "mixed mode" separations	no	L1	3, 5	12	17	2.0-7.5	tocopherols, fat-soluble vitamins, disinfectants	94	ODS-AL
		Polymer C18	96	polymethacrylate matrix, wide pH applicability	—	—	6	proprietary	C18 equivalent 10%	2.0-13.0	phenols, anilines, peptides in high pH, pharmaceuticals, quaternary amines	96	Polymer C18
	C8	Triart C8	13-40 216	multi-layer hybrid particle	yes	L1	1.9, 3, 5	12	17	1.0-12.0	acid, neutral, basic and chelating compounds, metabolites, "versatile" stationary phase	13 216	Triart C8
		Pro C8	72	C8, with very low residual non-specific interactions	yes	L7	3, 5	12	10	2.0-7.5	acidic, neutral, basic and chelating compounds, drugs and metabolites	72	Pro C8
		C8 (Octyl)	100	traditional C8	yes	L7	3, 5	12, 20, 30	10, 7, 4	2.0-7.5	proteins and peptides, estrogens, general purpose phase	100	C8 (Octyl)
		YMCbasic	98	monomeric bonded chains of C8 and smaller	—	L7	3, 5	proprietary	8	2.0-7.5	basic molecules w/o modifiers, anilines, alkaloids, antidepressants	98	YMCbasic
	C4	Ph (Phenyl)	102	monomeric bonded phenyl	yes	L11	3, 5	12, 30	9, 3	2.0-7.5	phenols, fullerenes, sweeteners	102	Ph (Phenyl)
		Pro C4	74	C4, with very low residual non-specific interactions	yes	L26	3, 5	12	8	2.0-7.5	polar acidic, neutral, basic and chelating compounds, polar peptides	74	Pro C4
		C4 (Butyl)	104	traditional C4	yes	L26	3, 5	12, 20, 30	7, 5, 3	2.0-7.5	biological separations, polar compounds	104	C4 (Butyl)
		PROTEIN-RP	158	high stability, good recovery rates	yes	L26	5	proprietary	—	1.5-7.5	proteins, peptides	158	PROTEIN-RP
		YMC-PAH	194	proprietary bonding chemistry	—	—	3, 5	—	—	2.0-8.0	polyaromatic hydrocarbons	194	YMC-PAH
Normal Phase / HILIC		TMS (C1)	106	trimethyl silane	—	L13	3, 5	12, 30	4, 3	2.0-7.5	water-soluble vitamins	106	TMS (C1)
		PVA-SIL	118	polyvinyl alcohol bonded on silica support	—	L24	5	12	—	2.0-9.5	proteins, phospholipids, retinoids, lipids	118	PVA-SIL
		Polyamine II (PBMN)	124	mixed secondary and tertiary amino derivative	—	—	5	12	—	2.0-7.5	malto-oligosaccharides, tocopherols, nucleotides, sugars	124	Polyamine II (PBMN)
		NH ₂ (Amino)	126	primary amino derivate	—	L8	3, 5	12	3	2.0-7.5	sugars, nucleotides, water-soluble vitamins	126	NH ₂ (Amino)
		CN (Cyano)	120	useful for SFC applications	yes	L10	3, 5	12, 30	7, 3	2.0-7.5	proteins, steroids, catechols	120	CN (Cyano)
		Triart Diol-HILIC	32	versatile HILIC column	—	L20	1.9; 3; 5	12	—	2.0-10.0	peptides, proteins, malto-oligosaccharides	32	Triart Diol-HILIC
		Diol	122 146	versatile alternative to silica for normal phase separations	—	L20	5	6, 12	—	2.0-7.5	peptides, proteins, malto-oligosaccharides	122 146	Diol
		SIL (Silica)	116 214	ultra high purity, high mechanical stability	—	L3	3, 5	6, 12, 20, 30	—	2.0-7.5	small organic molecules, fat-soluble vitamins, tocopherols	116 214	SIL (Silica)
IEX		BioPro QA / SP	136	high ion exchange capacity, porous hydrophilic polymer	—	—	5	100	—	2.0-12.0	proteins, peptides, nucleotides	136	BioPro QA / SP
		BioPro QA-F / SP-F	141	high ion exchange capacity, non-porous hydrophilic polymer	—	—	5	—	—	2.0-12.0	proteins, peptides, nucleotides	141	BioPro QA-F / SP-F
SEC		Diol-60, -120, -200, -300	122 146	versatile phase for gel filtration separations	—	L20	5	6, 12, 20, 30	—	5.0-7.5	peptides, proteins, malto-oligosaccharides	122 146	Diol-60, -120, -200, -300



ANALYTICAL SCIENCES

 **interchim**

211 bis, avenue JF Kennedy
BP 1140 - 03100 Montluçon
Fax 33 4 70 03 82 60



■ Hotline +33 4 70 03 73 09 ■



■ interchrom@interchim.com ■