

Stacker Parallel Purification System™

Rapid sequential and parallel purification in a 24 well micro titer plate footprint using standard filtration, phase separation and SPE columns...





The Stacker is designed to make your parallel and combinatorial chemistry workups and purification, quick and easy.

The principle of the Stacker allows efficient purification of samples from any standard 24 well format micro titer plate, as well as the GreenHouse reaction blocks, Lollipop, 24 x 7ml micro titer plates and other 24 position reactors/synthesisers.

The Stacker allows two levels of filtration, phase separation and/or SPE (Solid Phase Extraction) columns to be stacked, one above the other. This enables parallel and sequential filtration, phase separation and/or SPE in 24 samples.

Stacker Features & Benefits

- Stacked column design allows sequential purification e.g. phase separation and SPE or filtration and SPE
- Easily removable upper support plate allows unit to be used for single or sequential parallel purification
- Clear gas cover allows slight positive pressure to be applied to increase flow rate of difficult samples
- Design accommodates all standard 3ml and 6ml tabbed or tabless SPE or filtration columns
- Accepts all standard 24 well micro titer plates, as well as GreenHouse Reaction Blocks for fraction collection
- Unique recessed design of support plates allows positioning of standard format columns within micro titer plate footprint without column tabs obstructing each other
- RDT market a complete range of consumables dedicated to the Stacker and preparative chemistry. These include Filtration, Phase Separation and SPE Columns including Silica, SCX, Aminopropyl and SAX Acetate.



Solid Phase Extraction (SPE) Columns

RDT offer four standard types of SPE Columns in two sizes that will cover the vast majority of purifications for the Stacker. These are:-

| Cat No | Description | Volume | Pk Qty |
|---------|-------------------|--------|--------|
| RR99810 | Silica 500mg | 3ml | 50 |
| RR99811 | Silica 1.0g | 6ml | 30 |
| RR99812 | Aminopropyl 500mg | 3ml | 50 |
| RR99813 | Aminopropyl 1.0g | 6ml | 30 |
| RR99814 | SCX 500mg | 3ml | 50 |
| RR99815 | SCX 1.0g | 6ml | 30 |
| RR99816 | SAX Acetate 500mg | 3ml | 50 |
| RR99817 | SAX Acetate 1.0g | 6ml | 30 |

The following guidelines are intended to help in the initial selection of a column according to the main interactions of the solid phase, polar, non-polar or anion/cation exchange.

Silica

The main interaction is polar via the silanol functionality. Most useful for purifying moderately polar compounds where the impurities are either or both very non-polar (eluted first) or very polar (retained on column).

Aminopropyl





Normally used as a polar sorbent. The amino group can hydrogen bond with a variety of polar groups, e.g. hydroxyl, amino etc and will have polar interactions with polar groups. It will also function as a weak anion exchanger when protonated.

SCX

A strong cation exchanger used to extract basic compounds capable of holding a positive charge. Can be used for extracting basic impurities from non-basic products (which pass straight through), or for extracting basic products from non-basic impurities (which pass straight through). The basic products can then be eluted using a medium with a higher pH.

SAX Acetate

A strong anion exchanger used to extract compounds capable of carrying a negative charge (e.g. weak acids). Can be used for extracting acidic impurities from non-acidic products (which pass straight through), or for extracting acidic products from non-acidic impurities (which pass straight through). The acidic products can then be eluted using a medium with a lower pH.



20µm PE Filtration Columns

Typical applications include:

- Isolation of precipitates or crystallites from the reaction mixture
- Filtration of resins prior to cleavage
- Filtration of solid supported reagents post synthesis
- Removal of solid impurities prior to SPE

In the first two cases the Filtration Columns are used in the 'Lower Level' as subsequent SPE is not required. Pressure can be applied to the Columns using the Acrylic Gas Cover to increase the flow rate. If individual Columns are slow or become blocked pressure can be applied to individual columns using the appropriate tapered adapter. The crude product or resin is washed prior to isolation or cleavage etc as required.

Solid impurities, for example salts may be removed using the Filtration Columns in the 'Upper Level' prior to SPE in the 'Lower Level'. See the Stacker Instructions for further information.

Phase Separation Columns

The Phase Separation Columns contain a hydrophobic frit and are designed to separate organic phases that are denser than water (eg dichloromethane). They can be used in the 'Lower Level' if further purification is not required or in the 'Upper Level' if SPE purification is required after phase separation.

Following partition of a crude reaction mixture between for example dichloromethane and water, the mixture is transferred to a Phase Separation Column where the lower organic layer will run through the frit leaving the upper aqueous layer for disposal.

It should not be necessary to apply pressure, indeed doing so can allow the aqueous layer to break through the hydrophobic frit.

| Cat No | Description | Volume | Pk Qty |
|---------|---------------------------|--------|--------|
| RR99818 | 20µm PE Filtration Column | 3ml | 100 |
| RR99819 | 20µm PE Filtration Column | 6ml | 100 |
| RR99820 | Phase Separation Column | 3ml | 100 |
| RR99821 | Phase Separation Column | 6ml | 100 |

6 Channel Impact² Electronic Pipettor


The cordless, rechargeable Impact² pipettor accepts six extra-long, narrow 1250µl pipette tips, allowing you to reach to the very bottom of longer tubes or deep wells (including the GreenHouse Reaction Tubes). The Impact² allows programming and storage of five different protocols of up to 40 steps each, with one additional floating program. Plus 3 standard and 2 additional slow speed controls for working with viscous liquids. The Impact²'s unique expandable tip spacing system allows transfer of samples from 96 to 24 well formats with ease (freely expandable from 9mm to 19.81mm). [The Impact² pipettor is available in other volumes and with 8 or 12 channels - details on request.](#)



Stacker Parallel Purification System™

Rapid sequential and parallel purification
in a 24 well micro plate footprint using
standard filtration, phase separation
and SPE columns...

International patents pending



Stacker Accessories

Acrylic Gas Cover - This semi-disposable acrylic cover will eventually be attacked by solvents and will need replacing. Its simple design allow us to offer replacements at a very reasonable price.

Column Caps & Luer Tips - Disposable polypropylene red caps and clear Luer Tips are for sealing/plugging columns for storage or transportation.

Micro Titer Plates - Three popular 24 well micro titer plate styles ideal for fraction collection from the Stacker.



| Cat No | Description | Qty |
|---|--|-----|
| RR99800 | Stacker Parallel Purification System | 1 |
| SPE Columns | | |
| RR99810 | Silica 500mg/3ml | 50 |
| RR99811 | Silica 1.0g/6ml | 30 |
| RR99812 | Aminopropyl 500mg/3ml | 50 |
| RR99813 | Aminopropyl 1.0g/6ml | 30 |
| RR99814 | SCX 500mg/3ml | 50 |
| RR99815 | SCX 1.0g/6ml | 30 |
| RR99816 | SAX Acetate 500mg/3ml | 50 |
| RR99817 | SAX Acetate 1.0g/6ml | 30 |
| Filtration Columns | | |
| RR99818 | 20um PE Filtration 3ml | 100 |
| RR99819 | 20um PE Filtration 6ml | 100 |
| Phase Separation Columns | | |
| RR99820 | Phase Separation 3ml | 100 |
| RR99821 | Phase Separation 6ml | 100 |
| Accessories & Replacement Parts | | |
| RR99802 | Replacement Acrylic Gas Cover | 1 |
| RR99806 | Stepped Connector for 3ml Columns | 1 |
| RR99807 | Stepped Connector for 6ml Columns | 1 |
| RR99808 | Clear Gas Tubing, 1 metre | 1 |
| RR99809 | Spare Polypropylene Needle Tips | 50 |
| RR99822 | Column Caps 3ml | 100 |
| RR99823 | Column Caps 6ml | 100 |
| RR99824 | Luer Tip Caps | 100 |
| 24 Well Micro Titer Collection Plates | | |
| T-9000U | Titan 24 x 9ml Conical Bottomed PTFE Plate | 1 |
| RR98020 | Loose 7ml HDPE Pots | 100 |
| RR98021 | 24 Place Polycarbonate Pot Holder with Cover | 1 |
| RR99440 | 10ml x 24 Well Clear Polypropylene Micro Plate | 25 |
| Impact² Multi-channel Pipettors | | |
| RR70048 | Impact ² 6 x1250µl Multi-Channel Pipettor | 1 |
| RR70213 | 1250µl Extra-long Racked Pipette Tips | 720 |
| RR70051 | Impact ² Pipettor Stand | 1 |



A 24 position personal parallel synthesiser designed by medicinal chemists at GlaxoWellcome



Innovative micro titer plate solutions for compound storage, combi-chem, parallel synthesis and HTS

Radleys Discovery Technologies are specialists in combi-chem innovation...

As a dynamic organisation Radleys Discovery Technologies are best able to react to the needs of this constantly developing market sector. Radleys Discovery Technologies specific areas of expertise is focused on parallel synthesis and parallel purification consumables and apparatus. The essential products that are the cornerstone of your parallel chemistry program.

© Radleys Discovery Technologies 2001

INTERFINE CHEMICALS

interchim

211 bis, avenue JF Kennedy
BP 1140 - 03100 Montluçon
Fax +33 4 70 03 82 60

■ Hotline +33 4 70 03 73 01 • interfine@interchim.fr ■