

# Sample Prep - Solid Phase Extraction

## Introduction

### Example of publication with Interchim SPE product

#### **Upti-Clean™ CN**

Chadi Abbara and all, Laboratoire de Pharmacologie, Service de Pharmacie & Service d'addictologie Hôpital Paul BROUSSE 94800 Villejuif, France ; Journal of Pharmaceutical and Biomedical Analysis 41 (2006) 1011–1016 ; Development and validation of a method for the quantitation of 9tetrahydrocannabinol in human plasma by high performance liquid chromatography after solid-phase extraction.

#### **Upti-Clean™ SI**

F. Boudsocq, P. Benaim, Y. Canitrot, M. Knibiehler, F. Ausseil, J. P. Capp, A. Bieth, C. Long, B. David, I. Shevelev, E. Frierich-Heinecken, U. Hubscher, F. Amalric, G. Massiot, J. S. Hoffmann, and C. Cazaux, Equipe Instabilité Génétique et Cancer, Institut de Pharmacologie et de Biologie Structurale, Unité Mixte Recherche Centre National de la Recherche Scientifique 5089 (F.B., P.B., Y.C., J.P.C., A.B., J.S.H., C.C.), Centre de Recherche en Pharmacologie-Santé, Unité Mixte Recherche Centre National de la Recherche Scientifique/P. Fabre 2587 (M.K.), Centre de Criblage Pharmacologique, Unité Mixte Recherche Centre National de la Recherche Scientifique/P. Fabre 2646 (F.A.), and Chimie des Substances Naturelles Bioactives, Unité Mixte Recherche Centre National de la Recherche Scientifique/P. Fabre 2597 (C.L., B.D., G.M.), Institut de Sciences et Technologies du Médicament de Toulouse 3, Toulouse, France; and Institute of Veterinary Biochemistry and Molecular Biology, University of Zurich, Zurich, Switzerland (E.F.-H., U.H.) ; MOLECULAR PHARMACOLOGY Mol Pharmacol 67:1485–1492, 2005, Modulation of Cellular Response to Cisplatin by a Novel Inhibitor of DNA Polymerase B.

#### **Atoll™ XC**

R. Déporte and all, Department of Pharmacokinetic & Department of Biostatistic, Anticancer Centre René Gauducheau 44805 Nantes, France ; Journal of Chromatography B, xxx (2006) xxx–xxx ; High-performance liquid chromatographic assay with UV detection for measurement of dihydrouracil / uracil ratio in plasma.

#### **Upti-Clean™ SI/CN**

Roberto Alzaga and all, Environmental Chemistry Department, IIQAB-CSIC, Jordi Girona 18–26, E-08034 Barcelona, Spain ; Journal of Chromatography A, 1025 (2004) 133–138 ; Fast solid-phase extraction–gas chromatography–mass spectrometry procedure for oil fingerprinting Application to the Prestige oil spill.

#### **Upti-Clean™ C18U**

Sandrine Demanèche and all, Laboratoire de Biochimie et Biophysique des Systèmes Intégrés, Unité Mixte de Recherche CEA-CNRS-Université Joseph Fourier-UMR5092, Laboratoire de Chimie des Protéines, Département de Réponse et Dynamique Cellulaires, CEA-Grenoble, Grenoble, France ; Applied and Environmental Microbiology, 2004 November, p. 6714–6725 ; Identification and Functional Analysis of Two Aromatic-Ring-Hydroxylating Dioxygenases from a Sphingomonas Strain That Degrades Various Polycyclic Aromatic Hydrocarbons.

#### **Upti-Clean™ C18U**

Serge Krivobok and all, Laboratoire de Biochimie et Biophysique des Systèmes Intégrés and Laboratoire de Chimie des Protéines, Département de Réponse et Dynamique Cellulaires, CNRS UMR 5092, CEA-Grenoble, 38054 Grenoble Cedex 9, France ; JOURNAL OF BACTERIOLOGY, July 2003, p. 3828–3841 ; Identification of Pyrene-Induced Proteins in Mycobacterium sp. Strain 6PY1: Evidence for Two Ring-Hydroxylating Dioxygenases.

#### **Atoll™ XC & Upti-Clean™ Florisil**

J. Le Faouder, E. Bichon and all, LABERCA, Ecole Nationale Vétérinaire de Nantes, Route de Gachet, Atlanpôle La Chantrerie, BP 50707, 44087 Nantes Cedex 03, France ; Science Direct, Talanta 73 (2007) 710–717, Transfer assessment of fipronil residues from feed to cow milk.

#### **Atoll™ XWP**

A. Salvador and all, Université Claude Bernard UMR 5180 69622 Villeurbanne & CEPHAC Europe 86281 Saint-Benoit, France ; Chromatographia 2006, 63, 609–615 ; Simultaneous LC-MS-MS Analysis of Capecitabine and its Metabolites (5 $\beta$ -deoxy-5-fluorocytidine, 5 $\beta$ -deoxy-5-fluorouridine, 5-fluorouracil) After Off-Line SPE from Human Plasma.

#### **Atoll™ ATH**

Séverine Compain, Dimitri Schlemmer, Mikael Levi and all, CEA, Service de Pharmacologie et d'Immunologie, DSV/DRM, CEA/Saclay, 91191 Gif-sur-Yvette Cedex and SPIBio, Parc d'Activité du Pas du Lac, 10 bis avenue Ampère, F-78180 Montigny le Bretonneux, France ; JOURNAL OF MASS SPECTROMETRY, J. Mass Spectrom. 2005; 40: 9–18 ; Development and validation of a liquid chromatographic/tandem mass spectrometric assay for the quantitation of nucleoside HIV reverse transcriptase inhibitors in biological matrices.

# Sample Prep - Solid Phase Extraction

## Introduction

Solid Phase Extraction continues to be the fastest growing technique utilised for sample preparation. The ease of use and flexibility of SPE means that increasingly, this is the chosen pre-step adopted to clean and concentrate samples prior to analysis in HPLC, HPLC/MS, GC or GC/MS.

Advances in the analytical process are placing greater demands and expectations on sample cleaning and therefore increasing the quality required from SPE products. For this reason, polymeric sorbents with high loading capacities and spherical ultrapure silica have become widespread.

Recovery, capacity, selectivity & reproducibility are the principal sample prep. demands of today's analyst. We have developed a state-of-the-art SPE product range incorporating silica and polymer based technology. Upti-prep™ (silica) and Atoll™ (polymeric) push the boundaries of expectation from modern day sample preparation challenges.

Interchim's full range of SPE products is highlighted within this chapter. The typical procedure for the SPE process is described below and considerations toward sorbent and working format selection follow.

## Typical SPE procedure

### Bed volume definition :

The bed volume is defined as the minimum volume of solvent necessary to wet the defined quantity of sorbent within the column. This can vary depending on the nature of the sorbent.

e.g. : ~ 120 µl per 100 mg of silica gel sorbent 60 Å  
~ 600 µl per 500 mg of silica gel 60 Å

[Incomplete elution of compound of interest will occur if the sorbent mass is too large for the volume of solvent used.  
Incomplete retention of compounds of interest will occur if there is an inadequate sorbent mass leading to compound eluting in the fraction or in the washing solvent. Such cases lead to lower recovery rates].


## 1. Conditioning Step

Sorbent activation and functional group activation are achieved by passing a volume of an appropriate solvent or a mixture of solvent, through the column. Column frits are simultaneously solvated.

Methanol or acetonitrile are commonly used for activating hydrophobic sorbents, whilst hexane or dichloromethane activate hydrophilic sorbents. 2 to 4 bed volumes are typically recommended.

### Technical tips

- Check solvent miscibility before using.
- Keep the solvent level beyond the sorbent to maintain its activated status.
- For ion exchange bonded silica, activate with methanol, water and then water buffered to the pH of choice.
- Vacuum & positive pressure manifolds are recommended to achieve more reproducible extractions



# Sample Prep - Solid Phase Extraction

## Introduction

### 2. Sample loading step

Apply sample onto the upper part of the sorbent bed. Matrix contaminants may pass through the column unretained, and additionally, other matrix components may be more or less strongly retained on the sorbent surface. To get a maximum purification efficiency, the sample flow needs to be controlled.

To achieve faster flow of viscous sample through a column, 90 to 140  $\mu\text{m}$  sorbents can be used. The exchange capacity and selectivity are unaffected.

[It is necessary to analyze the unretained fraction to check if all compounds of interest have been retained]

### 3. Washing step

Passing solvents through columns washes away interfering compounds whilst leaving the analyte undisturbed on the sorbent bed. Different solvents or solvent mixtures may be used to improve the rinsing efficiency.

### 4. Drying step

Solvent traces are evaporated by circulating air through the column over a 2 to 10 minute time period. This improves the extraction yield.

### 5. Elution step

An appropriate solvent is passed through the column to disrupt the analyte-sorbent interaction and to elute 100% of compounds of interest.


The appropriate solvent must have maximum interaction with the compound of interest and a minimal interaction with the remaining impurities, leaving them undisturbed on the sorbent bed. In addition the volume of the elution solvent needs to be as small as possible to maximise the concentration factor.

[Sorbent with low particle size (e.g 30,50  $\mu\text{m}$ ) gives a lower elution volume than larger sorbent particle size (e.g 90, 140  $\mu\text{m}$ )].

### 6. Concentration step

Compounds of interest are concentrated by evaporating a part of the solvent. If necessary, dry the eluate with anhydrous sulfate to remove possible water traces. The concentrated sample is then ready for analysis.

Interchim recommends that all steps should be carefully optimized according to your specific extraction. This will improve the quality of the final analysis.



# Sample Prep - Solid Phase Extraction

## Selecting the appropriate sorbent

**Sorbent selection requires consideration of sample volume, the nature of the analyte, analyte concentration and the inherent properties of the sorbent itself. For environmental studies, a volume of several hundred millilitres might be necessary for a good pre-concentration (e.g., organic pollutants) whereas in the pharmaceutical industry, the sample volumes that require cleaning may only be millilitres.**

The selected sorbent needs to have an excellent affinity for the compounds of interest and at the same time a weak affinity for irrelevant compounds within the matrix. Choosing the correct sorbent results in a specific selectivity for the compounds of interest. A sufficient loading capacity also needs to be identified to optimise retention volumes of the desired compound.

There are four general modes used in Solid Phase Extraction : reversed phase, normal phase and ion exchange that require different sorbent types namely hydrophobic, hydrophilic, ion-exchange and mixed mode.

### Hydrophobic sorbents

In reversed phase, the non polar functional groups of the sorbent operate according to Van der Waals forces. Alkyl and aromatic chains are function groups that have affinity with non-polar and mid-polar compounds.

Free silanol groups left on the sorbent favour polar interactions.

For aromatic compounds, eg, pharma-based chemistries, we recommend to use phenyl selectivity (polystyrene divinylbenzene polymers).

### Hydrophilic sorbents

Normal phase provides an efficient cleaning of molecules with polar function groups.

Cyano (CN) functional groups can be used either in normal phase to extract polar compounds or in reversed phase for mid-polar compounds.

Diol functional groups can enhance polar compound extraction compared to virgin silica.

Amino sorbent ( $\text{NH}_2$ ) can be used either as weak anion exchangers (for strong acids), or as a polar sorbent that can interact with OH, NH, SH ...

### Ion exchange sorbents

Ion-exchange retention is based on ionic interaction. This sorbent creates a strong attraction with opposite functional groups of the sample compounds.

Ion exchange sorbent interactions depend essentially on counter-ion pH and ionic strength.

Strong anion exchange phases (SAX) possess a strong quaternary amine. They are used to extract weak acids which have one or more negative charges.

Strong cation exchange phases (SCX) contain sulfonic acid that are used to extract weak basic compounds which have one or more positive charges.

Weak anion exchange phases (DEAE,  $\text{NH}_2$ ) possess a diethyl amino ethyl and amino group. They are used to extract strong acids which have one or more negative charges.

Weak cation exchange phases (WCX) contain carboxylic acid that are used to extract strong basic compounds which have one or more positive charges.


### Mixed mode sorbents

Mixed mode sorbents exhibit the greatest selectivity. Ion exchange and hydrophobic chains are bonded onto the surface of silica providing unique selectivity.

This technique is used for basic compound extraction especially for drugs and metabolites within biologic fluids. Initially compounds that possess acid or basic functionality are retained by ion exchange functionality. A washing step with an appropriate pH, removes ionizable impurities. Passing an organic solvent through the column will then remove retained impurities that result from hydrophobic bonding.

# Sample Prep - Solid Phase Extraction

## Selecting the appropriate sorbent



### Ion-exchange capacity


The capacity of ion-exchange or mixed mode sorbent is measured by the mass of ionizable compounds that interact with a pre-defined sorbent mass.  
E.g : For 1g of bonded SCX sorbent with an exchange capacity of 0.6 meq/g, it is possible to retain 0.6 mmol of compound.

### Optimal pH

Acids and bases in the sample need to be in their ionized form to develop interactions with the sorbent.  
To maintain reproducible and repeatable recovery rates, sample and sorbent need to be buffered at an optimised pH.  
Eg : for an acid w/pKa 3.6 & a base w/pKa 10.6, the pH zone should be between pH :5.6 and 8.6.











### Relative selectivity of the counter-ions

A counter-ion is an ionic entity able to interact with an ion exchange sorbent.  
It improves the efficiency of cleaning steps including elution according to its concentration in solution and its affinity with the exchanger sorbent.



# Sample Prep - Solid Phase Extraction

## Sorbent selection guide

Upti-clean™ C18, C18U, LCC, C8, PH	
Recovery™ C18	
Atoll™ ATL	
Atoll™ XC	
Atoll™ XWP	
Atoll™ AEV	
Atoll™ ATH	
Upti-clean™ SI, CN, Diol	
Upti-clean™ NH₂, DEAE, SAX, WCX, SCX	
Upti-clean™ MM1, MM2, MM3, MM4	

Polarity of analytes




Polar

Mid-Polar

non-Polar

Capacity to work w/ aqueous matrix :

high



médium



low



Capacity to interact w/ analyte :





# Sample Prep - Solid Phase Extraction

## Sorbent Properties & Characteristics

### Sorbent Selection

Silicas & Polymers are the most popular sorbents used for Solid Phase Extraction. Polymer loading capacities are higher than silica sorbents, however, silica sorbents exhibit greater selectivity. Sorbent features are highlighted below and provide the analyst with an initial consideration for appropriate sorbent selection.

#### Polymers :

Polymer sorbents are very stable from pH 1 to 14, they exhibit high loading capacities allowing for the cleaning of a broad range of compounds through a variety of matrices (waters, oils, plasma, urines...). Interchim polymers have a very high specific surface area that maximises pi-pi interactions. The capacity of our polymers are typically 15% greater than competitive polymers and 25% higher than silicas. These polymers are particularly suited for polar compound cleaning. The polymer surface can be easily modified and facilitates a large selectivity range from hydrophobic to hydrophilic interactions.

#### Silicas :

Silica & bonded silica are rigid supports that do not shrink or swell with solvents. The silica surface can be easily modified, this creates a potential for a large selectivity for SPE from hydrophobic to hydrophilic interactions. The pH stability of bonded silica is limited, typically to within the range of 2 to 7.5, this is chemistry dependant. Interchim offers more than 30 different silica based selectivities. Our sorbents take advantage of our ultra pure spherical silica (Upti-prep), and this achieves greater reproducibility, and establishes repeatable extraction and optimized sample recoveries.

#### The Interchim Product range :

Atoll™ (polymer), Upti-Clean™ (silica) and Recovery™ (silica) are Interchim's sorbent technologies for Solid Phase Extraction.

Within any sample cleanup process, the analyst has to consider the characteristics of their sample relative to the features of the SPE products available. Interchim provides, on a custom basis, method development kits which combine sorbent type, quantity of media and suitable housing format.

Please contact Interchim's technical support center for assistance in your selection process.



# Sample Prep - Solid Phase Extraction

## Sorbent Properties & Characteristics

Description	Sorbent	Interchim type	Comment	Exchange capacity	End-capping	Porosity A	Specific surface area m <sup>2</sup> /g	Particle diameter µm	Rank	Purity
Atoll	PSDVB	30XC	High capacity		No	n.a.	1500	30	Spherical	Pure
Atoll	PSDVB	XC	High capacity		No	n.a.	1500	70	Spherical	Pure
Atoll	PSDVB	XWP	High capacity		No	Wide Pore	1200	90	Spherical	Pure
Atoll	n.a.	ATH	Hydrophilic		No	70	800	75	Spherical	Pure
Atoll	n.a.	30ATH	Hydrophilic		No	70	800	30	Spherical	Pure
Atoll	PSHEMA	30AEV	Hydrophilic / hydrophobic		No	70	800	30	Spherical	Pure
Atoll	PSHEMA	AEV	Hydrophilic / hydrophobic		No	70	800	75	Spherical	Pure
Atoll	PSDVB	ATL	Hydrophobic		No	70	800	100	Spherical	Pure
Recovery	Silica Upti-prep	REC	C18		n.a.	n.a.	n.a.	50	Spherical	Pure
Recovery	Silica Upti-prep	RSI	SI		No	n.a.	n.a.	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	C18-S	% C : 18		Yes	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	C18U-S	% C : 16		No	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	LCC	% C : 10		Yes	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	RPAQ	C18, 100% hydrophilic		No	n.a.	n.a.	75	Spherical	Pure
Upti-Clean	Silica Upti-prep	C18-S2F	High flow		Yes	60	500	140	Spherical	Pure
Upti-Clean	Silica Upti-prep	C18U-S2F	High flow		No	60	500	140	Spherical	Pure
Upti-Clean	Silica Upti-prep	C8-S	% C : 11		Yes	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	C8U-S	% C : 9		No	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	C8-S2F	High flow		Yes	60	500	140	Spherical	Pure
Upti-Clean	Silica Upti-prep	C2	% C : 6		Yes	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	PH-S	% C : 9		No	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	CH	% C : 10,5		No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	SI-S			No	60	500	50	Granular	Pure
Upti-Clean	Silica Upti-prep	NH2-S	% C : 5		No	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	CN-S	% C : 8		Yes	60	500	50	Spherical	Pure
Upti-Clean	Silica Upti-prep	OH	% C : 7,5		No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	SCX	Strong acid	0,7 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	WCX	Weak acid	0,22 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	DEAE	Weak base	0,33 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	SAX	Strong base	0,30 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	MM1	Mixed mode RP / SCX	0,09 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	MM2	Mixed mode RP / WCX	0,10 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	MM3	Mixed mode RP / SAX	0,14 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Silica Upti-prep	MM4	Mixed mode RP / WAX	0,15 meq/g	No	60	450	60	Granular	Pure
Upti-Clean	Florisil	FL	Standard grade		No	n.a.	n.a.	200	Granular	Std.
Upti-Clean	Florisil	FLPR	Pesticide Residue grade		No	n.a.	n.a.	200	Granular	Std.
Upti-Clean	Polyamide	P6	6		No	n.a.	n.a.	100	Granular	Std.
Upti-Clean	Amberlite	XAD2	Polystyrene		No	90	330	16-50 mesh	Spherical	Std.
Upti-Clean	Amberlite	XAD4	Polystyrene		No	50	750	16-50 mesh	Spherical	Std.
Upti-Clean	Amberlite	XAD7	Acrylic Ester		No	80	450	16-50 mesh	Spherical	Std.
Upti-Clean	Amberlite	XAD16	Polystyrene		No	100	825	16-50 mesh	Spherical	Std.
Upti-Clean	Alumina	ALA	Acid		No	n.a.	200	32-63	Granular	Std.
Upti-Clean	Alumina	ALN	Neutral		No	n.a.	200	32-63	Granular	Std.
Upti-Clean	Alumina	ALB	Basic		No	n.a.	200	32-63	Granular	Std.
Upti-Clean	Silica	BCD	β-Cyclodextrine		No	n.a.	n.a.	40	n.a.	Std.
Upti-Clean	Carbon	CG	Graphitized		No	non-porous	100	n.a.	Granular	Std.
Upti-Clean	Carbon	CA	Activated		No	n.a.	n.a.	40	Granular	Std.





# Sample Prep - Solid Phase Extraction

## Sorbent Properties & Characteristics









# Sample Prep - Solid Phase Extraction

## SPE format

Interchim has developed a comprehensive line of SPE products to support efficient and reproducible cleanup for a wide range of sample type and volume.

Our technological advances ensure our SPE range provide a margin of accuracy to within 1%. 96 well plates are delivered with a weighing certificate as a warranty of the real mass of sorbent within each well.

SPE format	Housing	Nature	Volume	Frits
Standard columns		Polypropylene MG	(1 - 3 - 6 - 15 -25 - 75 - 150) ml	(10 - 20 - 70) $\mu$ m Polyethylene or Teflon
LRC columns		Polypropylene MG	Robotic Large Capacity (LRC) 15ml	20 $\mu$ m Polyethylene or Teflon
Glass columns		Glass	6 ml	20 $\mu$ m Teflon
Cartridges		Polypropylene MG	300 - 600 - 900 mg	20 $\mu$ m Polyethylene
96 well plates		Polypropylene MG	2 ml square well	20 $\mu$ m Polyethylene
48 well plates		Polypropylene MG	5 ml square well	20 $\mu$ m Polyethylene

# Sample Prep - Solid Phase Extraction


## Upti-Clean™ - Silica based

Upti-Clean™ are Interchim's range of silica based SPE products. They have been developed from a spherical high purity (99.9%) silica called Upti-prep™. Upti-prep™ has strictly controlled porosity and particle size ( $\pm 10A$ ,  $\pm 5 \mu m$ ) working in a pH range of 1 to 7.5

Interchim bonding technology ensures greater batch to batch reproducibility for our bonded silica, there is no longer a requirement for batch reservations. Upti-Clean™ products subsequently achieve superior recovery rates relative to traditional irregular silicas, exhibiting excellent reproducibility & consistency

The Upti-Clean™ range achieves accuracy of  $\pm 1\%$  in mass per column thanks to a strictly controlled packing process. Upti-Clean™ products are thoroughly quality control tested in-house to guarantee traceability. Products are supplied with an individual certificate detailing the specific production number and sorbent batch.

Upti-Clean™ are supplied in a variety of formats relative to current extraction apparatus and automated systems. Packaging has been specifically designed to guarantee the integrity of the product.



### Upti-Clean™ cartridges

Upti-Clean™ cartridges are easy to use, efficient SPE cleaning devices with widespread applications throughout the pharmaceutical, toxicology and clinical fields

- Luer inlet & outlet fitting
- No specific SPE apparatus necessary
- Polypropylene MG housing
- C18, C8 bonded chemistries for non and mid-polar (aqueous) compounds
- Virgin silica for polar compound purification in non-polar solvents
- Useful for sample storage & transport



Weight	Qty	C18	C8	Silica
300 mg	25 u	C18-300/SC-25	C8-300/SC-25	SI-300/SC-25
300 mg	50 u	C18-300/SC-50	C8-300/SC-50	SI-300/SC-50
300 mg	100 u	C18-300/SC-100	C8-300/SC-100	SI-300/SC-100
600 mg	25 u	C18-600/SC-25	C8-600/SC-25	SI-600/SC-25
600 mg	50 u	C18-600/SC-50	C8-600/SC-50	SI-600/SC-50
600 mg	100 u	C18-600/SC-100	C8-600/SC-100	SI-600/SC-100
900 mg	25 u	C18-900/SC-25	C8-900/SC-25	SI-900/SC-25
900 mg	50 u	C18-900/SC-50	C8-900/SC-50	SI-900/SC-50
900 mg	100 u	C18-900/SC-100	C8-900/SC-100	SI-900/SC-100

# Sample Prep - Solid Phase Extraction

## Recovery™ - Silica based

Interchim Recovery™ columns address recovery and reproducibility problems, highlighted in recent studies, that are associated with only a part of the standard 60 Å silica's specific surface area being accessible in SPE silica based cleanup procedures.

Recovery™ columns utilise an optimised version of Upti-prep™ silica. They prevent physical phenomena associated with older generation silica sorbent and utilise 100% of their specific surface area. Recovery™ can be used in all solvent conditions (including 100% water) achieving greater reproducibility and consistency.

- **Recovery™ REC** : C18, fully end-capped for non-polar, mid-polar & polar compounds in aqueous environment.
- **Recovery™ RSI** : Virgin silica for polar and mid-polar compounds from organic matrices.

### Recovery™ columns

Weight	Volume	Qty	Recovery C18	Recovery Silica
<b>Std. columns / PE frits</b>				
50 mg	1 ml	50 u	REC-50/1	RSI-50/1
100 mg	1 ml	100 u	REC-100/1	RSI-100/1
100 mg	3 ml	50 u	REC-100/3	RSI-100/3
200 mg	3 ml	50 u	REC-200/3	RSI-200/3
500 mg	3 ml	50 u	REC-500/3	RSI-500/3
500 mg	6 ml	30 u	REC-500/6	RSI-500/6
1000 mg	6 ml	30 u	REC-1G/6	RSI-1G/6
2000 mg	6 ml	20 u	REC-2G/6	RSI-2G/6
2000 mg	15 ml	20 u	REC-2G/15	RSI-2G/15
2000 mg	25 ml	20 u	REC-2G/25	RSI-2G/25
<b>LRC columns / PE frits</b>				
100 mg	LRC 15	50 u	REC-100LRC	RSI-100LRC
200 mg	LRC 15	50 u	REC-200LRC	RSI-200LRC
500 mg	LRC 15	50 u	REC-500LRC	RSI-500LRC
<b>Std. columns / PTFE frits</b>				
50 mg	1 ml	50 u	REC-50/1T	RSI-50/1T
100 mg	1 ml	100 u	REC-100/1T	RSI-100/1T
100 mg	3 ml	50 u	REC-100/3T	RSI-100/3T
200 mg	3 ml	50 u	REC-200/3T	RSI-200/3T
500 mg	3 ml	50 u	REC-500/3T	RSI-500/3T
500 mg	6 ml	30 u	REC-500/6T	RSI-500/6T
1000 mg	6 ml	30 u	REC-1G/6T	RSI-1G/6T
<b>Glass columns / PTFE frits</b>				
200 mg	6 ml	30 u	REC-200/6G	RSI-200/6G
500 mg	6 ml	30 u	REC-500/6G	RSI-500/6G
1000 mg	6 ml	30 u	REC-1G/6G	RSI-1G/6G



Std. columns  
PP straight tube + 20µm PE frits



LRC columns  
PP tube + 20µm PE frits



Std. columns  
PP tube + 20µm PTFE frits



Glass columns  
Glass tube + 20µm PTFE frits

Glass tubes and PTFE frits are solvent resistant. These materials guarantee purifications without extractable.

# Sample Prep - Solid Phase Extraction

## Recovery™ - Silica based

### Recovery™ 96 well plates

Automated cleanup procedures are now an integral part of the modern laboratory.

Interchim Recovery™ is available in a 96 well format with a 2 ml volume per well. This high quality unit provides rapid throughput of sample, within an automated process, whilst establishing excellent well-to-well consistency.

Recovery™ 96 well plates utilise Interchim's state of the art automated weight machine. This packing technology guarantees unprecedented accuracy (+/- 1% mass per column) compared to the commonly utilised competitor volumetric systems. Recovery™ therefore raises the standards & client expectations for sample recovery /reproducibility.

Recovery™ 96 well plates are QC tested in-house to guarantee tracability. Products are supplied with an individual certificate detailing Mfg number, Sorbent batch number & specifications. An additional certificate is supplied that states real sorbent mass in each individual well.

Recovery™ 96 well plates are packed in a PEHD/Al bag for long term integral storage. This prevents potential damage from UV and moisture.

Recovery™ 96 well plates are manufactured from polypropylene, conform to a standard footprint (127,76 mm x 85,47 mm x 19,74 mm) and are compatible with existing systems on the market.



### Recovery® 96 Well plate - 2 ml

Weight	Qty	Recovery™ C18	Recovery™ Silica
30 mg	1 u	REC-30/WP20	RSI-30/WP20
50 mg	1 u	REC-50/WP20	RSI-50/WP20
60 mg	1 u	REC-60/WP20	RSI-60/WP20
100 mg	1 u	REC-100/WP20	RSI-100/WP20
150 mg	1 u	REC-150/WP20	RSI-150/WP20
200 mg	1 u	REC-200/WP20	RSI-200/WP20



# Sample Prep - Solid Phase Extraction

## Upti-Clean™ Series-S - Silica based

Upti-Clean™ Series-S columns are easy to use, efficient SPE cleaning devices with widespread applications throughout the pharmaceutical, toxicology and clinical fields.

- 30 selectivities available
- Polypropylene MG or Glass housing
- Std. PE or PTFE frits
- Luer outlet fitting

### Upti-Clean™ Series-S, columns

#### Hydrophobic bonding

Weight	Volume	Qty	C18-S	C18U-S	LCC	RPAQ	C8-S	C8U-S	C2	PH-S
<b>Std. columns / PE frits</b>										
50 mg	1 ml	50 u	C18-S-50/1	C18U-S-50/1	LCC-50/1	RPAQ-50/1	C8-S-50/1	C8U-S-50/1	C2-50/1	PH-S-50/1
100 mg	1 ml	100 u	C18-S-100/1	C18U-S-100/1	LCC-100/1	RPAQ-100/1	C8-S-100/1	C8U-S-100/1	C2-100/1	PH-S-100/1
100 mg	3 ml	50 u	C18-S-100/3	C18U-S-100/3	LCC-100/3	RPAQ-100/3	C8-S-100/3	C8U-S-100/3	C2-100/3	PH-S-100/3
200 mg	3 ml	50 u	C18-S-200/3	C18U-S-200/3	LCC-200/3	RPAQ-200/3	C8-S-200/3	C8U-S-200/3	C2-200/3	PH-S-200/3
500 mg	3 ml	50 u	C18-S-500/3	C18U-S-500/3	LCC-500/3	RPAQ-500/3	C8-S-500/3	C8U-S-500/3	C2-500/3	PH-S-500/3
500 mg	6 ml	30 u	C18-S-500/6	C18U-S-500/6	LCC-500/6	RPAQ-500/6	C8-S-500/6	C8U-S-500/6	C2-500/6	PH-S-500/6
1000 mg	6 ml	30 u	C18-S-1G/6	C18U-S-1G/6	LCC-1G/6	RPAQ-1G/6	C8-S-1G/6	C8U-S-1G/6	C2-1G/6	PH-S-1G/6
2000 mg	6 ml	20 u	C18-S-2G/6	C18U-S-2G/6	LCC-2G/6	RPAQ-2G/6	C8-S-2G/6	C8U-S-2G/6	C2-2G/6	PH-S-2G/6
2000 mg	15 ml	20 u	C18-S-2G/15	C18U-S-2G/15	LCC-2G/15	RPAQ-2G/15	C8-S-2G/15	C8U-S-2G/15	C2-2G/15	PH-S-2G/15
2000 mg	25 ml	20 u	C18-S-2G/25	C18U-S-2G/25	LCC-2G/25	RPAQ-2G/25	C8-S-2G/25	C8U-S-2G/25	C2-2G/25	PH-S-2G/25
<b>LRC columns / PE frits</b>										
100 mg	LRC 15	50 u	C18-S-100LRC	C18U-S-100LRC	LCC-100LRC	RPAQ-100LRC	C8-S-100LRC	C8U-S-100LRC	C2-100LRC	PH-S-100LRC
200 mg	LRC 15	50 u	C18-S-200LRC	C18U-S-200LRC	LCC-200LRC	RPAQ-200LRC	C8-S-200LRC	C8U-S-200LRC	C2-200LRC	PH-S-200LRC
500 mg	LRC 15	50 u	C18-S-500LRC	C18U-S-500LRC	LCC-500LRC	RPAQ-500LRC	C8-S-500LRC	C8U-S-500LRC	C2-500LRC	PH-S-500LRC
<b>Std. columns / PTFE frits</b>										
50 mg	1 ml	50 u	C18-S-50/1T	C18U-S-50/1T	LCC-50/1T	RPAQ-50/1T	C8-S-50/1T	C8U-S-50/1T	C2-50/1T	PH-S-50/1T
100 mg	1 ml	100 u	C18-S-100/1T	C18U-S-100/1T	LCC-100/1T	RPAQ-100/1T	C8-S-100/1T	C8U-S-100/1T	C2-100/1T	PH-S-100/1T
100 mg	3 ml	50 u	C18-S-100/3T	C18U-S-100/3T	LCC-100/3T	RPAQ-100/3T	C8-S-100/3T	C8U-S-100/3T	C2-100/3T	PH-S-100/3T
200 mg	3 ml	50 u	C18-S-200/3T	C18U-S-200/3T	LCC-200/3T	RPAQ-200/3T	C8-S-200/3T	C8U-S-200/3T	C2-200/3T	PH-S-200/3T
500 mg	3 ml	50 u	C18-S-500/3T	C18U-S-500/3T	LCC-500/3T	RPAQ-500/3T	C8-S-500/3T	C8U-S-500/3T	C2-500/3T	PH-S-500/3T
500 mg	6 ml	30 u	C18-S-500/6T	C18U-S-500/6T	LCC-500/6T	RPAQ-500/6T	C8-S-500/6T	C8U-S-500/6T	C2-500/6T	PH-S-500/6T
1000 mg	6 ml	30 u	C18-S-1G/6T	C18U-S-1G/6T	LCC-1G/6T	RPAQ-1G/6T	C8-S-1G/6T	C8U-S-1G/6T	C2-1G/6T	PH-S-1G/6T
<b>Glass columns / PTFE frits</b>										
200 mg	6 ml	30 u	C18-S-200/6G	C18U-S-200/6G	LCC-200/6G	RPAQ-200/6G	C8-S-200/6G	C8U-S-200/6G	C2-200/6G	PH-S-200/6G
500 mg	6 ml	30 u	C18-S-500/6G	C18U-S-500/6G	LCC-500/6G	RPAQ-500/6G	C8-S-500/6G	C8U-S-500/6G	C2-500/6G	PH-S-500/6G
1000 mg	6 ml	30 u	C18-S-1G/6G	C18U-S-1G/6G	LCC-1G/6G	RPAQ-1G/6G	C8-S-1G/6G	C8U-S-1G/6G	C2-1G/6G	PH-S-1G/6G



# Sample Prep - Solid Phase Extraction

## Upti-Clean™ Series-S - Silica based

### Upti-Clean™ Series-S 96 well plates

Automated cleanup procedures are now an integral part of the modern laboratory.

Interchim Upti-Clean™ is available in a 96 well format with a 2 ml volume per well. This high quality unit provides rapid throughput of sample, within an automated process, whilst establishing excellent well-to-well consistency.

Upti-Clean™ 96 well plates utilise Interchim's state of the art automated weight machine. This packing technology guarantees unprecedented accuracy (+/- 1% mass per column) compared to the commonly utilised competitor volumetric systems. Upti-Clean™ therefore raises the standards & client expectations for sample recovery /reproducibility.

Upti-Clean™ 96 well plates are QC tested in-house to guarantee tracability. Products are supplied with an individual certificate detailing Mfg number, Sorbent batch number & specifications. An additional certificate is supplied that states real sorbent mass in each individual well.

Upti-Clean™ 96 well plates are packed in a PEHD/Al bag for long term integral storage. This prevents potential damage from UV and moisture.

Upti-Clean™ 96 well plates are manufactured from polypropylene, conform to a standard footprint (127,76 mm x 85,47 mm x 19,74 mm) and are compatible with existing systems on the market.



### Upti-Clean™ Series-S, 96 well plates - 2 ml

#### Hydrophobic bonding

Weight	Qty	C18-S	C18U-S	LCC	RPAQ	C8-S
30 mg	1 u	C18-30/WP20	C18U-30/WP20	LCC-30/WP20	RPAQ-30/WP20	C8-30/WP20
50 mg	1 u	C18-50/WP20	C18U-50/WP20	LCC-50/WP20	RPAQ-50/WP20	C8-50/WP20
60 mg	1 u	C18-60/WP20	C18U-60/WP20	LCC-60/WP20	RPAQ-60/WP20	C8-60/WP20
100 mg	1 u	C18-100/WP20	C18U-100/WP20	LCC-100/WP20	RPAQ-100/WP20	C8-100/WP20
150 mg	1 u	C18-150/WP20	C18U-150/WP20	LCC-150/WP20	RPAQ-150/WP20	C8-150/WP20
200 mg	1 u	C18-200/WP20	C18U-200/WP20	LCC-200/WP20	RPAQ-200/WP20	C8-200/WP20
		C8U-S	C2	PH-S	CH	
30 mg	1 u	C8U-30/WP20	C2-30/WP20	PH-30/WP20	CH-30/WP20	
50 mg	1 u	C8U-50/WP20	C2-50/WP20	PH-50/WP20	CH-50/WP20	
60 mg	1 u	C8U-60/WP20	C2-60/WP20	PH-60/WP20	CH-60/WP20	
100 mg	1 u	C8U-100/WP20	C2-100/WP20	PH-100/WP20	CH-100/WP20	
150 mg	1 u	C8U-150/WP20	C2-150/WP20	PH-150/WP20	CH-150/WP20	
200 mg	1 u	C8U-200/WP20	C2-200/WP20	PH-200/WP20	CH-200/WP20	



# Sample Prep - Solid Phase Extraction

Upti-Clean™ Series-S - Silica based

## Hydrophilic columns

Weight	Volume	Qty	SI-S	NH2-S	CN-S	OH
<b>Std. columns / PE frits</b>						
50 mg	1 ml	50 u	SI-S-50/1	NH2-S-50/1	CN-S-50/1	OH-50/1
100 mg	1 ml	100 u	SI-S-100/1	NH2-S-100/1	CN-S-100/1	OH-100/1
100 mg	3 ml	50 u	SI-S-100/3	NH2-S-100/3	CN-S-100/3	OH-100/3
200 mg	3 ml	50 u	SI-S-200/3	NH2-S-200/3	CN-S-200/3	OH-200/3
500 mg	3 ml	50 u	SI-S-500/3	NH2-S-500/3	CN-S-500/3	OH-500/3
500 mg	6 ml	30 u	SI-S-500/6	NH2-S-500/6	CN-S-500/6	OH-500/6
1000 mg	6 ml	30 u	SI-S-1G/6	NH2-S-1G/6	CN-S-1G/6	OH-1G/6
2000 mg	6 ml	20 u	SI-S-2G/6	NH2-S-2G/6	CN-S-2G/6	OH-2G/6
2000 mg	15 ml	20 u	SI-S-2G/15	NH2-S-2G/15	CN-S-2G/15	OH-2G/15
2000 mg	25 ml	20 u	SI-S-2G/25	NH2-S-2G/25	CN-S-2G/25	OH-2G/25

### LRC columns / PE frits

100 mg	LRC 15	50 u	SI-S-100LRC	NH2-S-100LRC	CN-S-100LRC	OH-100LRC
200 mg	LRC 15	50 u	SI-S-200LRC	NH2-S-200LRC	CN-S-200LRC	OH-200LRC
500 mg	LRC 15	50 u	SI-S-500LRC	NH2-S-500LRC	CN-S-500LRC	OH-500LRC

### Std. columns / PTFE frits

50 mg	1 ml	50 u	SI-S-50/1T	NH2-S-50/1T	CN-S-50/1T	OH-50/1T
100 mg	1 ml	100 u	SI-S-100/1T	NH2-S-100/1T	CN-S-100/1T	OH-100/1T
100 mg	3 ml	50 u	SI-S-100/3T	NH2-S-100/3T	CN-S-100/3T	OH-100/3T
200 mg	3 ml	50 u	SI-S-200/3T	NH2-S-200/3T	CN-S-200/3T	OH-200/3T
500 mg	3 ml	50 u	SI-S-500/3T	NH2-S-500/3T	CN-S-500/3T	OH-500/3T
500 mg	6 ml	30 u	SI-S-500/6T	NH2-S-500/6T	CN-S-500/6T	OH-500/6T
1000 mg	6 ml	30 u	SI-S-1G/6T	NH2-S-1G/6T	CN-S-1G/6T	OH-1G/6T

### Glass columns / PTFE frits

200 mg	6 ml	30 u	SI-S-200/6G	NH2-S-200/6G	CN-S-200/6G	OH-200/6G
500 mg	6 ml	30 u	SI-S-500/6G	NH2-S-500/6G	CN-S-500/6G	OH-500/6G
1000 mg	6 ml	30 u	SI-S-1G/6G	NH2-S-1G/6G	CN-S-1G/6G	OH-1G/6G

## Hydrophilic 96 well plates - 2 ml

Weight	Qty	SI-S	NH2-S	CN-S	OH
30 mg	1 u	SI-30/WP20	NH2-30/WP20	CN-30/WP20	OH-30/WP20
50 mg	1 u	SI-50/WP20	NH2-50/WP20	CN-50/WP20	OH-50/WP20
60 mg	1 u	SI-60/WP20	NH2-60/WP20	CN-60/WP20	OH-60/WP20
100 mg	1 u	SI-100/WP20	NH2-100/WP20	CN-100/WP20	OH-100/WP20
150 mg	1 u	SI-150/WP20	NH2-150/WP20	CN-150/WP20	OH-150/WP20
200 mg	1 u	SI-200/WP20	NH2-200/WP20	CN-200/WP20	OH-200/WP20

## Ion exchange columns

Weight	Volume	Qty	SCX	WCX	DEAE	SAX
<b>Std. columns / PE frits</b>						
50 mg	1 ml	50 u	SCX-50/1	WCX-50/1	DEAE-50/1	SAX-50/1
100 mg	1 ml	100 u	SCX-100/1	WCX-100/1	DEAE-100/1	SAX-100/1
100 mg	3 ml	50 u	SCX-100/3	WCX-100/3	DEAE-100/3	SAX-100/3
200 mg	3 ml	50 u	SCX-200/3	WCX-200/3	DEAE-200/3	SAX-200/3
500 mg	3 ml	50 u	SCX-500/3	WCX-500/3	DEAE-500/3	SAX-500/3
500 mg	6 ml	30 u	SCX-500/6	WCX-500/6	DEAE-500/6	SAX-500/6
1000 mg	6 ml	30 u	SCX-1G/6	WCX-1G/6	DEAE-1G/6	SAX-1G/6
2000 mg	6 ml	20 u	SCX-2G/6	WCX-2G/6	DEAE-2G/6	SAX-2G/6
2000 mg	15 ml	20 u	SCX-2G/15	WCX-2G/15	DEAE-2G/15	SAX-2G/15
2000 mg	25 ml	20 u	SCX-2G/25	WCX-2G/25	DEAE-2G/25	SAX-2G/25

### LRC columns / PE frits

100 mg	LRC 15	50 u	SCX-100LRC	WCX-100LRC	DEAE-100LRC	SAX-100LRC
200 mg	LRC 15	50 u	SCX-200LRC	WCX-200LRC	DEAE-200LRC	SAX-200LRC
500 mg	LRC 15	50 u	SCX-500LRC	WCX-500LRC	DEAE-500LRC	SAX-500LRC

### Std. columns / PTFE frits

50 mg	1 ml	50 u	SCX-50/1T	WCX-50/1T	DEAE-50/1T	SAX-50/1T
100 mg	1 ml	100 u	SCX-100/1T	WCX-100/1T	DEAE-100/1T	SAX-100/1T
100 mg	3 ml	50 u	SCX-100/3T	WCX-100/3T	DEAE-100/3T	SAX-100/3T
200 mg	3 ml	50 u	SCX-200/3T	WCX-200/3T	DEAE-200/3T	SAX-200/3T
500 mg	3 ml	50 u	SCX-500/3T	WCX-500/3T	DEAE-500/3T	SAX-500/3T
500 mg	6 ml	30 u	SCX-500/6T	WCX-500/6T	DEAE-500/6T	SAX-500/6T
1000 mg	6 ml	30 u	SCX-1G/6T	WCX-1G/6T	DEAE-1G/6T	SAX-1G/6T

### Glass columns / PTFE frits

200 mg	6 ml	30 u	SCX-200/6G	WCX-200/6G	DEAE-200/6G	SAX-200/6G
500 mg	6 ml	30 u	SCX-500/6G	WCX-500/6G	DEAE-500/6G	SAX-500/6G
1000 mg	6 ml	30 u	SCX-1G/6G	WCX-1G/6G	DEAE-1G/6G	SAX-1G/6G

## Ion exchange :es - 2 ml

Weight	Qty	SCX	WCX	DEAE	SAX
30 mg	1 u	SCX-30/WP20	WCX-30/WP20	DEAE-30/WP20	SAX-30/WP20
50 mg	1 u	SCX-50/WP20	WCX-50/WP20	DEAE-50/WP20	SAX-50/WP20
60 mg	1 u	SCX-60/WP20	WCX-60/WP20	DEAE-60/WP20	SAX-60/WP20
100 mg	1 u	SCX-100/WP20	WCX-100/WP20	DEAE-100/WP20	SAX-100/WP20
150 mg	1 u	SCX-150/WP20	WCX-150/WP20	DEAE-150/WP20	SAX-150/WP20
200 mg	1 u	SCX-200/WP20	WCX-200/WP20	DEAE-200/WP20	SAX-200/WP20

# Sample Prep - Solid Phase Extraction

Upti-Clean™ Series-S - Silica based

## Mixed mode columns

Weight	Volume	Qty	MM1	MM2	MM3	MM4
<b>Std. columns / PE frits</b>						
50 mg	1 ml	50 u	MM1-50/1	MM2-50/1	MM3-50/1	MM4-50/1
100 mg	1 ml	100 u	MM1-100/1	MM2-100/1	MM3-100/1	MM4-100/1
100 mg	3 ml	50 u	MM1-100/3	MM2-100/3	MM3-100/3	MM4-100/3
200 mg	3 ml	50 u	MM1-200/3	MM2-200/3	MM3-200/3	MM4-200/3
500 mg	3 ml	50 u	MM1-500/3	MM2-500/3	MM3-500/3	MM4-500/3
500 mg	6 ml	30 u	MM1-500/6	MM2-500/6	MM3-500/6	MM4-500/6
1000 mg	6 ml	30 u	MM1-1G/6	MM2-1G/6	MM3-1G/6	MM4-1G/6
2000 mg	6 ml	20 u	MM1-2G/6	MM2-2G/6	MM3-2G/6	MM4-2G/6
2000 mg	15 ml	20 u	MM1-2G/15	MM2-2G/15	MM3-2G/15	MM4-2G/15
2000 mg	25 ml	20 u	MM1-2G/25	MM2-2G/25	MM3-2G/25	MM4-2G/25

<b>LRC columns / PE frits</b>						
100 mg	LRC 15	50 u	MM1-100LRC	MM2-100LRC	MM3-100LRC	MM4-100LRC
200 mg	LRC 15	50 u	MM1-200LRC	MM2-200LRC	MM3-200LRC	MM4-200LRC
500 mg	LRC 15	50 u	MM1-500LRC	MM2-500LRC	MM3-500LRC	MM4-500LRC

<b>Std. columns / PTFE frits</b>						
50 mg	1 ml	50 u	MM1-50/1T	MM2-50/1T	MM3-50/1T	MM4-50/1T
100 mg	1 ml	100 u	MM1-100/1T	MM2-100/1T	MM3-100/1T	MM4-100/1T
100 mg	3 ml	50 u	MM1-100/3T	MM2-100/3T	MM3-100/3T	MM4-100/3T
200 mg	3 ml	50 u	MM1-200/3T	MM2-200/3T	MM3-200/3T	MM4-200/3T
500 mg	3 ml	50 u	MM1-500/3T	MM2-500/3T	MM3-500/3T	MM4-500/3T
500 mg	6 ml	30 u	MM1-500/6T	MM2-500/6T	MM3-500/6T	MM4-500/6T
1000 mg	6 ml	30 u	MM1-1G/6T	MM2-1G/6T	MM3-1G/6T	MM4-1G/6T

<b>Glass columns - PTFE frits</b>						
200 mg	6 ml	30 u	MM1-200/6G	MM2-200/6G	MM3-200/6G	MM4-200/6G
500 mg	6 ml	30 u	MM1-500/6G	MM2-500/6G	MM3-500/6G	MM4-500/6G
1000 mg	6 ml	30 u	MM1-1G/6G	MM2-1G/6G	MM3-1G/6G	MM4-1G/6G



Std. columns  
PP straight tube + 20µm PE frits



LRC columns  
PP tube + 20µm PE frits



Std. columns  
PP straight tube + 20µm PTFE  
frits



Glass columns  
Glass tube + 20µm PTFE frits



## Mixed mode 96 well plates - 2 ml

Weight	Qty	MM1	MM2	MM3	MM4
30 mg	1 u	MM1-30/WP20	MM2-30/WP20	MM3-30/WP20	MM4-30/WP20
50 mg	1 u	MM1-50/WP20	MM2-50/WP20	MM3-50/WP20	MM4-50/WP20
60 mg	1 u	MM1-60/WP20	MM2-60/WP20	MM3-60/WP20	MM4-60/WP20
100 mg	1 u	MM1-100/WP20	MM2-100/WP20	MM3-100/WP20	MM4-100/WP20
150 mg	1 u	MM1-150/WP20	MM2-150/WP20	MM3-150/WP20	MM4-150/WP20
200 mg	1 u	MM1-200/WP20	MM2-200/WP20	MM3-200/WP20	MM4-200/WP20

# Sample Prep - Solid Phase Extraction

## Upti-Clean™ Series-S2F - Silica based

Upti-Clean™ Series-S2F columns are easy to use, efficient SPE cleaning devices. 140 µm particle size achieves faster flow and helps to achieve non-polar and mid-polar compound extraction from crude matrices without clogging.

- 3 selectivities available
- Polypropylene MG or Glass housing
- Std. PE or PTFE frits
- Luer outlet fitting

Applications : urines, plasmas, oils...

### Upti-Clean™ Series-S2F, columns

Weight	Volume	Qty	C18-S2F	C18U-S2F	C8-S2F
<b>Std. columns / PE frits</b>					
50 mg	1 ml	50 u	C18-S2F-50/1	C18U-S2F-50/1	C8-S2F-50/1
100 mg	1 ml	100 u	C18-S2F-100/1	C18U-S2F-100/1	C8-S2F-100/1
100 mg	3 ml	50 u	C18-S2F-100/3	C18U-S2F-100/3	C8-S2F-100/3
200 mg	3 ml	50 u	C18-S2F-200/3	C18U-S2F-200/3	C8-S2F-200/3
500 mg	3 ml	50 u	C18-S2F-500/3	C18U-S2F-500/3	C8-S2F-500/3
500 mg	6 ml	30 u	C18-S2F-500/6	C18U-S2F-500/6	C8-S2F-500/6
1000 mg	6 ml	30 u	C18-S2F-1G/6	C18U-S2F-1G/6	C8-S2F-1G/6
2000 mg	6 ml	20 u	C18-S2F-2G/6	C18U-S2F-2G/6	C8-S2F-2G/6
2000 mg	15 ml	20 u	C18-S2F-2G/15	C18U-S2F-2G/15	C8-S2F-2G/15
2000 mg	25 ml	20 u	C18-S2F-2G/25	C18U-S2F-2G/25	C8-S2F-2G/25



Std. columns  
PP straight tube + 20µm PE frits

<b>LRC columns / PE frits</b>					
100 mg	LRC 15	50 u	C18-S2F-100LRC	C18U-S2F-100LRC	C8-S2F-100LRC
200 mg	LRC 15	50 u	C18-S2F-200LRC	C18U-S2F-200LRC	C8-S2F-200LRC
500 mg	LRC 15	50 u	C18-S2F-500LRC	C18U-S2F-500LRC	C8-S2F-500LRC



LRC columns  
PP tube + 20µm PE frits

<b>Std. columns / PTFE frits</b>					
50 mg	1 ml	50 u	C18-S2F-50/1T	C18U-S2F-50/1T	C8-S2F-50/1T
100 mg	1 ml	100 u	C18-S2F-100/1T	C18U-S2F-100/1T	C8-S2F-100/1T
100 mg	3 ml	50 u	C18-S2F-100/3T	C18U-S2F-100/3T	C8-S2F-100/3T
200 mg	3 ml	50 u	C18-S2F-200/3T	C18U-S2F-200/3T	C8-S2F-200/3T
500 mg	3 ml	50 u	C18-S2F-500/3T	C18U-S2F-500/3T	C8-S2F-500/3T
500 mg	6 ml	30 u	C18-S2F-500/6T	C18U-S2F-500/6T	C8-S2F-500/6T
1000 mg	6 ml	30 u	C18-S2F-1G/6T	C18U-S2F-1G/6T	C8-S2F-1G/6T



Std. columns  
PP straight tube + 20µm PTFE frits

<b>Glass columns / PTFE frits</b>					
200 mg	6 ml	30 u	C18-S2F-200/6G	C18U-S2F-200/6G	C8-S2F-200/6G
500 mg	6 ml	30 u	C18-S2F-500/6G	C18U-S2F-500/6G	C8-S2F-500/6G
1000 mg	6 ml	30 u	C18-S2F-1G/6G	C18U-S2F-1G/6G	C8-S2F-1G/6G



Glass columns  
Glass tube + 20µm PTFE frits

# Sample Prep - Solid Phase Extraction

## Upti-Clean™ Specialty

The Upti-Clean™ Specialty Series columns complete the Series S range. The Specialty series are designed for specific applications highlighted below.

- Polypropylene MG housing
- Std. PE or PTFE frits
- Luer outlet fitting

### Florisil® columns

Florisil is a magnesia-loaded silica gel used to trap polar impurities within non-polar matrices. Florisil can also be used as an alternative to traditional virgin silica when working with viscous solvents.

Florisil PR is a grade for Pesticide applications

Weight	Volume	Qty	FL	FLPR
<b>Std. columns / PE frits</b>				
500 mg	3 ml	50 u	FL-500/3	FLPR-500/3
500 mg	6 ml	30 u	FL-500/6	FLPR-500/6
1000 mg	6 ml	30 u	FL-1G/6	FLPR-1G/6
2000 mg	6 ml	20 u	FL-2G/6	FLPR-2G/6
2000 mg	15 ml	20 u	FL-2G/15	FLPR-2G/15
2000 mg	25 ml	20 u	FL-2G/25	FLPR-2G/25

<b>Std. columns / PTFE frits</b>				
500 mg	3 ml	50 u	FL-500/3T	FLPR-500/3T
500 mg	6 ml	30 u	FL-500/6T	FLPR-500/6T
1000 mg	6 ml	30 u	FL-1G/6T	FLPR-1G/6T

### Alumina columns

The aluminium atom lacks two electrons within its center that are responsible for ion pair interaction.

The acidic treatment of Alumina favors the retention of cationic species whilst a basic treatment of Alumina leads to the retention of anionic species.

Neutral Alumina is suitable to clean non ionisable compounds with polar function.

Applications : Environmental (dioxines, PCB,...)

Weight	Volume	Qty	ALA	ALB	ALN
<b>Std. columns / PE frits</b>					
500 mg	3 ml	50 u	ALA-500/3	ALB-500/3	ALN-500/3
500 mg	6 ml	30 u	ALA-500/6	ALB-500/6	ALN-500/6
1000 mg	6 ml	30 u	ALA-1G/6	ALB-1G/6	ALN-1G/6
2000 mg	6 ml	20 u	ALA-2G/6	ALB-2G/6	ALN-2G/6
2000 mg	15 ml	20 u	ALA-2G/15	ALB-2G/15	ALN-2G/15
2000 mg	25 ml	20 u	ALA-2G/25	ALB-2G/25	ALN-2G/25

<b>Std. columns / PTFE frits</b>					
500 mg	3 ml	50 u	ALA-500/3T	ALB-500/3T	ALN-500/3T
500 mg	6 ml	30 u	ALA-500/6T	ALB-500/6T	ALN-500/6T
1000 mg	6 ml	30 u	ALA-1G/6T	ALB-1G/6T	ALN-1G/6T

### Upti-Clean™ DRY columns

These columns remove trace water found in organic solvents or non aqueous samples. They contain a treated, high purity sodium sulfate group within the structure.

Description	Weight	Volume	Qty	P/N
Upti-Clean DRY	300 mg	1 ml	50 u	DRY-300/1
Upti-Clean DRY	500 mg	3 ml	50 u	DRY-500/3
Upti-Clean DRY	2000 mg	3 ml	50 u	DRY-2G/3



# Sample Prep - Solid Phase Extraction

Upti-Clean™ Specialty

## Amberlite™ columns

Amberlite® is the first generation of polymer resins. They are used for fast separation of a variety of compounds from biological fluids. Amberlite® suffer from weak selectivity.



Weight	Volume	Qty	XAD-2	XAD-4	XAD-7	XAD-16
<b>Std. columns / PE frits</b>						
100 mg	1 ml	100 u	XAD2-100/1	XAD4-100/1	XAD7-100/1	XAD16-100/1
200 mg	3 ml	50 u	XAD2-200/3	XAD4-200/3	XAD7-200/3	XAD16-200/3
500 mg	3 ml	50 u	XAD2-500/3	XAD4-500/3	XAD7-500/3	XAD16-500/3
500 mg	6 ml	30 u	XAD2-500/6	XAD4-500/6	XAD7-500/6	XAD16-500/6
1000 mg	6 ml	30 u	XAD2-1G/6	XAD4-1G/6	XAD7-1G/6	XAD16-1G/6
1000 mg	12 ml	20 u	XAD2-1G/12	XAD4-1G/12	XAD7-1G/12	XAD16-1G/12
2000 mg	6 ml	30 u	XAD2-2G/6	XAD4-2G/6	XAD7-2G/6	XAD16-2G/6
2000 mg	12 ml	20 u	XAD2-2G/12	XAD4-2G/12	XAD7-2G/12	XAD16-2G/12
5000 mg	35 ml	20 u	XAD2-5G/35	XAD4-5G/35	XAD7-5G/35	XAD16-5G/35
10000 mg	60 ml	12 u	XAD2-10G/60	XAD4-10G/60	XAD7-10G/60	XAD16-10G/60
20000 mg	60 ml	12 u	XAD2-20G/60	XAD4-20G/60	XAD7-20G/60	XAD16-20G/60

## Polyamide columns

Amide functionality upon a Nylon support. Polyamide columns are typically used for aromatic and natural compound extraction such as PAH or flavanoids

Weight	Volume	Qty	P6
<b>Std. columns / PE frits</b>			
100 mg	1 ml	100 u	P6-100/1
100 mg	3 ml	50 u	P6-100/3
200 mg	3 ml	50 u	P6-200/3
500 mg	3 ml	50 u	P6-500/3
500 mg	6 ml	30 u	P6-500/6
1000 mg	6 ml	30 u	P6-1G/6
2000 mg	6 ml	20 u	P6-2G/6
2000 mg	15 ml	20 u	P6-2G/15
2000 mg	25 ml	20 u	P6-2G/25

## β-Cyclodextrine columns

β-Cyclodextrine is useful for cleaning small molecular mass oligosaccharides.

Weight	Volume	Qty	BCD
<b>Std. columns / PE frits</b>			
100 mg	1 ml	100 u	B-CD-100/1
100 mg	3 ml	50 u	B-CD-100/3
200 mg	3 ml	50 u	B-CD-200/3
500 mg	3 ml	50 u	B-CD-500/3
500 mg	6 ml	30 u	B-CD-500/6
1000 mg	6 ml	30 u	B-CD-1G/6



# Sample Prep - Solid Phase Extraction

Upti-Clean™ Specialty

## Graphitized carbon columns

The electronic interactions of graphitized carbon assure the retention of a large variety of compounds. This sorbent is especially suited to polar compound extraction within an aqueous matrix.

Weight	Volume	Qty	CG-ENV	CG
<b>Std. columns / PE frits</b>				
50 mg	1 ml	100 u	CG-ENV-50/1	CG-50/1
100 mg	1 ml	100 u	CG-ENV-100/1	CG-100/1
150 mg	3 ml	50 u		CG-150/3
200 mg	3 ml	50 u	CG-ENV-200/3	CG-200/3
250 mg	3 ml	50 u		CG-250/3
250 mg	6 ml	30 u		CG-250/6
500 mg	3 ml	50 u	CG-ENV-500/3	CG-500/3
500 mg	6 ml	30 u	CG-ENV-500/6	CG-500/6
1000 mg	6 ml	-	CG-ENV-1G/6	CG-1G/6
1000 mg	15 ml	20 u	CG-ENV-1G/15	CG-1G/1

## Florisil & Alumina 96 well plate

Weight	Qty	FL	FLPR	ALA	ALB	ALN
50 mg	1 u	FL-50/WP20	FLPR-50/WP20	ALA-50/WP20	ALB-50/WP20	ALN-50/WP20
100 mg	1 u	FL-100/WP20	FLPR-100/WP20	ALA-100/WP20	ALB-100/WP20	ALN-100/WP20
150 mg	1 u	FL-150/WP20	FLPR-150/WP20	ALA-150/WP20	ALB-150/WP20	ALN-150/WP20



# Sample Prep - Solid Phase Extraction

## Upti-Clean™ 48 well plates

The Upti-Clean® 48 well-plate is a flash purification device for chemists working in Drug Discovery Laboratories and Medicinal Chemistry Departments. The 48 well format maximises sample purification throughput and is an essential tool for library purification. The large volume wells (5 & 7 ml) allow high sample volume load.

Interchim silica and bonded silica are rigid supports that do not shrink or swell with solvents. The pH stability of bonded silica is limited, typically to within the range of 2 to 7.5. This is dependant on the chemistry and Interchim offers 8 different selectivities in this format. Our sorbents take advantage of our ultra pure spherical silica (Upti-prep™), and this achieves greater purification, and establishes optimized sample recovery.

Upti-Clean™ 48 well plates utilise Interchim's state of the art automated weight machine. This packing technology guarantees unprecedented accuracy (+/- 1% mass per column) compared to the commonly utilised competitor volumetric systems. Upti-Clean™ therefore raise the standards & client expectations for sample recovery /reproducibility.

Upti-Clean™ 48 well plates are QC tested in-house to guarantee tracability. Products are supplied with an individual certificate detailing Mfg number, Sorbent batch number & specifications. An additional certificate is supplied that states real sorbent mass in each individual well.

Upti-Clean™ 48 well plates are packed in a PEHD/Al bag for long term integral storage. This prevents potential damage from UV and moisture.

Upti-Clean® 48 well plates are manufactured from polypropylene MG, conform to a standard footprint (127,76 mm x 85,47 mm) and are compatible with existing systems on the market.

Sorbent	Mass / well 100 mg / 5 ml	Mass / well 250 mg / 5 ml	Mass / well 500 mg / 7 ml	Mass / well 1000 mg / 7 ml
Silica	SI-100/WP50	SI-250/WP50	SI-500/WP70	SI-1G/WP70
C18	C18-100/WP50	C18-250/WP50	C18-500/WP70	C18-1G/WP70
Penyl	PH-100/WP50	PH-250/WP50	PH-500/WP70	PH-1G/WP70
Amino	NH2-100/WP50	NH2-250/WP50	NH2-500/WP70	NH2-1G/WP70
Strong Anion Exchanger	SAX-100/WP50	SAX-250/WP50	SAX-500/WP70	SAX-1G/WP70
Strong Cation Exchanger	SCX-100/WP50	SCX-250/WP50	SCX-500/WP70	SCX-1G/WP70
Alumina basic	ALB-100/WP50	ALB-250/WP50	ALB-500/WP70	ALB-1G/WP70
Florisil	FL-100/WP50	FL-250/WP50	FL-500/WP70	FL-1G/WP70



# Sample Prep - Solid Phase Extraction

## Atoll™ - Polymer based

Atoll™ is a comprehensive product range based upon porous polymers consisting of ultra pure spherical particles with strictly controlled particle & pore size distribution. Atoll™ is the latest generation of polymer that does not shrink or swell with a range of standard solvents.

Atoll™ can be used through a wide range of pH (1-14) and is compatible with all solvents and existing samples. It has a very high surface area, Atoll™ XC exhibiting the highest surface area currently available on market (1500 sq M /g) approximately three times greater loading capacity than traditional silica sorbents. The higher capacity of Atoll™ allows for smaller bed masses.

Interchim bonding technology ensures greater batch to batch reproducibility for the Atoll™ family of products, therefore there is no longer a need for batch reservations. Our products subsequently achieve superior recovery rates relative to traditional media, exhibiting excellent reproducibility & consistency.

The Atoll™ family of products : Atoll Xtrem™ Capacity (XC) - Atoll™ Xtrem Capacity Wide Pore (XWP) - Atoll™ Hydrophilic (ATH) - Atoll™ Environment (AEV) - Atoll™ Lipophilic (ATL)

The Atoll™ range achieves accuracy of +/- 1% in mass per column thanks to a strictly controlled packing process. Atoll™ products are thoroughly quality control tested in-house to guarantee tracability. Products are supplied with an individual certificate detailing the specific production number and sorbent batch.

Atoll™ columns are packed in a PEHD/Al bag for long term integral storage. This prevents potential damage from UV and moisture.

Type	Product code	Sorbent	Particle size	Porosity	Surface area
Atoll Xtrem Capacity	30XC	PSDVB	30 µm	n.a.	1500 m <sup>2</sup> /g
	XC	PSDVB	70 µm	n.a.	1500 m <sup>2</sup> /g
Atoll Xtrem Capacity Wide Pore	XWP	PSDVB	90 µm	n.a.	1200 m <sup>2</sup> /g
Atoll Hydrophile	30ATH	n.a.	30 µm	70 Å	800 m <sup>2</sup> /g
	ATH	n.a.	75 µm	70 Å	800 m <sup>2</sup> /g
Atoll Environment	AEV	PSHEMA	75 µm	70 Å	800 m <sup>2</sup> /g
Atoll Lipophile	ATL	PSDVB	100 µm	70 Å	800 m <sup>2</sup> /g

### Typical Atoll procedure

#### Conditioning

100 % MeOH : 1 ml  
100 % H<sub>2</sub>O : 1 ml

#### Sample loading

0,2 - 1000 ml

#### Washing

H<sub>2</sub>O/MeOH (95/5) : 1 ml

#### Drying

2 min

#### Elution

100 % MeOH : 1 ml

#### Analysis

# Sample Prep - Solid Phase Extraction

## Atoll™ XC columns - Polymer based

### Atoll™ Xtrem Capacity (XC)

Atoll™ XC exhibits the highest surface area currently available on the market (1500 sq M /g), approximately three times greater loading capacity than traditional silica sorbents. Atoll XC improves, 2 to 3 fold, the retention of a compound of interest, compared to traditional silica.

Atoll XC polymer is a versatile sorbent for the extraction and pre-concentration of non-polar, polar, acidic, basic or neutral compounds.

- Std. 70 µm and 30 µm
- pH stability : 1 to 14
- Pharmaceutical applications : drugs & metabolites in biological fluids
- Environmental applications : PAH, PCB, carbamates, phenyl-ureas, acrylamide, glyphosate

Weight	Volume	Qty	Atoll XC	Atoll 30 µm XC
<b>Std. columns / PE frits</b>				
30 mg	1 ml	50 u	XC-30/1	30XC-30/1
50 mg	1 ml	50 u	XC-50/1	30XC-50/1
60 mg	1 ml	50 u	XC-60/1	30XC-60/1
75 mg	1 ml	50 u	XC-75/1	30XC-75/1
100 mg	1 ml	50 u	XC-100/1	30XC-100/1
100 mg	3 ml	50 u	XC-100/3	On request
150 mg	3 ml	50 u	XC-150/3	On request
200 mg	3 ml	50 u	XC-200/3	On request
250 mg	3 ml	50 u	XC-250/3	On request
500 mg	6 ml	30 u	XC-500/6	On request

<b>LRC columns / PE frits</b>				
50 mg	LRC 15	50 u	XC-50LRC	30XC-50LRC
75 mg	LRC 15	50 u	XC-75LRC	30XC-75LRC
100 mg	LRC 15	50 u	XC-100LRC	30XC-100LRC
200 mg	LRC 15	50 u	XC-200LRC	On request

<b>Std. columns / PTFE frits</b>				
30 mg	1 ml	50 u	XC-30/1T	30XC-30/1T
50 mg	1 ml	50 u	XC-50/1T	30XC-50/1T
60 mg	1 ml	50 u	XC-60/1T	30XC-60/1T
75 mg	1 ml	50 u	XC-75/1T	30XC-75/1T
100 mg	1 ml	50 u	XC-100/1T	30XC-100/1T
100 mg	3 ml	50 u	XC-100/3T	On request
150 mg	3 ml	50 u	XC-150/3T	On request
200 mg	3 ml	50 u	XC-200/3T	On request
250 mg	3 ml	50 u	XC-250/3T	On request

<b>Glass columns / PTFE frits</b>				
200 mg	6 ml	30 u	XC-200/6G	On request
500 mg	6 ml	30 u	XC-500/6G	On request



Std. columns  
PP straight tube + 20µm PE frits



LRC columns  
PP tube + 20µm PE frits



Std. columns  
PP straight tube + 20µm PTFE frits



Glass columns  
Glass tube + 20µm PTFE frits

# Sample Prep - Solid Phase Extraction

## Atoll™ ATH columns - Polymer based

### Atoll™ Hydrophilic (ATH)

Atoll™ Hydrophilic contains proprietary chemical modifications that achieve unsurpassed compound extraction within 100% aqueous matrices. This technology facilitates mixed hydrophilic /hydrophobic interactions.

- Std. 70 µm
- pH stability : 2 to 12
- Pharmaceutical applications : drugs & metabolites in biological fluids
- Environmental applications : PAH, PCB, carbamates, phenyl-ureas, acrylamide, glyphosate

Weight	Volume	Qty	Atoll ATH	Atoll 30 µm ATH
<b>Std. columns / PE frits</b>				
30 mg	1 ml	50 u	ATH-30/1	30ATH-30/1
50 mg	1 ml	50 u	ATH-50/1	30ATH-50/1
60 mg	1 ml	50 u	ATH-60/1	30ATH-60/1
75 mg	1 ml	50 u	ATH-75/1	30ATH-75/1
100 mg	1 ml	50 u	ATH-100/1	30ATH-100/1
100 mg	3 ml	50 u	ATH-100/3	On request
150 mg	3 ml	50 u	ATH-150/3	On request
200 mg	3 ml	50 u	ATH-200/3	On request
250 mg	3 ml	50 u	ATH-250/3	On request
500 mg	6 ml	30 u	ATH-500/6	On request
1000 mg	6 ml	30 u	ATH-1G/6	On request

<b>LRC columns / PE frits</b>				
50 mg	LRC 15	50 u	ATH-50LRC	30ATH-50LRC
75 mg	LRC 15	50 u	ATH-75LRC	30ATH-75LRC
100 mg	LRC 15	50 u	ATH-100LRC	30ATH-100LRC
200 mg	LRC 15	50 u	ATH-200LRC	On request

<b>Std. columns / PTFE frits</b>				
30 mg	1 ml	50 u	ATH-30/1T	30ATH-30/1T
50 mg	1 ml	50 u	ATH-50/1T	30ATH-50/1T
60 mg	1 ml	50 u	ATH-60/1T	30ATH-60/1T
75 mg	1 ml	50 u	ATH-75/1T	30ATH-75/1T
100 mg	1 ml	50 u	ATH-100/1T	30ATH-100/1T
100 mg	3 ml	50 u	ATH-100/3T	On request
150 mg	3 ml	50 u	ATH-150/3T	On request
200 mg	3 ml	50 u	ATH-200/3T	On request
250 mg	3 ml	50 u	ATH-250/3T	On request
500 mg	6 ml	30 u	ATH-500/6T	On request
1000 mg	6 ml	30 u	ATH-1G/6T	On request

<b>Glass columns / PTFE frits</b>				
200 mg	6 ml	30 u	ATH-200/6G	On request
500 mg	6 ml	30 u	ATH-500/6G	On request

# Sample Prep - Solid Phase Extraction

## Atoll™ AEV columns - Polymer based

### Atoll™ Environment (AEV)

Atoll™ AEV is a hydroxyethylmethacrylate - polystyrene polymer that supports the extraction of compounds weakly retained on C18. Atoll™ AEV exhibits high extraction flow of non-polar and mid-polar compounds that allows to clean and collect multi-residue samples in a single fraction prior to analysis. These characteristics identify this support as suitable to perform environmental extractions.

- 75 µm
- pH stability : 1 to 12
- Environmental applications : PAH, PCB, carbamates, phenyl-ureas

Weight	Volume	Qty	AEV
<b>Std. columns / PE frits</b>			
30 mg	1 ml	50 u	AEV-30/1
50 mg	1 ml	50 u	AEV-50/1
60 mg	1 ml	50 u	AEV-60/1
75 mg	1 ml	50 u	AEV-75/1
100 mg	1 ml	50 u	AEV-100/1
100 mg	3 ml	50 u	AEV-100/3
150 mg	3 ml	50 u	AEV-150/3
200 mg	3 ml	50 u	AEV-200/3
250 mg	3 ml	50 u	AEV-250/3
500 mg	6 ml	30 u	AEV-500/6
1000 mg	6 ml	30 u	AEV-1G/6

<b>LRC columns / PE frits</b>			
50 mg	LRC 15	50 u	AEV-50LRC
75 mg	LRC 15	50 u	AEV-75LRC
100 mg	LRC 15	50 u	AEV-100LRC
200 mg	LRC 15	50 u	AEV-200LRC

<b>Std. columns / PTFE frits</b>			
30 mg	1 ml	50 u	AEV-30/1T
50 mg	1 ml	50 u	AEV-50/1T
60 mg	1 ml	50 u	AEV-60/1T
75 mg	1 ml	50 u	AEV-75/1T
100 mg	1 ml	50 u	AEV-100/1T
100 mg	3 ml	50 u	AEV-100/3T
150 mg	3 ml	50 u	AEV-150/3T
200 mg	3 ml	50 u	AEV-200/3T
250 mg	3 ml	50 u	AEV-250/3T
500 mg	6 ml	30 u	AEV-500/6T
1000 mg	6 ml	30 u	AEV-1G/6T

<b>Glass columns / PTFE frits</b>			
200 mg	6 ml	30 u	AEV-200/6G
500 mg	6 ml	30 u	AEV-500/6G



Std. columns  
PP straight tube + 20µm PE frits



LRC columns  
PP tube + 20µm PE frits



Std. columns  
PP straight tube + 20µm PTFE frits



Glass columns  
Glass tube + 20µm PTFE frits



# Sample Prep - Solid Phase Extraction

## Atoll™ ATL columns - Polymer based

### Atoll™ Lipophilic (ATL)

Atoll™ ATL is a 100 µm PSDVB optimized for the extraction of hydrophobic compounds in a large variety of matrices and is an excellent alternative to high flow silica for crude samples.

- pH stability : 1 to 14
- Applications : oils, muds, reactional environment after synthesis, ...

Weight	Volume	Qty	Atoll ATL
<b>Std. columns / PE frits</b>			
30 mg	1 ml	50 u	ATL-30/1
50 mg	1 ml	50 u	ATL-50/1
60 mg	1 ml	50 u	ATL-60/1
75 mg	1 ml	50 u	ATL-75/1
100 mg	1 ml	50 u	ATL-100/1
100 mg	3 ml	50 u	ATL-100/3
150 mg	3 ml	50 u	ATL-150/3
200 mg	3 ml	50 u	ATL-200/3
250 mg	3 ml	50 u	ATL-250/3
500 mg	6 ml	30 u	ATL-500/6
1000 mg	6 ml	30 u	ATL-1G/6

<b>LRC columns / PE frits</b>			
50 mg	LRC 15	50 u	ATL-50LRC
75 mg	LRC 15	50 u	ATL-75LRC
100 mg	LRC 15	50 u	ATL-100LRC
200 mg	LRC 15	50 u	ATL-200LRC

<b>Std. columns / PTFE frits</b>			
30 mg	1 ml	50 u	ATL-30/1T
50 mg	1 ml	50 u	ATL-50/1T
60 mg	1 ml	50 u	ATL-60/1T
75 mg	1 ml	50 u	ATL-75/1T
100 mg	1 ml	50 u	ATL-100/1T
100 mg	3 ml	50 u	ATL-100/3T
150 mg	3 ml	50 u	ATL-150/3T
200 mg	3 ml	50 u	ATL-200/3T
250 mg	3 ml	50 u	ATL-250/3T
500 mg	6 ml	30 u	ATL-500/6T
1000 mg	6 ml	30 u	ATL-1G/6T

<b>Glass columns / PTFE frits</b>			
200 mg	6 ml	30 u	ATL-200/6G
500 mg	6 ml	30 u	ATL-500/6G

# Sample Prep - Solid Phase Extraction

## Atoll™ XWP columns - Polymer based

### Atoll™ Xtrem Capacity Wide Pore (XWP)

Atoll™ XWP is a highly cross linked PSDVB exhibiting a very high surface area (1200 m<sup>2</sup>/g) and a 90 µm particle size. The apparent porosity of Atoll™ XWP makes it particularly suitable for protein & peptide cleaning from biologic fluids within a cut-off limit of ~ 500 KD. Its characteristics guarantee an excellent throughput and achieve efficient extractions without clogging.

- pH stability : 1 to 14
- Applications : proteins and peptides in polar matrices, molecules with high molecular weight in all solvent types.

Weight	Volume	Qty	Atoll XWP
<b>Std. columns / PE frits</b>			
30 mg	1 ml	50 u	XWP-30/1
50 mg	1 ml	50 u	XWP-50/1
60 mg	1 ml	50 u	XWP-60/1
75 mg	1 ml	50 u	XWP-75/1
100 mg	1 ml	50 u	XWP-100/1
100 mg	3 ml	50 u	XWP-100/3
150 mg	3 ml	50 u	XWP-150/3
200 mg	3 ml	50 u	XWP-200/3
250 mg	3 ml	50 u	XWP-250/3
500 mg	6 ml	30 u	XWP-500/6

<b>LRC columns / PE frits</b>			
50 mg	LRC 15	50 u	XWP-50LRC
75 mg	LRC 15	50 u	XWP-75LRC
100 mg	LRC 15	50 u	XWP-100LRC
200 mg	LRC 15	50 u	XWP-200LRC

<b>Std. columns / PTFE frits</b>			
30 mg	1 ml	50 u	XWP-30/1T
50 mg	1 ml	50 u	XWP-50/1T
60 mg	1 ml	50 u	XWP-60/1T
75 mg	1 ml	50 u	XWP-75/1T
100 mg	1 ml	50 u	XWP-100/1T
100 mg	3 ml	50 u	XWP-100/3T
150 mg	3 ml	50 u	XWP-150/3T
200 mg	3 ml	50 u	XWP-200/3T
250 mg	3 ml	50 u	XWP-250/3T

<b>Glass columns / PTFE frits</b>			
200 mg	6 ml	30 u	XWP-200/6G
500 mg	6 ml	30 u	XWP-500/6G



Std. columns  
PP straight tube + 20µm PE frits



LRC columns  
PP tube + 20µm PE frits



Std. columns  
PP straight tube + 20µm PTFE frits



Glass columns  
Glass tube + 20µm PTFE frits



# Sample Prep - Solid Phase Extraction

## Atoll™ Bio-Ion Exchange columns - Polymer based

### Atoll™ Bio-Ion Exchange

Atoll™ Bio-Ion Exchange bonded chemistry exhibits Ion exchange interaction and offers distinct selectivity to support the extraction and purification of proteins in aqueous matrices, such as biological fluids.

Adsorbent vol.	Volume	Qty	SPCE	WPCE	WPAE	SPAE
0,5 ml	3 ml	50 u	SPCE-0.5/3	WPCE-0.5/3	WPAE-0.5/3	SPAE-0.5/3
1 ml	3 ml	50 u	SPCE-1X/3	WPCE-1X/3	WPAE-1X/3	SPAE-1X/3
2 ml	6 ml	50 u	SPCE-2X/6	WPCE-2X/6	WPAE-2X/6	SPAE-2X/6

# Sample Prep - Solid Phase Extraction

## Atoll™ 96 well plates - Polymer based

Atoll™ is a comprehensive product range based upon porous polymers consisting of ultra pure spherical particles with strictly controlled particle & pore size distribution. Atoll™ is the latest generation of polymer that does not shrink or swell within a range of standard solvents. The Atoll™ family of products : Atoll Xtrem™ Capacity (XC) - Atoll™ Xtrem Capacity Wide Pore (XWP) - Atoll™ Hydrophilic (ATH) - Atoll™ Environment (AEV) - Atoll™ Lipophilic (ATL)  
See page A.43 to A.49 for relevant information about Atoll™ characteristics.

Automated cleanup procedures are now an integral part of the modern laboratory.

Atoll™ is available in a 96 well format with a 2 ml volume per well. This high quality unit provides rapid throughput of sample, within an automated process, whilst establishing excellent well-to-well consistency.

Atoll™ 96 well plates utilise Interchim's state of the art automated weight machine.

This packing technology guarantees unprecedented accuracy (+/- 1% mass per column) compared to the commonly utilised competitor volumetric systems. Atoll™ therefore raises the standards & client expectations for sample recovery /reproducibility.

Atoll™ 96 well plates are QC tested in-house to guarantee tracibility. Products are supplied with an individual certificate detailing Mfg number, Sorbent batch number and specifications. An additional certificate is supplied that states real sorbent mass in each individual well.

Atoll™ 96 well plates are packed in a PEHD/Al bag for long term integral storage. This prevents potential damage from UV and moisture.

Atoll™ 96 well plates are manufactured from polypropylene, conform to a standard footprint (127,76 mm x 85,47 mm x 19,74 mm) and are compatible with existing systems on the market.

Weight	Qty	Atoll XC	Atoll 30XC	Atoll ATH	Atoll 30ATH
30 mg	1 u	XC-30/WP20	30XC-30/WP20	ATH-30/WP20	30ATH-30/WP20
50 mg	1 u	XC-50/WP20	30XC-50/WP20	ATH-50/WP20	30ATH-50/WP20
60 mg	1 u	XC-60/WP20	30XC-60/WP20	ATH-60/WP20	30ATH-60/WP20
75 mg	1 u	XC-75/WP20	30XC-75/WP20	ATH-75/WP20	30ATH-75/WP20
100 mg	1 u	XC-100/WP20	30XC-100/WP20	ATH-100/WP20	30ATH-100/WP20
150 mg	1 u	XC-150/WP20	n.a.	n.a.	n.a.
200 mg	1 u	XC-200/WP20	n.a.	n.a.	n.a.

Weight	Qty	Atoll AEV	Atoll ATL	Atoll XWP
30 mg	1 u	AEV-30/WP20	ATL-30/WP20	XWP-30/WP20
50 mg	1 u	AEV-50/WP20	ATL-50/WP20	XWP-50/WP20
60 mg	1 u	AEV-60/WP20	ATL-60/WP20	XWP-60/WP20
75 mg	1 u	AEV-75/WP20	ATL-75/WP20	XWP-75/WP20
100 mg	1 u	AEV-100/WP20	ATL-100/WP20	XWP-100/WP20





# Sample Prep - Solid Phase Extraction

## Polymer/ Silica kits - method development & optimization

The rapid development of robust, reproducible & repeatable Solid Phase Extraction methods is a primary challenge for today's analyst. To assist with the achievement of this goal Interchim has developed a complete range of kits. The analyst can combine, within the kits, sorbent nature (silicas or polymers), mass and container volume. The combination of these parameters is critical to obtain optimum extraction yield as well as achieving good reproducibility.

Method development kits provide quick and efficient selection of the most suitable column. The analyst should know a basic background of information regarding the nature of the compound of interest, the impurities and the matrix. Upon identification of these parameters the optimization kit should optimize extraction yield and elution volume.

Interchim is happy to assist you in building an appropriate kit for method development and optimization.

### Pharmaceutical 1

Qty	P/N SPE-D9 Kit 100 mg / 1 mL	P/N SPE-D10 Kit 200 mg / 3 mL
6 columns C18-S	C18-S-100/1	C18-S-200/3
6 columns C8-S	C8-S-100/1	C8-S-200/3
6 columns CN-S	CN-S-100/1	CN-S-200/3
6 columns MM1	MM1-100/1	MM1-200/3
6 columns RPAQ	RPAQ-100/1	RPAQ-200/3
6 columns SCX	SCX-100/1	SCX-200/3
6 columns SAX	SAX-100/1	SAX-200/3

### Pharmaceutical 2

Qty	P/N SPE-D11 Kit 100 mg / 1 mL	P/N SPE-D12 Kit 200 mg / 3 mL
6 columns C18-S	C18-S-100/1	C18-S-200/3
6 columns C8-S	C8-S-100/1	C8-S-200/3
6 columns RPAQ	RPAQ-100/1	RPAQ-200/3
6 columns MM1	MM1-100/1	MM1-200/3
6 columns SCX	SCX-100/1	SCX-200/3
6 columns XC	XC-100/1	XC-200/3
6 columns AEV	AEV-100/1	AEV-200/3

# Sample Prep - Solid Phase Extraction

## Polymer/ Silica kits - method development & optimization

### Non-polar kit

P/N SPE-D54 100 mg / 1 mL	P/N SPE-D55 200 mg / 3 mL	P/N SPE-D56 500 mg / 6 mL
6 columns C18-S	6 columns C18-S	6 columns C18-S
6 columns C18U-S	6 columns C18U-S	6 columns C18U-S
6 columns RPAQ	6 columns RPAQ	6 columns RPAQ
6 columns C8-S	6 columns C8-S	6 columns C8-S
6 columns C2	6 columns C2	6 columns C2
6 columns PH-S	6 columns PH-S	6 columns PH-S
6 columns CH	6 columns CH	6 columns CH
6 columns ATH	6 columns ATH	6 columns ATH
6 columns XC	6 columns XC	6 columns XC

### Ion exchange & mixed mode kit

P/N SPE-D57 200 mg / 3 mL	P/N SPE-D58 500 mg / 6 mL
6 columns SCX	6 columns SCX
6 columns SAX	6 columns SAX
6 columns WCX	6 columns WCX
6 columns DEAE	6 columns DEAE
6 columns NH2	6 columns NH2
6 columns MM1	6 columns MM1
6 columns MM2	6 columns MM2
6 columns MM3	6 columns MM3
6 columns MM4	6 columns MM4

### Polymer kit

P/N SPE-D15 100 mg / 1 mL	P/N SPE-D59 50 mg / 1 mL
6 columns XC	6 columns XC
6 columns XWP	6 columns XWP
6 columns AEV	6 columns AEV
6 columns ATH	6 columns ATH
6 columns ATL	6 columns ATL
	6 columns 30ATH
	6 columns 30XC

### Optimization kit

QtyP/N SPE-DO60	P/N SPE-DO61
6 columns C18-S, 30mg/1ml	6 columns XC, 20mg/1ml
6 columns C18-S, 40mg/1ml	6 columns XC, 30mg/1ml
6 columns C18-S, 50mg/1ml	6 columns XC, 40mg/1ml
6 columns C18-S, 60mg/1ml	6 columns XC, 50mg/1ml
6 columns C18-S, 70mg/1ml	6 columns XC, 60mg/1ml
6 columns C18-S, 80mg/1ml	6 columns XC, 70mg/1ml
6 columns C18-S, 90mg/1ml	6 columns XC, 80mg/1ml
6 columns C18-S, 100mg/1ml	



# Sample Prep - Solid Phase Extraction

## Polymer/ Silica - 96 well plate custom kits


The 96 well format allows for a faster, easier process of SPE method development. Kit selection should be based upon selection of appropriate sorbent and the identification of appropriate sorbent mass per well (20 to 200 mg relative to sorbent density).

Description	P/N
96 well plate development kit	96WP-DO63191

Interchims technical center can assist with identification of your specific requirements. Plate proposals are provided within 2 working days of receiving your enquiry.

Your specified plate is supplied with a weight certificate that specifies exact sorbent mass in every well to a tolerance of +/- 1% .

E.g. 96 well-plate pharmaceutical method development



### Upti-clean™ silica media

Hydrophobic & Hydrophilic	Ion Exchange & Mixed
C18-S end capped	NH <sub>2</sub> weak anion exchanger
C18U-S non-end capped	DEAE medium anion exchanger
LCC end capped	SAX strong anion exchanger
RP.AQ hydrophilic	MM4 hydrophobic /weak anion exchanger
C8-S end capped	MM2 hydrophobic /strong anion exchanger
C8U-S non-end capped	WCX weak cation exchanger
CN-S Cyano	SCX strong cation exchanger
PH-S Phenyl	MM3 hydrophobic /weak cation exchanger
SI-S Silica	MM1 hydrophobic /strong cation exchanger

### Atoll® Polymer

Hydrophobic, Hydrophilic & Mixed
ATL Medium hydrophobic
XC Strong hydrophobic, 70 µm
30XC Strong hydrophobic, 30 µm
XWP Strong hydrophobic, wide pore
AEV Hydrophobic /hydrophilic
ATH Medium hydrophilic, 70 µm
30ATH Medium hydrophilic, 30 µm

# Sample Prep - Solid Phase Extraction

## Interchim custom packing capabilities

Interchim's state of the art automated weight machine facilitates the custom packing of columns and multi-well plates that guarantees unprecedented accuracy ( $\pm 1\%$  mass per column) compared to the commonly utilised competitor volumetric systems.

Our systems therefore raise the standards & client expectations for sample recovery /reproducibility.

Interchim is able to assist clients in identification of the necessary parameters for custom supply i.e.

- Sorbent type
- Sorbent mass
- Column, plate or container format
- Column, plate or container volume
- Frit type and porosity
- Quantity of units
- Weight certification

Considerations toward production parameters are outlined on the following page.

### Proposals

Proposals are provided within 2 working days of receiving your enquiry.

This will incorporate a confidentiality agreement to maintain the integrity of both parties.



# Sample Prep - Solid Phase Extraction

## Interchim custom packing capabilities

### Sorbent considerations

A number of client enquiries are for the packing of their own sorbent or from a third party manufacturer. In such instances, sorbent type, physical characteristics and a safety data sheet need to be specified.

The Sorbent Mass is, in part, governed by the parameters of the column/ plate to be packed. The current range for this service is between 15 mg and 70 g with a +/- 1% guaranteed accuracy (according to our standard 1 ml to 150 ml containers).

Interchim's sorbent selection guide on page A.28 can assist your selection of appropriate sorbent. The guide highlights a range of ~ 50 silica or polymer based selectivities.

### Sorbent columns or plates

#### Column

- Polypropylene straight reservoir 1 - 3 - 6 - 15 - 25 - 75 - 150 ml
- Polypropylene LRC (Large Reservoir Capacity) 15 ml
- Glass straight reservoir 6 ml

#### Plate

- 96 Well plates 2 ml
- 48 Well plates 5 - 7 ml

### Frit type and porosity

- PTFE, Teflon® or glass fiber frits for polypropylene straight tubes & LRC reservoirs
- Teflon® frits for straight glass tubes
- PTFE frits for 48 & 96 plates

### Miscellaneous

The packing of specialty containers are considered relative to the parameters of our packing system.



# Sample Prep - Solid Phase Extraction


## Vacuum manifold for SPE columns

Vacuum manifold systems are time saving devices for sample preparation that establish greater consistency and reproducibility.

Interchim manifolds are compatible with all Luer fitting SPE columns existing on the market and allow the simultaneous processing of 12-to-24 samples.

Units are easy to set up, with a simple vacuum adjustment. Stop cocks are installed on each path.

Vacuum pumps minimize flow problems associated with crude sample or insufficient vacuum.



### Vacuum manifold protocol P/N 518100

1. Screw the black lid of the upper part of the manifold
2. Introduce the plastic waste in the Glass basin
3. Insert collection needles on male Luer outlets under the lid
4. Set taps on the lids female Luer inlets
5. Place lid cover on the glass basin
6. Install SPE columns on closed taps
7. Close unused lids on female Luer inlets
8. Open vacuum on the Glass basin outlet valve
9. Establish SPE steps at constant flow and depression (the vacuum may be monitored with the Glass basin valve)
10. Remove waste after the rinsing step & insert collection rack with glass tubes
11. Replace lid cover, add vacuum to Glass basin, elute to collect sample
12. Insert dry cover instead of the lid. re-apply vacuum
13. Attach dry cover (with heated nitrogen outlet) for evaporation and sample concentration.

Description	12 position P/N	Qty	16 position P/N	Qty	24 position P/N	Qty
Vacuum manifold Set-Complete	518100	u	336570	u	Q72030	u
Glass chamber	Q71530	u	Q71900	u	Q72230	u
Cover, gasket, & 12 stopcocks	Q71540	u	Q71910	u	Q72240	u
Gaskets	Q71550	2 u	Q71920	2 u	Q72250	2 u
Vacuum gauge, valve, & glass chamber	Q71560	u	Q71930	u	Q72260	u
Needles - Polypropylene	Q57820	12 u	Q71940	16 u	Q57830	24 u
Needles - Stainless Steel	Q71570	12 u	Q71950	16 u	Q72270	24 u
Collection Rack-shelves, legs, clips, & posts	Q71580	u	Q71960	u	Q72280	u
Plate - 13 mm	Q71590	u	Q71970	u	Q72290	u
Plate - volumetric flask	Q71600	u				
Plate - 16 mm test tube	Q71610	u	Q71980	u	Q72300	u
Plate - autosampler vial	Q71640	u				
Plate - dimple	Q71660	u	Q71990	u	Q72310	u
Plate - base	Q71670	u	Q72000	u	Q72320	u
Stopcocks	Q71680	12 u	Q72010	16 u	Q72330	24 u
Drying attachment	Q71690	u	Q72020	u	Q72340	u
PP vacuum waste container	Q71700	10 u				
Support posts for rack	BI4330	u	BI4340	u	BI4340	u

# Sample Prep - Solid Phase Extraction

## Vacuum manifold

### Accessories

Description	P/N	Qty
Female Luer Fittings	Q72360	2 u
Male Luer Fittings	Q72370	2 u
Support posts for rack	Q72380	3 u
Legs for cover - black	Q72390	4 u
Vacuum gauge & valve assembly	Q72400	u
Valve assembly only	Q72420	u
Vacuum gauge	Q72440	u
Retaining clips	Q72450	12 u
Vacuum manifold plugs	Q72460	50 u
Glass collection tube 12 x 75 mm	CD9520	1000 u
Glass collection tube 16 x 100 mm	CD9530	1000 u

### Needles and control valve

Description	P/N	Qty
Control Valve	Q72470	25 u
Control Valve	Q72471	50 u
Teflon Needles	Q72500	100 u
Teflon Needles	Q72501	500 u
Teflon Flow Needle	CD7022	100 u
Teflon Flow Needle	CD7023	500 u

### Flash Purification

#### Vacuum manifold system - 10 position

accepts large capacity columns (25, 75 & 150 ml).

Description	P/N	Qty
Flash Vacuum Manifold Set-Complete	BU3010	u
Glass chamber	BU3020	u
Cover, gasket, & 10 stopcocks	BU3030	u
Gaskets	BU3040	2 u
Vacuum gauge, valve, & glass chamber	BU3050	u
Needles - Polypropylene	BU3060	10 u
Needles - Stainless Steel	BU3070	10 u
Collection Rack-shelves, legs, clips, & posts	BU3080	u
Plate - 19 mm	BU3090	u
Plate - 25 mm	BU3100	u
Plate - dimple	BU3110	u
Plate - base	BU3120	u
Stopcocks	BU3140	10 u
Drying attachment	BU3160	u
Support posts for rack	BI4340	u




# Sample Prep - Solid Phase Extraction

## Positive pressure manifold

The positive pressure manifold provides an ultra efficient, reproducible, simple robust way to simultaneously process up to 48 samples. Individually regulated positions provide even pressure to each column. Manifolds can accommodate a large range of column sizes whilst providing an even flow for each position.

- Individually regulated positions to provide even pressure to each column.
- Four banks - 12 positions each. Flow is always provided to bank one. User can control flow to each of the other three banks.
- Capacity - 1-to-48 columns of either 1 mL, 3 mL, 6 mL, 10 mL or 15 mL. Switch between column sizes by using appropriate adaptor plate.
- Air Supply : N<sub>2</sub> or compressed air regulated to 75 psi & filtered to 10 µm.
- Dual pressure regulators allow users different pressure setting for extraction and column drying.
- Waste reservoir included that can be emptied between each waste step if required.
- Single switch raises & lowers sample racks and creates an airtight seal.

Description	P/N	Qty
Completed manifold system with 10 mL / 15 mL rack, collection rack, waste rack and waste bin	BY9950	1 u
10 mL & 15 mL SPE rack	CD3550	1 u
Adapter extraction plates for 1 mL columns	CD3520	1 u
Adapter extraction plates for 3 mL columns	CD3530	1 u
Adapter extraction plates for 6 mL columns	CD3540	1 u
Installation kit	CK4570	1 u
Waste rack	CD3570	1 u
Collection rack 13 x 100 mm glass tubes	CH6220	1 u
Collection rack 16 x 100 mm glass tubes	CD3560	1 u
Replacement in-line air filter	CD3580	1 u






# Sample Prep - Solid Phase Extraction

## Universal vacuum manifold for 48 & 96 well plates

The universal well plate vacuum manifold is a versatile and effective vacuum manifold system suited to both 48 & 96 well plates for fast and efficient automated sample processing techniques.

- Filtration
- Protein crash filtration
- Solid Phase Extraction (SPE)

Description	P/N	Qty
Well plate manifold (48 & 96), complete system	AN1530	1 u
24 well collection plates 10 ml	BN9650	25 u
24 well pre-pierced lid cover	CA2320	100 u
24 well collection holder	BM4880	1 u
24 well extraction holder	BM4890	1 u
24 well rack collection tubes 16 x 60 mm	CE3800	24 u
48 well collection plate 5 ml	BG7400	25 u
48 well pre-pierced lid cover	BG7410	100 u
48 well collection holder	BM4900	1 u
48 well extraction holder	BM4910	1 u
48 well rack collection tubes 11 x 60 mm	CE3790	48 u
96 well collection plates 2 ml	U90380	50 u
96 well pre-pierced lid cover	BD7730	100 u
Standard base	BM4930	1 u
Up cover	CD3250	1 u
Universal Tube	CD3270	1 u
EPDM Gasket	CD3300	1 u
Neoprene Gasket	CD3330	1 u
Wedge 12.7 mm	CD3360	1 u
Wedge 25.4 mm	CD3390	1 u
Gauge (Manometer) / Liquid trap	CD3420	1 u
48 wells manifold for 3 ml tubes without flange	CD3500	1 u
24 wells manifold for 3 ml tubes without flange	CD3510	1 u
24 well plate caps	CH6230	24 u
48 well plate caps	CH6240	48 u



[Circular flanges are necessary for 3 & 6 ml formats. A wedge can be added to tighten the upper holder higher and allows to collect sample in approp. glass tubes. 24 & 48 position collecting plates can also be used.]

# Sample Prep - Solid Phase Extraction

## Vacuum manifold for 47 & 90 mm disc

These vacuum manifold systems allow 47 and 90 mm membrane filtration. The KEL-F membrane holder is suited for many solvents without filtrate contamination (80°C limit).

Manifolds are available with 1, 3 or 6 independent stations with a stop valve on each path. Aluminium clamps ensure a very tight seal between loading reservoir, the base and membrane holder.

Description	P/N	Qty
1 station manifold assembly (47 mm). Includes : 1 station manifold, KEL-F screen, funnel, base, clamp	<b>BX4370</b>	1 u
1 station manifold assembly (90 mm). Includes : 1 station manifold, KEL-F screen, funnel, base, clamp	<b>BX4380</b>	1 u
3 station manifold assembly (47 mm). Includes : 3 station manifold, KEL-F screen, funnel, base, clamp	<b>BX4350</b>	1 u
3 station manifold assembly (90 mm). Includes : 3 station manifold, KEL-F screen, funnel, base, clamp	<b>BX4360</b>	1 u
6 station manifold assembly (47 mm). Includes : 6 station manifold, KEL-F screen, funnel, base, clamp	<b>BX2030</b>	1 u
6 station manifold assembly (90 mm). Includes : 6 station manifold, KEL-F screen, funnel, base, clamp	<b>BX4340</b>	1 u
1 station manifold	<b>BX4400</b>	1 u
3 station manifold	<b>BX4410</b>	1 u
6 station manifold	<b>BX4420</b>	1 u
47 mm aluminum clamp	<b>BX4430</b>	1 u
90 mm aluminum clamp	<b>BX4440</b>	1 u
47 mm support base	<b>BX4450</b>	1 u
47 mm 300 ml funnel	<b>BX4460</b>	1 u
90 mm support base	<b>BX4470</b>	1 u
90 mm 1000 ml funnel	<b>BX4480</b>	1 u
47 mm KEL-F screen	<b>BX4490</b>	1 u
90 mm KEL-F screen	<b>BX4500</b>	1 u
Cartridge adaptor	<b>BX4510</b>	1 u
Teflon stopcocks for block manifolds	<b>CH6260</b>	6 u



# Sample Prep - Solid Phase Extraction

## SPE Accessories

### Polypropylene tubes

Volume	Qty	P/N
<b>Empty column</b>		
1 ml	100 u	541410
4 ml	100 u	541420
8 ml	100 u	541430
15 ml	100 u	541440
25 ml	100 u	541450
75 ml	50 u	823370
150 ml	25 u	S28581
<b>Column + one polyethylene 20 µm frits</b>		
1 ml	100 u	F97660
4 ml	100 u	F97710
8 ml	100 u	F97730
15 ml	100 u	F97750
25 ml	100 u	F97760
75 ml	50 u	F97780
<b>Column + one polyethylene 20 µm frits + caps</b>		
1 ml	100 u	F97800
4 ml	100 u	F97810
8 ml	100 u	F97820
15 ml	100 u	F97830
25 ml	100 u	F97840
75 ml	50 u	F97860
<b>Column + one PTFE 20 µm frits + caps</b>		
1 ml	100 u	F97870
4 ml	100 u	F97890
8 ml	100 u	F97900
15 ml	100 u	F97910
25 ml	100 u	F97920
75 ml	50 u	F97940
<b>Column + two polyethylene 20 µm frits</b>		
1 ml	100 u	F97960
4 ml	100 u	F97970
8 ml	100 u	F97980
15 ml	100 u	F97990
25 ml	100 u	F98000
75 ml	50 u	F98020




# Sample Prep - Solid Phase Extraction

## SPE Accessories

### Polypropylene tubes

Volume	Porosity	Membranes	Qty	P/N
<b>Column + two polyethylene 20 µm frits</b>				
1 ml			100 u	F98040
4 ml			100 u	F98050
8 ml			100 u	F98060
15 ml			100 u	F98070
25 ml			100 u	F98080
75 ml			50 u	F98100
<b>1/16" PE frits</b>				
1 ml	20 µm	PE	100 u	779530
4 ml	20 µm	PE	100 u	841880
8 ml	20 µm	PE	100 u	858750
15 ml	20 µm	PE	100 u	823280
25 ml	20 µm	PE	100 u	885460
75 ml	20 µm	PE	50 u	823380
<b>1/8" PE frits</b>				
3 ml	20 µm	PE	100 u	AM2340
8 ml	20 µm	PE	100 u	AZ3340
15 ml	20 µm	PE	100 u	S08600
<b>1/8" PE frits</b>				
25 ml	20 µm	PE	100 u	S08610
75 ml	20 µm	PE	50 u	S08620
150 ml	20 µm	PE	50 u	S28600
<b>Caps</b>				
1 ml			100 u	F97350
4 ml			100 u	F97360
8 ml			100 u	F97370
15 ml			100 u	F97440
25 ml			100 u	F97470
75 ml			50 u	F97490
<b>End caps</b>			100 u	F97510



# Sample Prep - Solid Phase Extraction

## SPE Accessories

### Glass columns

- Length : 300 mm
- Internal Ø : 20.4 mm
- Fitting : male Luer
- Option : PTFE or PE taps




### SPE column adapter

Suitable for 1, 3 and 6 ml SPE column tops with female Luer fitting

Applications :

- Increase reservoir volume that can be applied to the SPE column.
- Allow multiple selectivities when stack packed with different sorbent columns.



### Uptiplate Collection Plate & Tools for Chemistry

Uptiplate collection plates are compatible with all vacuum manifolds and SPE /filtration automated systems for the collection of solvents & samples. Manufactured from high quality polymer, they are available in round or square well formats and the well bottoms have a "U" or "V" shape design for sample management.

Interchim generally recommends the square "U" shape format. The "U" format enhances sample collection whilst the square well design assists the introduction of a flexible sealing lid.



Collection plate	Well	Well bottom	Max. volume	Qty	P/N
24 well 10 ml	Square	V	10 ml	25 u	BN9650
24 well 10 ml	Square	U	9 ml	10 u	CA2260
48 well 5 ml	Square	V	5 ml	25 u	BG7400
48 well 7,5 ml	Square	V	7,3 ml	30 u	BG7430
96 well 2 ml	Square	V	2,1 ml	50 u	U90380
96 well 1 ml	Square	U	1 ml	50 u	U90370
96 well 350 µl	Square	U	350 µl	50 u	U90360
96 well 1,3 ml	Square	V	1,35 ml	25 u	BN1940
96 well 1 ml	Round	U	1 ml	50 u	U90350
384 well 150 µl	Square	V	160 µl	25 u	BD3390
384 well 50 µl	Square	U	55 µl	60 u	BM4970
384 well 50 µl	Square	U	55 µl	30 u	BM4971

# Sample Prep - Solid Phase Extraction

## SPE Accessories

### Uptiplate soft cover

Uptiplate soft covers maintain sample integrity during delivery and storage. They are compatible with all collection plates on the market.

Pre-pierced covers are available to avoid needle clogging during automated analysis.

Santoprene is currently the soft cover of choice exhibiting a chemical resistance similar to Neoprene. They are not compatible with chlorinated or aromatic solvents such as Benzene,  $\text{CCl}_4$ ,  $\text{CHCl}_3$ ,  $\text{CHCl}_2$ ,  $\text{CH}_3\text{Cl}$ , Chlorobenzene, Cyclohexane, Toluene, Xylene. In such instances Teflon /Silicone lids are usually more suitable.

Collection plate	Well bottom	Pierceable	Nature	Qty	P/N
24 well	Square	Yes	Santoprene	100 u	CA2320
48 well	Square	Yes	Santoprene	100 u	BG7410
48 well	Square	No	Santoprene	100 u	BG7420
96 well	Square	Yes	Santoprene	100 u	BD7730
96 well	Square	Yes	Santoprene	50 u	BD7731
96 well	Square	Yes	Teflon /Silicone	20 u	CD9840
96 well	Square	No	Santoprene	100 u	AP3131
96 well	Square	No	Santoprene	50 u	AP3130
96 well	Round	No	Elastomer	50 u	BB2740
96 well	Round	Yes	Silicone encased Teflon	50 u	CD9490
96 well	Round	Yes	Teflon /Silicone	20 u	CD9510
384 well	Square	Yes	Silicone encased Teflon	30 u	CD9500



### 96 Vial plates

96 Vial plates are polymer based holders for type 1 borosilicate glass collection vials. They are compatible with all solvents and facilitate the collection of samples after filtration or extraction to 380  $\mu\text{m}$ .

An elastomer coated PTFE soft cover is inserted on every vial for sample delivery & storage. This is easily achieved with a Mat Capper. Vials can be directly accessed with a syringe. A knife can be used to access one or several closed vials.

Description	Qty	P/N
Plate + 96 vials	5 u	863610
Glass vials	500 u	863620
Soft cover	5 u	863640
Cutter	1 u	863660
Mat Capper	1 u	AL5880





# Sample Prep - Solid Phase Extraction

## Tools for Chemistry

### Tools for Chemistry

#### Uptiplate solvent storage plates

Polypropylene reagent reservoirs are generally used for reagent & solvent storage. They are suitable for a number of liquid handling configurations, from 8 or 12 channel pipettes, to 96 or 384 pipette tip heads.

Plate wells have a pyramid construction to reduce dead volume. Suited to any robotic system.


Description	Qty	P/N
<b>Uptiplate solvent storage</b>		
24 square well - 7,5 ml PP	30 u	<b>BG7430</b>
96 square well - 0,7 ml PP - Flat Bottom	50 u	<b>CD2210</b>
96 square well - 1,5 ml PP - Flat Bottom	50 u	<b>CD9810</b>
Universal lid - Robotic - Polystyrene - 127,8 x 85,5 mm	100 u	<b>BU4410</b>
Universal lid - Polystyrene - 127 x 84,7 mm	100 u	<b>BV0090</b>
Reservoir Baffle 105,7 x 69,7 mm	10 u	<b>CA2280</b>
Low Profile - PP - Flat Bottom	25 u	<b>BM4550</b>
<b>Uptiplate solvent storage - partitioned</b>		
PP - 12 columns - Pyramid Bottom - 4,4 cm height	25 u	<b>BM8790</b>
PP - 8 Row - Pyramid Bottom - 4,4 cm height	25 u	<b>BM8730</b>
PP - 8 Row - Pyramid Bottom - Low Profile - 1,9 cm height	25 u	<b>BM8780</b>
PP - 12 columns - Pyramid Bottom - Low Profile - 1,9 cm height	25 u	<b>BJ9450</b>
PP - 12 columns - Pyramid Bottom - Low Profile - 4,4 cm height	25 u	<b>BJ0870</b>
PP - 2 x 12 columns - Pyramid Bottom - Low Profile - 1,9 cm height - 3,5 ml max vol.	25 u	<b>CA2170</b>
PP - 16 Row - Pyramid Bottom - Low Profile - 1,9 cm height	25 u	<b>BM8830</b>
PP - 24 columns - Pyramid Bottom - Low Profile - 1,9 cm height	25 u	<b>BM8860</b>
PP - 4 columns + Reservoir equivalent to 20 columns - Pyramid Bottom	25 u	<b>BM8760</b>
Low Profile - 1,9 cm height		
PP - 4 columns - Pyramid Bottom - 73 ml max vol. per well	25 u	<b>CG7170</b>



BJ0870



BJ9450



BM8760



BM8780



BM8790

# Sample Prep - Solid Phase Extraction

## UptiTip™ - Microvolume preparation

### Micro Extraction - UptiTip™

UptiTip™ pipette tips are for the management of very small (0.1 µl) sample volumes and are an excellent alternative to packed cartridges.

#### UptiTip™ Coated

UptiTip™ coated pipette tips are chemically activated with a range of SPE media for ultra-efficient clean-up of small sample volumes.

##### Features :

- Faster sample preparation with reduced sample loss
- Maximised surface area-to-sample for ultra-efficient clean-up
- Direct coating removes potential sample contamination from the support
- Sample volumes as small as 0.1 µl
- Available in volumes of : 0.1-10 µl and 10-200 µl

##### Applications :

- Desalting
- MALDI
- Mass spectroscopy
- Electrophoresis
- Protein purification
- HPCE, HPLC, CEC

#### UptiTip™ Packed


UptiTip™ capped tips are filled with a range of activated sorbents for a variety of selectivities. The fine slit on the tip bottom (1-2 µm) allows liquid to pass through whilst chromatographic material (20-30 µm diameter) is retained.

No filter is required so dead volume is minimal.

The UptiTip™ can also be used as centrifugal spin column. Centrifuge adaptors are provided.

##### Features :

- Faster sample preparation with minimal sample loss
- Sample volumes as small as 0.1 µl
- Available in volumes of : 0.1-10 µl and 10-200 µl



Media	UptiTip™ -Coated*		UptiTip™-Packed*	
	1-10 µl	10-200 µl	1-10 µl	10-200 µl
C-18	BI5010	BI5020	BI5270	BI5280
C-08	BI5030	BI5040	BI5290	BI5300
C-04	BI5050	BI5060	BI5310	BI5320
Carbon	BV7460	BV7470	BU3190	BU3210
HILIC	CC6880	CC6890	CH7060	CH7070
HILIC SDS Removal	BI5100	BI5110	BI5390	BI5400
PolyCAT A	BI5120	BI5130	BI5410	BI5420
SDS-Removal	BI5150	BI1130	BI5440	BI5450
TiO <sub>2</sub>	BH3750	BH3760	BT3530	BU3630
ZrO <sub>2</sub>	BH3730	BH3740	CA8260	BX5810
<b>Affinity Media</b>				
Silica IMAC	BI5170	BI5180	BI5460	BI5470
Ni IMAC	BI5190	BI5200	BI5480	BI5500
Fe	CA8080	CH7580	CA8100	CH7490
Protein A	BI5210	BI5220	BI5510	BI5520
Protein G			BI5540	BI5560
Lectin ConA			BJ3650	BJ3770
Lectin WGA			BJ3780	BJ3790
Trypsin	BH3770	BI5230	BI5570	BI5580
Streptavidin	CH5900			

# Sample Prep - Solid Phase Extraction

## On-line extraction

### Upti-trap™ On-line extraction

Upti-trap™ ensures extraction and /or sample pre-concentration, prior to HPLC injection.


Upti-trap™ units are available either in 15 x 4.0 mm or 15 x 2.1 mm configurations. The loading capacity for a 15 x 4.0mm unit is approximately 8 to 12 mg whilst the 15 x 2.1 mm unit is approximately 1.5 to 3 mg.

The speed and efficiency of Upti-trap™ make this an excellent device for medicinal research & DMPK Labs whilst Upti-traps™ pre-concentration characteristics make it suitable to environmental applications.

This technique allows for a fast, robust and reproducible method development process that can easily be automated.


#### 1. Extraction mode

The compound of interest is trapped on the sorbent. Impurities are washed by the mobile phase A.



#### 2. Elution mode

The compound of interest is eluted by the mobile phase B.



(This example requires a 6 port - two position valve and two HPLC pumps).

A 10 port, 2 position valve increases productivity as one sample can be extracted while the other is analyzed.

Sorbent	Functionality /comment	Particle size	P/N 15 x 4,0 mm	P/N 15 x 2,1 mm
Upti-clean™ C18	C18 end-capped	50 µm	U25750	U25840
Upti-clean™ C8	C18 end-capped	50 µm	BG6090	CE0580
Upti-clean™ RPAQ	C18 100% hydrophile	75 µm	U25760	U25850
Upti-clean™ NH <sub>2</sub>	Amino	50 µm	BV3700	BV3720
Upti-clean™ SCX	Strong cation exchange	60 µm	U25770	U25860
Upti-clean™ SAX	Strong anion exchange	60 µm	BB8650	BG6950
Upti-clean™ MM1	RP / Strong cation exchange	60 µm	U25780	U25870
Atoll™ XC	Hydrophobic	30 µm	U70480	U70500
Atoll™ ATH	Hydrophilic	30 µm	CE0560	CE0590
Atoll™ AEV	Hydrophilic /Hydrophobic	75 µm	U25820	U25910
Atoll™ ATL	Hydrophobic	75 µm	U25800	U25890
Atoll™ XWP	Hydrophobic	90 µm	BU5550	BU5560



# Sample Prep - Solid Phase Extraction

United Chemical Technologies (UCT) product guide



## CLEAN SCREEN®

Drugs of Abuse Columns

- DAU = Acidic, Basic & Neutral Drugs
- THC = Carboxy THC
- THCA = THC  $\Delta^9$  Carboxylic Acid
- GHB = Gamma-Hydroxybutyrate
- ETG = Ethyl Glucuronide
- BNZ = Benzodiazepine
- CLEAN-THRU® Tips



## STYRE SCREEN®

Polymeric Based Columns

- DBX = Copolymeric
- DVB = Divinylbenzene
- BCX = Benzenesulfonic Acid
- C18 = Octadecyl C18
- CCX = Carboxylic Acid
- QAX = Quaternary Amine



## CLEAN-UP®

Solid Phase Extraction Columns

- Ion Exchange
- Hydrophobic
- Hydrophilic
- Copolymeric
- Covalent



## ENVIRO-CLEAN®

- Products
- Environmental SPE Cartridge
- SPE Inert Glass Syringe Barrels

# Sample Prep - Solid Phase Extraction


## UCT product guide

### First of its kind in the field of forensic and clinical toxicology !

The recently published *Forensic and Clinical Applications of Solid Phase Extraction*, by Michael J. Telepchak, Thomas August and Glynn Chaney, has been met with enthusiasm by those in the SPE field, and has been recommended as a valuable laboratory reference. Dr. Terry Danielson, Ph.D., who reviewed the book for the American Society of Crime Laboratory Directors, calls attention to the "extensive details of many currently available SPE separation procedures", and describes the book as a convenient compendium of SPE technology, and is relevant to the development, implementation and practice of modern SPE appropriate to students, and experienced practitioners."


In the *Canadian Society of Forensic Science Journal*, Dr. Karen Woodall, Ph.D., of the Toronto Centre of Forensic Sciences, calls the book a "must read" for anyone interested in SPE, especially helpful in resolving the day-to-day problems that can occur when using [SPE] and gives some excellent examples of how to deal with some of these occurrences such as recovery variability, contamination, flow, and non-extraction problems."

BOOK "FORENSIC AND CLINICAL APPLICATIONS OF SOLID PHASE EXTRACTION"  
P/N : ZZ3801



### Reservoirs for Bonded Phase Extractions

Stated Volume (mL)	Tube Configuration	Bed Diameter (mm)	Sorbent Mass (mg)
1	Cylindrical	5.5	50-200
3	Cylindrical	8.5	50-1000
6	Cylindrical	12.5	200-2000
10	Expanded	8.5	50-1000
15	Cylindrical	15.5	500-2000
25	Cylindrical	20	500-5000
75	Cylindrical	27.5	1000-10000
150	Cylindrical	38.0	10000-70000



1 mL



3 mL




6 mL




10 mL




15 mL



25 mL



75 mL



150 mL



# Sample Prep - Solid Phase Extraction

UCT Clean Screen® columns

## Copolymeric bonded phases for drug abuse testing

Analytical demand for more efficient, robust and clean extraction of drugs from biological matrices led to the development of Clean Screen® sorbents. Since 1986, Clean Screen® has led the industry with dependable and reproducible Solid Phase Extraction products and applications. Clean Screen phases are true copolymeric sorbents that contain hydrophobic and ion exchange functional groups uniquely polymerized to a silica substrate. The design and quality of Clean Screen provides superior sample clean up, recovery and reproducibility.

Mixed mode separations allow maximum selectivity for extraction of acids, neutrals and bases. This selectivity makes Clean Screen ideal for both screening and confirmation analysis for virtually all drug categories. Clean Screen DAU, THC, and GHB columns are used extensively by forensic and clinical chemists including :

- Post Mortem Investigations
- Criminal Investigations
- Urine Drug Testing
- Athletic Drug Testing
- Racing Laboratories
- Therapeutic Drug Monitoring
- Medical Drug Screening



Recent additions to this product line include Clean Screen® Ethyl Glucuronide and Clean Screen® Benzodiazepines.

## Mechanism of Clean Screen®

Carboxylic acid functionalities present in the sample are ionized. This creates a repulsion between the column and many sample borne interferences, thereby reducing the likelihood of their adsorbing onto the column.

At this pH, ibuprofen & barbiturates are not ionized and are hydrophobically adsorbed onto the column (figure 1). At the same time, drugs with amine functionalities such as cocaine and phencyclidine adsorb onto the column via both hydrophobic and ionic attraction (fig. 1).

The column can then be washed with water or weak aqueous buffers at or below pH 6 without risking loss of the analytes. After drying the column, it is possible to elute the hydrophobically bound analytes using solvents of minimal polarity such as methylene chloride or a hexane/ethyl acetate mixture (fig. 2). Cationic analytes will remain bound to the column. Many compounds of intermediate polarity and potential interferences will also remain on the column. The majority of these potential interferences can be removed by using a methanol wash.







# Sample Prep - Solid Phase Extraction

## UCT Clean Screen® columns

Cationic analytes bound to the column can be eluted after another drying step. The drying steps are necessary to remove water which would have prevented the water-immiscible elution solvents from optimally interacting with the analytes (fig. 3).



To elute the cationic analytes, an organic solution with a high pH (between 11 & 12) should be used. A methylene chloride isopropanol-ammonium hydroxide mixture will simultaneously disrupt these ionic interactions and successfully elute the desired compound (fig. 4).



### Clean Screen® DAU

This column is copolymerized on a rigid, purified silica gel support. The two functional groups include a reverse phase, and an ion exchanger, benzenesulfonic acid. This column is commonly used for analyzing a wide range of drugs of abuse, including acidic, basic & neutral drugs.

**Application :** Dual functionality for weak bases and hydrophobic compounds.

P/N	Weight /Vol.	Qty
CSDAU131	130 mg/1 mL	100
CSDAU133	130 mg/3 mL	50
CSDAU203	200 mg/3 mL	50
CSDAU303	300 mg/3 mL	50
CSDAU503	500 mg/3 mL	50
CSDAU206	200 mg/6 mL	50
CSDAU506	500 mg/6 mL	50
CSDAU1M6	1 g/6 mL	30
ZSDAU005	50 mg/10 mL	50
ZSDAU013	130 mg/10 mL	50
ZSDAU020	200 mg/10 mL	50
CSDAU515	500 mg/15 mL	50

### Clean Screen® THC

This column is copolymerized on a rigid, purified silica gel support. The two functional groups include a reverse phase, and an ion exchanger, primary amine. This column is commonly used for analyzing THC and its metabolites.


**Application :** Dual functionality for acids and hydrophobic compounds.

P/N	Weight /Vol.	Qty
CSTHC131	130 mg/1 mL	100
CSTHC203	200 mg/3 mL	50
CSTHC303	300 mg/3 mL	50
CSTHC503	500 mg/3 mL	50
CSTHC206	200 mg/6 mL	50
CSTHC506	500 mg/6 mL	50
CSTHC1M6	1 g/6 mL	30
ZSTHC013	130 mg/10 mL	50
ZSTHC020	200 mg/10 mL	50
CSTHC515	500 mg/15 mL	50

### Clean Screen® GHB

The small polar nature of the molecule and the lack of a UV chromophore complicate the chromatographic and spectrophotometric analysis of GHB. Chemically, GHB is unstable and readily forms Gamma-butyrolactone when heated in acid conditions. Most analytical methods are based upon the interconversion to the lactone and chemical derivatization to form the TMS derivative. This column is for the extraction of free GHB.

P/N	Weight /Vol.	Qty
CSGHB203	200 mg/3 mL	50
ZSGHB020	200 mg/10 mL	50
ZCGHB020	200 mg/10 mL	50



## UCT Clean-Up® columns

## Hydrophobic Extraction Columns


This sorbent is composed of a silica backbone bonded with hydrocarbon chains. It is used to extract compounds which exhibit non-polar or neutral characteristics out of complex matrices. The C18 phase is the most widely used for non-polar interactions because of its nonselective nature; C18 will extract a large number of compounds with differing chemical properties. To enhance selectivity, UCT offers a wide range of hydrophobic sorbents, from C2 to C20. Multiple chain configurations are available for some sorbents. Endcapped or unendcapped sorbents are available for all chain lengths.



## Mechanism of hydrophobic bonding



Compounds are retained by non-polar interactions from polar solvents or matrix environments. They are bound by dispersion forces / van der Waals forces. Elution, or disruption of the non-polar interactions is achieved by solvents or solvent mixtures with sufficient non-polar character. Some polar solvents, such as acetonitrile have enough non-polar characteristics to disrupt non-polar binding to cause elution of a compound from the sorbent. Methanol can be used as well, although it should be noted that it will take off both polar & non-polar analytes of interest & interferences.

### Example of Hydrophobic Bonding



## Unendcapped vs. Endcapped

Bonded phases are manufactured by the reaction of organosilanes with activated silica. During the polymerization reaction of carbon chains to the silica backbone, a very stable silyl ether linkage forms. Our unendcapped columns allow hydroxyl sites to remain, thus making these columns slightly hydrophilic. In order to decrease this slight polarity, these hydroxyl sites are deactivated. Proprietary bonding techniques ensure that these sites are 100% reacted, leading to a complete endcapping. Because there are no hydroxyl sites left, our endcapped columns are more hydrophobic than our unendcapped columns.



## Functionalized hydrophobic silica based phases

Sorbent	Product code	Structure	% Carbon
C2 ethyl	C02	-SiCH <sub>2</sub> CH <sub>3</sub>	6.60
C3 propyl	C03	-Si-(CH <sub>2</sub> ) <sub>2</sub> CH <sub>3</sub>	7.60
C4 n-butyl	Cn4	-Si-(CH <sub>2</sub> ) <sub>3</sub> CH <sub>3</sub>	8.50
Ci4 isobutyl	Ci4	-Si-CH <sub>2</sub> CH(CH <sub>3</sub> ) <sub>2</sub>	8.80
Ct4 tertiary butyl	Ct4	-Si-C(CH <sub>3</sub> ) <sub>3</sub>	8.50
C5 pentyl	C05	-Si-(CH <sub>2</sub> ) <sub>4</sub> CH <sub>3</sub>	9.50
C6 hexyl	C06	-Si-(CH <sub>2</sub> ) <sub>5</sub> CH <sub>3</sub>	11.00
C7 heptyl	C07	-Si-(CH <sub>2</sub> ) <sub>6</sub> CH <sub>3</sub>	11.00
C8 octyl	C08	-Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	11.10
C10 decyl	C10	-Si-(CH <sub>2</sub> ) <sub>9</sub> CH <sub>3</sub>	15.70
C12 dodecyl	C12	-Si-(CH <sub>2</sub> ) <sub>11</sub> CH <sub>3</sub>	not tested
C18 octadecyl	C18	-Si-(CH <sub>2</sub> ) <sub>17</sub> CH <sub>3</sub>	21.70
C20 eicosyl	C20	-Si-(CH <sub>2</sub> ) <sub>19</sub> CH <sub>3</sub>	24.30
C30 tricontyl	C30	-Si-(CH <sub>2</sub> ) <sub>29</sub> CH <sub>3</sub>	26.00
Cyclohexyl	CYH1	-Si-(CH)	11.60
Phenyl	PHY1	-Si-(PH)	11.00

# Sample Prep - Solid Phase Extraction

## UCT Clean-Up® columns

Weight/Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50mg/1mL	C2, Ethyl	CEC021L1	CUC021L1	C4, n-Butyl*	CECN41L1	CUCN41L1	C6, Hexyl	CEC061L1	CUC061L1	100
100mg/1mL		CEC02111	CUC02111		CECN4111	CUCN4111		CEC06111	CUC06111	100
100mg/3mL		CEC02113	CUC02113		CECN4113	CUCN4113		CEC06113	CUC06113	50
200mg/3mL		CEC02123	CUC02123		CECN4123	CUCN4123		CEC06123	CUC06123	50
500mg/3mL		CEC02153	CUC02153		CECN4153	CUCN4153		CEC06153	CUC06153	50
500mg/6mL		CEC02156	CUC02156		CECN4156	CUCN4156		CEC06156	CUC06156	50
1g/6mL		CEC021M6	CUC021M6		CECN41M6	CUCN41M6		CEC061M6	CUC061M6	30
100mg/10mL		CEC0211Z	CUC0211Z		CECN411Z	CUCN411Z		CEC0611Z	CUC0611Z	50
200mg/10mL		CEC0212Z	CUC0212Z		CECN412Z	CUCN412Z		CEC0612Z	CUC0612Z	50
500mg/10mL		CEC0215Z	CUC0215Z		CECN415Z	CUCN415Z		CEC0615Z	CUC0615Z	50
50mg/1mL	C3, Propyl	CECN31L1	CUCN31L1	C5, Pentyl	CEC051L1	CUC051L1	C7, Heptyl	CEC071L1	CUC071L1	100
100mg/1mL		CECN3111	CUCN3111		CEC05111	CUC05111		CEC07111	CUC07111	100
100mg/3mL		CECN3113	CUCN3113		CEC05113	CUC05113		CEC07113	CUC07113	50
200mg/3mL		CECN3123	CUCN3123		CEC05123	CUC05123		CEC07123	CUC07123	50
500mg/3mL		CECN3153	CUCN3153		CEC05153	CUC05153		CEC07153	CUC07153	50
500mg/6mL		CECN3156	CUCN3156		CEC05156	CUC05156		CEC07156	CUC07156	50
1g/6mL		CECN31M6	CUCN31M6		CEC051M6	CUC051M6		CEC071M6	CUC071M6	30
100mg/10mL		CECN311Z	CUCN311Z		CEC0511Z	CUC0511Z		CEC0711Z	CUC0711Z	50
200mg/10mL		CECN312Z	CUCN312Z		CEC0512Z	CUC0512Z		CEC0712Z	CUC0712Z	50
500mg/10mL		CECN315Z	CUCN315Z		CEC0515Z	CUC0515Z		CEC0715Z	CUC0715Z	50

\*Available on request C4 Isobutyl and C4 Tertiary Butyl

Weight/Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50mg/1mL	C8, Octyl	CEC081L1	CUC081L1	C12, nDodecyl	CEC121L1	CUC121L1	Cyclohexyl	CECYH1L1	CUCYH1L1	100
100mg/1mL		CEC08111	CUC08111		CEC12111	CUC12111		CECYH111	CUCYH111	100
100mg/3mL		CEC08113	CUC08113		CEC12113	CUC12113		CECYH113	CUCYH113	50
200mg/3mL		CEC08123	CUC08123		CEC12123	CUC12123		CECYH123	CUCYH123	50
500mg/3mL		CEC08153	CUC08153		CEC12153	CUC12153		CECYH153	CUCYH153	50
500mg/6mL		CEC08156	CUC08156		CEC12156	CUC12156		CECYH156	CUCYH156	50
1g/6mL		CEC081M6	CUC081M6		CEC121M6	CUC121M6		CECYH1M6	CUCYH1M6	30
100mg/10mL		CEC0811Z	CUC0811Z		CEC1211Z	CUC1211Z		CECYH11Z	CUCYH11Z	50
200mg/10mL		CEC0812Z	CUC0812Z		CEC1212Z	CUC1212Z		CECYH12Z	CUCYH12Z	50
500mg/10mL		CEC0815Z	CUC0815Z		CEC1215Z	CUC1215Z		CECYH15Z	CUCYH15Z	50
50mg/1mL	C10, nDecyl	CEC101L1	CUC101L1	C18, Octadecyl	CEC181L1	CUC181L1	Phenyl	CEPHY1L1	CUPHY1L1	100
100mg/1mL		CEC10111	CUC10111		CEC18111	CUC18111		CEPHY111	CUPHY111	100
100mg/3mL		CEC10113	CUC10113		CEC18113	CUC18113		n.a.	n.a.	50
200mg/3mL		CEC10123	CUC10123		CEC18123	CUC18123		CEPHY123	CUPHY123	50
500mg/3mL		CEC10153	CUC10153		CEC18153	CUC18153		CEPHY153	CUPHY153	50
500mg/6mL		CEC10156	CUC10156		CEC18156	CUC18156		CEPHY156	CUPHY156	50
1g/6mL		CEC101M6	CUC101M6		CEC181M6	CUC181M6		CEPHY1M6	CUPHY1M6	30
100mg/10mL		CEC1011Z	CUC1011Z		CEC1811Z	CUC1811Z		CEPHY11Z	CUPHY11Z	50
200mg/10mL		CEC1012Z	CUC1012Z		CEC1812Z	CUC1812Z		CEPHY12Z	CUPHY12Z	50
500mg/10mL		CEC1015Z	CUC1015Z		CEC1815Z	CUC1815Z		CEPHY15Z	CUPHY15Z	50

# Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

## Hydrophilic Normal Phase Columns

This sorbent is composed of a silica backbone bonded with carbon chains containing polar functional groups. Groups which will possess such polarity include amines, hydroxyls and carbonyls.


## Functionalized hydrophilic silica based phases

Sorbent	Product code	Structure	% Carbon
Silica	SIL1	-SiOH	N/A
Diol	DOL1	-Si-(CH <sub>2</sub> ) <sub>3</sub> OCH <sub>2</sub> CHOHCH <sub>2</sub> OH	8.00
Cyanopropyl	CNP1	-Si-(CH <sub>2</sub> ) <sub>3</sub> CN	6.90
Florisil®	FLS		N/A
Alumina, Acidic	ALA		N/A
Alumina, Neutral	ALN		N/A
Alumina, Basic	ALB		N/A
Carbon	CARB		N/A

## Mechanism of hydrophilic bonding

Compounds are retained on hydrophilic sorbents through polar interactions including hydrogen bonding, pi-pi or dipole-dipole interaction. These types of interactions occur when a distribution of electrons between individual atoms in functional groups is unequal, causing negative and positive polarity. Compounds typically extracted on a hydrophilic column include analytes which have polar groups, including amines, hydroxyls and carbonyls. Elution is performed by strong polar solvents.

Example of  
Hydrophilic  
Bonding




Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50mg/1mL	Unbonded Silica	CUSIL1L1	Florisil®	CUFLS1L1	Alumina, Acidic	CUALA1L1	Alumina, Basic	CUALB1L1	100
100mg/1mL		CUSIL111		CUFLS111		CUALA111		CUALB111	100
200mg/3mL		CUSIL123		CUFLS123		CUALA123		CUALB123	50
500mg/3mL		CUSIL153		CUFLS153		CUALA153		CUALB153	50
500mg/6mL		CUSIL156		CUFLS156		CUALA156		CUALB156	50
1g/6mL		CUSIL1M6		CUFLS1M6		CUALA1M6		CUALB1M6	30
100mg/10mL		CUSIL11Z		CUFLS11Z		CUALA11Z		CUALB11Z	50
200mg/10mL		CUSIL12Z		CUFLS12Z		CUALA12Z		CUALB12Z	50
500mg/10mL		CUSIL15Z		CUFLS15Z		CUALA15Z		CUALB15Z	50
50mg/1mL	Alumina, neutral	CUALN1L1	CN, Cyanopropyl Endcapped	CECNP1L1	CN, Cyanopropyl Unendcapped	CUCNP1L1	Diol	CUDOL1L1	100
100mg/1mL		CUALN111		CECNP111		CUCNP111		CUDOL111	100
200mg/3mL		CUALN123		CECNP123		CUCNP123		CUDOL123	50
500mg/3mL		CUALN153		CECNP153		CUCNP153		CUDOL153	50
500mg/6mL		CUALN156		CECNP156		CUCNP156		CUDOL156	50
1g/6mL		CUALN1M6		CECNP1M6		CUCNP1M6		CUDOL1M6	30
100mg/10mL		CUALN11Z		CECNP11Z		CUCNP11Z		CUDOL11Z	50
200mg/10mL		CUALN12Z		CECNP12Z		CUCNP12Z		CUDOL12Z	50
500mg/10mL		CUALN15Z		CECNP15Z		CUCNP15Z		CUDOL15Z	50

# Sample Prep - Solid Phase Extraction

UCT Clean-Up® columns

## Ion Exchange extraction columns

This sorbent is composed of a silica backbone bonded with a carbon chain terminated by a negatively or positively charged functional group. Ion exchange interactions occur between a sorbent that carries a charge and a compound of opposite charge. This electrostatic interaction is reversible by neutralizing the sorbent and/or analyte. Ion exchange bonds can also be disrupted by introduction of a "counter ion" to compete with the analyte for binding sites on the sorbent.



Sorbent	Code	Structure	Pka	% Carbon	meq /g
Aminopropyl (1° amine)	NAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> NH <sub>2</sub>	9.8	6.65	0.310
N-2 Aminoethyl (1° & 2° amine)	PSA1	-Si-(CH <sub>2</sub> ) <sub>3</sub> NH(CH <sub>2</sub> ) <sub>2</sub> NH <sub>2</sub>	10.1, 10.9	9.70	0.320
Diethylamino (3° amine)	DAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> N(CH <sub>2</sub> CH <sub>3</sub> ) <sub>2</sub>	10.6	8.40	0.280
Quaternary Amine Chloride	QAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> N <sup>+</sup> -(CH <sub>3</sub> ) <sub>3</sub> Cl <sup>-</sup>	always charged	8.40	0.250
Quaternary Amine Hydroxide	CHQAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> N <sup>+</sup> -(CH <sub>3</sub> ) <sub>3</sub> CH <sub>3</sub> CO <sub>2</sub> <sup>-</sup>	always charged	8.40	0.250
Quaternary Amine Acetate	CAQAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> N <sup>+</sup> -(CH <sub>3</sub> ) <sub>3</sub> OH <sup>-</sup>	always charged	8.40	0.250
Quaternary Amine Formate	CFQAX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> N <sup>+</sup> -(CH <sub>3</sub> ) <sub>3</sub> CHO <sub>2</sub> <sup>-</sup>	always charged	8.40	0.250
Cation Exchange					
Carboxylic Acid	CCX1	-Si-CH <sub>2</sub> COOH	4.8	9.10	0.170
Propylsulfonic Acid	PCX1	-Si-(CH <sub>2</sub> ) <sub>3</sub> SO <sub>3</sub> H	<1	7.10	0.180
Benzenesulfonic Acid	BCX1	-Si-(CH <sub>2</sub> ) <sub>2</sub> -(PH)-SO <sub>3</sub> H	always charged	11.00	0.320
Benzenesulfonic Acid High Load	BCXH1	-Si-(CH <sub>2</sub> ) <sub>2</sub> -(PH)-SO <sub>3</sub> H	always charged	15.00	0.650

## Mechanism of Ion Exchange bonding

Compounds are retained on the sorbent through ionic bonds. Therefore, it is essential that the sorbent and the analyte to be extracted are charged. Generally, the number of molecules with charged cationic groups increases at pH values below the molecule's pKa value. The number of molecules with charged anionic groups decreases at pH values below the molecule's pKa value. To ensure 99% or more ionization, the pH should be at least two pH units below the pKa of the cation and two pH units above the pKa of the anion. Elution occurs by using a solvent to raise the pH above the pKa of the cationic group or to lower the pH below the pKa of the anion to disrupt retention. At this point, the sorbent or compound will be neutralized.

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50mg/1mL	Aminopropyl	CUNAX1L1	"PSA	CUPSA1L1	Diethylamino	CUDAX1L1	"Quaternary	CUQAX1L1	100
100mg/1mL		CUNAX111	N-2 Aminoethyl"	CUPSA111		CUDAX111	Amine	CUQAX111	100
200mg/3mL		CUNAX123		CUPSA123		CUDAX123	with Chloride	CUQAX123	50
500mg/3mL		CUNAX153		CUPSA153		CUDAX153	Counter Ion**	CUQAX153	50
500mg/6mL		CUNAX156		CUPSA156		CUDAX156		CUQAX156	50
1g/6mL		CUNAX1M6		CUPSA1M6		CUDAX1M6		CUQAX1M6	30
100mg/10mL		CUNAX11Z		CUPSA11Z		CUDAX11Z		CUQAX11Z	50
200mg/10mL		CUNAX12Z		CUPSA12Z		CUDAX12Z		CUQAX12Z	50
500mg/10mL		CUNAX15Z		CUPSA15Z		CUDAX15Z		CUQAX15Z	50
50mg/1mL	Carboxylic Acid	CUCCX1L1	Propyl	CUPCX1L1	Benzene-	CUBCX1L1	Benzene-	CUBCX1H1L1	100
100mg/1mL		CUCCX111	sulfonic	CUPCX111	-sulfonic	CUBCX111	-sulfonic	CUBCX1HL11	100
200mg/3mL		CUCCX123	Acid	CUPCX123	Acid	CUBCX123	Acid High	CUBCX1HL23	50
500mg/3mL		CUCCX153		CUPCX153		CUBCX153	Load	CUBCX1HL53	50
500mg/6mL		CUCCX156		CUPCX156		CUBCX156		CUBCX1HL56	50
1g/6mL		CUCCX1M6		CUPCX1M6		CUBCX1M6		CUBCX1HLM6	30
100mg/10mL		CUCCX11Z		CUPCX11Z		CUBCX11Z		CUBCX1HL1Z	50
200mg/10mL		CUCCX12Z		CUPCX12Z		CUBCX12Z		CUBCX1HL2Z	50
500mg/10mL		CUCCX15Z		CUPCX15Z		CUBCX15Z		CUBCX1HL5Z	50

\*\*Available with Acetate , Hydroxide & Formate Counter Ion



## UCT Clean-Up® columns

### Copolymeric Extraction columns (ion Exchange + hydrophobicity)


This sorbent is composed of a silica backbone with two types of functional chains attached - an ion exchanger or polar chain and a hydrophobic carbon chain. Our copolymeric phases are produced in a way to allow for equal parts of each functional group to attach to the silica backbone. This copolymerization technique yields reproducible bonded phases and unique copolymeric chemistries which allow the controlled use of mixed mode separation mechanisms. This type of dual chemistry is beneficial especially when one is looking for both a neutral & charged compound. This is common when a neutral parent drug metabolizes and becomes a charged compound.



### Functionalized mixed mode silica based phases



Sorbent	Product code	Structure	% Carbon	meq/g
Aminopropyl + C8	NAX2	-Si-(CH <sub>2</sub> ) <sub>3</sub> NH <sub>2</sub> & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	12.3	0.163
Quaternary Amine + C8	QAX2	-Si-(CH <sub>2</sub> ) <sub>3</sub> N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub> & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	13.60	0.160
Carboxylic Acid + C8	CCX2	-Si-CH <sub>2</sub> COOH & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	12.50	0.105
Propylsulfonic Acid + C8	PCX2	-Si-(CH <sub>2</sub> ) <sub>3</sub> SO <sub>3</sub> H & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	14.62	0.114
Benzenesulfonic Acid + C8	BCX2	-Si-(CH <sub>2</sub> ) <sub>2</sub> -(PH)-SO <sub>3</sub> H & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	12.30	0.072
Cyanopropyl + C8	CNP2	-Si-(CH <sub>2</sub> ) <sub>3</sub> CN & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	14.60	0.163
Cyclohexyl + C8	CYH2	-Si-(PH) & -Si-(CH <sub>2</sub> ) <sub>7</sub> CH <sub>3</sub>	N/A	N/A

### Example of Copolymeric Bonding



## Mechanism of mixed mode bonding

Using a sample composed of a theoretical neutral parent drug and its charged (acidic) metabolite, it is applied at a pH of 6 (figure 1). At this pH, many amine groups are positively charged. Since the column is also positively charged, compounds with this chemistry (cations) are repelled. Depending on the pKa of the metabolite, carboxylic acid groups may be negatively charged, allowing the metabolite to bond to the positively charged sorbent. Since the column also possesses a hydrophobic chain, the neutral parent drug also bonds to the column. Water or a weak aqueous buffer (pH6) washes away hydrophilically bound interferences (figure 2). The column is then dried, careful to free the column of any residual aqueous phase that would interfere with elution.






# Sample Prep - Solid Phase Extraction

## UCT Clean-Up® columns


Elution 1  
(Figure 3)

The hydrophobically bound neutral parent drug is eluted with a solvent of minimal polarity, such as hexane/ethyl acetate - 80:20.



Elution 2  
(Figure 4)

The final elution employs an acid to neutralize the charge of acidic analytes. Ionic interaction is released, and analytes are eluted in an appropriate solvent mixture.



Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Hydrophobic	CUPCX2L1	Hydrophobic	CUCCX2L1	Hydrophobic	CUBCX2L1	Octadecyl	CUBCX3L1	100
100 mg/1 mL	plus	CUPCX211	plus	CUCCX211	plus	CUBCX211	plus	CUBCX311	100
200 mg/3 mL	Propylsulfonic	CUPCX223	Carboxylic	CUCCX223	Benzene-	CUBCX223	Benzene-	CUBCX323	50
500 mg/3 mL	Acid	CUPCX253	Acid	CUCCX253	-sulfonic	CUBCX253	-sulfonic	CUBCX353	50
500 mg/6 mL		CUPCX256		CUCCX256	Acid	CUBCX256	Acid	CUBCX356	50
1 g/6 mL		CUPCX2M6		CUCCX2M6		CUBCX2M6		CUBCX3M6	30
100 mg/10 mL		CUPCX21Z		CUCCX21Z		CUBCX21Z		CUBCX31Z	50
200 mg/10 mL		CUPCX22Z		CUCCX22Z		CUBCX22Z		CUBCX32Z	50
500 mg/10 mL		CUPCX25Z		CUCCX25Z		CUBCX25Z		CUBCX35Z	50
50 mg/1 mL	Hydrophobic	CUPSA2L1	Octadecyl	CUPSA3L1	Hydrophobic	CUQAX2L1	Hydrophobic	CUNAX2L1	100
100 mg/1 mL	plus	CUPSA211	plus	CUPSA311	plus	CUQAX211	plus	CUNAX211	100
200 mg/3 mL	N-2	CUPSA223	N-2	CUPSA323	Quaternary	CUQAX223	Aminopropyl	CUNAX223	50
500 mg/3 mL	Aminoethyl	CUPSA253	Aminoethyl	CUPSA353	Amine	CUQAX253		CUNAX253	50
500 mg/6 mL		CUPSA256		CUPSA356		CUQAX256		CUNAX256	50
1 g/6 mL		CUPSA2M6		CUPSA3M6		CUQAX2M6		CUNAX2M6	30
100 mg/10 mL		CUPSA21Z		CUPSA31Z		CUQAX21Z		CUNAX21Z	50
200 mg/10 mL		CUPSA22Z		CUPSA32Z		CUQAX22Z		CUNAX22Z	50
500 mg/10 mL		CUPSA25Z		CUPSA35Z		CUQAX25Z		CUNAX25Z	50

# Sample Prep - Solid Phase Extraction

## UCT Styre Screen® Polymeric Resin columns

Styre Screen® extraction columns contain an ultra clean, highly cross-linked styrene and divinylbenzene copolymer sorbent that is functionalized with both a reverse phase, hydrophobic component and a strong cation exchanger. High & reproducible recoveries for acidic, neutral and basic compounds are achievable with a single column. The Styre Screen® particles have an average particle size of 30 microns and a very high analyte capacity making them ideal for standard solid phase extraction applications. The increased analyte capacity means that less sorbent bed mass is needed which results in faster flow rates and less solvent use. Higher throughput and less solvent waste disposal translate into significant savings in both time and money. In addition, no conditioning steps are required for most drugs of abuse applications.



### Advantage

- No conditioning steps
- Copolymer allows for extraction of acids, neutrals and bases
- High and reproducible recoveries
- Clean extractions
- Highly cross-linked styrene/divinylbenzene polymer
- Reduction in sorbent mass
- Faster flow rates
- pH stable (1 to 14)
- Reduction in solvent use
- High sorbent capacity
- Methods for NIDA/SAMHSA 5 Drugs

Weight /Vol.	Sorbent	P/N	Application	Qty
30 mg/1 mL	DBX - Benzenesulfonic Acid + C18	SSDBX031	Dual functionality for weak acids and hydrophobic compounds.	100
30 mg/3 mL		SSDBX033		50
50 mg/6 mL		SSDBX056		50
30 mg/1 mL	DVB - Polystyrene Divinylbenzene	SSDBX031	n.a.	100
30 mg/3 mL		SSDBX033	n.a.	50
50 mg/6 mL		SSDBX056	n.a.	50
30 mg/1 mL	BCX - Reverse Phase	SSBCX031	Scavenger for amines, alcohols and other nucleophiles.	100
30 mg/3 mL		SSBCX033		50
50 mg/6 mL		SSBCX056		50
30 mg/1 mL	C18 -Reverse Phase	SSC18031	Removes hydrophobic impurities, de-salting and purification of hydrophobic compounds.	100
30 mg/3 mL		SSC18033		50
50 mg/6 mL		SSC18056		50
30 mg/1 mL	CCX - Carboxylic Acid	SSCCX031	Scavenger for strong anions (Quaternary amines or metals)	100
30 mg/3 mL		SSCCX033		50
50 mg/6 mL		SSCCX056		50
30 mg/1 mL	QAX - Quaternary Amine	SSQAX031	"Removes large or more hydrophobic compounds."	100
30 mg/3 mL		SSQAX033		50
50 mg/6 mL		SSQAX056		50



# Sample Prep - Solid Phase Extraction

## UCT Enviro-Clean® - Universal cartridges



The Enviro-Clean® Universal Cartridge is the choice of modern contract labs. This inexpensive, easy to use cartridge provides consistent extractions with clean blanks. Built in flow control allows for consistent flow rates. Enviro-Clean® sorbents UCT polypropylene, and PTFE frits offer a clean blank with every batch. Designed for the environmental lab, the cartridge is made to handle large volumes of waste water. An optional bottle holder is available for continuous feed from Boston Round and wide mouth bottles.



SPE-DEX® is a registered trademark of Horizon Technology, Inc.

Product Name	P/N	Weight /Vol.	Description / Application	Qty
UNIVERSAL C18	<b>ECUNIC18</b>	1100 mg/83 mL	1100 mg of endcapped C18 for pesticides, PCBs and a large assortment of applications.	8
UNIVERSAL 525	<b>ECUNI525</b>	1500 mg/83 mL	1500 mg of our special C18 blend. This cartridge is specifically designed for EPA Method 525.	8
UNIVERSAL PAH / DRO	<b>ECUNIPAH</b>	2000 mg/83 mL	2000 mg of C18 specifically designed for PAH extraction.	8
UNIVERSAL OIL & GREASE	<b>ECUNIOAG</b>	4000 mg/83 mL	4000 mg of large particle C18 with an assortment of PE frit filters. No more liquid/liquid emulsions or clogged disks.	15

The cartridge will fit all standard manifolds and disk manifolds with adapter.

# Sample Prep - Solid Phase Extraction

## UCT dispersive SPE – “QuEChERS”

**Quick, Easy, Cheap, Effective, Rugged and Safe Approach for determining pesticide residues in fruits, vegetables and other foods.**

The QuEChERS method is gaining in popularity around the world as the method of choice for food testing. The QuEChERS method offers the advantages of high recoveries, accurate results, high sample throughput, low solvent and glassware usage (no chlorinated solvents), less labor and bench space and lower reagent costs. Organic acids and other potential contaminants are removed during the cleanup process.

UCT provides a variety of QuEChERS products containing primary secondary amine (PSA), C18, magnesium sulfate anhydrous and graphitized carbon black. These products are used in the method's clean-up step. Bulk, pre-cleaned magnesium sulfate anhydrous is available for the extraction part of the method.

PSA is used to remove various polar organic acids, polar pigments, some sugars and fatty acid co-extractables from QuEChERS extracts. Combined with C18, samples containing less than 1% lipids can be cleaned of most lipids and sterols. Graphitized carbon is used to remove sterols and pigments such as chlorophyll. The downside to carbon is its ability to retain planar molecules. Schenck and Vega (April 2001) reported that 3/1 acetone : toluene performed well at eluting many compounds from carbon.



P/N	Description	Qty
CUMPSCB2CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA & 50mg Carbon	100
CUMPS2CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA	100
CUMPS18CT	2 mL micro-centrifuge tubes with 150 mg Anhydrous Magnesium Sulfate, 50 mg PSA & 50 mg endcapped C18	100
ECMPSCB15CT	15 mL centrifuge tubes with 900 mg Anhydrous Magnesium Sulfate, 300 mg PSA & 50 mg endcapped C18	50
ECMPSC1815CT	15 mL centrifuge tubes with 900 mg Anhydrous Magnesium Sulfate, 300 mg PSA & 150 mg endcapped C18	50
ECPSACB6	6 mL columns with 400 mg PSA on bottom, 200 mg Graphitized Carbon-Black on top, separated by a Teflon frit*	30
ECPSACB256	6 mL columns with 250 mg Graphitized Carbon on top, 500 mg PSA on the bottom, separated with a Teflon frit*	30
ECPSACB506	6 mL columns with 500 mg Graphitized Carbon on top, 500 mg PSA on the bottom, separated with a Teflon frit*	30
ECMSSA50CT	50 mL PP centrifuge tube with 6 g Anhydrous Magnesium Sulfate, 1.5 g Anhydrous Sodium Acetate	250
ECMSSC50CT	50 mL PP centrifuge tube with 4 g Anhydrous Magnesium Sulfate, 1 g NaCl	250

### Example of procedure

1. Transfer 15 g of homogenized sample into a 50 ml FEP centrifuge tube.
2. Add 15 ml of 1% acetic acid in acetonitrile, 1.5 g sodium acetate anhydrous, 6 g of UCT magnesium sulfate anhydrous and an internal standard.
3. Shake vigorously for 1 minute.
4. Centrifuge for 3 minutes at 3700 rpm.
5. Transfer an aliquot of the supernatant to the UCT product.
6. Shake for 1 minute.
7. Centrifuge for 3 minutes at 3700 rpm.
8. Analyze.

\* Products available with Polyethylene or Teflon frits. Choice depends application and price requirements.

# Sample Prep - Solid Phase Extraction

## UCT Enviro-Clean® columns

### Polypropylene and Inert Glass Extraction columns

Enviro-Clean® solid phase extraction columns are designed especially for the isolation and separation of environmental analytes such as pesticides, herbicides, polyaromatic hydrocarbons, polychlorinated biphenyls and other environmentally related compounds.

Enviro-Clean® offers a selection of high quality solid phase extraction columns geared to support the environmental chemist with a very broad range of analytical applications. The most important function of the solid phase extraction column for the environmental chemist is the clean separation of an analyte from a variety of compounds. An important function of the extraction column is that it will concentrate a low level of analyte from large samples for accurate analysis. When evaluating analyte extraction or separation, Enviro-Clean® offers nonpolar, polar, ion-exchange and copolymeric phases for application in the environmental laboratory.

### Hydrophobic Extraction Columns

Non-polar phases are often referred to as hydrophobic and function by the interactions of the carbon-hydrogen bond of the analyte and the sorbent. C18 is the most widely used of these phases. EPA approved methods for analyzing organics in drinking water specify the C18 hydrophobic phase. This method requires that large sample volumes (liters) be analyzed which utilizes the compound concentration function of the hydrophobic sorbent.

Weight /Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50 mg/1 mL	C2, Ethyl	EEC021L1	EUC021L1	C4, n-Butyl	EECN41L1	EUCN41L1	C6, Hexyl	EEC061L1	EUC061L1	100
100 mg/1 mL		EEC02111	EUC02111		EECN4111	EUCN4111		EEC06111	EUC06111	100
200 mg/3 mL		EEC02123	EUC02123		EECN4123	EUCN4123		EEC06123	EUC06123	50
500 mg/3 mL		EEC02153	EUC02153		EECN4153	EUCN4153		EEC06153	EUC06153	50
500 mg/6 mL		EEC02156	EUC02156		EECN4156	EUCN4156		EEC06156	EUC06156	50
1000 mg/6 mL		EEC021M6	EUC021M6		EECN41M6	EUCN41M6		EEC061M6	EUC061M6	30
100 mg/10 mL		EEC0211Z	EUC0211Z		EECN411Z	EUCN411Z		EEC0611Z	EUC0611Z	50
200 mg/10 mL		EEC0212Z	EUC0212Z		EECN412Z	EUCN412Z		EEC0612Z	EUC0612Z	50
500 mg/10 mL		EEC0215Z	EUC0215Z		EECN415Z	EUCN415Z		EEC0615Z	EUC0615Z	50
2000 mg/15 mL		EEC0212M15	EUC0212M15		EECN412M15	EUCN412M15		EEC0612M15	EUC0612M15	20
50 mg/1 mL	C3, Propyl	EECN31L1	EUCN31L1	C5, Pentyl	EEC051L1	EUC051L1	C7, Heptyl	EEC071L1	EUC071L1	100
100 mg/1 mL		EECN3111	EUCN3111		EEC0511	EUC05111		EEC07111	EUC07111	100
200 mg/3 mL		EECN3123	EUCN3123		EEC05123	EUC05123		EEC07123	EUC07123	50
500 mg/3 mL		EECN3153	EUCN3153		EEC05153	EUC05153		EEC07153	EUC07153	50
500 mg/6 mL		EECN3156	EUCN3156		EEC05156	EUC05156		EEC07156	EUC07156	50
1000 mg/6 mL		EECN31M6	EUCN31M6		EEC051M6	EUC051M6		EEC071M6	EUC071M6	30
100 mg/10 mL		EECN311Z	EUCN311Z		EEC0511Z	EUC0511Z		EEC0711Z	EUC0711Z	50
200 mg/10 mL		EECN312Z	EUCN312Z		EEC0512Z	EUC0512Z		EEC0712Z	EUC0712Z	50
500 mg/10 mL		EECN315Z	EUCN315Z		EEC0515Z	EUC0515Z		EEC0715Z	EUC0715Z	50
2000 mg/15 mL		EECN312M15	EUCN312M15		EEC0512M15	EUC0512M15		EEC0712M15	EUC0712M15	20



# Sample Prep - Solid Phase Extraction

UCT Enviro-Clean® columns



## Hydrophobic Extraction Columns

Weight /Vol.	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Sorbent	Endcapped	Unendcapped	Qty
50 mg/1 mL	C8, Octyl	EEC081L1	EUC081L1	C12, nDodecyl	EEC121L1	EUC121L1	Cyclohexyl	EECYH1L1	EUCYH1L1	100
100 mg/1 mL		EEC08111	EUC08111		EEC12111	EUC12111		EECYH111	EUCYH111	100
200 mg/3 mL		EEC08123	EUC08123		EEC12123	EUC12123		EECYH123	EUCYH123	50
500 mg/3 mL		EEC08153	EUC08153		EEC12153	EUC12153		EECYH153	EUCYH153	50
500 mg/6 mL		EEC08156	EUC08156		EEC12156	EUC12156		EECYH156	EUCYH156	50
1000 mg/6 mL		EEC081M6	EUC081M6		EEC121M6	EUC121M6		EECYH1M6	EUCYH1M6	30
100 mg/10 mL		EEC0811Z	EUC0811Z		EEC1211Z	EUC1211Z		EECYH11Z	EUCYH11Z	50
200 mg/10 mL		EEC0812Z	EUC0812Z		EEC1212Z	EUC1212Z		EECYH12Z	EUCYH12Z	50
500 mg/10 mL	C10, nDecyl	EEC0815Z	EUC0815Z	C18, Octadecyl	EEC1215Z	EUC1215Z	Phenyl	EECYH15Z	EUCYH15Z	50
2000 mg/15 mL		EEC0812M15	EUC0812M15		EEC1212M15	EUC1212M15		EECYH12M15	EUCYH12M15	20
50 mg/1 mL	C10, nDecyl	EEC101L1	EUC101L1	C18, Octadecyl	EEC18111	EUC18111	Phenyl	EEPHY1L1	EUPHY1L1	100
100 mg/1 mL		EEC10111	EUC10111		EEC18123	EUC18123		EEPHY111	EUPHY111	100
200 mg/3 mL		EEC10123	EUC10123		EEC18153	EUC18153		EEPHY123	EUPHY123	50
500 mg/3 mL		EEC10153	EUC10153		EEC18156	EUC18156		EEPHY153	EUPHY153	50
500 mg/6 mL		EEC10156	EUC10156		EEC181M6	EUC181M6		EEPHY156	EUPHY156	50
1000 mg/6 mL		EEC101M6	EUC101M6		EEC1811Z	EUC1811Z		EEPHY1M6	EUPHY1M6	30
100 mg/10 mL		EEC1011Z	EUC1011Z		EEC1812Z	EUC1812Z		EEPHY11Z	EUPHY11Z	50
200 mg/10 mL		EEC1012Z	EUC1012Z		EEC1815Z	EUC1815Z		EEPHY12Z	EUPHY12Z	50
500 mg/10 mL	C10, nDecyl	EEC1015Z	EUC1015Z	C18, Octadecyl	EEC1812M15	EUC1812M15	Phenyl	EEPHY15Z	EUPHY15Z	50
2000 mg/15 mL		EEC1012M15	EUC1012M15					EEPHY12M15	EUPHY12M15	20

## Hydrophilic Extraction Columns

**Polar or hydrophilic** phases function by hydrogen bonding, pi-pi and dipole-dipole interaction. Ion exchange interactions occur between the sorbent and the analyte of opposite charge.

Enviro-Clean® sorbents are available in both cation or anion exchangers exhibiting both weak and strong characteristics.

**Copolymeric phases** offer a new approach to the environmental analyst by providing very clean extracts and high compound recovery. Dual functionalities, hydrophobic plus ion-exchange or polar allow a higher degree of selectivity than was previously possible. Analytes retained by multiple mechanisms can be washed by disrupting only one mechanism. Careful selection of the solvent strength results in a greater removal of chromatographic contamination.

Weight /Vol.	Qty	P/N
50 mg/1 mL	100	EUCARB1L1
100 mg/1 mL	100	EUCARB111
200 mg/3 mL	50	EUCARB123
500 mg/3 mL	50	EUCARB153
200 mg/6 mL	50	EUCARB126
500 mg/6 mL	50	EUCARB156
1000 mg/6 mL	30	EUCARB1M6
100 mg/10 mL	50	EUCARB11Z
200 mg/10 mL	50	EUCARB12Z
500 mg/10 mL	50	EUCARB15Z
1000 mg/15 mL	20	EUCARB1M15
2000 mg/15 mL	20	EUCARB12M15

## Carbon-Graphitized

**Application :** Carbon supports have been used to isolate extremely polar organic compounds. They work by a hydrophobic mechanism with a high surface area and ion exchange. This interaction can happen in a wide range of polar and non-polar solvents.



# Sample Prep - Solid Phase Extraction

## UCT Enviro-Clean® columns

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Unbonded Silica	EUSIL1L1	Florisil®	EUFLS1L1	Alumina, Acidic	EUALA1L1	Alumina, Basic	EUALB1L1	100
100 mg/1 mL		EUSIL111		EUFLS111		EUALA111		EUALB111	100
200 mg/3 mL		EUSIL123		EUFLS123		EUALA123		EUALB123	50
500 mg/3 mL		EUSIL153		EUFLS153		EUALA153		EUALB153	50
500 mg/6 mL		EUSIL156		EUFLS156		EUALA156		EUALB156	50
1000 mg/6 mL		EUSIL1M6		EUFLS1M6		EUALA1M6		EUALB1M6	30
100 mg/10 mL		EUSIL11Z		EUFLS11Z		EUALA11Z		EUALB11Z	50
200 mg/10 mL		EUSIL12Z		EUFLS12Z		EUALA12Z		EUALB12Z	50
500 mg/10 mL		EUSIL15Z		EUFLS15Z		EUALA15Z		EUALB15Z	50
2000 mg/15 mL		EUSIL12M15		EUFLS12M15		EUALA12M15		EUALB12M15	20
50 mg/1 mL	Alumina, neutral	EUALN1L1	CN, Cyanopropyl Endcapped	EECNP1L1	CN, Cyanopropyl Unendcapped	EUCNP1L1	Diol	EUDOL1L1	100
100 mg/1 mL		EUALN111		EECNP111		EUCNP111		EUDOL111	100
200 mg/3 mL		EUALN123		EECNP123		EUCNP123		EUDOL123	50
500 mg/3 mL		EUALN153		EECNP153		EUCNP153		EUDOL153	50
500 mg/6 mL		EUALN156		EECNP156		EUCNP156		EUDOL156	50
1000 mg/6 mL		EUALN1M6		EECNP1M6		EUCNP1M6		EUDOL1M6	30
100 mg/10 mL		EUALN11Z		EECNP11Z		EUCNP11Z		EUDOL11Z	50
200 mg/10 mL		EUALN12Z		EECNP12Z		EUCNP12Z		EUDOL12Z	50
500 mg/10 mL		EUALN15Z		EECNP15Z		EUCNP15Z		EUDOL15Z	50
2000 mg/15 mL		EUALN12M15		EECNP12M15		EUCNP12M15		EUDOL12M15	20

## Ion exchange and mixed mode Extraction Columns

Weight /Vol.	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Sorbent	P/N	Qty
50 mg/1 mL	Diethylamino	EUDAX1L1	Aminopropyl	EUNAX1L1	Quaternary Amine with Chloride Counter Ion	EUQAX1L1	Carboxylic Acid	EUCCX1L1	100
100 mg/1 mL		EUDAX111		EUNAX111		EUQAX111		EUCCX111	100
200 mg/3 mL		EUDAX123		EUNAX123		EUQAX123		EUCCX123	50
500 mg/3 mL		EUDAX153		EUNAX153		EUQAX153		EUCCX153	50
500 mg/6 mL		EUDAX156		EUNAX156		EUQAX156		EUCCX156	50
1000 mg/6 mL		EUDAX1M6		EUNAX1M6		EUQAX1M6		EUCCX1M6	30
100 mg/10 mL		EUDAX11Z		EUNAX11Z		EUQAX11Z		EUCCX11Z	50
200 mg/10 mL		EUDAX12Z		EUNAX12Z		EUQAX12Z		EUCCX12Z	50
500 mg/10 mL		EUDAX15Z		EUNAX15Z		EUQAX15Z		EUCCX15Z	50
2000 mg/15 mL		EUDAX12M15		EUNAX12M15		EUQAX12M15		EUCCX12M15	20
50 mg/1 mL	Benzene- -sulfonic Acid	EUBCX1L1	Hydrophobic plus Carboxylic Acid	EUCCX2L1	Hydrophobic plus Benzene- -sulfonic Acid	EUBCX2L1	Hydrophobic plus Quaternary Amine	EUQAX2L1	100
100 mg/1 mL		EUBCX111		EUCCX211		EUBCX211		EUQAX211	100
200 mg/3 mL		EUBCX123		EUCCX223		EUBCX223		EUQAX223	50
500 mg/3 mL		EUBCX153		EUCCX253		EUBCX253		EUQAX253	50
500 mg/6 mL		EUBCX156		EUCCX256		EUBCX256		EUQAX256	50
1000 mg/6 mL		EUBCX1M6		EUCCX2M6		EUBCX2M6		EUQAX2M6	30
100 mg/10 mL		EUBCX11Z		EUCCX21Z		EUBCX21Z		EUQAX21Z	50
200 mg/10 mL		EUBCX12Z		EUCCX22Z		EUBCX22Z		EUQAX22Z	50
500 mg/10 mL		EUBCX15Z		EUCCX25Z		EUBCX25Z		EUQAX25Z	50
2000 mg/15 mL		EUBCX12M15		EUCCX22M15		EUBCX22M15		EUQAX22M15	50