



interchim®

(U)HPLC Analysis : C18AQ stationary phases

www.interchim.com

www.blog.interchim.com

www.forum.interchim.com

www.facebook.com/interchim



(U)HPLC, prep-LC & Flash columns

Uptisphere® CS Evolution

Core Shell columns for fast & highly efficient identification & quantification of **small molecules**.

Uptisphere® 120Å

HPLC & prep LC columns for the identification, quantification & purification of **small molecules & pharma compounds**.

Uptisphere® Strategy™

(U)HPLC, Analytical & prep LC columns with **high loadability** for identification, quantification & purification of **small molecules & pharma compounds**.

Uptisphere® X-serie™

HPLC & prep LC columns for the identification, quantification & purification of **small molecules & bio-drugs at high & low pH**.

Uptisphere® 300Å

HPLC & prep LC columns for identification, quantification & purification of **Proteins, Peptides & Polypeptides**.

puriFlash® Prep

prep LC columns for sophisticated purification of **small & bio-molecules & pharma compounds**.

puriFlash®

Flash columns for routine purification of **small & bio-molecules & pharma compounds**.



(U)HPLC Analysis



Small molecules

Oligonucleotides, Peptides
& Polypeptides

70% C18
30% Hilic, C8, C4,
Ph/Hex, PFP, ...

50% C18
35% C8, C4
15% Hilic

Purification

Small molecules


Oligonucleotides,
Peptides & Polypeptides

25% Reverse phase
75% Normal phase

80% Reverse phase
20% Normal phase/Hilic



C18AQ stationary phases Versatility



Using polar or hydrophilic endcapping along with bonding of longer alkyl chains such as C18 is a successful development approach for stationary phases that can retain polar analytes reproducibly under highly aqueous conditions.

For solubility reasons, many polar compounds prefer a highly aqueous mobile phase and can be retained only with a minimal concentration of organic modifier, sometimes less than 5%....

... we demonstrated and discussed the phenomenon of phase collapse and mentioned approaches to solve the problem and successfully separate polar analytes using mobile-phase systems with a very high percentage of water, even as great as 100% water to separate polar compounds in highly aqueous environments.

LCGC Europe dec 2002 ; Columns for Reversed-Phase LC Separations in Highly Aqueous Mobile Phases


Ronald E. Majors, Agilent Technologies, Wilmington, Delaware, USA, and Matthew Przybyciel, ES Industries, West Berlin, New Jersey, USA

Benefits of C18AQ stationary phases


- ✓ Stable under 100% H₂O conditions
- ✓ Retentive and Selective for polar compounds which represent the main molecules to be analyzed and purified today
- ✓ Good peak shape with basic compounds
- ✓ Available from analytical scale to purification, Core-Shell particles up to preparative => Interchim, unique offer on the market!



History & Evolution of Interchim C18 AQ technologies




Interchim C18-AQ : stable under 100% H2O conditions



Interchim C18-AQ Phases

=> Perfectly usable with 100% aqueous mobile phases
Repeatability of analysis times under 100% H₂O conditions



C18 Silica



=> The C18 bonding condenses towards the silica surface using mobile phases containing > 95% H₂O
=> No or loss of retention & separation



Interchim C18-AQ : retention with polar compounds

Interchim C18-AQ Phases

=> Better retention of polar compounds



Conventional C18 silica



=> Less retention of polar compounds



Interchim C18-AQ : selective for polar compounds

Coelution

Acidic Drugs Separation



W..... S..... C18 5 μ m 250 x 4.6mm

Uptisphere C18-HDO 5 μ m 250 x 4.6mm

1) Barbituric Acid, 2) Atropine, 3) Codeine, 4) Pavaverine, 5) Noscapine
30/70 ACN-buffer pH: 2.3 @ 1 ml/min – 30 $^{\circ}$ T, UV: 220 nm

Interchim C18-AQ : selectivity with basic polar compounds

	α 2,6 diMPyridine/Pyridine	
CS Evolution 2.6 μ m C18-AQ	1,18	😊
CS Evolution 2.6 μ m C18-RP	1,12	😊
AMT Halo 2.7 μ m C18	1,00	😞
Uptisphere 3 μ m C18-HDO	1,08	😊
Uptisphere 5 μ m C18-HDO	1,09	😊
Uptisphere Strategy 3 μ m C18-RP	1,06	😊
Y.. 3 μ m ODS-AQ	1,05	😊
PuriFlash 5 μ m C18AQ 70A	1,43	😊


Test Conditions

Mobile Phase : Acetonitrile/Buffer pH=5.5 - (20/80)

Flow Rate (mL/min) : 1.0




Interchim C18-AQ : good peak shape with basic compounds




Interchim C18-AQ

- ⇒ Low activity of surface silanols
- ⇒ Good peak shape for basic compounds



Core Shell Competitor
2.7µm 50x4.6mm

⇒



Uptisphere CS Evolution C18-AQ
2.6µm 50x4.6mm

1) Pyridine, 2) 2,6 Dimethylpyridine

Mobile phase 20/80 ACN-buffer pH: 5.7 @ 1 ml/min – UV: 254nm

CS Evolution 2.6µm C18-AQ

CS Evolution 2.6µm C18-RP

Uptisphere 2.2µm C18-HDO



UHPLC > 650 bars

CS Evolution 2.6µm C18-AQ

CS Evolution 2.6µm C18-RP

Uptisphere 2.2µm C18-HDO



HPLC > 400 bars

PuriFlash 5µm C18-AQ

Uptisphere 3 & 5µm C18-HDO

Uptisphere Strategy 3 & 5µm C18-RP

Uptisphere X-Series 3 & 5µm C18-AQ



HPLC < 400 bars

PuriFlash 5-10-15µm C18-AQ

Uptisphere Strategy 5-10-15µm C18-RP




Prep LC




Availability of Interchim C18AQ stationary phases

	UHPLC	HPLC	Prep LC	Flash
Uptisphere C18-HDO	☺	☺		
Uptisphere Strategy C18-RP	☺	☺	☺	
Uptisphere X-Series C18-AQ		☺	☺	
puriFlash PF C18-AQ		☺	☺	☺
Uptisphere CS Evolution C18-RP	☺	☺		
Uptisphere CS Evolution C18-AQ	☺	☺		


Uptisphere C18-HDO




CS Evolution C18-AQ




puriFlash PF C18-AQ




X-Series C18-AQ



CS Evolution C18-RP



Strategy C18-RP





C18-RP

Selectivity

Uptisphere® CS Evolution™

Capacity

Productivity



Uptisphere® CS Evolution™ C18-RP

85Å - 130m²/g

2.6 μm

Bonding: C18 monofunctional

%C: 6,0

End-capping: multi step mixte

pH stability: 1.5 to 8.0

Suitable for mid & non polar compounds separation. RP shows excellent mechanical stability that provides long-life & make it an excellent tool for analysis under acidic or basic conditions.

USP code: L1

Application:

mid-polar organic compounds



C18-AQ

Selectivity

Uptisphere® CS Evolution™

Capacity

Productivity



Uptisphere® CS Evolution™ C18-AQ

85Å - 130m²/g

2.6 μm

Bonding: C18 monofunctional

%C: 6.5%

End-capping: Mixte

pH stability: 1.5 to 7.0

*Suitable for mid & non polar compounds separation.
RP shows excellent mechanical stability under 100%
aqueous mobile phase conditions.*

USP code: L1

Application:

mid-polar organic compounds

100% water compatible


C18-RP

Selectivity

Uptisphere® Strategy™

Capacity

Productivity



Strategy™ C18-RP

100Å - 425m²/g

2.2, 3, 5, 10 & 15 μm

Bonding: C18 monofunctional

%C: 16

End-capping: multi step mixte

pH stability: 1.5 to 8.0

Suitable for mid & non polar compounds separation. RP shows excellent mechanical stability that make it an excellent tool for purification under acidic or basic conditions.

USP code: L1

Application:

mid-polar organic compounds


C18-HDO

Selectivity

Uptisphere® 120Å

Capacity

Productivity



Uptisphere® C18-HDO

120Å - 320m²/g

2.2, 3, 5µm

Bonding: C18 monofunctional

%C: 17

End-capping: Mixte

pH stability: 1.5 to 7.0

Suitable for mid & non polar compound separation. Shows excellent stability under 100% aqueous mobile phase condition.

USP code: L1

Application:

mid-polar organic compounds

100% water compatible


C18-AQ

Selectivity

Uptisphere® X-serie™

Capacity

Productivity



Uptisphere® X-serie™ C18-AQ

220Å - 200m²/g

3, 5 & 15 μm

Bonding: C18

%C: 14

End-capping: mixte

pH stability: 1.0 to 10.0

*mid-polar BioDrugs & Peptides with medium molecular weight.
100% water compatible*

Reusable columns

Application:

Bio-molecules, Bio-Drugs

C18-AQ

Selectivity

puriFlash® Prep

Capacity

Productivity



puriFlash® Prep C18-AQ

Spherical silica

60Å - 500m²/g

5, 10, 15 & 30 µm

Bonding: C18

%C: 14

End-capping: one-step hydrophilic

pH stability: 2.0 to 7.5


*The bonding chemistry allow to start gradient with 100% of water.
Suitable for the purification of mid and non polar compounds.*

Reusable columns

Application:

mid-polar organic compounds

100% water compatible



interchim

211 bis avenue JF Kennedy – BP1140 – 03100 Montluçon France

Tel. +33 4 70 03 88 55 – Fax +33 4 70 03 82 60

interchim@interchim.com – www.interchim.com