



TOSOH



HYDROPHOBIC INTERACTION CHROMATOGRAPHY

HIC
HYDROPHOBIC
INTERACTION
CHROMATO
GRAPHY

TOSOH BIOSCIENCE

<p>1 TOSOH BIOSCIENCE GMBH ZETTACHRING 6 70567 STUTTGART GERMANY</p> <p>T · + 49 (0) 711 · 13257-0 F · + 49 (0) 711 · 13257-89 INFO.SEP.EU@TOSOH.COM WWW.TOSOHBIOSCIENCE.COM</p>	<p>2 TOSOH BIOSCIENCE LLC 156 KEYSTONE DRIVE MONTGOMERYVILLE, PA 18936 - 9637, USA</p> <p>T · +1 215-283-5000 F · +1 215-283-5035 INFO.TBL@TOSOH.COM WWW.TOSOHBIOSCIENCE.COM</p>	<p>3 TOSOH CORPORATION 3-8-2 SHIBA, MINATO-KU TOKYO 105-8623 JAPAN</p> <p>T · +81 3-5427-5180 F · +81 3-5427-5220 INFO@TOSOH.CO.JP WWW.TOSOHBIOSCIENCE.COM</p>	<p>4 TOSOH SHANGHAI CO. LTD ROOM 2618A, INTERN. TRADE CENTER, NO. 2201 YAN-AN WEST ROAD SHANGHAI, 200336 CHINA</p> <p>T · +86 21-6270-2810 F · +86 21-6270-2820 PAN@TOSOH.COM.CN WWW.TOSOHBIOSCIENCE.COM</p>
--	--	--	--



➤ TOSOH HISTORY

- 1935 FOUNDING OF TOYO SODA MANUFACTURING CO., LTD.
- 1936 OPERATION OF NANYO MANUFACTURING COMPLEX BEGINS
- 1971 SCIENTIFIC INSTRUMENTS DIVISION FORMED, FIRST GPC COLUMN USING TSK-GEL DEVELOPED BY TOSOH
- 1974 HIGH PERFORMANCE LIQUID CHROMATOGRAPHY COLUMN PLANT IS COMPLETED
- 1979 TOSOH DEVELOPS TOYOPEARL MEDIA
- 1983 TOSOH DEVELOPS HYDROPHOBIC INTERACTION MEDIA
- 1987 TOSOHAAS US OPERATIONS FORMED IN MONTGOMERYVILLE
- 1989 TOSOHAAS GMBH OPERATIONS FORMED IN STUTTGART
- 1995 TOSOH NANYO GEL FACILITY RECEIVES ISO 9001
- 2000 IN NOVEMBER FORMER TOSOHAAS US OPERATIONS BECOMES TOSOH BIOSEP LLC, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2001 IN JANUARY FORMER TOSOHAAS GMBH EUROPEAN OPERATIONS BECOMES TOSOH BIOSEP GMBH, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2002/2003 TOSOH CORPORATION ANNOUNCES THAT ALL TOSOH AFFILIATED SCIENTIFIC AND DIAGNOSTIC SYSTEM RELATED COMPANIES IN EUROPE, WILL BE UNIFIED UNDER THE NEW NAME TOSOH BIOSCIENCE.
- 2008 ECOSEC , THE 7TH GENERATION GPC SYSTEM IS INTRODUCED GLOBALLY
- 2009 20TH ANNIVERSARY OF TOSOH BIOSCIENCE GMBH, STUTTGART

HYDROPHOBIC INTERACTION CHROMATOGRAPHY




Hydrophobic Interaction Chromatography (HIC) is a widely-used technique for separation and purification of proteins and peptides. HIC sorts biomolecules by degree of their surface hydrophobicity. Samples are adsorbed to the resin at relatively high salt concentrations and eluted by applying a decreasing salt gradient. The mild conditions used in HIC separation of peptides and proteins typically maintain protein structure and biologic activity. This makes HIC a powerful tool for the process purification of biomolecules.

An optimum HIC process step will balance high dynamic binding capacity (DBC), adequate selectivity, good mass recovery and retention of biological activity. The key parameter is selecting the best resin for the given separation problem. Proteins show varying degrees of hydrophobicity depending on their amino acid composition, structure and size. Separation can therefore be optimized either by varying the mobile phase or by using different HIC packings. Matching the hydrophobicity of the target compound to the resin hydrophobicity is critical for the best overall purification performance. This is the reason why Tosoh Bioscience offers seven product lines of ToyoPearl HIC resins using five different ligands. The different degrees of hydrophobicity and selectivity support the user in selecting the best solution for a given target.

The hydrophobicity increases through the ligand series: Ether, Polypropylenglycol (PPG), Phenyl, Butyl, Hexyl. ToyoPearl HIC resins are available in three different average particle sizes (35 µm (S), 65 µm (M) & 100 µm (C)) for intermediate purification or capture chromatography. For high resolution HIC Tosoh Bioscience offers TSK-GEL resins with particle sizes of 20 and 30 µm.

FIGURE 1
HIC Ligand Candidates



HIC

HOW IT WORKS

Many theories and models have been proposed to describe the HIC retention mechanism but none of them has gained universal acceptance. HIC is based upon interactions between hydrophobic patches on the surface of biomolecules and the hydrophobic ligands of the stationary phase. It is commonly believed that the driving force of interaction is the entropy gain arising from changes in the order of the water molecules surrounding the interacting hydrophobic groups. Protein binding to HIC adsorbents is promoted by moderately high concentrations of anti-chaotropic salts. Elution is achieved by a linear or stepwise decrease of salt in the mobile phase.

Selectivity

The hydrophobicity of a target with known structure can be roughly estimated as it often increases with the size of the protein surface. Nevertheless, practical screening experiments under standard buffer conditions are essential to select the optimum resin. The hydrophobicity of the resin determines the salt concentration necessary to adsorb the target. With low-hydrophobic ligands the difference between adsorption and precipitation might be so small that certain proteins may partially precipitate under binding conditions. On the other hand a high-hydrophobic stationary phase might cause irreversible binding of hydrophobic proteins.

HIC Method Development

The goal in purification method development is optimizing conditions for maximum capacity and recovery of the target molecules. There are several parameters which affect HIC separations in addition to the hydrophobicity of the ligand:

- Salt type
- pH
- Buffer concentration
- Temperature
- Gradient type, slope
- Particle and pore size
- Column dimensions





FIGURE 2
Binding Capacities of TSKgel Phenyl-5PW



Optimizing Salt Type and Concentration

Besides the hydrophobicity of the resin, the eluent salts make a major impact on a HIC separation. Ammonium sulfate and sodium chloride are most commonly used for HIC applications. Sometimes citrate-buffers or dual salt systems are used to improve resolution. While the type of salt affects retention and selectivity the initial salt concentration is the key to maximize binding capacity for the target. The salt concentration required for binding is related to the size of the surface area of the protein. Small, hydrophilic proteins will need high salt, e.g. up to 3 M ammonium sulfate, for efficient binding but it can decrease below 1 M for very large proteins. Figure 2 shows the influence of salt concentration on binding capacity of TSKgel Phenyl-5PW for various proteins.

Other Parameters

pH can be used for fine tuning. A good starting pH is 7.0, irrespective of the component's isoelectric point. The pH can influence not only retention but also DBC (Figure 3).


Most HIC applications are performed at room temperature or at 4°C. A higher temperature might be used to influence binding strength and selectivity.

Elution is typically performed by gradient elution. The sample is applied at a salt concentration high enough for adsorption of the targets. As the salt concentration is lowered, proteins become increasingly desorbed and move down the column. Resolution can be increased by decreasing gradient slope. In manufacturing scale processes step gradients are more common than linear gradients.

Resolution in HIC can be improved by increasing the column lengths, since the full length of the column bed interacts continuously with sample components.

Organic modifiers can speed up a HIC separation or alter the selectivity. For purification of small molecules up to 20% ethanol might be used.

FIGURE 3
Influence of pH





HIC TSK-GEL AND TOYOPEARL HIC RESINS

The **particle size** depends on the sample and the required resolution. Capturing steps from a crude feedstock are usually performed with coarse particles (Toyopearl C). In intermediate purification steps medium size particles (Toyopearl S or M) are used, whereas for polishing the even smaller TSK-GEL materials with 20 µm or 30 µm particles are ideal. TSK-GEL columns with 10 µm beads are best suited for analytical purposes or for small scale purifications (Figure 4).

Toyopearl HIC Material

Toyopearl and TSK-GEL HIC resins are specifically designed for use in biopharmaceutical production. Their rigid methacrylic polymer structure shows excellent pressure/flow properties enabling high process throughput. Large pore diameters and narrow particle size distribution allow rapid adsorption kinetics and exceptional resolution. For seamless scale-up Tosoh Bioscience offers a complete HIC toolbox, ranging from analytical TSK-GEL HPLC columns up to bulk media used for pilot and production scale.

HIC Ligands

The wide range of Toyopearl and TSK-GEL HIC selectivities enables a developer to optimize protein separations at the extremes of the hydrophobic spectrum. The hydrophobicity of the resins increases through the series:

Ether < PPG < Phenyl < Butyl < Hexyl


Highly retentive Hexyl and Butyl resins are used to separate hydrophilic proteins and should be considered for separations requiring a low ionic strength. Toyopearl Ether resin is used for the purification of very hydrophobic targets such as certain monoclonal antibodies and membrane proteins. PPG and Phenyl complement the other HIC ligands and offer alternatives for mid-range hydrophobic proteins.

Today high-resolution HIC applications gain more and more interest. TSK-GEL 5PW media with small particle sizes are ideally suited if high resolution is an issue. TSK-GEL 5PW bulk material is available with the ligands Ether and Phenyl. TSK-GEL columns in various dimensions are available with Ether, Phenyl and Butyl chemistry.

Regulatory Support

Pharmaceutical industry all over the world successfully uses Toyopearl HIC resins in the downstream processing of a variety of biologically active proteins, including several FDA-approved therapeutic drugs. For Toyopearl HIC resins 'Regulatory Support Files', describing the specifications, the manufacturing and the QA/QC of the product are registered at the FDA. In addition, Tosoh Bioscience's application specialists are available for discussion of your specific separation challenge or process validation issues.

➤ FIGURE 4
HIC Resins



HIC

TSK-GEL AND TOYOPEARL HIC RESINS

Dynamic Binding Capacity

In downstream processing steps, the dynamic binding capacity (DBC) of the resin for the target is even more important than selectivity. Selecting media with a different pore size is an option, if DBCs are not satisfying. Tosoh Bioscience provides resins designed for maximum dynamic binding capacity for dedicated proteins. The standard Toyopearl resins have an average pore size of 1000 Å, suitable for most targets.


Smaller Pore Size HIC Resins

The accessible surface area of a porous bead increases by decreasing the mean pore diameter and so does the dynamic binding capacity. This led to the development of two specialty lines of HIC materials with smaller pores. For monoclonal antibodies a pore size of 750 Å is sometimes favorable. Toyopearl resins exhibiting this pore size are available with three ligands: PPG-600, Phenyl-600 and Butyl-600. For smaller molecules such as peptides Toyopearl resins with even narrower pore diameter (500 Å) are used to create the SuperButyl-550C resin.

The variety of HIC phases increases the probability of matching a resin best to the given target, at the same time making the screening procedure more complex. Figure 6 shows all available Toyopearl resins sorted according to their pore size and relative hydrophobicities.

FIGURE 5

Typical Dynamic Binding Capacities of Lysozyme



Binding capacity (mg/mL - gel)
(10% Breakthrough)

Ether-650M	12.5
Phenyl-650M	27.5
Butyl-650M	32.2
Hexyl-650C	33.2
PPG-600M	44.2
SuperButyl-550C	54.3

Conditions

Column size: 7.8 mm ID x 20 cm L
 Samples: 1 mg lysozyme in 0.1 M phosphate buffer + 1.8 M sodium sulfate (pH 7.0)
 Linear velocity: 100 cm/h
 Detection: UV @280 nm

The Toyopearl Phenyl-600M resin was designed as a high-sub type. The higher ligand density results in a higher hydrophobicity than Toyopearl Phenyl-650 resins.

FIGURE 6

Hydrophobicity and Average Pore Size of Toyopearl HIC Resins





Figure 5 and 7 show the dynamic binding capacities of Toyopearl resins for Lysozyme and a monoclonal antibody. For a small protein such as Lysozyme the SuperButyl-550C is the best choice (Figure 5). Figure 7 demonstrates the superior DBC of the Butyl-600M and Phenyl-600M resins for large proteins.

FIGURE 7

Breakthrough Curves of a Polyclonal IgG on Various HIC Resins



Binding capacity (mg/mL - gel)
(10% breakthrough)

Toyopearl Phenyl-600M	40.0
Agarose Type Phenyl	32.0
Toyopearl Butyl-600M	29.0
Toyopearl Phenyl-650M	16.0
Toyopearl PPG-600M	3.0

Column size:

7.8 mm ID x 20 cm

Sample:

polyclonal human IgG

Binding buffer:

1 g/l IgG in 0.8 mol/l (NH₄)₂SO₄ + 0.1 mol/l sodium phosphate (pH 7.0)

Linear velocity:

300 cm/h

Temperature:

25°C



HIC SCREENING

ToyoScreen® for Easy Resin Scouting


In order to simplify the screening process, Tosoh Bioscience offers sets of prepacked columns with different resins. They provide a convenient way to screen different resins effectively for both, target retention and recovery. ToyoScreen is available with 1 and 5 ml bed volumes for most Toyopearl resins and can be connected to common laboratory liquid chromatography instruments. If the LC system is equipped with automated solvent and column switching valves, screening of resins at various buffer conditions can be easily performed in overnight runs.

The effect of the different hydrophobicities of Toyopearl resins on retention and resolution of standard proteins are illustrated in Figure 8. A standard mixture of proteins was separated using ToyoScreen columns. Fast screening of a larger number of resins under various conditions can be realized by applying robotic fluid handling systems and high throughput screening tools in 96 well plate formats.

Comparison of HIC Resins

Non-specific binding effects from the base material of the resin can alter resolution and selectivity. The matrix of Toyopearl and TSK-GEL HIC resins is a uniform, hydrophilic polymer. HIC resins from other manufacturers, based on different base resins, might exhibit different properties regarding hydrophobicity, selectivity and resolution even if they are functionalized with the same ligand. This is important to consider when screening resins of various manufacturers.

FIGURE 8 Screening of Toyopearl HIC Resins - Standard Proteins



Column: ToyoScreen (1 ml)
 Eluent A: 0,1 M Phosphate Buffer + 1.8 M Sodium Sulfate (pH 7.0)
 Eluent B: 0.1 M Phosphate Buffer (pH 7.0)
 Flow Rate: 1 ml/min
 Gradient: 30 min linear
 Inj.Vol.: 50 l
 Samples: Ribonuclease A, Lysozyme, α -Chymotrypsinogen, 1 mg/ml



HIC SCALE UP

Seamless scale up


In terms of cost efficiency a production step should deliver maximum yield of the active product in short time. It will always be a compromise between throughput, resolution and recovery. The capacity of the column must fit to the yield of the upstream process or of the previous purification steps respectively. The target capacity determines the column dimensions, while the nature of the sample and the approached resolution determine the particle size.

The chemistry of the resins is very similar from the pre-packed TSK-GEL PW HPLC columns to the TSK-GEL-5PW and Toyopearl bulk resins. This offers the opportunity to find the ideal particle size for the intended use regardless of whether it is laboratory scale purification, a process polishing, intermediate or capture step. Figure 10 shows the separation of four standard proteins on various Phenyl media. Increasing the bead size from 10 µm (TSKgel Phenyl-5PW) over 35 µm and 65 µm up to 100 µm only reduces resolution but does not impair selectivity.

Superior pressure/flow characteristics


High flow rates reduce process cycle time and increase productivity. The rigid polymeric backbone of Toyopearl and TSK-GEL HIC resins assures superior pressure/flow characteristics over a wide range of flow rates. Figure 9 shows the excellent pressure flow/curves for all grades of Toyopearl Butyl-650, determined on a production size column with 40 cm ID and 20 cm length.

FIGURE 9
Pressure Flow Curve



Resin: Toyopearl Butyl-650
 Column Size: 40 cm ID x 20 cm L
 Eluent: Water
 Temp.: Room Temperature

FIGURE 10
Improvement of Performance by Reducing Particle Size



Column: 7.5 mm ID x 7.5 cm L
 Sample: 1. Myoglobin, 2. Ribonuclease A, 3. Lysozyme, 4.α-Chymotrypsinogen
 Injection: 200 µl
 Elution: 60 min. linear gradient from 2.0 M to 0 M of (NH₄)₂SO₄ in 0.1 M phosphate buffer, pH 7.0
 Flow rate: 1.0 ml/min.
 Detection: UV @ 280 nm



HIC APPLICATIONS

Applications

Toyopearl and TSK-GEL HIC resins are used in downstream purification of a variety of biopharmaceuticals. HIC is often used in capture steps following an ammonium sulphate precipitation. It is decreasing the salt concentration at the same time as conducting a purification step. HIC is a common intermediate process step for the purification of monoclonal antibodies. It is typically used to remove leached Protein A and aggregates subsequent to an affinity step. A typical industrial purification scheme for the isolation of mAbs from a cell culture supernatant is shown in Figure 12.

Monoclonal Antibodies

The diverse hydrophobic nature of mAbs is shown in Figure 11. The retention time as an indicator of hydrophobicity was measured for 51 different mouse IgGs on a TSKgel Phenyl-5PW analytical column. The elution time differs by a factor of 2-3 indicating very different hydrophobicities. The Toyopearl series of HIC ligands with different hydrophobicities offers a range of options for finding the right resin for the target molecule. For the highly hydrophobic mouse anti-chicken 14 kDa lectin the hydrophilic Ether ligand works well. Figure 13 shows the purification of this antibody from ascites fluid with Toyopearl Ether-650M material.

Aggregate Removal


HIC in flow through mode is often used to remove aggregates generated in Protein A purification steps for mAbs. These impurities have chemical properties very similar to the target but they will generally be more hydrophobic than the native protein. Therefore they bind at relatively low salt concentrations to Butyl or Phenyl resins allowing the target to flow through the column.

In addition to the mentioned examples HIC is used successfully for a variety of other applications such as plasmid purification and endotoxin removal.

FIGURE 11

Hydrophobic Diversity of Mouse mAbs


Plot of chromatographic elution times for 51 different mouse mAbs



Column: TSKgel Phenyl-5PW
 Eluent: (A) 0.1 M phosphate buffer containing 1.8 M ammonium sulfate (pH 7.0)
 (B) 0.1 M phosphate buffer (pH 7.0)
 Flow rate: 1 mL/min
 Gradient: (B) 0% (0 min)--0% (5 min)--100% (65 min) linear
 Samples: 51 kinds of mouse monoclonal antibodies

FIGURE 12

Example of Industrial mAb Purification



Even glycoproteins, which often bind irreversibly to saccharide based media, can be purified by HIC on polymer based resins.


Regeneration of the Column

The type and frequency of regeneration of a column naturally depends on the samples applied. Standard cleaning procedures involve washing with high pH (e.g. 0.5 N NaOH). Toyopearl and TSK-GEL HIC resins are recommended for use from pH 2.0 to 12.0, although short exposures to higher pH for cleaning in place are possible.

FIGURE 13

Purification of mAbs from Ascites Fluid

65 µm Toyopearl Ether-650M



Column: Toyopearl Ether-650M, 7.5 mm ID x 7.5 cm L
 Sample: anti-chicken 14 kDa lectin, diluted ascites fluid, 0.76 mg in 50 µl
 Elution: 60 min. linear gradient from 1.5 M to 0 M (NH₄)₂SO₄ in 0.1 M phosphate buffer (pH 7.0)
 Flow rate: 136 cm/h
 Detection: UV @ 280 nm

ORDERING INFORMATION

► ORDERING INFORMATION

Toyopearl HIC resins:

Hydrophobicity	Chemical Structure	Product description	Container size (mL)	Part #	Particle size (µm)	Pore Size (Å)				
weak	HW65-(OCH ₂ CH ₂) _n -OH	Ether-650S	25	43151	20-50	1000				
			100	16172						
			1,000	16174						
			5,000	16176						
			25	19805						
	Ether-650M	100	16173	40-90	1000					
		1,000	16175							
		5,000	16177							
		medium	HW60-(OCH(CH ₃)-CH ₂) _n -OH			PPG-600M	25	21301	40-90	750
							100	21302		
1,000	21303									
5,000	21304									
25	43152			20-50	1000					
HW65-OC ₆ H ₅	Phenyl-650S		100			14477				
			1,000			14784				
			5,000			14935				
			25			19818	40-90	1000		
			Phenyl-650M	100	14478					
1,000	14783									
5,000	14943									
25	43126	50-150		1000						
Phenyl-650C	100				14479					
	1,000		14785							
	5,000		21887		40-90	750				
	HW60-OC ₆ H ₅		Phenyl-600M				100	21888		
		1,000		21889						
5,000		21890								
strong		HW65-O-(CH ₂) ₃ -CH ₃		Butyl-650S			25	43153	20-50	1000
					100	07476				
	1,000		14701							
	Butyl-650M		5,000	07975	40-90	1000				
			25	19802						
			100	07477						
	Butyl-650C		1,000	14702	50-150	1000				
			5,000	07976						
			25	43127						
	HW60-O-(CH ₂) ₃ -CH ₃	Butyl-600M	100	07478	40-90	750				
			1,000	14703						
			5,000	07977						
		SuperButyl-550C	25	21448			50-150	500		
			100	19955						
			1,000	19956						
	HW65-O-(CH ₂) ₅ -CH ₃	Hexyl-650C	5,000	19957	50-150	1000				
			25	44465						
			100	19026						
Hexyl-650C		1,000	19027	50-150			1000			
		5,000	19028							
		25	44465							

Toyopearl LABPAK:

Part #	Product description	Container size (mL)	Particle size (µm)
43150	HICPAK HP (Ether, Phenyl, Butyl-650S)	3 x 25 mL	35
19806	HICPAK (Ether, Phenyl, Butyl-650M)	3 x 25 mL	65
43125	HICPAK-C (Phenyl, Butyl, Hexyl-650C)	3 x 25 mL	100

ORDERING INFORMATION



ORDERING INFORMATION

TSK-GEL 5PW HIC resins for high resolution:

Hydrophobicity	Chemical Structure	Product description	Container size (mL)	Part #	Particle size (µm)	Pore Size (Å)
weak	5PW-(OCH ₂ CH ₂) _n -OH	Ether-5PW (20)	25	43276	10-30	1000
			250	16052		
			1,000	16053		
		Ether-5PW (30)	5,000	18437	20-40	1000
			25	43176		
			250	16050		
		Phenyl-5PW (20)	1,000	16051	10-30	1000
			5,000	18439		
			25	43277		
medium	5PW-OC ₆ H ₅	Phenyl-5PW (20)	25	43277	10-30	1000
			250	14718		
			1,000	14719		
		Phenyl-5PW (30)	5,000	18438	20-40	1000
			25	43177		
			250	14720		
		Phenyl-5PW (30)	1,000	14721	20-40	1000
			5,000	17210		

ToyoScreen process development columns for HIC:

Part #	Product description	Package
21372	ToyoScreen Ether-650M, 1 mL	1 mL x 6 each
21373	ToyoScreen Ether-650M, 5 mL	5 mL x 6 each
21374	ToyoScreen Phenyl-650M, 1 mL	1 mL x 6 each
21375	ToyoScreen Phenyl-650M, 5 mL	5 mL x 6 each
21376	ToyoScreen Butyl-650M, 1 mL	1 mL x 6 each
21377	ToyoScreen Butyl-650M, 5 mL	5 mL x 6 each
21378	ToyoScreen Hexyl-650C, 1 mL	1 mL x 6 each
21379	ToyoScreen Hexyl-650C, 5 mL	5 mL x 6 each
21380	ToyoScreen PPG-600M, 1 mL	1 mL x 6 each
21381	ToyoScreen PPG-600M, 5 mL	5 mL x 6 each
21495	ToyoScreen Butyl-600M, 1 mL	1 mL x 6 each
21494	ToyoScreen Butyl-600M, 5 mL	5 mL x 6 each
21892	ToyoScreen Phenyl-600M, 1 mL	1 mL x 6 each
21893	ToyoScreen Phenyl-600M, 5 mL	5 mL x 6 each
21382	ToyoScreen SuperButyl-550C, 1 mL	1 mL x 6 each
21383	ToyoScreen SuperButyl-550C, 5 mL	5 mL x 6 each
21398	ToyoScreen HIC Mix Pack, 1 mL	1 mL x 6 Grades x 1 each
21399	ToyoScreen HIC Mix Pack, 5 mL	5 mL x 6 Grades x 1 each

ToyoScreen column accessories

21400	ToyoScreen Column Holder
-------	--------------------------

TSK-GEL LABPAK:

Part #	Product description	Container size (mL)	Particle size (µm)
43278	HICPAK PW (20) (Ether-5PW, Phenyl-5PW)	2 x 25 mL	10-30
43175	HICPAK PW (30) (Ether-5PW, Phenyl-5PW)	2 x 25 mL	20-40