



CY_{anine} Alkyne

CY_{anine} fluorophores functionalized for labeling biomolecules by click chemistry

Products Description

The table below gives main physical and fluorescence characteristics of the activated dyes.

| Product name cat.number/gty* | MW (g·mol ⁻¹) | λ abs./em. | mol. abs. | Comment, structure |
|---|------------------------------|---------------|---------------------|--|
| 1.5 | +added MW | (nm) | $(M^{-1}cm^{-1})$ | |
| CY _{anine} 3 – Alkyne | 530.14 | 555 / 570 | 150 000 OY: 0.31 | |
| (M) $FP-1A6320, Img$ | | | | a |
| Red powder Soluble in DMF, DMSO, dichloromethane), insoluble in water | | | | HN |
| Mono-Sulfo-CY _{anine} 3 – | 573.75 | 550 / 567 | 96 800 | 0 |
| Alkyne | | | QY: 0.15 | |
| ^(M) FP-1C4620, 1mg | | | | - The second sec |
| Soluble in DMF, DMSO, MetOH, | | | | 2 |
| EtOH and acetonitrile. Poor in water. | | | | |
| | | | | \mathcal{F}^{N_3} (bears 1 sulfo, 1 methyl) |
| DiSulfo-CY _{anine} 3 – Alkyne, | 675.79 | 548 / 567 | 162 000 | [©] o₃s so ⁰ ₃ |
| Na ⁺ salt | | | QY: 0.15 | |
| ^(M) FP-1C4630, 1mg | | | | N N N N N N N N N N N N N N N N N N N |
| Soluble in water and polar organic | | | | (B) Na |
| solvent (DMF, DMSO), insoluble in organic solvents (DCM, chloroform) | | | | |
| DiSulfo-CY _{anine} 3 – Alkyne, | 691.9 | 548/563 | | EC:162000; QY:0.1; |
| K^+ salt | | | | CF260:0.03; CF280:0.06 |
| ^(M) FP-0B8390, 1mg | | | | Solubility is good in water, DMF, DMSO |

interchim

211 bis, av BP 1140 -Fax +33 4

Kennedy

03100 Mont 70 03 82 60



| FT-1A6320 | | | | |
|---|---|-----------------------|---------------------|--|
| Product name cat.number/qty* | MW (g·mol ⁻¹) +added MW | λ abs./em. (nm) | mol. abs. | Comment, structure |
| $\begin{array}{ll} \textbf{DiSulfo-CY}_{anine}\textbf{3}-\textbf{Alkyne},\\ \text{TFA}^{}\text{salt}\\ \text{(M)} & \text{FP-LQV030, 1mg} \end{array}$ | 781.86 | 555 / 565 | | $\begin{array}{c} HO_{3}S \\ H_{3}C \\ CF_{3}CO_{2}^{-} \\ CF_{2}H_{5} \end{array} \xrightarrow{O} \\ N^{H} \\ N^{H$ |
| Tri-Sulfo-CY _{anine} 3-Alkyne ^(M) FP-1C8831, 5mg | 761.93 | 555 / 580 | 151 000 | (3 Sulfo, 1 butylSulfo |
| CY _{anine} 5 – Alkyne ^(M) FP-OO5590, 1mg dark blue powder / solution Solubility is good in dichloromethane, DMF, DMSO, alcohols, low in water | 556.18 (520.3) | 646 / 662 | 250 000 QY: 0.2 | CI- |
| Di-Sulfo-CY _{anine} 5 – Alkyne, Na salt ^(M) FP-1C4640, 1mg Soluble in water and polar organic solvent (DMF, DMSO), insoluble in organic solvents (DCM, chloroform) | 717.94 (701.83) | 646 / 662 | 271 000 QY: 0.28 | eo ₃ s h h h h h h h h h h h h h |
| Di-Sulfo-CY _{anine} 5 – Alkyne, K salt FP-SJI060, 1mg Very high solubility in water | 717.94 | 649/662 | 271 000 QY: 0.2 | (bears 2 Sulfo) QY:0.28; CF260:0.04; CF280:0.04 – Solubility is very good in Water, DMF, DMSO |
| Di-Sulfo-CY _{anine} 5 – Alkyne , TFA salt ^(M) FP-LQV090, 1mg | 807.90 | 649 / 665 | 5 | HO ₃ S H ₃ C CH ₃ H ₃ C CH ₃ H ₅ C CH ₃ SO ₃ H CF ₃ CO ₂ $\overset{N^+}{C_2H_5}$ N (bears 2 Sulfo, 1 Ethyl) |
| Tri-Sulfo-CY anine 5-Alkyne(M)FP-1C8841, 5mgSoluble in DMSO, DMF, MeOH, WaterAbs/Em = 646/663 nm; EC: 250000 | 787.98 | 646 / 663 | 271 000 | (bears 3 Sulfo, 1 butylSulfo) |
| CY _{anine} 5.5 – Alkyne, Cl ⁻ ^(M) FP-SJH910, 1mg | 656.30 | 673 / 707 | 209 000 QY : 0.2 | |



| Product name | MW | λ | mol. abs. | Comment, structure |
|--|-------------------------------------|------------------|---|---|
| cat.number/qty* | (g·mol ⁻¹) +added MW | abs./em. (nm) | (M ⁻¹ cm ⁻¹) | |
| Tri-Sulfo-CY _{anine} 5.5 – | 1257.69 | 678 / 694 | 251 000 | о с сн. и с сн. и он |
| Alkyne, TFA salt | (TEA salt) | | EC: | |
| ^(M) FP-MRV061 , 5mg | 934.12(protonated) | | CF_{260} : 0.05, CF_{280} : 0.18 | |
| ; Soluble in DMSO, DMF, MeOH, Water Abs/Em = 678/694 nm, | | | | |
| | 10(0.14 | (70 / 701 | | (bears 3 Sulfo including 1 butylSulfo) |
| Tetra-Sulfo-CY _{anine} 5.5 – | 1068.14 | 6/8///01 | | |
| Alkyne, TFA salt | | | | · QX XII. |
| ^(M) FP-LQV320, 1mg | | | | |
| | | | | (bears 4 Sulfo 1 Ethyl) |
| $CV \cdot 7 - Alkyne Cl salt$ | 622.28 | 749 / 776 | " | Solubility is good in DMSO, DMF, |
| $^{(M)}$ FP-WZE120, 1mg | [+586.4] | | | alcohols |
| DiSulfo-CY _{anine} 7 – Alkyne. | 818.91 | 749 / 776 | | HO ₃ S H ₃ C CH ₂ H ₂ C CH ₃ SO ₃ H |
| TFA salt | | | | A a a a a a a a a a a a a a a a a a a a |
| ^(M) FP-LQV250, 1mg | | | | CF ₃ CO ₂ · N^{*} $CF_{3}CO_{2}$ · $C_{2}H_{5}$ |
| Solubility in DMSO | | | | NH NH |
| | | | | (bears 2 Sulfo, |
| Tri-Sulfo-CY _{anine} 7 – Alkvne. | 1313.80 | 753/775 | 255 000 | |
| TFA salt | (TEA salt), | | | |
| ^(M) FP-111220. 1mg | d | | | (3Sulfo, 1 Buty) |
| FP-111221, 5mg | | | | 17 |
| Soluble in DMSO, DMF, MeOH, Water Abs/Em = $753/775$ nm, EC: 255000 ; CF ₂₆₀ : 0.05, CF ₂₈₀ ; 0.05 | | | | |

Storage: -20° C for long term, protected from light ^(M)

Cyanine fluorophores

LIFE SCIENCES

A variety of **Cyanine dyes** has been used to label proteins, nucleic acids and other biomolecules for fluorescence techniques (imaging, biochemical analysis). They replace advantageously the conventional fluorochromes such as Fluorescein(FITC) and rhodamines (TRITC, RRX).

CY_{anine}3 can replace orange-fluorescent dyes, like Tetramethylrhodamine (TRITC).

70 03 82 60

 $Cy_{anine}3$ is one of the most broadly used fluorophores which can be detected by various fluorometers, imagers, and microscopes. Due to inherently high extinction coefficient, this dye is also easily detected by naked eye on gels, and in solution. See also alternative superior dye: <u>FluoProbes547H</u>.

CY_{anine}3.5 can replace SulfoRhodamine 101.

See also alternative superior dye: FluoProbes594.

 $CY_{anine}\mathbf{5}$ can replace far red red fluorescent dyes.

During last years, CY_{anine}5 flurophore has become an incredibly popular label in life science research and diagnostics.



FT-1A6320

Fluorophore emission has maximum in red region, where many CCD detectors have maximum sensitivity, and biological objects have low background. Dye color is very intense, therefore quantity as small as 1 nanomol can be detected in gel electrophoresis by naked eye. See also alternative superior dye: <u>FluoProbes647H</u>

 CY_{anine} 5.5 can replace near infrared fluorescent dyes.

See also alternative superior dye: <u>FluoProbes682</u>.

CY_{anine}7 is a near infrared red fluorophores used in in vivo imaging applications.

See also alternative superior dye: FluoProbes752.

CY_{anine}7.5 is a near infrared red fluorophores used for *in vivo* imaging applications.

See also alternative superior dye: FluoProbes800.

 $Sulfo - CY_{anine}$ dyes are water-soluble derivates of the CY_{anine} dyes. DiSulfonated compounds are the most classic, but some tri- and quadri-sulfonated forms are available as well, for even higher hydrosolubility.

Click Chemistry

Alkyne reacts with Azide groups (standard Click reaction using Cu catalyzer). Among the vast variety of organic reactions, Click Chemistry has been selected as a conjugation chemistry reaction because of several advantages:

► It is very selective. Click Chemistry reaction takes place only between azide and alkyne components. It is does not interfere with most any other organic groups present in DNA and proteins being labeled, such as amino and carboxy groups.

► There are no azides and alkynes in native biomolecules. These groups should be specially introduced into DNA and proteins. Alkyne-containing DNA can be prepared with alkyne phosphoramidite during standard oligo synthesis. Proteins labeled with azide and alkyne can be made using azide activated ester and alkyne activated ester.

► Click Chemistry takes place in water. Aqueous DMSO, DMF, acetonitrile, alcohols, or pure water and buffers can be used for the reaction. The reaction is biocompatible and can take place in living cells.

► Reaction is quick and quantitative. Click Chemistry is a tool that allows preparation of nanomoles of conjugates in diluted solutions.

► The reaction is pH-insensitive. Unlike reaction of NHS esters with amines, and some other conjugation chemistries, there is no need to control pH in reaction mixture. There is no need to add any special buffer, acid or base - Click Chemistry works well in pH interval of 4-11.

Click Chemistry thus became a tool for universal modification of DNA, proteins, conjugate preparation, and fluorescent labeling. This is where our reagents come to help: you can perform easy preparation of conjugates in your lab.

• Click Chemistry reaction – CuAAC

Click Chemistry is a reaction between azide and alkyne yielding covalent product - 1,5-disubstituted 1,2,3-triazole. This process is also known as CuAAC - Cu catalyzed alkyne azide cycloaddition.



The catalyst is often introduced as Cu-TBTA (Copper(II)-TBTA complex FY2780).

General protocols for standard click chemistry (Cu-catalyzed)

See the technical sheet FT-FY2780 for more information on Click Chemistry and protocols.



FT-1A6320 References



- Ngo et al. (2012) State-selective metabolic labeling of cellular proteins. ACS Chemical Biology. 7(8):1326.
- Rubino et al. (2012) Chemoselective Modification of Viral Surfaces via Bioorthogonal Click Chemistry. J. Vis. Exp. 66 :4246. Yao et al. (2012) Fluorophore Targeting to Cellular Proteins via Enzyme-Mediated Azide Ligation and Strain-Promoted Cycloaddition. J. Am. Chem. Soc. 134 :4246.
- Best et al. (2009) Click chemistry and bioorthogonal reactions: unprecedented selectivity in the labeling of biological molecules. Biochemistry. 48(28):6571.
- Ning et al. (2008) Visualizing metabolically labeled glycoconjugates of living cells by copper-free and fast Huisgen cycloaddition. Anweg. Chem. Int. Ed. 47(12) :2253.
- Baskin et al. (2007) Copper-free click chemistry for dynamic in vivo imaging. PNAS. 104(43) :16793.
- Dieterich et al. (2007) Labeling, detection and identification of newly synthesized proteomes with bioorthogonal noncanonical amino acid tagging. Nature Protocols. 2(3):532

Related products

* Other CY_{anine} dyes functionalized by NHS (<u>FT-BB7493</u>), Maleimide (<u>FT-JO6660</u>), Azide (<u>FT-HO7250</u>), Alkyne (<u>FT-1A6320</u>), DBCO (<u>FT-1A6320</u>: CycloAlkynes, for strain-promoted Click reactions), Tetrazine (<u>FT-WXS720</u>), Hydrazide (<u>LQV050</u>), Amino group (<u>CY3AM0</u>), Carboxyl group (<u>CY3CA0</u>). 3Dye 2D DIGE (CY2/CY3/CY5) labeling kit (<u>EV0870</u>)

* Related labels • Superior FluoProbes fluorescent dyes,

- activated by –NHS (<u>list</u>), i.e. FP488-NHS #<u>BA68000</u> - activated by –Azide, i.e. FP488-Azide #<u>YE4970</u> • Classic dyes such as **FAM**, **R110**, **JOE TAMRA**, and **ROX**.

Ordering information

Catalog size quantities and prices may be found at <u>www.interchim.com/</u>. Please inquire for bigger quantities and for any information, please ask : FluoProbes[®] / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from FluoProbes[®] are sold for research use only, and are not intended for food, drug, household, or cosmetic use. FluoProbes[®] is not liable for any damage resulting from handling or contact with this product.